

PhilHealth

eClaims

Implementation Guide

APPROVAL SHEET

This document along with its attachments have been approved as the Official Implementation Guide for the Electronic Claims (eClaims) Project.

Following approval of this document, the documents will be distributed to Health Facilities and Health Information System / Electronic Medical Record (HIS/EMR) Service Providers for compliance to the said project

Prepared by:

MANUEL S. MANLANGIT
System Developer

GUILLERMO I. RAZON, JR.
System Developer

Recommended for Release by:

SENECA L. BELCHEZ
OIC-Division Chief
Information Systems Management Division
IT Management Department

Approved by:

ALLAN M. DE VILLA
Acting Senior Manager
IT Management Department

JOVITA V. ARAGONA
Senior Vice President/Chief Information Officer
Office of the Senior Vice President Chief Information
Officer, Information Management Sector

REVISION HISTORY

Revision Number	Remarks
20240215	<ul style="list-style-type: none"> Initial Release for the SearchCaseRate, GenerateToken, ValidateEsoa and ValidateCF5 methods
20240216	<ul style="list-style-type: none"> Changed the parameter of SearchCaseRate method from object input parameters to comma separated list of input parameters
20240228	<ul style="list-style-type: none"> Added the requestQrAuthorization and inquireQrTrackingNo methods for the implementation of a use case for QR code of the eGov super app in the PhilHealth claims processing
20240306	<ul style="list-style-type: none"> Added the uploadClaims method; Added Annex B, C, D, E, F
20240418	<ul style="list-style-type: none"> Added softwareCertId in the header of uploadClaims method
20240423	<ul style="list-style-type: none"> Added the following methods: <ul style="list-style-type: none"> addRequiredDocument eClaimsFileCheck getClaimStatus getDoctorPAN getMemberPIN getUploadedClaimsMap getVoucherDetails isDoctorAccredited searchEmployer getDBServerDateTime getServerDateTime getServerVersion generatePBEFPDF
20240823	<ul style="list-style-type: none"> Changed the attribute pActualCharges to pChargesNetOfApplicableVat in eSOA XML, Updated the expected value of patientIs parameter in isClaimEligible method to the following: <ul style="list-style-type: none"> 'M' – patient is member (Self) 'S' – patient is spouse 'C' – patient is child 'P' – patient is parent, Added Data Dictionary of CF5, Added eSOA Libraries
20240910	<ul style="list-style-type: none"> Added softwareCertificateId in the header of Upload eClaims method (to be removed on the next version) Updated the attribute memberpPIN to memberPIN in isClaimEligible method
20240911	<ul style="list-style-type: none"> Added Medicine Library in the eSOA Library (Annex F)
20241111	<ul style="list-style-type: none"> Updated following items: <ul style="list-style-type: none"> Upload eClaims Method - The correct value of the attribute pUserName must be colon plus software certificate ID; Removed the softwareCertificateId header Add Required Document Method - The body should be a raw JSON; Updated the documentation Get Uploaded Claims Map Method - query parameter name "ReceiptTicketNo" is now "receiptTicketNumber"

	<ul style="list-style-type: none"> ○ Is Doctor Accredited Method - updated the key names documentation ○ Search Employer Method - updated the key names documentation ○ Is Claim Eligible Method - Added a definition for hospitalCode; Updated the sample JSON result documentation; ○ Generate PBEF PDF Method - updated the documentation ● Moved the Annexes B to F to a different file.
20250217	<ul style="list-style-type: none"> ● Revised the Data Dictionary and Document Type Definition (DTD) of eSOA. ● Changed the Drug and Medicine Library for eSOA. ● Updated the Data Dictionary of CF5 to set "N" (None) as the valid value for Laterality. ● Removed the requestQrAuthorization and inquireQrTrackingNo methods as this requirement has been deferred.

TABLE OF CONTENTS

APPROVAL SHEET	2
REVISION HISTORY	3
TABLE OF CONTENTS	5
INTRODUCTION	1
APPLICATION PROGRAM INTERFACE (API)	2
WEB SERVICE METHODS	3
1. Get Token Method	3
2. Validate eSOA Method	4
3. Validate CF5 Method	11
4. Upload eClaims Method	14
5. Search Case Rates Method	32
6. Add Required Document Method	37
7. EClaims File Check Method	39
8. Get Claim Status Method	41
9. Get Doctor PAN Method	44
10. Get Member PIN Method	46
11. Get Uploaded Claims Map Method	48
12. Get Voucher Details Method	50
13. Is Doctor Accredited Method	55
14. Search Employer Method	57
15. Get DB Server Date Time Method	59
16. Get Server Date Time Method	61
17. Get Server Version Method	63
18. Is Claim Eligible Method	64
19. Generate PBEF PDF Method	68
Annex A	70
Guidelines for the Data Encryption Using the Cipher Key of the Health Facility	70
Annex B	72
Document Type Code and Description	72
Annex C	74
Data Dictionary eClaimsUpload	74
Annex D	81
Data Dictionary ValidateEsoa	81
Annex E	83
Data Dictionary validateCF5	83
Annex E	84
eSOA Item Library	84

INTRODUCTION

The Philippine Health Insurance Corporation (PhilHealth) has been implementing the Electronic Claims (eClaims) System since 2016 for the All Case Rates (ACR) provider payment mechanism. This system was developed to automate the claims process, enabling improved turnaround times for claims processing and reimbursement of providers. Despite notable improvements in the efficiency of claims processing, there exists continued delays in turnaround time, and high numbers of “Return-to-Hospital” claims. These can be mainly attributed to the existence of manual procedures in the claims process, including the adjudication of claims, and the lack of front-end claims data validation mechanisms.

The Universal Health Care Act of 2019 aims to address these issues by mandating PhilHealth to shift to paying providers prospectively using Diagnosis-Related Groups (DRG). A key component of this reform is an enhanced version of the current eClaims system which accommodates the collection of data required for DRG grouping.

This guide is designed to provide in-house hospital developers and service providers essential insights and instructions on how to successfully navigate the electronic claims submission process within the context of the PhilHealth DRG system.

Disclaimer:

All information and content in this material is provided in good faith by Philippine Health Insurance Corporation-Universal Health Care Surge Team (PhilHealth-UHCST), and is based on sources believed to be reliable and accurate at the time of development. The PhilHealth-UHCST, and their respective officers, employees and supervisors, do not accept legal liability or responsibility for the Material, or any consequences arising from its use.

Philippine Health Insurance Corporation (PhilHealth) is committed to ensuring ease of availability at the point of care for all its beneficiaries.

APPLICATION PROGRAM INTERFACE (API)

The Application Programming Interface (API) of the PhilHealth e-Claims Web Service (PECWS) serves as the primary interface for interactions with PhilHealth's electronic claims processing system. This existing facade API is the gateway through which healthcare institutions and service providers can connect to and communicate with PhilHealth's infrastructure for the purpose of submitting electronic claims.

It acts as the bridge, allowing these entities to send claim-related data and information to PhilHealth, enabling the processing and evaluation of healthcare claims. Through the PECWS API, users can submit data regarding patient information, medical procedures, diagnoses, and other relevant details necessary for the billing and reimbursement of healthcare services.

WEB SERVICE METHODS

1. Get Token Method

This method will generate a token as authorization key to access API methods using the encrypted Accreditation Number and Software Certification ID of the health facility.

Endpoint

https://{pecws.domain}/PHIC/Claims3.0/getToken

Method

GET

Inputs

Header

Key	Value/Remarks
accreditationNo	The PhilHealth Accreditation Number (PAN) of the Health Facility
softwareCertificateId	The software certification ID for PECWS 3.0 of the Health Facility

Output

JSON object containing the following key-value pairs:

Key	Value/Remarks
success	A value of 'true' that data has been retrieved successfully
message	If an error was encountered during the execution of this method, this will contain the error message
result	The generated Token

Sample Output Payload

```
{
  "message": "Token is valid for 20 seconds",
  "result": "eyJhbGciOiJIUzI1NiBeyJhbGciOiJIUzI1Ni",
  "success": true
}
```


2. Validate eSOA Method

This method validates an encrypted eSOA XML file against the Document Type Definition (DTD), ensuring compliance with the required data format and valid values set by PhilHealth.

Important Note:

After successful validation, the eSOA XML must be encrypted using the PhilHealth Public Key and submitted as an attachment to electronic claims through the **EclaimsUpload** method, with the Document Type set to **ESA**.

Endpoint

`https://{pecws.domain}/PHIC/Claims3.0/validateeSOA`

Method

POST

Inputs

Header

Key	Values
token	PECWS authentication token

Body

The eSOA XML must be encrypted using the cipher key issued by PhilHealth to the Health Facility. The encrypted text contains the following key-value pairs.

Key	Value/Remarks
docMimeType	"text/xml"
hash	The computed hash value of the unencrypted text.
key1	"" (Empty string)
key2	"" (Empty string)
iv	The initialization vector used in the AES encryption
doc	The result of the encryption of eSOA XML text using the cipher key of the HF

Output

JSON object containing the following key-value pairs:

Key	Value/Remarks														
success	A value of 'true' that data has been retrieved successfully														
message	If an error was encountered during the execution of this method, this will contain the error message														
result	<div>Encrypted of result which contains the following key-value pairs</div> <table><tr><th>Key</th><th>Value/Remarks</th></tr><tr><td>docMimeType</td><td>"text/xml"</td></tr><tr><td>hash</td><td>The computed hash value of the unencrypted text.</td></tr><tr><td>key1</td><td>"" (Empty string)</td></tr><tr><td>key2</td><td>"" (Empty string)</td></tr><tr><td>iv</td><td>The initialization vector used in the AES encryption</td></tr><tr><td>doc</td><td>The result of the encryption of eSOA XML text using the cipher key of the HF</td></tr></table> <div>When decrypted, it may contains a response JSON object containing the error details. An array of errors as shown below:</div> <pre>{ "errors": ["string1", "string2", "string3"] }</pre>	Key	Value/Remarks	docMimeType	"text/xml"	hash	The computed hash value of the unencrypted text.	key1	"" (Empty string)	key2	"" (Empty string)	iv	The initialization vector used in the AES encryption	doc	The result of the encryption of eSOA XML text using the cipher key of the HF
Key	Value/Remarks														
docMimeType	"text/xml"														
hash	The computed hash value of the unencrypted text.														
key1	"" (Empty string)														
key2	"" (Empty string)														
iv	The initialization vector used in the AES encryption														
doc	The result of the encryption of eSOA XML text using the cipher key of the HF														

eSOA DTD

```
<!--
    Philippine Health Insurance Corporation
    Electronic Statement of Account Data Type Definition Version
    0.1 : 2022-08-19 12:24pm: Initial
    0.2 : 2023-02-13 03:51pm: Revised
    0.3 : 2023-05-30 09:14am: Added pCategory in ItemizedBillingItem
    0.4 : 2023-08-14 05:17pm: Changed pActualCharges to pChargesNetOfApplicableVat
    0.5 : 2025-02-17 01:08pm: Added the specs for the "Others" element and "Others"
    category
-->

<!ELEMENT eSOA (SummaryOfFees, ProfessionalFees, ItemizedBillingItems)>

<!ATTLIST eSOA
    pHciPan CDATA #REQUIRED
    pHciTransmittalId CDATA #REQUIRED>

<!ELEMENT SummaryOfFees (RoomAndBoard, DrugsAndMedicine, LaboratoryAndDiagnostic,
    OperatingRoomFees, MedicalSupplies, Others, PhilHealth, Balance)>

<!ELEMENT ProfessionalFees (ProfessionalFee*, PhilHealth, Balance)>

<!ELEMENT ItemizedBillingItems (ItemizedBillingItem+)>

<!ELEMENT PhilHealth EMPTY>
<!ATTLIST PhilHealth
    pTotalCaseRateAmount CDATA #REQUIRED
>

<!ELEMENT Balance EMPTY>
<!ATTLIST Balance
    pAmount CDATA #REQUIRED
>

<!ELEMENT OtherFundSource EMPTY>
<!ATTLIST OtherFundSource
    pDescription CDATA #REQUIRED
    pAmount CDATA #REQUIRED>

<!ELEMENT SummaryOfFee EMPTY>
<!ATTLIST SummaryOfFee
    pChargesNetOfApplicableVat CDATA #REQUIRED
    pSeniorCitizenDiscount CDATA #REQUIRED
    pPWDDiscount CDATA #REQUIRED
    pPCSO CDATA #REQUIRED
    pDSWD CDATA #REQUIRED
    pDOHMAP CDATA #REQUIRED
    pHMO CDATA #REQUIRED
>

<!ELEMENT RoomAndBoard (SummaryOfFee, OtherFundSource*)>

<!ELEMENT DrugsAndMedicine (SummaryOfFee, OtherFundSource*)>

<!ELEMENT LaboratoryAndDiagnostic (SummaryOfFee, OtherFundSource*)>

<!ELEMENT OperatingRoomFees (SummaryOfFee, OtherFundSource*)>
```

```

<!ELEMENT MedicalSupplies (SummaryOfFee, OtherFundSource*)>

<!ELEMENT Others (SummaryOfFee, OtherFundSource*)>

<!ELEMENT ProfessionalFee (ProfessionalInfo, SummaryOfFee)>

<!ELEMENT ProfessionalInfo EMPTY>
<!--ATTLIST ProfessionalInfo
  pPAN CDATA #REQUIRED
  pFirstName CDATA #REQUIRED
  pMiddleName CDATA #REQUIRED
  pLastName CDATA #REQUIRED
  pSuffixName CDATA #REQUIRED
-->

<!ELEMENT ItemizedBillingItem EMPTY>
<!--ATTLIST ItemizedBillingItem
  pServiceDate CDATA #REQUIRED
  pItemCode CDATA #REQUIRED
  pItemName CDATA #REQUIRED
  pUnitOfMeasurement CDATA #REQUIRED
  pUnitPrice CDATA #REQUIRED
  pQuantity CDATA #REQUIRED
  pTotalAmount CDATA #REQUIRED
  pCategory
(RoomAndBoard|DrugsAndMedicine|LaboratoryAndDiagnostic|OperatingRoomFees|MedicalSupplies|Others) #REQUIRED
-->

```

Sample Unencrypted eSOA XML

```

<eSOA
  pHciPan="HXXXXX678"
  pHciTransmittalId="HXXXXXXXXXX">
  <SummaryOfFees>
    <RoomAndBoard>
      <SummaryOfFee
        pChargesNetOfApplicableVat="2000.00"
        pSeniorCitizenDiscount="0"
        pPWDDiscount="0"
        pPCSO="2000.00"
        pDSWD="0"
        pDOHMAP="0"
        pHMO="0"
      />
    <OtherFundSource
      pDescription="Hello World Foundation"
      pAmount="3000.00"
    />
  </RoomAndBoard>
  <DrugsAndMedicine>
    <SummaryOfFee
      pChargesNetOfApplicableVat="3800"
      pSeniorCitizenDiscount="0"
      pPWDDiscount="0"
      pPCSO="0"
      pDSWD="1000.00"
      pDOHMAP="0"
      pHMO="0"
    />
  </DrugsAndMedicine>
</eSOA>

```

```

</DrugsAndMedicine>
<LaboratoryAndDiagnostic>
  <SummaryOfFee
    pChargesNetOfApplicableVat ="4000.00"
    pSeniorCitizenDiscount="0"
    pPWDDiscount="0"
    pPCSO="0"
    pDSWD="0"
    pDOHMAP="0"
    pHMO="0"
  />
</LaboratoryAndDiagnostic>
<OperatingRoomFees>
  <SummaryOfFee
    pChargesNetOfApplicableVat ="12000.00"
    pSeniorCitizenDiscount="0"
    pPWDDiscount="0"
    pPCSO="0"
    pDSWD="0"
    pDOHMAP="0"
    pHMO="0"
  />
</OperatingRoomFees>
<MedicalSupplies>
  <SummaryOfFee
    pChargesNetOfApplicableVat ="2000.00"
    pSeniorCitizenDiscount="1000.00"
    pPWDDiscount="0"
    pPCSO="0"
    pDSWD="0"
    pDOHMAP="0"
    pHMO="0"
  />
</MedicalSupplies>
<Others>
  <SummaryOfFee
    pChargesNetOfApplicableVat ="100.00"
    pSeniorCitizenDiscount="0"
    pPWDDiscount="0"
    pPCSO="0"
    pDSWD="0"
    pDOHMAP="0"
    pHMO="0"
  />
</Others>
<PhilHealth
  pTotalCaseRateAmount="15000.00"
/>
<Balance
  pAmount="3500.00"
/>
</SummaryOfFees>
<ProfessionalFees>
  <ProfessionalFee>
    <ProfessionalInfo
      pPAN="1234-4567890-1"
      pFirstName="JUAN"
      pMiddleName="TAMAD"
      pLastName="DELA CRUZ"
      pSuffixName=""
    />
    <SummaryOfFee
      pChargesNetOfApplicableVat ="0"
      pSeniorCitizenDiscount="0"
      pPWDDiscount="0"

```

```

        pPCSO="0"
        pDSWD="0"
        pDOHMAP="0"
        pHMO="0"
    />
</ProfessionalFee>
<PhilHealth
    pTotalCaseRateAmount="2000.00"
/>
<Balance
    pAmount="27180.00"
/>
</ProfessionalFees>
<ItemizedBillingItems>
    <ItemizedBillingItem
        pServiceDate="06-30-2021"
        pItemCode="PARAC0000000015SUPP40000000000"
        pCategory="DrugsAndMedicine"
        pItemName="PARACETAMOL 250 mg SUPPOSITORY"
        pUnitOfMeasurement=""
        pUnitPrice="700.00"
        pQuantity="5"
        pTotalAmount="3500.00"
    />
    <ItemizedBillingItem
        pServiceDate="06-30-2021"
        pItemCode=""
        pCategory="DrugsAndMedicine"
        pItemName="SAMPLE MEDICINE NOT INCLUDED IN THE LIBRARY"
        pUnitOfMeasurement="SAMPLE UNIT"
        pUnitPrice="300.00"
        pQuantity="1"
        pTotalAmount="300.00"
    />
    <ItemizedBillingItem
        pServiceDate="06-30-2021"
        pItemCode="1545"
        pCategory="LaboratoryAndDiagnostic"
        pItemName="HEMATOLOGY: CBC"
        pUnitOfMeasurement=""
        pUnitPrice="500.00"
        pQuantity="3"
        pTotalAmount="4000.00"
    />
    <ItemizedBillingItem
        pServiceDate="06-30-2021"
        pItemCode="1898"
        pCategory="OperatingRoomFees"
        pItemName="OPERATING ROOM"
        pUnitOfMeasurement=""
        pUnitPrice="12000"
        pQuantity="1"
        pTotalAmount="12000"
    />
    <ItemizedBillingItem
        pServiceDate="06-30-2021"
        pItemCode="175"
        pCategory="MedicalSupplies"
        pItemName="GLOVES STERILE, S-6.0 NON LATEX"
        pUnitOfMeasurement="PAIR"
        pUnitPrice="500.00"
        pQuantity="1"
        pTotalAmount="500"
    />
    <ItemizedBillingItem

```

```

        pServiceDate="06-30-2021"
        pItemCode="1288"
        pCategory="MedicalSupplies"
        pItemName="FACE MASK (EAR LOOP)DISPOSABLE"
        pUnitOfMeasurement="Box"
        pUnitPrice="500.00"
        pQuantity="3"
        pTotalAmount="1500.00"
    />
    <ItemizedBillingItem
        pServiceDate="06-30-2021"
        pItemCode=""
        pCategory="Others"
        pItemName="SAMPLE ITEM FOR OTHERS CATEGORY"
        pUnitOfMeasurement="SAMPLE UNIT"
        pUnitPrice="25.00"
        pQuantity="4"
        pTotalAmount="100.00"
    />
</ItemizedBillingItems>
</eSOA>

```

Sample Encrypted XML as Input Payload

```

{
  "docMimeType": "text/xml",
  "hash": "dfe701c0c9bbca0678300540591fac928f7dc8f4d74d9771a3841317e7a99aec",
  "key1": "",
  "key2": "",
  "iv": "y1jPMxvQE2aJPVnqqn1pDQ==",
  "doc": "Pms1FWFZT+odAp0qf2zMmroSUr3lYgDFnhYeJqBkuhJNMJU5geSEN=="
}

```

Sample Output Payload

```

{
  "message": "",
  "result": {
    "doc": "t1b2b2NF4PeAS1VoVOX3xIk/GMyB0s/bR6C56VUG2f2Xhka==",
    "docMimeType": "text/xml",
    "hash": "91cb1c3361f70f56ad5df326a000f5c",
    "iv": "oW/4PLliOnwgPNax1d6QWQ==",
    "key1": "",
    "key2": ""
  },
  "success": true
}

```

3. Validate CF5 Method

This method validates an encrypted CF5 XML file against the Document Type Definition (DTD), ensuring compliance with the required data format and values set by PhilHealth.

Important Note:

After successful validation, CF5 XML should be encrypted using PhilHealth Public Key and submitted as attachment to electronic claims using the **EclaimsUpload Method**, with the **Document Type** set to **CF5**.

Endpoint

`https://{pecws.domain}/PHIC/Claims3.0/validateCF5`

Method

POST

Inputs

Headers

Key	Values
token	PECWS authentication token

Body

The JSON object contains the following key-value pairs:

Name	Values
cf5	The CF5 XML text encrypted using cipher key of the HF <pre>{ "docMimeType": "text/xml", "hash": "dc8f4d74d977dfe701c0c9bbca000540591fac928f71a3841317e7a99aec", "key1": "", "key2": "", "iv": "y1jPMxvQE2aJPVnqqn1pDQ==", "doc": "PMs1FWFZT+odAp0qf2zMmroSUR3lYgDFnhYeJqBkuhJNMJU5geEN==" }</pre>
eclaims	The e-Claims XML text encrypted using cipher key of the HF <pre>{ "docMimeType": "text/xml", "hash": "dc8f4d74d977dfe701c0c9bbca0678300540591fa28f71a3841317e7a99aec", "key1": "", "key2": "", "iv": "y1jPMxvQE2aJPVnqqn1pDQ==", "doc": "PMs1FWFZT+ag2dbj3gs123fsfq3ahqasfwnosfaqwrgssrjlsrqarewr==" }</pre>

CF5 DTD

```
<?xml version='1.0' encoding='UTF-8'?>
<!--
  DRG E-Claims DTD History
  DIAGNOSTIC RELATED GROUP Version 1.0
  Version History
    1.0 2023-05-16 by JHEONARD
      : Testing part since e form not in thier final form

  DRG E-Claims DTD History
  DIAGNOSTIC RELATED GROUP Version 1.1
  Version History
    1.1 2023-06-7 by JHEONARD
      : Final Forms of DRG E-Forms
      : Add possible required for claims but not yet final

  DRG E-Claims DTD History
  DIAGNOSTIC RELATED GROUP Version 1.3
  Version History
    1.2 2023-07-12 by JHEONARD
      : Remove duplicate data

  DRG E-Claims DTD History
  DIAGNOSTIC RELATED GROUP Version 1.3
  Version History
    1.3 2024-02-15 by JHEONARD
      : Remove duplicate data
      : Remove Attributes (Series,Lhio,Admission Time)
      : Rename Tag (DRG - CF5)
-->

<!ELEMENT CF5 (DRGCLAIM)>
<!ATTLIST CF5 pHospitalCode CDATA #REQUIRED
>

<!ELEMENT DRGCLAIM (SECONDARYDIAGS,PROCEDURES)>
<!ATTLIST DRGCLAIM
  PrimaryCode CDATA #REQUIRED
  NewBornAdmWeight CDATA #REQUIRED
  Remarks CDATA #REQUIRED
  ClaimNumber CDATA #REQUIRED
>

<!ELEMENT SECONDARYDIAGS (SECONDARYDIAG)>
<!ELEMENT SECONDARYDIAG EMPTY>
<!ATTLIST SECONDARYDIAG
  SecondaryCode CDATA #REQUIRED
  Remarks CDATA #REQUIRED
>

<!ELEMENT PROCEDURES (PROCEDURE)+>
<!ELEMENT PROCEDURE EMPTY>
<!ATTLIST PROCEDURE
  RvsCode CDATA #REQUIRED
  Laterality CDATA #REQUIRED
  Ext1 CDATA #REQUIRED
  Ext2 CDATA #REQUIRED
  Remarks CDATA #REQUIRED
>
```

Sample Unencrypted CF5 XML

```
<CF5
  pHospitalCode="300806">
<DRGCLAIM
  ClaimNumber="300806-20221216-1-1"
  PrimaryCode="A00.0"
  NewBornAdmWeight=""
  Remarks="">
  <SECONDARYDIAGS>
    <SECONDARYDIAG
      SecondaryCode="A00.1"
      Remarks=""/>
    </SECONDARYDIAGS>
  <PROCEDURES>
    <PROCEDURE
      RvsCode=""
      Laterality=""
      Ext1=""
      Ext2=""
      Remarks=""/>
    </PROCEDURES>
  </DRGCLAIM>
</CF5>
```

Sample Input Payload

```
{
  "cf5": {
    "docMimeType": "text/xml",
    "hash": "dc8f4d74d977dfe701c0c9bbca0678300540591fac928f71a3841317e7a99aec",
    "key1": "",
    "key2": "",
    "iv": "y1jPMxvQE2aJPVnqqn1pDQ==",
    "doc": "PMs1FWFZT+odAp0qf2zMmroSUR3lYgDFnhYeJqBkuhJNMJU5geEN=="
  },
  "eclaims": {
    "docMimeType": "text/xml",
    "hash": "dc8f4d74d977dfe701c0c9bbca0678300540591fac928f71a3841317e7a99aec",
    "key1": "",
    "key2": "",
    "iv": "y1jPMxvQE2aJPVnqqn1pDQ==",
    "doc": "PMs1FWFZT+odAp0qf2zMmroSUR3lYgDFnhYeJqBkuhJNMJU5geEN=="
  }
}
```

4. Upload eClaims Method

This method allows health facilities to transmit electronic claims (eClaims) files to PhilHealth. It ensures Document Type Definition (DTD) compliance and validates XML element attributes based on the eClaims XML Elements Attribute Definition table.

The transmission date serves as the official date received for the uploaded claims, which will be used to measure the Turnaround Time (TAT).

Endpoint

`https://{pecws.domain}/PHIC/Claims3.0/uploadeClaims`

Method

POST

Inputs

Header

Key	Value
token	PECWS authentication token

Body

The body is the eclaims XML in JSON format as encrypted using the cipher key of the Health Facility. The JSON object contains the following key-value pairs:

Key	Value/Remarks
docMimeType	"text/xml"
hash	The computed hash value of the unencrypted text.
key1	"" (Empty string)
key2	"" (Empty string)
iv	The initialization vector used in the AES encryption
doc	The encrypted e-claim XML text.

Output

This method returns a JSON object with the key-value pairs described below:

Key	Value/Remarks
-----	---------------

success	A value of 'true' indicates a successful operation														
message	If an error was encountered during the execution of this method, this will contain the error message														
result	<p>The JSON object as the result of the encryption (using the cipher key of the health facility) of the XML text containing the Receipt Ticket Number and other data about the processing of the submitted e-claim data.</p> <table border="1"> <thead> <tr> <th>Key</th><th>Value/Remarks</th></tr> </thead> <tbody> <tr> <td>docMimeType</td><td>"text/xml"</td></tr> <tr> <td>hash</td><td>The computed hash value of the unencrypted text.</td></tr> <tr> <td>key1</td><td>"" (Empty string)</td></tr> <tr> <td>key2</td><td>"" (Empty string)</td></tr> <tr> <td>iv</td><td>The initialization vector used in the AES encryption</td></tr> <tr> <td>doc</td><td>The encrypted text of the JSON object that contains the records of the matching benefit packages. Sample XML text and the DTD of the XML text is shown below.</td></tr> </tbody> </table>	Key	Value/Remarks	docMimeType	"text/xml"	hash	The computed hash value of the unencrypted text.	key1	"" (Empty string)	key2	"" (Empty string)	iv	The initialization vector used in the AES encryption	doc	The encrypted text of the JSON object that contains the records of the matching benefit packages. Sample XML text and the DTD of the XML text is shown below.
Key	Value/Remarks														
docMimeType	"text/xml"														
hash	The computed hash value of the unencrypted text.														
key1	"" (Empty string)														
key2	"" (Empty string)														
iv	The initialization vector used in the AES encryption														
doc	The encrypted text of the JSON object that contains the records of the matching benefit packages. Sample XML text and the DTD of the XML text is shown below.														

eClaims DTD

<!--

Philippine Health Insurance Corporation
eClaims Document Type Definition Version 1.9

Version History

1.9.0 06-23-2017 08:16PM (MSM)

- : removed the enumerated list values for the pDocumentType attribute for flexibility of adding new elements
- : added pServiceProvider attribute to the eCLAIMS element to indicate the provider of the system used to encode and submit the e-claims XML file

1.8.0 09-19-2016 05:13PM (MSM)

- : added "IMRT" element as additional to the list of supported repetitive procedures
- : added "APR" element to CF2 for saving the data of the Consent to Access Patient Record section of RF2
- : added "pHasAttachedSOA" attribute to CF2 element
- : added new "CATARACTINFO" element to hold info about data like IOL sticker number for cataract operation. The "CATARACT" element will be deprecated later
- : changed the comments for the pThumbMarkedBy attribute

1.7.6 07-31-2015 10:38PM

- : remove "()" in this part (BENEFITS)

1.7.5 revised: 06-23-2015 1:19PM

- : added "P" (Lifetime Member) in the pMembershipType attribute acceptable values
- : added "ANR" (Anesthesia Record) & "HDR" (Hemodialysis Record) in the pDocumentType attribute acceptable values

1.7.4 revised: 11-03-2014 9:52am

- : added "MRF" (PhilHealth Member Registration Form) in the pDocumentType

1.7.3 revised: 01-28-2014 9:58am

- : added pPreAuthDate in ZBENEFIT element
- : added pCaseRateAmount in CASERATE element
- : added pDoctorSignDate in PROFESSIONALS element
- : added pPatientType & pIsEmergency in CLAIM element
- : transfer pCataractPreAuth to CATARACT element

1.7.2 revised: 12-16-2013 5:03pm

- : abstracted the Particulars from CF3 outside
- : added the CF3_OLD and CF3_NEW elements

1.7.1 revised: 11-04-2013 3:31pm

- : updated the elements based on the claim forms version 11_04_2013

1.7 revised: 09-24-2013 08:56am

- : major revision to cater all case rate policy

1.6.1 revised: 01-29-2013 01:29pm

- : added CSF in list of document types

1.6 revised: 01-25-2013 01:11pm

- : updated the new coding of documents

1.5 revised: 10-18-2012 10:01am

- : abstracted the Observation Codes from ZBenefit outside
- : added Official Receipt details

1.4 revised: 07-16-2012 01:40pm

- : added support for observation for z-benefits

1.3 revised: 06-01-2012 09:04am

- : added elements for document urls
- : added elements for Z Benefit

1.2 revised: 02-24-2012 06:16pm

- : added code for packages
- : added codes for case-rates

1.1 revised: 08-19-2011 11:02am

- : added requirements for case rate

1.0 revised: 11-12-2010 04:42pm

- : initial version

-->

```

<!ENTITY Ntilde "&#209;">
<!ENTITY ntilde "&#241;">
<!ENTITY nbsp "&#160;">

<!ELEMENT eCLAIMS (eTRANSMITTAL)>
<!ATTLIST eCLAIMS
  pUserName CDATA #REQUIRED
  pUserPassword CDATA #REQUIRED
  pHospitalCode CDATA #REQUIRED
  pHospitalEmail CDATA #REQUIRED
  pServiceProvider CDATA #IMPLIED>

<!ELEMENT eTRANSMITTAL (CLAIM+)>
<!ATTLIST eTRANSMITTAL
  pHospitalTransmittalNo CDATA #REQUIRED
  pTotalClaims CDATA #REQUIRED>

<!ELEMENT CLAIM (CF1, CF2, (ALLCASERATE | ZBENEFIT), CF3?, PARTICULARS?, RECEIPTS?, DOCUMENTS)>
<!ATTLIST CLAIM
  pClaimNumber CDATA #REQUIRED
  pTrackingNumber CDATA #REQUIRED
  pPhilhealthClaimType (ALL-CASE-RATE|Z-BENEFIT) #REQUIRED
  pPatientType (I|O) #REQUIRED
  pIsEmergency (Y|N) #REQUIRED>

<!ELEMENT CF1 EMPTY>
<!ATTLIST CF1
  pMemberPIN CDATA #REQUIRED
  pMemberLastName CDATA #REQUIRED
  pMemberFirstName CDATA #REQUIRED
  pMemberSuffix CDATA #REQUIRED
  pMemberMiddleName CDATA #REQUIRED
  pMemberBirthDate CDATA #REQUIRED
  pMemberShipType (S|G|I|NS|NO|PS|PG|P) #REQUIRED
  pMailingAddress CDATA #REQUIRED
  pZipCode CDATA #REQUIRED
  pMemberSex (M|F) #REQUIRED
  pLandlineNo CDATA #REQUIRED
  pMobileNo CDATA #REQUIRED
  pEmailAddress CDATA #REQUIRED
  pPatientIs (M|S|C|P) #REQUIRED
  pPatientPIN CDATA #REQUIRED
  pPatientLastName CDATA #REQUIRED
  pPatientFirstName CDATA #REQUIRED
  pPatientSuffix CDATA #REQUIRED
  pPatientMiddleName CDATA #REQUIRED
  pPatientBirthDate CDATA #REQUIRED
  pPatientSex (M|F) #REQUIRED
  pPEN CDATA #REQUIRED
  pEmployerName CDATA #REQUIRED>

<!ELEMENT CF2 (DIAGNOSIS, SPECIAL, PROFESSIONALS+, CONSUMPTION, APR?)>
<!ATTLIST CF2
  pPatientReferred (Y|N) #REQUIRED
  pReferredIHCPAccreCode CDATA #REQUIRED
  pAdmissionDate CDATA #REQUIRED
  pAdmissionTime CDATA #REQUIRED
  pDischargeDate CDATA #REQUIRED
  pDischargeTime CDATA #REQUIRED
  pDisposition (I|R|H|A|E|T) #REQUIRED
  pExpiredDate CDATA #REQUIRED
  pExpiredTime CDATA #REQUIRED
  pReferralIHCPAccreCode CDATA #REQUIRED
  pReferralReasons CDATA #REQUIRED
  pAccommodationType (P|N) #REQUIRED

```

```

    pHasAttachedSOA (Y|N) #IMPLIED>

<!ELEMENT DIAGNOSIS (DISCHARGE+)>
<!ATTLIST DIAGNOSIS
    pAdmissionDiagnosis CDATA #REQUIRED>

<!ELEMENT DISCHARGE ((ICDCODE+, RVSCODES*)|(ICDCODE*, RVSCODES+))>
<!ATTLIST DISCHARGE
    pDischargeDiagnosis CDATA #REQUIRED>

<!ELEMENT ICDCODE EMPTY>
<!ATTLIST ICDCODE
    pICDCode CDATA #REQUIRED>

<!ELEMENT RVSCODES EMPTY>
<!ATTLIST RVSCODES
    pRelatedProcedure CDATA #REQUIRED
    pRVSCode CDATA #REQUIRED
    pProcedureDate CDATA #REQUIRED
    pLaterality (L|R|B|N) #REQUIRED>

<!ELEMENT SPECIAL (PROCEDURES?, MCP?, TBDOTS?, ABP?, NCP?, HIVAIDS?, CATARACTINFO?)>
<!ELEMENT PROCEDURES ((HEMODIALYSIS?, PERITONEAL?, LINAC?, COBALT?, TRANSFUSION?,
BRACHYTHERAPY?, CHEMOTHERAPY?, DEBRIDEMENT?, IMRT?),
(HEMODIALYSIS|PERITONEAL|LINAC|COBALT|TRANSFUSION|BRACHYTHERAPY|CHEMOTHERAPY|DEBRIDEMENT|IM
RT))>
<!ELEMENT HEMODIALYSIS (SESSIONS*)>
<!ELEMENT PERITONEAL (SESSIONS*)>
<!ELEMENT LINAC (SESSIONS*)>
<!ELEMENT COBALT (SESSIONS*)>
<!ELEMENT TRANSFUSION (SESSIONS*)>
<!ELEMENT BRACHYTHERAPY (SESSIONS*)>
<!ELEMENT CHEMOTHERAPY (SESSIONS*)>
<!ELEMENT DEBRIDEMENT (SESSIONS*)>
<!ELEMENT IMRT (SESSIONS*)>

<!ELEMENT SESSIONS EMPTY>
<!ATTLIST SESSIONS
    pSessionDate CDATA #REQUIRED>

<!ELEMENT MCP EMPTY>
<!ATTLIST MCP
    pCheckUpDate1 CDATA #REQUIRED
    pCheckUpDate2 CDATA #REQUIRED
    pCheckUpDate3 CDATA #REQUIRED
    pCheckUpDate4 CDATA #REQUIRED>

<!ELEMENT TBDOTS EMPTY>
<!ATTLIST TBDOTS
    pTBType (I|M) #REQUIRED
    pNTPCardNo CDATA #REQUIRED>

<!ELEMENT ABP EMPTY>
<!ATTLIST ABP
    pDay0ARV CDATA #REQUIRED
    pDay3ARV CDATA #REQUIRED
    pDay7ARV CDATA #REQUIRED
    pRIG CDATA #REQUIRED
    pABPOthers CDATA #REQUIRED
    pABPSpecify CDATA #REQUIRED>

<!ELEMENT NCP (ESSENTIAL?)>
<!ATTLIST NCP
    pEssentialNewbornCare (Y|N) #REQUIRED
    pNewbornHearingScreeningTest (Y|N) #REQUIRED

```

```

    pNewbornScreeningTest (Y|N) #REQUIRED
    pFilterCardNo CDATA #REQUIRED>

<!ELEMENT ESSENTIAL EMPTY>
<!ATTLIST ESSENTIAL
    pDrying (Y|N) #REQUIRED
    pSkinToSkin (Y|N) #REQUIRED
    pCordClamping (Y|N) #REQUIRED
    pProphylaxis (Y|N) #REQUIRED
    pWeighing (Y|N) #REQUIRED
    pVitaminK (Y|N) #REQUIRED
    pBCG (Y|N) #REQUIRED
    pNonSeparation (Y|N) #REQUIRED
    pHepatitisB (Y|N) #REQUIRED>

<!ELEMENT HIVAIDS EMPTY>
<!ATTLIST HIVAIDS
    pLaboratoryNumber CDATA #REQUIRED>

<!ELEMENT CATARACTINFO EMPTY>
<!ATTLIST CATARACTINFO
    pCataractPreAuth CDATA #REQUIRED
    pLeftEyeIOLStickerNumber CDATA #REQUIRED
    pLeftEyeIOLExpiryDate CDATA #REQUIRED
    pRightEyeIOLStickerNumber CDATA #REQUIRED
    pRightEyeIOLExpiryDate CDATA #REQUIRED>

<!ELEMENT PROFESSIONALS EMPTY>
<!ATTLIST PROFESSIONALS
    pDoctorAccreCode CDATA #REQUIRED
    pDoctorLastName CDATA #REQUIRED
    pDoctorFirstName CDATA #REQUIRED
    pDoctorMiddleName CDATA #REQUIRED
    pDoctorSuffix CDATA #REQUIRED
    pWithCoPay (Y|N) #REQUIRED
    pDoctorCoPay CDATA #REQUIRED
    pDoctorSignDate CDATA #REQUIRED>

<!ELEMENT CONSUMPTION (BENEFITS|(HCIFEES, PROFFEES, PURCHASES)))>
<!ATTLIST CONSUMPTION
    pEnoughBenefits (Y|N) #REQUIRED>

<!ELEMENT APR (APRBYPATSIG|APRBYPATREPSIG|APRBYTHUMBMARK)>

<!ELEMENT APRBYPATSIG EMPTY>
<!ATTLIST APRBYPATSIG
    pDateSigned CDATA #REQUIRED>

<!ELEMENT APRBYPATREPSIG ((DEFINEDPATREPREL|OTHERPATREPREL),
(DEFINEDREASONFORSIGNING|OTHERREASONFORSIGNING)))>
<!ATTLIST APRBYPATREPSIG
    pDateSigned CDATA #REQUIRED>

<!ELEMENT DEFINEDPATREPREL EMPTY>
<!--      pRelCode:          One of
                                S: Spouse
                                C: Child
                                P: Parent
                                I: Siblings
                                O: Others
                                -->

<!ATTLIST DEFINEDPATREPREL
    pRelCode (S|C|P|I) #REQUIRED>
<!ELEMENT OTHERPATREPREL EMPTY>
<!ATTLIST OTHERPATREPREL

```



```

        pRelCode CDATA #FIXED "0"
        pRelDesc CDATA #REQUIRED>

<!ELEMENT DEFINEDREASONFORSIGNING EMPTY>
<!--      pReasonCode:      One of
                                I: Patient is incapacitated
                                O: Other reasons. Should be specified in pReasonDesc
                                -->
<!ATTLIST DEFINEDREASONFORSIGNING
        pReasonCode (I) #REQUIRED>
<!ELEMENT OTHERREASONFORSIGNING EMPTY>
<!ATTLIST OTHERREASONFORSIGNING
        pReasonCode CDATA #FIXED "0"
        pReasonDesc CDATA #REQUIRED>

<!ELEMENT APRBYTHUMBMARK EMPTY>
<!--      pThumbmarkedBy:      One of
                                P: by patient/member
                                R: by representative
                                -->
<!ATTLIST APRBYTHUMBMARK
        pThumbmarkedBy (P|R) #REQUIRED>

<!ELEMENT BENEFITS EMPTY>
<!ATTLIST BENEFITS
        pTotalHCIFees CDATA #REQUIRED
        pTotalProfFees CDATA #REQUIRED
        pGrandTotal CDATA #REQUIRED>

<!ELEMENT HCIFEES EMPTY>
<!ATTLIST HCIFEES
        pTotalActualCharges CDATA #REQUIRED
        pDiscount CDATA #REQUIRED
        pPhilhealthBenefit CDATA #REQUIRED
        pTotalAmount CDATA #REQUIRED
        pMemberPatient (Y|N) #REQUIRED
        pHMO (Y|N) #REQUIRED
        pOthers (Y|N) #REQUIRED>

<!ELEMENT PROFFEES EMPTY>
<!ATTLIST PROFFEES
        pTotalActualCharges CDATA #REQUIRED
        pDiscount CDATA #REQUIRED
        pPhilhealthBenefit CDATA #REQUIRED
        pTotalAmount CDATA #REQUIRED
        pMemberPatient (Y|N) #REQUIRED
        pHMO (Y|N) #REQUIRED
        pOthers (Y|N) #REQUIRED>

<!ELEMENT PURCHASES EMPTY>
<!ATTLIST PURCHASES
        pDrugsMedicinesSupplies (Y|N) #REQUIRED
        pDMSTotalAmount CDATA #REQUIRED
        pExaminations (Y|N) #REQUIRED
        pExamTotalAmount CDATA #REQUIRED>

<!ELEMENT ALLCASERATE (CASERATE+)>
<!ELEMENT CASERATE (CATARACT?)>
<!ATTLIST CASERATE
        pCaseRateCode CDATA #REQUIRED
        pICDCode CDATA #REQUIRED
        pRVSCode CDATA #REQUIRED
        pCaseRateAmount CDATA #REQUIRED>

```

```

<!ELEMENT CATARACT EMPTY>
<!ATTLIST CATARACT
    pCataractPreAuth CDATA #REQUIRED>

<!ELEMENT ZBENEFIT EMPTY>
<!ATTLIST ZBENEFIT
    pZBenefitCode
        (Z0011|Z0012|Z0013|Z0021|Z0022|Z003|Z0041|Z0042|Z0051|Z0052|Z0061|Z0062|Z0071|Z0072|Z0081|Z0
082|Z0091|Z0092) #REQUIRED
    pPreAuthDate CDATA #REQUIRED>

<!ELEMENT CF3 (CF3_OLD?, CF3_NEW?)>
<!ELEMENT CF3_OLD (PHEX, MATERNITY?)>
<!ATTLIST CF3_OLD
    pChiefComplaint CDATA #REQUIRED
    pBriefHistory CDATA #REQUIRED
    pCourseWard CDATA #REQUIRED
    pPertinentFindings CDATA #REQUIRED>

<!ELEMENT MATERNITY (PRENATAL, DELIVERY, POSTPARTUM)>
<!ELEMENT PRENATAL (CLINICALHIST, OBSTETRIC, MEDISURG, CONSULTATION+)>
<!ATTLIST PRENATAL
    pPrenatalConsultation CDATA #REQUIRED
    pMCPOrientation (Y|N) #REQUIRED
    pExpectedDeliveryDate CDATA #REQUIRED>

<!ELEMENT CLINICALHIST EMPTY>
<!ATTLIST CLINICALHIST
    pVitalSigns (Y|N) #REQUIRED
    pPregnancyLowRisk (Y|N) #REQUIRED
    pLMP CDATA #REQUIRED
    pMenarcheAge CDATA #REQUIRED
    pObstetricG CDATA #REQUIRED
    pObstetricP CDATA #REQUIRED
    pObstetric_T CDATA #REQUIRED
    pObstetric_P CDATA #REQUIRED
    pObstetric_A CDATA #REQUIRED
    pObstetric_L CDATA #REQUIRED>

<!ELEMENT OBSTETRIC EMPTY>
<!ATTLIST OBSTETRIC
    pMultiplePregnancy (Y|N) #REQUIRED
    pOvarianCyst (Y|N) #REQUIRED
    pMyomaUteri (Y|N) #REQUIRED
    pPlacentaPrevia (Y|N) #REQUIRED
    pMiscarriages (Y|N) #REQUIRED
    pStillBirth (Y|N) #REQUIRED
    pPreEclampsia (Y|N) #REQUIRED
    pEclampsia (Y|N) #REQUIRED
    pPrematureContraction (Y|N) #REQUIRED>

<!ELEMENT MEDISURG EMPTY>
<!ATTLIST MEDISURG
    pHypertension (Y|N) #REQUIRED
    pHeartDisease (Y|N) #REQUIRED
    pDiabetes (Y|N) #REQUIRED
    pThyroidDisaster (Y|N) #REQUIRED
    pObesity (Y|N) #REQUIRED
    pAsthma (Y|N) #REQUIRED
    pEpilepsy (Y|N) #REQUIRED
    pRenalDisease (Y|N) #REQUIRED
    pBleedingDisorders (Y|N) #REQUIRED
    pPreviousCS (Y|N) #REQUIRED
    pUterineMyomectomy (Y|N) #REQUIRED>

```

```

<!ELEMENT CONSULTATION EMPTY>
<!ATTLIST CONSULTATION
  pVisitDate CDATA #REQUIRED
  pAOGWeeks CDATA #REQUIRED
  pWeight CDATA #REQUIRED
  pCardiacRate CDATA #REQUIRED
  pRespiratoryRate CDATA #REQUIRED
  pBloodPressure CDATA #REQUIRED
  pTemperature CDATA #REQUIRED>

<!ELEMENT DELIVERY EMPTY>
<!ATTLIST DELIVERY
  pDeliveryDate CDATA #REQUIRED
  pDeliveryTime CDATA #REQUIRED
  pObstetricIndex CDATA #REQUIRED
  pAOGLMP CDATA #REQUIRED
  pDeliveryManner CDATA #REQUIRED
  pPresentation CDATA #REQUIRED
  pFetalOutcome CDATA #REQUIRED
  pSex CDATA #REQUIRED
  pBirthWeight CDATA #REQUIRED
  pAPGARScore CDATA #REQUIRED
  pPostpartum CDATA #REQUIRED>

<!ELEMENT POSTPARTUM EMPTY>
<!ATTLIST POSTPARTUM
  pPerinealWoundCare (Y|N) #REQUIRED
  pPerinealRemarks CDATA #REQUIRED
  pMaternalComplications (Y|N) #REQUIRED
  pMaternalRemarks CDATA #REQUIRED
  pBreastFeeding (Y|N) #REQUIRED
  pBreastFeedingRemarks CDATA #REQUIRED
  pFamilyPlanning (Y|N) #REQUIRED
  pFamilyPlanningRemarks CDATA #REQUIRED
  pPlanningService (Y|N) #REQUIRED
  pPlanningServiceRemarks CDATA #REQUIRED
  pSurgicalSterilization (Y|N) #REQUIRED
  pSterilizationRemarks CDATA #REQUIRED
  pFollowupSchedule (Y|N) #REQUIRED
  pFollowupScheduleRemarks CDATA #REQUIRED>

<!ELEMENT CF3_NEW (ADMITREASON?, COURSE?)>

<!ELEMENT ADMITREASON (CLINICAL+, LABDIAG+, PHEX)>
<!ATTLIST ADMITREASON
  pBriefHistory CDATA #REQUIRED
  pReferredReason CDATA #REQUIRED
  pIntensive (Y|N) #REQUIRED
  pMaintenance (Y|N) #REQUIRED>

<!ELEMENT CLINICAL EMPTY>
<!ATTLIST CLINICAL
  pCriteria CDATA #REQUIRED>

<!ELEMENT LABDIAG EMPTY>
<!ATTLIST LABDIAG
  pCriteria CDATA #REQUIRED>

<!ELEMENT PHEX EMPTY>
<!ATTLIST PHEX
  pBP CDATA #REQUIRED
  pCR CDATA #REQUIRED
  pRR CDATA #REQUIRED
  pTemp CDATA #REQUIRED
  pHEENT CDATA #REQUIRED

```

```

    pChestLungs CDATA #REQUIRED
    pCVS CDATA #REQUIRED
    pAbdomen CDATA #REQUIRED
    pGUIE CDATA #REQUIRED
    pSkinExtremities CDATA #REQUIRED
    pNeuroExam CDATA #REQUIRED>

<!ELEMENT COURSE (WARD+)>
<!ELEMENT WARD EMPTY>
<!ATTLIST WARD
    pCourseDate CDATA #REQUIRED
    pFindings CDATA #REQUIRED
    pAction CDATA #REQUIRED>

<!ELEMENT PARTICULARS ((DRGMED+|XLSO+), (DRGMED*|XLSO*))>
<!ELEMENT DRGMED EMPTY>
<!ATTLIST DRGMED
    pPurchaseDate CDATA #REQUIRED
    pDrugCode CDATA #REQUIRED
    pPNDFCode CDATA #REQUIRED
    pGenericName CDATA #REQUIRED
    pBrandName CDATA #REQUIRED
    pPreparation CDATA #REQUIRED
    pQuantity CDATA #REQUIRED>
<!ELEMENT XLSO EMPTY>
<!ATTLIST XLSO
    pDiagnosticDate CDATA #REQUIRED
    pDiagnosticType (IMAGING|LABORATORY|SUPPLIES|OTHERS) #REQUIRED
    pDiagnosticName CDATA #REQUIRED
    pQuantity CDATA #REQUIRED>

<!ELEMENT RECEIPTS (RECEIPT+)>
<!ELEMENT RECEIPT (ITEM+)>
<!ATTLIST RECEIPT
    pCompanyName CDATA #REQUIRED
    pCompanyTIN CDATA #REQUIRED
    pBIRPermitNumber CDATA #REQUIRED
    pReceiptNumber CDATA #REQUIRED
    pReceiptDate CDATA #REQUIRED
    pVATExemptSale CDATA #REQUIRED
    pVAT CDATA #REQUIRED
    pTotal CDATA #REQUIRED>
<!ELEMENT ITEM EMPTY>
<!ATTLIST ITEM
    pQuantity CDATA #REQUIRED
    pUnitPrice CDATA #REQUIRED
    pDescription CDATA #REQUIRED
    pAmount CDATA #REQUIRED>

<!ELEMENT DOCUMENTS (DOCUMENT+)>
<!ELEMENT DOCUMENT EMPTY>
<!ATTLIST DOCUMENT
    pDocumentType CDATA #REQUIRED
    pDocumentURL CDATA #REQUIRED>

```

Sample eClaims XML

```
<eCLAIMS
  pUserName=":SOFTWARE-CERTIFICATE-ID-HERE"
  pUserPassword=""
  pHospitalCode="123456"
  pHospitalEmail=email@yahoo.com
  pServiceProvider="SAMPLE SERVICE PROVIDER">
  <eTRANSMITTAL
    pHospitalTransmittalNo="20160901"
    pTotalClaims="1">
    <CLAIM
      pClaimNumber="123456-20160930-2"
      pTrackingNumber=""
      pPhilhealthClaimType="ALL-CASE-RATE"
      pPatientType="I"
      pIsEmergency="N">
      <CF1
        pMemberPIN="072007271094"
        pMemberLastName="DELA CRUZ"
        pMemberFirstName="JUAN"
        pMemberSuffix=""
        pMemberMiddleName="OCAMPO"
        pMemberBirthDate="09-19-1973"
        pMemberShipType="G"
        pMailingAddress="PHILIPPINES"
        pZipCode="1234"
        pMemberSex="M"
        pLandlineNo=""
        pMobileNo=""
        pEmailAddress="delacruzjuan@yahoo.com.ph"
        pPatientIs="M"
        pPatientPIN="072007271094"
        pPatientLastName="DELA CRUZ"
        pPatientFirstName="JUAN"
        pPatientSuffix=""
        pPatientMiddleName="OCAMPO"
        pPatientBirthDate="09-19-1973"
        pPatientSex="M"
        pPEN="11-047400000-2"
        pEmployerName="PHILHEALTH"/>
      <CF2
        pPatientReferred="Y"
        pReferredIHCPAccreCode="H12345678"
        pAdmissionDate="09-01-2016"
        pAdmissionTime="01:00:00PM"
        pDischargeDate="09-03-2016"
        pDischargeTime="03:00:00PM"
        pDisposition="I"
        pExpiredDate=""
        pExpiredTime=""
        pReferralIHCPAccreCode=""
        pReferralReasons=""
        pAccommodationType="N"
        pHasAttachedSOA="Y">
      <DIAGNOSIS
        pAdmissionDiagnosis="PNEUMONIA">
      <DISCHARGE
        pDischargeDiagnosis="END STAGE RENAL DISEASE">
        <ICDCODE pICDCode="013.012"/>
        <ICDCODE pICDCode="013.013"/>
        <RVSCODES
          pRelatedProcedure="HEMODIALYSIS"
          pRVSCode="90935"
          pProcedureDate="08-26-2009"
```

```

        pLaterality="L"/>
</DISCHARGE>
<DISCHARGE
    pDischargeDiagnosis="DENGUE">
        <ICDCODE pICDCode="A90.0"/>
    </DISCHARGE>
</DIAGNOSIS>
<SPECIAL>
    <!-- For Repetitive Procedures -->
    <PROCEDURES>
        <HEMODIALYSIS>
            <SESSIONS pSessionDate="08-25-2009"/>
            <SESSIONS pSessionDate="08-26-2009"/>
        </HEMODIALYSIS>
        <CHEMOTHERAPY>
            <SESSIONS pSessionDate="08-27-2009"/>
        </CHEMOTHERAPY>
    </PROCEDURES>
    <!-- For MCP Package -->
    <!--MCP
        pCheckUpDate1="08-25-2009"
        pCheckUpDate2="08-26-2009"
        pCheckUpDate3="08-27-2009"
        pCheckUpDate4="08-28-2009"/-->
    <!-- For TB DOTS Package -->
    <!--TBDOTS
        pTBType="I"
        pNTPCardNo=""/-->
    <!-- For Animal Bite Package -->
    <!--ABP
        pDay0ARV="08-25-2009"
        pDay3ARV="08-25-2009"
        pDay7ARV="08-25-2009"
        pRIG="08-25-2009"
        pABPOthers=""
        pABPSpecify=""/-->
    <!-- For Newborn Care Package -->
    <!--NCP
        pEssentialNewbornCare="Y"
        pNewbornHearingScreeningTest="N"
        pNewbornScreeningTest="N"
        pFilterCardNo="">
        <ESSENTIAL
            pDrying="Y"
            pSkinToSkin="Y"
            pCordClamping="Y"
            pProphylaxis="Y"
            pWeighing="Y"
            pVitaminK="Y"
            pBCG="Y"
            pNonSeparation="Y"
            pHepatitisB="Y"/>
    </NCP-->
    <!-- For Outpatient HIV/AIDS Treatment Package -->
    <!--HIVAIDS
        pLaboratoryNumber=""/-->
    <!--CATARACTINFO
        pCataractPreAuth=""
        pLeftEyeIOLStickerNumber=""
        pLeftEyeIOLExpiryDate=""
        pRightEyeIOLStickerNumber=""
        pRightEyeIOLExpiryDate=""/-->
</SPECIAL>
<PROFESSIONALS
    pDoctorAccreCode="1234-1527066-1"

```

```

pDoctorLastName="TEST2"
pDoctorFirstName="TEST2"
pDoctorMiddleName="TEST2"
pDoctorSuffix=""
pWithCoPay="Y"
pDoctorCoPay="1000"
pDoctorSignDate="08-25-2009"/>
<CONSUMPTION
  pEnoughBenefits="Y">
    <!-- if pEnoughBenefits="Y"-->
    <BENEFITS
      pTotalHCIFees="1000"
      pTotalProfFees="1000"
      pGrandTotal="2000"/>
    <!-- if pEnoughBenefits="N"-->
    <!--
    <HCIFEES
      pTotalActualCharges="2000"
      pDiscount="1800"
      pPhilhealthBenefit="1500"
      pTotalAmount="300"
      pMemberPatient="Y"
      pHMO="N"
      pOthers="N"/>
    <PROFFEES
      pTotalActualCharges="3000"
      pDiscount="2500"
      pPhilhealthBenefit="1500"
      pTotalAmount="1000"
      pMemberPatient="Y"
      pHMO="Y"
      pOthers="N"/>
    <PURCHASES
      pDrugsMedicinesSupplies="Y"
      pDMSTotalAmount="1000"
      pExaminations="N"
      pExamTotalAmount=""/>
    -->
  </CONSUMPTION>
</APR>
  <APRBYPATREPSIG pDateSigned="08-26-2009">
    <DEFINEDPATREPREL pRelCode="S" />
    <OTHERREASONFORSIGNING pReasonDesc="MEMBER IS MISSING" />
  </APRBYPATREPSIG>
</APR>
</CF2>
<!-- pPhilhealthClaimType="ALL-CASE-RATE" -->
<ALLCASERATE>
  <CASERATE
    pCaseRateCode="CR0001"
    pICDCode="A90"
    pRVSCode=""
    pCaseRateAmount="10000"/>
  <CASERATE
    pCaseRateCode="CR0002"
    pICDCode=""
    pRVSCode="90935"
    pCaseRateAmount="2600"/>
</ALLCASERATE>
<!-- pPhilhealthClaimType="Z-BENEFIT" -->
<!--ZBENEFIT
  pZBenefitCode="Z0011">
</ZBENEFIT-->
<CF3>
  <CF3_OLD

```

```

pChiefComplaint=""
pBriefHistory=""
pCourseWard=""
pPertinentFindings="">
<PHEX
  pBP=""
  pCR=""
  pRR=""
  pTemp=""
  pHEENT=""
  pChestLungs=""
  pCVS=""
  pAbdomen=""
  pGUIE=""
  pSkinExtremities=""
  pNeuroExam=""/>
<MATERNITY>
  <PRENATAL
    pPrenatalConsultation="08-26-2009"
    pMCPOrientation="Y"
    pExpectedDeliveryDate="08-26-2009">
  <CLINICALHIST
    pVitalSigns="Y"
    pPregnancyLowRisk="Y"
    pLMP="08-26-2009"
    pMenarcheAge="21"
    pObstetricG=""
    pObstetricP=""
    pObstetric_T=""
    pObstetric_P=""
    pObstetric_A=""
    pObstetric_L=""/>
  <OBSTETRIC
    pMultiplePregnancy="N"
    pOvarianCyst="N"
    pMyomaUteri="N"
    pPlacentaPrevia="N"
    pMiscarriages="N"
    pStillBirth="N"
    pPreEclampsia="N"
    pEclampsia="N"
    pPrematureContraction="N"/>
  <MEDISURG
    pHypertension="N"
    pHeartDisease="N"
    pDiabetes="N"
    pThyroidDisaster="N"
    pObesity="N"
    pAsthma="N"
    pEpilepsy="N"
    pRenalDisease="N"
    pBleedingDisorders="N"
    pPreviousCS="N"
    pUterineMyomectomy="Y"/>
  <CONSULTATION
    pVisitDate="08-20-2009"
    pAOGWeeks=""
    pWeight="50"
    pCardiacRate=""
    pRespiratoryRate="26"
    pBloodPressure="160/100"
    pTemperature="38.5 C"/>
  <CONSULTATION
    pVisitDate="08-25-2009"
    pAOGWeeks=""

```



```

        pWeight="60"
        pCardiacRate=""
        pRespiratoryRate="26"
        pBloodPressure="160/100"
        pTemperature="38.5 C"/>
    <CONSULTATION
        pVisitDate="08-30-2009"
        pAOGWeeks=""
        pWeight="65"
        pCardiacRate=""
        pRespiratoryRate="26"
        pBloodPressure="160/100"
        pTemperature="38.5 C"/>
</PRENATAL>
<DELIVERY
    pDeliveryDate="09-01-2009"
    pDeliveryTime="12:00AM"
    pObstetricIndex=""
    pAOGLMP=""
    pDeliveryManner=""
    pPresentation=""
    pFetalOutcome=""
    pSex="M"
    pBirthWeight="5000"
    pAPGARScore=""
    pPostpartum=""/>
<POSTPARTUM
    pPerinealWoundCare="Y"
    pPerinealRemarks=""
    pMaternalComplications="Y"
    pMaternalRemarks=""
    pBreastFeeding="Y"
    pBreastFeedingRemarks=""
    pFamilyPlanning="Y"
    pFamilyPlanningRemarks=""
    pPlanningService="Y"
    pPlanningServiceRemarks=""
    pSurgicalSterilization="Y"
    pSterilizationRemarks=""
    pFollowupSchedule="Y"
    pFollowupScheduleRemarks=""/>
</MATERNITY>
</CF3_OLD>
<!--CF3_NEW>
    <ADMITREASON
        pBriefHistory=""
        pReferredReason=""
        pIntensive="N"
        pMaintenance="N">
        <CLINICAL pCriteria="COUGH"/>
        <CLINICAL pCriteria="COLDS"/>
        <CLINICAL pCriteria="FEVER"/>
        <CLINICAL pCriteria="RR=26"/>
        <CLINICAL pCriteria="T= 38.5 C"/>
        <CLINICAL pCriteria="BP = 160/100"/>
        <LABDIAG pCriteria="CHEST X-RAY- PNEUMONIA"/>
        <LABDIAG pCriteria="CBC-INCREASE WBC"/>
        <LABDIAG pCriteria="URINALYSIS"/>
    <PHEX
        pBP=""
        pCR=""
        pRR=""
        pTemp=""
        pHEENT=""
        pChestLungs=""

```

```

        pCVS=""
        pAbdomen=""
        pGUIE=""
        pSkinExtremities=""
        pNeuroExam=""/>
    </ADMITREASON>
    <COURSE>
        <WARD
            pCourseDate=""
            pFindings="CHEST X-RAY"
            pAction="FOR ADMISSION"/>
        <WARD
            pCourseDate=""
            pFindings="URINALYSIS"
            pAction="START PENICILLIN IV EVERY 6 HRS."/>
    </COURSE>
</CF3_NEW-->
</CF3>
<PARTICULARS>
    <DRGMED
        pPurchaseDate="08-26-2009"
        pDrugCode="X0001234"
        pPNDFCode=""
        pGenericName="PARACETAMOL"
        pBrandName="GAYAGESIC"
        pPreparation="TABLET 250MG"
        pQuantity="3"/>
    <DRGMED
        pPurchaseDate="08-26-2009"
        pDrugCode="X0001235"
        pPNDFCode=""
        pGenericName="PARACETAMOL"
        pBrandName="GAYAGESIC"
        pPreparation="TABLET 250MG"
        pQuantity="3"/>
    <XLSO
        pDiagnosticDate="08-26-2009"
        pDiagnosticType="IMAGING"
        pDiagnosticName="ULTRASOUND"
        pQuantity="2"/>
    <XLSO
        pDiagnosticDate="08-26-2009"
        pDiagnosticType="IMAGING"
        pDiagnosticName="ULTRASOUND"
        pQuantity="2"/>
</PARTICULARS>
<RECEIPTS>
    <RECEIPT
        pCompanyName="COMPANY"
        pCompanyTIN="123-456-789"
        pBIRPermitNumber="12345"
        pReceiptNumber="00001"
        pReceiptDate="08-25-2009"
        pVATExemptSale="0.00"
        pVAT="12.00"
        pTotal="100.00">
        <ITEM
            pQuantity="10"
            pUnitPrice="5"
            pDescription="BIOGESIC"
            pAmount="50.00"/>
        <ITEM
            pQuantity="5"
            pUnitPrice="10"
            pDescription="NEOZEP"

```

```

                pAmount="50.00"/>
            </RECEIPT>
        </RECEIPTS>
        <DOCUMENTS>
            <DOCUMENT
                pDocumentType="CSF"
                pDocumentURL=
                "https://hospitalwebserver/eclaims/claimnumber/yyyymmdd000001.pdf"/>
            <DOCUMENT
                pDocumentType="OPR"
                pDocumentURL=
                "https://hospitalwebserver/eclaims/claimnumber/yyyymmdd000002.pdf"/>
            <DOCUMENT
                pDocumentType="SOA"
                pDocumentURL=
                "https://hospitalwebserver/eclaims/claimnumber/yyyymmdd000003.pdf"/>
        </DOCUMENTS>
    </CLAIM>
</eTRANSMITTAL>
</eCLAIMS>

```

Sample encrypted XML:

```

{
  "docMimeType": "text/xml",
  "hash": " fac928f71a3841317e7a99aec dc8f4d74d977dfe701c0c9bbca0678300540591",
  "key1": "",
  "key2": "",
  "iv": "nhYeJqBkuhJNM y1jPMxvQE2aJPVnqqn1pDQs==",
  "doc": "aJPVnqqn1pDQ Pms1FWFZT+odAp0qf2zManmroSUr3lYgDFnhYeJqBkuhJNMJU5geEN=="
}

```

Sample encrypted output value:

```

{
  "result": {
    "docMimeType": "text/xml",
    "hash": "dc8f4d74d977dfe701c0c9bbca0678300540591fac928f71a3841317e7a99aec",
    "key1": "",
    "key2": "",
    "iv": "y1jPMxvQE2aJPVnqqn1pDQ==",
    "doc": "Pms1FWFZT+odAp0qf2zManmroSUr3lYgDFnhYeJqBkuhJNMJU5geEN=="
  },
  "success": true,
  "message": ""
}

```

A. Successfully Received

```

<eRECEIPT
  pUserName=""
  pUserPassword=""
  pHospitalCode="123456"
  pHospitalTransmittalNo="001"
  pTotalClaims="1"
  pTransmissionControlNumber="1234-5601-1234-1253"
  pTransmissionDate="08-26-2009"
  pTransmissionTime="00:00:00AM"
  pReceiptTicketNumber="1234-5601-1234"

```

</eRECEIPT>

B. Unsuccessfully Received

<eRECEIPT

pUserName=""

pUserPassword=""

pHospitalCode="123456"

pHospitalTransmittalNo="001"

pTotalClaims="1"

pTransmissionControlNumber=""

pTransmissionDate="08-26-2009"

pTransmissionTime="00:00:00AM" >

<REMARKS pErrCode="T01" pErrDescription="Invalid parameter value: pAmtActual" />

<REMARKS pErrCode="T02" pErrDescription="Invalid parameter value: pOperationDate" />

</eRECEIPT>

5. Search Case Rates Method

This method allows the caller to retrieve case rate amounts for target benefit packages, including different applicable periods.

Endpoint

`https://{pecws.domain}/PHIC/Claims3.0/searchCaseRates`

Method

POST

Inputs

Header

Key	Value
token	PECWS authentication token

Body

The JSON object contains the following key-value pairs:

Name	Values
icdcode	The ICD code of the target benefit package
rvscode	The RVS code of the target benefit package
description	Substring or full text of the name or description of the target benefit packages
targetdate	Format: mm-dd-yyyy. If this parameter has a valid date value (which, normally the admission date of a target claim), this method will return only the record for the applicable period that covers the target date. If an empty string is passed as the value is passed of this parameter, all the applicable periods of the target benefit packages will be returned

Input Payload

{

```

"icdcode": "",
"rvscode": "",
"description": "DENGUE",
"targetdate": "02-14-2024"
}

```

Output

This method returns a JSON object with the key-value pairs described below:

Key	Value/Remarks														
success	A value of 'true' indicates a successful operation														
message	If an error was encountered during the execution of this method, this will contain the error message														
result	<p>The JSON object of the result of the encryption of the XML text containing the records of the matching benefit packages using the cipher key of the health facility. The JSON object contains the following key value pairs:</p> <table> <tr> <th>Key</th><th>Value/Remarks</th></tr> <tr> <td>docMimeType</td><td>"text/xml"</td></tr> <tr> <td>hash</td><td>The computed hash value of the unencrypted text.</td></tr> <tr> <td>key1</td><td>"" (Empty string)</td></tr> <tr> <td>key2</td><td>"" (Empty string)</td></tr> <tr> <td>iv</td><td>The initialization vector used in the AES encryption</td></tr> <tr> <td>doc</td><td>The encrypted text of the JSON object containing the records of the matching benefit packages. Sample XML text and the DTD of the XML text is shown below.</td></tr> </table>	Key	Value/Remarks	docMimeType	"text/xml"	hash	The computed hash value of the unencrypted text.	key1	"" (Empty string)	key2	"" (Empty string)	iv	The initialization vector used in the AES encryption	doc	The encrypted text of the JSON object containing the records of the matching benefit packages. Sample XML text and the DTD of the XML text is shown below.
Key	Value/Remarks														
docMimeType	"text/xml"														
hash	The computed hash value of the unencrypted text.														
key1	"" (Empty string)														
key2	"" (Empty string)														
iv	The initialization vector used in the AES encryption														
doc	The encrypted text of the JSON object containing the records of the matching benefit packages. Sample XML text and the DTD of the XML text is shown below.														

Sample Encrypted Output Payload

```

{
  "success": true,
  "message": "",
  "result": {
    "docMimeType": "text/xml",
    "hash": "dc8f4d74d977dfe701c0c9bbca0678300540591fac928f71a3841317e7c",
    "key1": "",
    "key2": "",
    "iv": "y1jPMxvQE2aJPVnqqn1pDQ==",
    "doc": "PMs1FWFZT+odArTp0qf2zMmroSUR31YgDFnhYeJqBkuhJNMJU5geEN=="
  }
}

```

Sample Decrypted XML Text of the Output Payload

```
{ "eCASERATES": {  
  "CASERATES": [  
    { "AMOUNT": {  
      "pPrimaryProfFee": "350.0",  
      "pPrimaryHCIFee": "2250.0",  
      "pCheckFacilityMAT": "F",  
      "pPrimaryCaseRate": "2600.0",  
      "pSecondaryHCIFee": "0.00",  
      "pSecondaryProfFee": "0.00",  
      "pSecondaryCaseRate": "0.00",  
      "pCheckFacilityH1": "F",  
      "pCheckFacilityH2": "F",  
      "pCheckFacilityFSDC": "F",  
      "pCheckFacilityH3": "F",  
      "pCheckFacilityASC": "F",  
      "pCheckFacilityPCF": "T",  
      "pCheckFacilityTSEKAP": "F",  
      "pCheckFacilityABTC": "F",  
      "pCheckFacilityTBDOTSC": "F",  
      "pCheckFacilityOPMC": "F",  
      "pCheckFacilityRHU": "F",  
      "pCheckFacilityDATRC": "F",  
      "pCheckFacilityHIVTH": "F",  
      "pCheckFacilityFPC": "F",  
      "pCheckFacilityCIU": "F",  
      "pCheckFacilityDSP": "F",  
      "pCheckFacilityPCB": "F",  
    },  
    "pCaseRateCode": "CR0389",  
    "pCaseRateDescription": "DIALYSIS PROCEDURE OTHER THAN HEMODIALYSIS (E.G.  
PERITONEAL, HEMOFILTRATION)",  
    "pItemCode": "90945",  
    "pItemDescription": "DIALYSIS PROCEDURE OTHER THAN HEMODIALYSIS (E.G. PERITONEAL,  
HEMOFILTRATION)",  
    "pEffectivityDate": "09-15-2015",  
    "pEffectivityEndDate": "02-13-2024"  
  },  
  { "AMOUNT": {  
    "pPrimaryProfFee": "500.0",  
    "pPrimaryHCIFee": "3500.0",  
    "pCheckFacilityMAT": "F",  
    "pPrimaryCaseRate": "4000.0",  
    "pSecondaryHCIFee": "0.00",  
    "pSecondaryProfFee": "0.00",  
    "pSecondaryCaseRate": "0.00",  
    "pCheckFacilityH1": "F",  
    "pCheckFacilityH2": "F",  
    "pCheckFacilityFSDC": "F",  
    "pCheckFacilityH3": "F",  
    "pCheckFacilityASC": "F",  
    "pCheckFacilityPCF": "T",  
    "pCheckFacilityTSEKAP": "F",  
    "pCheckFacilityABTC": "F",  
    "pCheckFacilityTBDOTSC": "F",  
    "pCheckFacilityOPMC": "F",  
    "pCheckFacilityRHU": "F",  
    "pCheckFacilityDATRC": "F",  
    "pCheckFacilityHIVTH": "F",  
    "pCheckFacilityFPC": "F",  
    "pCheckFacilityCIU": "F",  
    "pCheckFacilityDSP": "F",  
  }  
}
```

```

        "pCheckFacilityPCB": "F"
    },
    "pCaseRateCode": "CR0389",
    "pCaseRateDescription": "DIALYSIS PROCEDURE OTHER THAN HEMODIALYSIS (E.G.
PERITONEAL, HEMOFILTRATION)",
    "pItemCode": "90945",
    "pItemDescription": "DIALYSIS PROCEDURE OTHER THAN HEMODIALYSIS (E.G. PERITONEAL,
HEMOFILTRATION)",
    "pEffectivityDate": "12-02-2014",
    "pEffectivityEndDate": "09-14-2015"
},
{ "AMOUNT": {
    "pPrimaryProfFee": "500.0",
    "pPrimaryHCIFee": "3500.0",
    "pCheckFacilityMAT": "F",
    "pPrimaryCaseRate": "4000.0",
    "pSecondaryHCIFee": "0.00",
    "pSecondaryProfFee": "0.00",
    "pSecondaryCaseRate": "0.00",
    "pCheckFacilityH1": "F",
    "pCheckFacilityH2": "F",
    "pCheckFacilityFSDC": "F",
    "pCheckFacilityH3": "F",
    "pCheckFacilityASC": "F",
    "pCheckFacilityPCF": "T",
    "pCheckFacilityTSEKAP": "F",
    "pCheckFacilityABTC": "F",
    "pCheckFacilityTBDOTSC": "F",
    "pCheckFacilityOPMC": "F",
    "pCheckFacilityRHU": "F",
    "pCheckFacilityDATRC": "F",
    "pCheckFacilityHIVTH": "F",
    "pCheckFacilityFPC": "F",
    "pCheckFacilityCIU": "F",
    "pCheckFacilityDSP": "F",
    "pCheckFacilityPCB": "F"
},
    "pCaseRateCode": "CR0389",
    "pCaseRateDescription": "DIALYSIS PROCEDURE OTHER THAN HEMODIALYSIS (E.G.
PERITONEAL, HEMOFILTRATION)",
    "pItemCode": "90945",
    "pItemDescription": "DIALYSIS PROCEDURE OTHER THAN HEMODIALYSIS (E.G. PERITONEAL,
HEMOFILTRATION)",
    "pEffectivityDate": "01-01-2014",
    "pEffectivityEndDate": "12-01-2014"
},
{ "AMOUNT": {
    "pPrimaryProfFee": "455.0",
    "pPrimaryHCIFee": "2925.0",
    "pCheckFacilityMAT": "F",
    "pPrimaryCaseRate": "3380.0",
    "pSecondaryHCIFee": "0.00",
    "pSecondaryProfFee": "0.00",
    "pSecondaryCaseRate": "0.00",
    "pCheckFacilityH1": "F",
    "pCheckFacilityH2": "F",
    "pCheckFacilityFSDC": "F",
    "pCheckFacilityH3": "F",
    "pCheckFacilityASC": "F",
    "pCheckFacilityPCF": "T",
    "pCheckFacilityTSEKAP": "F",
    "pCheckFacilityABTC": "F",
    "pCheckFacilityTBDOTSC": "F",

```



```

        "pCheckFacilityOPMC": "F",
        "pCheckFacilityRHU": "F",
        "pCheckFacilityDATRC": "F",
        "pCheckFacilityHIVTH": "F",
        "pCheckFacilityFPC": "F",
        "pCheckFacilityCIU": "F",
        "pCheckFacilityDSP": "F",
        "pCheckFacilityPCB": "F"
    },
    "pCaseRateCode": "CR0389",
    "pCaseRateDescription": "DIALYSIS PROCEDURE OTHER THAN HEMODIALYSIS (E.G.
PERITONEAL, HEMOFILTRATION)",
    "pItemCode": "90945",
    "pItemDescription": "DIALYSIS PROCEDURE OTHER THAN HEMODIALYSIS (E.G. PERITONEAL,
HEMOFILTRATION)",
    "pEffectivityDate": "02-14-2024",
    "pEffectivityEndDate": "12-31-9999"
}
]
}

```

6. Add Required Document Method

This method allows the addition of required documents in compliance with RTH claims.

Endpoint

`https://{pecws.domain}/PHIC/Claims3.0/addRequiredDocument`

Method

POST

Header

Key	Value
token	PECWS authentication token

Body

JSON object containing the following key-value pairs

Key	Value/Remarks																				
pSeriesLhioNo	Series Lhio Number																				
pXML	<div>JSON Object list of documents to be added to the claim submitted<table><tr><th>Key</th><th>Value/Remarks</th></tr><tr><td>docMimeType</td><td>"text/xml"</td></tr><tr><td>hash</td><td>The computed hash value of the unencrypted text.</td></tr><tr><td>key1</td><td>"" (Empty string)</td></tr><tr><td>key2</td><td>"" (Empty string)</td></tr><tr><td>iv</td><td>The initialization vector used in the AES encryption</td></tr><tr><td>doc</td><td>The encrypted text of the XML object containing the document type and document URL. Sample XML text is shown below.</td></tr></table><div>XML Document attributes<table><tr><th>Key</th><th>Value/Remarks</th></tr><tr><td>pDocumentType</td><td>Document Type</td></tr><tr><td>pDocumentURL</td><td>Document URL</td></tr></table></div></div>	Key	Value/Remarks	docMimeType	"text/xml"	hash	The computed hash value of the unencrypted text.	key1	"" (Empty string)	key2	"" (Empty string)	iv	The initialization vector used in the AES encryption	doc	The encrypted text of the XML object containing the document type and document URL. Sample XML text is shown below.	Key	Value/Remarks	pDocumentType	Document Type	pDocumentURL	Document URL
Key	Value/Remarks																				
docMimeType	"text/xml"																				
hash	The computed hash value of the unencrypted text.																				
key1	"" (Empty string)																				
key2	"" (Empty string)																				
iv	The initialization vector used in the AES encryption																				
doc	The encrypted text of the XML object containing the document type and document URL. Sample XML text is shown below.																				
Key	Value/Remarks																				
pDocumentType	Document Type																				
pDocumentURL	Document URL																				

Sample JSON:

```
{"pSeriesLhioNo": "",
```

```

"pXML": {
    "docMimeType": "",
    "hash": "",
    "key1": "",
    "key2": "",
    "iv": "",
    "doc": ""
}

```

Sample XML to be encrypted as pXML value:

```

<DOCUMENTS>
  <DOCUMENT pDocumentType="CF1" pDocumentURL="http://sample/file/other/cf1.pdf"/>
  <DOCUMENT pDocumentType="CF2" pDocumentURL="http://sample/file/other/cf2.pdf"/>
  <DOCUMENT pDocumentType="OPR" pDocumentURL="http://sample/file/other/opr.pdf"/>
</DOCUMENTS>

```

Output

This method returns a JSON object with the key-value pairs described below:

Key	Value/Remarks
success	A value of 'true' indicates a successful operation
message	If an error was encountered during the execution of this method, this will contain the error message
result	Possible results: "Claims has already been paid" "Claims has already been denied"

7. EClaims File Check Method

This method allows the caller to validate the eClaims XML File.

Endpoint

`https://{pecws.domain}/PHIC/Claims3.0/eClaimsFileCheck`

Method

POST

Header

Key	Value
token	PECWS authentication token

Body

The body is the eclaims XML in JSON format as encrypted using the cipher key of the Health Facility. The JSON object contains the following key-value pairs:

Key	Value/Remarks
docMimeType	"text/xml"
hash	The computed hash value of the unencrypted text.
key1	"" (Empty string)
key2	"" (Empty string)
iv	The initialization vector used in the AES encryption
doc	The encrypted e-claim XML text.

Output

This method returns a JSON object with the key-value pairs described below:

Key	Value/Remarks
success	A value of 'true' indicates a successful operation
message	If an error was encountered during the execution of this method, this will contain the error message
result	The JSON object as the result of the encryption (using the cipher key of the health facility) of the XML text containing the Receipt Ticket Number and

	other data about the processing of the submitted e-claim data.														
	<table> <tr> <th>Key</th><th>Value/Remarks</th></tr> <tr> <td>docMimeType</td><td>"text/xml"</td></tr> <tr> <td>hash</td><td>The computed hash value of the unencrypted text.</td></tr> <tr> <td>key1</td><td>"" (Empty string)</td></tr> <tr> <td>key2</td><td>"" (Empty string)</td></tr> <tr> <td>iv</td><td>The initialization vector used in the AES encryption</td></tr> <tr> <td>doc</td><td>The encrypted text of the JSON object containing the records of the matching benefit packages. Sample XML text and the DTD of the XML text is shown below.</td></tr> </table>	Key	Value/Remarks	docMimeType	"text/xml"	hash	The computed hash value of the unencrypted text.	key1	"" (Empty string)	key2	"" (Empty string)	iv	The initialization vector used in the AES encryption	doc	The encrypted text of the JSON object containing the records of the matching benefit packages. Sample XML text and the DTD of the XML text is shown below.
Key	Value/Remarks														
docMimeType	"text/xml"														
hash	The computed hash value of the unencrypted text.														
key1	"" (Empty string)														
key2	"" (Empty string)														
iv	The initialization vector used in the AES encryption														
doc	The encrypted text of the JSON object containing the records of the matching benefit packages. Sample XML text and the DTD of the XML text is shown below.														

Sample encrypted output value:

```
{
  "result": {
    "docMimeType": "text/xml",
    "hash": "dc8f4d74d977dfe701c0c9bbca0678300540591fac928f71a3841317e7a99aec",
    "key1": "",
    "key2": "",
    "iv": "y1jPMxvQE2aJPVnqqn1pDQ==",
    "doc": "PMs1FWFZT+odAp0qf2zManmroSUr3lYgDFnhYeJqBkuhJNMJU5geEN=="
  },
  "success": true,
  "message": ""
}
```

8. Get Claim Status Method

This method allows the retrieval of the processing status of one or more claims using claim series numbers.

Endpoint

`https://{pecws.domain}/PHIC/Claims3.0/getClaimStatus`

Method

GET

Header

Key	Value
token	PECWS authentication token

Body

The JSON object contains the following key-value pairs:

Key	Value
serieslhionos	JSON Object list of Strings containing multiple claims Series Nos

Output

This method returns a JSON object with the key-value pairs described below:

Key	Value/Remarks												
success	A value of 'true' indicates a successful operation												
message	If an error was encountered during the execution of this method, this will contain the error message												
result	<div>The JSON object as the result of the encryption (using the cipher key of the health facility) of the XML text containing the Receipt Ticket Number and other data about the processing of the submitted e-claim data.<table><tr><th>Key</th><th>Value/Remarks</th></tr><tr><td>docMimeType</td><td>"text/xml"</td></tr><tr><td>hash</td><td>The computed hash value of the unencrypted text.</td></tr><tr><td>key1</td><td>"" (Empty string)</td></tr><tr><td>key2</td><td>"" (Empty string)</td></tr><tr><td>iv</td><td>The initialization vector used in the AES encryption</td></tr></table></div>	Key	Value/Remarks	docMimeType	"text/xml"	hash	The computed hash value of the unencrypted text.	key1	"" (Empty string)	key2	"" (Empty string)	iv	The initialization vector used in the AES encryption
Key	Value/Remarks												
docMimeType	"text/xml"												
hash	The computed hash value of the unencrypted text.												
key1	"" (Empty string)												
key2	"" (Empty string)												
iv	The initialization vector used in the AES encryption												

	<div> <div>doc</div> <div>The encrypted text of the JSON object containing the records of the matching benefit packages. Sample XML text and the DTD of the XML text is shown below.</div> </div>
--	---

Sample JSON input payload:

```
{ "serieslhionos": ["1","2"] }
```

Sample encrypted output value:

```
{
  "result": {
    "docMimeType": "text/xml",
    "hash": "dc8f4d74d977dfe701c0c9bbca0678300540591fac928f71a3841317e7a99aec",
    "key1": "",
    "key2": "",
    "iv": "y1jPMxvQE2aJPVnqqn1pDQ==",
    "doc": "PMs1FWFZT+odAp0qf2zManmroSUr3lYgDFnhYeJqBkuhJNMJU5geEN=="
  },
  "success": true,
  "message": ""
}
```

Sample decrypted result value:

```
{
  "CLAIMS": [
    {
      "STATUS": {
        "CLAIM": {
          "TRAIL": {
            "PROCESS": [
              {
                "pProcessStage": "VALIDATION",
                "pProcessDate": "07-25-2012"
              },
              {
                "pProcessStage": "EDITING",
                "pProcessDate": "07-25-2012"
              },
              {
                "pProcessStage": "VALIDATION",
                "pProcessDate": "07-23-2012"
              },
              {
                "pProcessStage": "EDITING (RECEIVING)",
                "pProcessDate": "07-23-2012"
              },
              {
                "pProcessStage": "ENCODING",
                "pProcessDate": "07-23-2012"
              },
              {
                "pProcessStage": "RECEIVING",
                "pProcessDate": "07-23-2012"
              }
            ]
          }
        }
      },
      "pClaimSeriesLhio": "120723190000119",
      "pPin": "190892937993",
    }
  ]
}
```

```

        "pPatientLastName": "LASTNAME",
        "pPatientFirstName": "FIRSTNAME",
        "pPatientMiddleName": "MIDDLENAME",
        "pPatientSuffix": "",
        "pAdmissionDate": "05-02-2012",
        "pDischargeDate": "05-06-2012",
        "pClaimDateReceived": "05-15-2012",
        "pClaimDateRefile": "",
        "pStatus": "IN PROCESS"
    },
    "pAsOf": "07-25-2012",
    "pAsOfTime": "04:46:23PM"
}
},
{
    "STATUS": {
        "CLAIM": {
            "TRAIL": {
                "PROCESS": [
                    {
                        "pProcessStage": "VALIDATION",
                        "pProcessDate": "07-25-2012"
                    },
                    {
                        "pProcessStage": "EDITING",
                        "pProcessDate": "07-25-2012"
                    },
                    {
                        "pProcessStage": "VALIDATION",
                        "pProcessDate": "07-23-2012"
                    },
                    {
                        "pProcessStage": "EDITING (RECEIVING)",
                        "pProcessDate": "07-23-2012"
                    },
                    {
                        "pProcessStage": "ENCODING",
                        "pProcessDate": "07-23-2012"
                    },
                    {
                        "pProcessStage": "RECEIVING",
                        "pProcessDate": "07-23-2012"
                    }
                ]
            },
            "pClaimSeriesLhio": "120723190000119",
            "pPin": "190592937994",
            "pPatientLastName": "LASTNAME",
            "pPatientFirstName": "FIRSTNAME",
            "pPatientMiddleName": "MIDDLENAME",
            "pPatientSuffix": "III",
            "pAdmissionDate": "05-02-2012",
            "pDischargeDate": "05-06-2012",
            "pClaimDateReceived": "05-15-2012",
            "pClaimDateRefile": "",
            "pStatus": "IN PROCESS"
        },
        "pAsOf": "07-25-2012",
        "pAsOfTime": "04:46:23PM"
    }
}
]
}

```


9. Get Doctor PAN Method

This method retrieves the PhilHealth Accreditation Number (PAN) of a health care professional (HCP) using the professional's complete name and date of birth.

Endpoint

`https://{pecws.domain}/PHIC/Claims3.0/getDoctorPAN`

Method

POST

Header

Key	Value
token	PECWS authentication token

Body

The body is in JSON format as encrypted using the cipher key of the Health Facility. The JSON object contains the following key-value pairs:

Key	Value
lastname	Health Care Professional Last Name
firstname	Health Care Professional First Name
middlename	Health Care Professional Middle Name
suffix	Health Care Professional Suffix Name
birthdate	Health Care Professional Birth Date (Date Format should be : 'MM-DD-YYYY')

Output

This method returns a JSON object with the key-value pairs described below:

Key	Value/Remarks				
success	A value of 'true' indicates a successful operation				
message	If an error was encountered during the execution of this method, this will contain the error message				
result	The JSON object as the result of the encryption (using the cipher key of the health facility) of the JSON object containing the PAN of the submitted encrypted JSON input payload. <table><tr><th>Key</th><th>Value/Remarks</th></tr><tr><td>docMimeType</td><td>"text/xml"</td></tr></table>	Key	Value/Remarks	docMimeType	"text/xml"
Key	Value/Remarks				
docMimeType	"text/xml"				

	hash	The computed hash value of the unencrypted text.
	key1	"" (Empty string)
	key2	"" (Empty string)
	iv	The initialization vector used in the AES encryption
	doc	The encrypted text of the JSON object containing the PAN and the decrypted result is shown below.

Sample encrypted JSON input payload:

```
{
  "docMimeType": "text/xml",
  "hash": "dc8f4d74d977dfe701c0c9bbca067830054.....",
  "key1": "",
  "key2": "",
  "iv": "y1jPMxvQE2aJPV.....",
  "doc": "PMs1FWFZT+odAp0qf2zManmroSUr3lYgDF..... .=="
}
```

Sample decrypted JSON input payload:

```
{ "lastname": "LASTNAME",
  "firstname": "FIRSTNAME",
  "middlename": "MIDDLENAME",
  "suffix": "III",
  "birthdate": "01-01-1990" }
```

Sample encrypted JSON output payload:

```
{
  "message": "",
  "result": {
    "doc": "TxUkrmwTSKWuP.....",
    "docMimeType": "text/xml",
    "hash": "dac3ccd4a4b726a3305c17c4df5efdda2a9273724ac.....",
    "iv": "I/yCJcJG3ZAe2NNmR.....",
    "key1": "",
    "key2": ""
  },
  "success": true
}
```

Sample decrypted result JSON output payload:

```
{ "pan": "0000-00000000-0" }
```

10. Get Member PIN Method

This method retrieves the PhilHealth Identification Number (PIN) of a PhilHealth member using the member's complete name and date of birth.

Endpoint

`https://{pecws.domain}/PHIC/Claims3.0/getMemberPIN`

Method

POST

Header

Key	Value
token	PECWS authentication token

Body

The body is in JSON format as encrypted using the cipher key of the Health Facility. The JSON object contains the following key-value pairs:

Key	Value
lastname	Member Last Name
firstname	Member First Name
middlename	Member Middle Name
suffix	Member Suffix Name
birthdate	Member Birth Date (Date Format should be : 'MM-DD-YYYY')

Output

This method returns a JSON object with the key-value pairs described below:

Key	Value/Remarks		
success	A value of 'true' indicates a successful operation		
message	If an error was encountered during the execution of this method, this will contain the error message		
result	The JSON object as the result of the encryption (using the cipher key of the health facility) of the JSON object containing the PAN of the submitted encrypted JSON input payload. <table><tr><th>Key</th><th>Value/Remarks</th></tr></table>	Key	Value/Remarks
Key	Value/Remarks		

	docMimeType	"text/xml"
	hash	The computed hash value of the unencrypted text.
	key1	"" (Empty string)
	key2	"" (Empty string)
	iv	The initialization vector used in the AES encryption
	doc	The encrypted text of the JSON object containing the PIN and the decrypted result is shown below.

Sample encrypted JSON input payload:

```
{
  "docMimeType": "text/xml",
  "hash": "dc8f4d74d977dfe701c0c9bbca0678300540591...",
  "key1": "",
  "key2": "",
  "iv": "y1jPMxvQE2aJPV.....",
  "doc": "Pms1FWFZT+odAp0qf2zManmroSUR3lYgDFnhYeJ....."
}
```

Sample decrypted JSON input payload:

```
{
  "lastname": "LASTNAME",
  "firstname": "FIRSTNAME",
  "middlename": "MIDDLENAME",
  "suffix": "III",
  "birthdate": "01-01-1990"
}
```

Sample encrypted JSON output payload:

```
{
  "message": "",
  "result": {
    "doc": "TxUkrmwTSKWuP.....",
    "docMimeType": "text/xml",
    "hash": "dac3ccd4a4b726a3305c17c4df5efdda2a9273724ac.....",
    "iv": "I/yCJcJG3ZAe2NNmR.....",
    "key1": "",
    "key2": ""
  },
  "success": true
}
```

Sample decrypted result JSON output payload:

```
{ "pin": "000000000000" }
```

11. Get Uploaded Claims Map Method

This method facilitates claim reconciliation and verification by mapping the health facility's system-generated ID to the PhilHealth Claim Series Number. The Claim Series Number serves as a common reference between the health facility and PhilHealth for a given claim.

The method processes claims included in the **e-claim XML file** submitted through the **eClaimsUpload** method. The **receiptTicketNo** parameter must match the **pReceiptTicketNumber** attribute returned by the eClaimsUpload method for the corresponding set of claims.

Endpoint

`https://{pecws.domain}/PHIC/Claims3.0/getUploadedClaimsMap`

Method

GET

Header

Key	Value
token	PECWS authentication token

Parameter

Key	Value
receiptTicketNumber	The receipt ticket number is generated after calling the eClaimsUpload method

Output

This method returns a JSON object with the key-value pairs described below:

Key	Value/Remarks				
success	A value of 'true' indicates a successful operation				
message	If an error was encountered during the execution of this method, this will contain the error message				
result	The JSON object as the result of the encryption (using the cipher key of the health facility) of the XML text containing the Receipt Ticket Number and other data about the processing of the submitted e-claim data. <table><tr><th>Key</th><th>Value/Remarks</th></tr><tr><td>docMimeType</td><td>"text/xml"</td></tr></table>	Key	Value/Remarks	docMimeType	"text/xml"
Key	Value/Remarks				
docMimeType	"text/xml"				

	hash	The computed hash value of the unencrypted text.
	key1	"" (Empty string)
	key2	"" (Empty string)
	iv	The initialization vector used in the AES encryption
	doc	The encrypted text of the JSON object containing the records of the matching benefit packages. Sample XML text and the DTD of the XML text is shown below.

Sample encrypted result:

```
{
  "result": {
    "docMimeType": "text/xml",
    "hash": "dc8f4d74d977dfe701c0c9bbca0678300540591fac928f71a3841317e7a99aec",
    "key1": "",
    "key2": "",
    "iv": "y1jPMxvQE2aJPVnqqn1pDQ==",
    "doc": "PMs1FWFZT+odAp0qf2zManmroSUR3lYgDFnhYeJqBkuhJNMJU5geEN=="
  },
  "success": true,
  "message": ""
}
```

Sample decrypted result:

```
{
  "eCONFIRMATION": {
    "MAPPING": {
      "pClaimNumber": "09-08-01-006",
      "pPatientLastName": "LASTNAME",
      "pPatientFirstName": "MARIA",
      "pPatientMiddleName": "C",
      "pPatientSuffix": "",
      "pAdmissionDate": "08-25-2009",
      "pDischargeDate": "08-25-2009",
      "pClaimSeriesLhio": "090801990000199"
    },
    "pReceiptTicketNumber": "071311000005",
    "pHospitalCode": "300832",
    "pHospitalTransmittalNo": "3008321107000008",
    "pTotalClaims": "1",
    "pReceivedDate": "09-13-2009"
  }
}
```

12. Get Voucher Details Method

This method facilitates the reconciliation of paid claims by returning the voucher details and other payment information associated with claims.

The **voucherNo** parameter should be set to the **pVoucherNo** attribute returned by the **getClaimStatus** method if a voucher has already been prepared for the claim's payment. Note that a single voucher may contain multiple claims.

Endpoint

`https://{pecws.domain}/PHIC/Claims3.0/getVoucherDetails`

Method

GET

Header

Key	Value
token	PECWS authentication token

Parameter

Key	Value
voucherNo	Valid voucher number

Output

This method returns a JSON object with the key-value pairs described below:

Key	Value/Remarks								
success	A value of 'true' indicates a successful operation								
message	If an error was encountered during the execution of this method, this will contain the error message								
result	<div>The JSON object as the result of the encryption (using the cipher key of the health facility) of the XML text containing the Receipt Ticket Number and other data about the processing of the submitted e-claim data.<table><tr><th>Key</th><th>Value/Remarks</th></tr><tr><td>docMimeType</td><td>"text/xml"</td></tr><tr><td>hash</td><td>The computed hash value of the unencrypted text.</td></tr><tr><td>key1</td><td>"" (Empty string)</td></tr></table></div>	Key	Value/Remarks	docMimeType	"text/xml"	hash	The computed hash value of the unencrypted text.	key1	"" (Empty string)
Key	Value/Remarks								
docMimeType	"text/xml"								
hash	The computed hash value of the unencrypted text.								
key1	"" (Empty string)								

	key2	"" (Empty string)
	iv	The initialization vector used in the AES encryption
	doc	The encrypted text of the JSON object containing the records of the matching benefit packages. Sample XML text and the DTD of the XML text is shown below.

Sample encrypted result:

```
{
  "result": {
    "docMimeType": "text/xml",
    "hash": "dc8f4d74d977dfe701c0c9bbca0678300540591fac928f71a3841317e7a99aec",
    "key1": "",
    "key2": "",
    "iv": "y1jPMxvQE2aJPVnqqn1pDQ==",
    "doc": "PMS1FWFZT+odAp0qf2zManmroSUR3lYgDFnhYeJqBkuhJNMJU5geEN=="
  },
  "success": true,
  "message": ""
}
```

Sample decrypted result:

```
{
  "VOUCHER": {
    "CLAIM": [
      {
        "CHARGE": [
          {
            "pPayeeType": "C",
            "pPayeeCode": "2XX25",
            "pPayeeName": "DR. DJXXXXX XXX X. SXXX",
            "pRMBD": "0.00",
            "pDRUGS": "0.00",
            "pXRAY": "0.00",
            "pOPRM": "0.00",
            "pSPFee": "0.00",
            "pGPFee": "0.00",
            "pSURFee": "4000.00",
            "pANESFee": "0.00",
            "pGrossAmount": "4000.00",
            "pTaxAmount": "0.00",
            "pNetAmount": "4000.00"
          },
          {
            "pPayeeType": "C",
            "pPayeeCode": "3XX25",
            "pPayeeName": "DR. IXX OLXXXX A. CANXXXX",
            "pRMBD": "0.00",
            "pDRUGS": "0.00",
            "pXRAY": "0.00",
            "pOPRM": "0.00",
            "pSPFee": "0.00",

```



```

        "pGPFee": "0.00",
        "pSURFee": "0.00",
        "pANESFee": "1200.00",
        "pGrossAmount": "1200.00",
        "pTaxAmount": "0.00",
        "pNetAmount": "1200.00"
    },
    {
        "pPayeeType": "H",
        "pPayeeCode": "30XX04",
        "pPayeeName": "XXXX CITY XXXXXX HOSPITAL",
        "pRMBD": "800.00",
        "pDRUGS": "507.50",
        "pXRAY": "994.40",
        "pOPRM": "3490.00",
        "pSPFee": "0.00",
        "pGPFee": "0.00",
        "pSURFee": "0.00",
        "pANESFee": "0.00",
        "pGrossAmount": "5791.90",
        "pTaxAmount": "0.00",
        "pNetAmount": "5791.90"
    },
    {
        "pPayeeType": "M",
        "pPayeeCode": "P192003617072",
        "pPayeeName": "ALOXXX , BERNAXXX X",
        "pRMBD": "0.00",
        "pDRUGS": "2544.00",
        "pXRAY": "141.00",
        "pOPRM": "0.00",
        "pSPFee": "0.00",
        "pGPFee": "0.00",
        "pSURFee": "0.00",
        "pANESFee": "0.00",
        "pGrossAmount": "2685.00",
        "pTaxAmount": "0.00",
        "pNetAmount": "2685.00"
    }
],
    "pClaimSeriesLhio": "060516030019903",
    "pPin": "192003610605",
    "pPatientLastName": "ALOXXX",
    "pPatientFirstName": "BERNAXXX",
    "pPatientMiddleName": "B",
    "pPatientSuffix": "",
    "pAdmissionDate": "04-18-2006",
    "pDischargeDate": "04-20-2006",
    "pClaimDateReceived": "05-12-2006",
    "pClaimDateRefile": "",
    "pIsAdjustment": "F"
},
{
    "CHARGE": [
        {
            "pPayeeType": "C",
            "pPayeeCode": "27XX2",
            "pPayeeName": "DR. VIRXXX XX. X. DXXXO",
            "pRMBD": "0.00",
            "pDRUGS": "0.00",
            "pXRAY": "0.00",
            "pOPRM": "0.00",
            "pSPFee": "0.00",
            "pGPFee": "600.00",
            "pSURFee": "0.00",

```

```

        "pANESFee": "0.00",
        "pGrossAmount": "600.00",
        "pTaxAmount": "0.00",
        "pNetAmount": "600.00"
    },
    {
        "pPayeeType": "H",
        "pPayeeCode": "3XXX04",
        "pPayeeName": "XXXX CITY XXXXX HOSPITAL",
        "pRMBD": "2400.00",
        "pDRUGS": "240.00",
        "pXRAY": "671.25",
        "pOPRM": "0.00",
        "pSPFee": "0.00",
        "pGPFee": "0.00",
        "pSURFee": "0.00",
        "pANESFee": "0.00",
        "pGrossAmount": "3311.25",
        "pTaxAmount": "0.00",
        "pNetAmount": "3311.25"
    },
    {
        "pPayeeType": "M",
        "pPayeeCode": "P19200XX31034",
        "pPayeeName": "ANDXXX , CONCHXXX X",
        "pRMBD": "0.00",
        "pDRUGS": "1157.50",
        "pXRAY": "863.00",
        "pOPRM": "0.00",
        "pSPFee": "0.00",
        "pGPFee": "0.00",
        "pSURFee": "0.00",
        "pANESFee": "0.00",
        "pGrossAmount": "2020.50",
        "pTaxAmount": "0.00",
        "pNetAmount": "2020.50"
    }
},
    "pClaimSeriesLhio": "060516030031234",
    "pPin": "192005981034",
    "pPatientLastName": "ANDXXX",
    "pPatientFirstName": "CONCHXXX",
    "pPatientMiddleName": "XXXX",
    "pPatientSuffix": "",
    "pAdmissionDate": "04-21-2006",
    "pDischargeDate": "04-27-2006",
    "pClaimDateReceived": "05-12-2006",
    "pClaimDateRefile": "",
    "pIsAdjustment": "F"
},
    "SUMMARY": {
        "PAYEE": [
            {
                "pPayeeType": "C",
                "pPayeeCode": "30XX04",
                "pPayeeName": "HC- XXXX CITY XXXXX HOSPITAL ",
                "pRMBD": "0.00",
                "pDRUGS": "0.00",
                "pXRAY": "0.00",
                "pOPRM": "0.00",
                "pSPFee": "0.00",
                "pGPFee": "600.00",
                "pSURFee": "4000.00",
                "pANESFee": "1200.00",
            }
        ]
    }
}

```

```

        "pGrossAmount": "5800.00",
        "pTaxAmount": "0.00",
        "pNetAmount": "5800.00",
        "pCheckNo": "0000XXX429",
        "pCheckDate": "06-19-2006"
    },
    {
        "pPayeeType": "H",
        "pPayeeCode": "30XX04",
        "pPayeeName": "XXXX CITY XXXXX HOSPITAL",
        "pRMBD": "3200.00",
        "pDRUGS": "747.50",
        "pXRAY": "1665.65",
        "pOPRM": "3490.00",
        "pSPFee": "0.00",
        "pGPFee": "0.00",
        "pSURFee": "0.00",
        "pANESFee": "0.00",
        "pGrossAmount": "9103.15",
        "pTaxAmount": "0.00",
        "pNetAmount": "9103.15",
        "pCheckNo": "0000XXX430",
        "pCheckDate": "06-19-2006"
    },
    {
        "pPayeeType": "M",
        "pPayeeCode": "P1920XX987072",
        "pPayeeName": "ALOXXX , BERNXXXX X",
        "pRMBD": "0.00",
        "pDRUGS": "2544.00",
        "pXRAY": "141.00",
        "pOPRM": "0.00",
        "pSPFee": "0.00",
        "pGPFee": "0.00",
        "pSURFee": "0.00",
        "pANESFee": "0.00",
        "pGrossAmount": "2685.00",
        "pTaxAmount": "0.00",
        "pNetAmount": "2685.00",
        "pCheckNo": "0000XXX431",
        "pCheckDate": "06-19-2006"
    },
    {
        "pPayeeType": "M",
        "pPayeeCode": "P1920XX731034",
        "pPayeeName": "ANXXXX , CONXXXXX X",
        "pRMBD": "0.00",
        "pDRUGS": "1157.50",
        "pXRAY": "863.00",
        "pOPRM": "0.00",
        "pSPFee": "0.00",
        "pGPFee": "0.00",
        "pSURFee": "0.00",
        "pANESFee": "0.00",
        "pGrossAmount": "2020.50",
        "pTaxAmount": "0.00",
        "pNetAmount": "2020.50",
        "pCheckNo": "0000XXX432",
        "pCheckDate": "06-19-2006"
    }
},
    "pTotalAmount": "19608.65",
    "pNumberOfClaims": "2"
},
    "pVoucherNo": "201-062001-06I03",

```

```
}
    "pVoucherDate": "06-14-2006"
}
```

13. Is Doctor Accredited Method

This method verifies whether a Health Care Professional (HCP) is accredited based on the Accreditation Number, Admission Date, and Discharge Date.

Endpoint

`https://{pecws.domain}/PHIC/Claims3.0/isDoctorAccredited`

Method

POST

Header

Key	Value
token	PECWS authentication token

JSON Payload

JSON object containing the following key-value pairs

Key	Value
accrecode	Health Care Professional Accreditation Number
admissiondate	Date of Admission (Date Format should be : 'MM-DD-YYYY')
dischargedate	Date of Discharge (Date Format should be : 'MM-DD-YYYY')

Output

This method returns a JSON object with the key-value pairs described below:

Key	Value/Remarks		
success	A value of 'true' indicates a successful operation		
message	If an error was encountered during the execution of this method, this will contain the error message		
result	The JSON object as the result of the encryption (using the cipher key of the health facility) of the XML text containing the essential information for verifying accreditation status and validity. <table><tr><th>Key</th><th>Value/Remarks</th></tr></table>	Key	Value/Remarks
Key	Value/Remarks		

	<table> <tr> <td>isaccredited</td><td></td></tr> <tr> <td>accrecode</td><td>Health Care Professional Accreditation Number</td></tr> <tr> <td>admissiondate</td><td>Date of Admission (Date Format : 'MM-DD-YYYY')</td></tr> <tr> <td>dischargedate</td><td>Date of Discharge (Date Format : 'MM-DD-YYYY')</td></tr> <tr> <td>accreditation start</td><td>Date of Start of Accreditation (Date Format : 'MM-DD-YYYY')</td></tr> <tr> <td>accreditation end</td><td>Date of End of Accreditation (Date Format : 'MM-DD-YYYY')</td></tr> </table>	isaccredited		accrecode	Health Care Professional Accreditation Number	admissiondate	Date of Admission (Date Format : 'MM-DD-YYYY')	dischargedate	Date of Discharge (Date Format : 'MM-DD-YYYY')	accreditation start	Date of Start of Accreditation (Date Format : 'MM-DD-YYYY')	accreditation end	Date of End of Accreditation (Date Format : 'MM-DD-YYYY')
isaccredited													
accrecode	Health Care Professional Accreditation Number												
admissiondate	Date of Admission (Date Format : 'MM-DD-YYYY')												
dischargedate	Date of Discharge (Date Format : 'MM-DD-YYYY')												
accreditation start	Date of Start of Accreditation (Date Format : 'MM-DD-YYYY')												
accreditation end	Date of End of Accreditation (Date Format : 'MM-DD-YYYY')												

Sample JSON input payload:

```
{
  "accrecode": "1234567890",
  "admissiondate": "01-01-1990",
  "dischargedate": "01-01-1990"
}
```

14. Search Employer Method

This method allows searching for an employed member's employer using the PhilHealth Employer Number (PEN) and/or employer name.

Endpoint

`https://{pecws.domain}/PHIC/Claims3.0/searchEmployer`

Method

POST

Header

Key	Value
token	PECWS authentication token

Body

JSON object containing the following key-value pairs

Key	Value
philhealthno	PhilHealth Employer Number
employername	PhilHealth Employer Name

Output

This method returns a JSON object with the key-value pairs described below:

Key	Value/Remarks										
success	A value of 'true' indicates a successful operation										
message	If an error was encountered during the execution of this method, this will contain the error message										
result	<div>The JSON object as the result of the encryption (using the cipher key of the health facility) of the XML text containing the Receipt Ticket Number and other data about the processing of the submitted e-claim data.<table><tr><th>Key</th><th>Value/Remarks</th></tr><tr><td>docMimeType</td><td>"text/xml"</td></tr><tr><td>hash</td><td>The computed hash value of the unencrypted text.</td></tr><tr><td>key1</td><td>"" (Empty string)</td></tr><tr><td>key2</td><td>"" (Empty string)</td></tr></table></div>	Key	Value/Remarks	docMimeType	"text/xml"	hash	The computed hash value of the unencrypted text.	key1	"" (Empty string)	key2	"" (Empty string)
Key	Value/Remarks										
docMimeType	"text/xml"										
hash	The computed hash value of the unencrypted text.										
key1	"" (Empty string)										
key2	"" (Empty string)										

	iv	The initialization vector used in the AES encryption
	doc	The encrypted text of the JSON object containing the records of the matching benefit packages. Sample XML text and the DTD of the XML text is shown below.

Sample JSON input:

```
{
  "PEN": "123456789012",
  "employerName": "EMP1"
}
```

Sample JSON output:

```
{
  "result": {
    "docMimeType": "text/xml",
    "hash": "dc8f4d74d977dfe701c0c9bbca0678300540591fac928f71a3841317e7a99aec",
    "key1": "",
    "key2": "",
    "iv": "y1jPMxvQE2aJPVnqqn1pDQ==",
    "doc": "Pms1FWFZT+odAp0qf2zManmroSUr3lYgDFnhYeJqBkuhJNMJU5geEN=="
  },
  "success": true,
  "message": ""
}
```

Sample decrypted result:

```
{
  "eEMPLOYERS": {
    "employer": [{
      "pPEN": "123456789012",
      "pEmployerName": "EMP1",
      "pEmployerAddress": "EMP1"
    },
    {
      "pPEN": "123456789013",
      "pEmployerName": "EMP2",
      "pEmployerAddress": ""
    }
  ]},
  "ASOF": "01-01-2024"
}
```

15. Get DB Server Date Time Method

This method retrieves the current database server date and time, allowing the caller to verify the availability of one of PhilHealth's databases and determine if it is down or inaccessible.

Endpoint

`https://{pecws.domain}/PHIC/Claims3.0/getDBServerDateTime`

Method

GET

Header

Key	Value
token	PECWS authentication token

Output

This method returns a JSON object with the key-value pairs described below:

Key	Value/Remarks								
success	A value of 'true' indicates a successful operation								
message	If an error was encountered during the execution of this method, this will contain the error message								
result	The JSON Array containing the following key-value pairs <table><tr><th>Key</th><th>Value/Remarks</th></tr><tr><td>server</td><td>Database Server</td></tr><tr><td>dateTime</td><td>Database current date and time</td></tr><tr><td>remarks</td><td></td></tr></table>	Key	Value/Remarks	server	Database Server	dateTime	Database current date and time	remarks	
Key	Value/Remarks								
server	Database Server								
dateTime	Database current date and time								
remarks									

Sample JSON output:

```
{
  "result": [{
    "server": "DB Server 1",
    "datetime": "01-01-2024 01:20:20 PM",
    "remarks": ""
  },
  {
    "server": "DB Server 2",
    "datetime": "01-01-2024 01:20:20 PM",
    "remarks": ""
  }
],
  "success": true,
  "message": ""
}
```

16. Get Server Date Time Method

This method retrieves the current server date and time, allowing the caller to monitor the system's operational status and determine whether the API or server is active or experiencing downtime.

Endpoint

`https://{pecws.domain}/PHIC/Claims3.0/getServerDateTime`

Method

GET

Header

Key	Value
token	PECWS authentication token

Output

This method returns a JSON object with the key-value pairs described below:

Key	Value/Remarks								
success	A value of 'true' indicates a successful operation								
message	If an error was encountered during the execution of this method, this will contain the error message								
result	The JSON Object containing the following key-value pairs <table><tr><th>Key</th><th>Value/Remarks</th></tr><tr><td>server</td><td>Database Server</td></tr><tr><td>dateTime</td><td>Database current date and time</td></tr><tr><td>remarks</td><td></td></tr></table>	Key	Value/Remarks	server	Database Server	dateTime	Database current date and time	remarks	
Key	Value/Remarks								
server	Database Server								
dateTime	Database current date and time								
remarks									

Sample JSON output:

```
{
  "result": {
    "server": "DB Server 1",
    "datetime": "01-01-2024 01:20:20 PM",
    "remarks": ""
  },
  "success": true,
  "message": ""
}
```

}

17. Get Server Version Method

This method retrieves the current server version.

Endpoint

`https://{pecws.domain}/PHIC/Claims3.0/getServerVersion`

Method

GET

Header

Key	Value
token	PECWS authentication token

Output

This method returns a JSON object with the key-value pairs described below:

Key	Value/Remarks
success	A value of 'true' indicates a successful operation
message	If an error was encountered during the execution of this method, this will contain the error message
result	PECWS Version

Sample JSON output:

```
{
  "result": "PECWS 3.0",
  "success": true,
  "message": ""
}
```

18. Is Claim Eligible Method

This method verifies the eligibility of a member and their valid dependents for claims availment.

Eligibility should be checked upon admission (initial call), while the final call generates a Tracking Number to confirm the member's eligibility for their claims availment.

Endpoint

`https://{pecws.domain}/PHIC/Claims3.0/isClaimEligible`

Method

POST

Header

Key	Value
token	PECWS authentication token

Body

The body is in JSON format as encrypted using the cipher key of the Health Facility. The JSON object contains the following key-value pairs:

Key	Value
hospitalCode	Accreditation number of the health facility
isForOPDHemodialysisClaim	If the purpose of checking eligibility is for hemodialysis claim. Indicate "Y" for Yes and "N" for No.
memberPIN	PIN or PhilHealth Identification Number is an identification number issued by PhilHealth to the member
memberBasicInformation	Basic Information
patientIs	Flag whether patient is the member or if dependent the relationship of patient with the member. <ul style="list-style-type: none">• 'M' – patient is member (Self)• 'S' – patient is spouse• 'C' – patient is child• 'P' – patient is parent
admissionDate	Refer to the admission date. The format should be "MM-DD-YYYY"
patientPIN	The PIN to the patient. It is either member or dependent PIN.
PatientBasicInformation	Basic Information
membershipType	PhilHealth membership type of the member
pEN	The PhilHealth issued number to the employer of the member, but not limited to the following: <ul style="list-style-type: none">• 'S' – Employed Private• 'G' – Employer Government• 'I' – Indigent• 'NS' – Individually Paying

	<ul style="list-style-type: none"> • ‘NO’ – OFW ‘ • PS’ – Non Paying Private ‘ • PG’ – Non Paying Government • ‘P’ – Lifetime member
employerName	Applicable only if the member is employed. Refer to the employer's name of the member.
isFinal	An indicator if the call for eligibility checking is final. Indicate “1” for final and “0” for not final.

Basic Information

Key	Value
lastname	The last name of the member or patient
firstname	The first name of the member or patient
middlename	The middle name of the member or patient
maidenname	The maiden name of the member or patient
suffix	The suffix or extension name of the member or patient
sex	The sex of the member or patient (M=for Male and F=for Female)
dateOfBirth	The date of birth of the member or patient. Format should be “MM-DD-YYYY”

Output

This method returns a JSON object with the key-value pairs described below:

Key	Value/Remarks																
success	A value of ‘true’ indicates a successful operation																
message	If an error was encountered during the execution of this method, this will contain the error message																
result	<p>An encrypted JSON object which contains the following key-value pairs</p> <table> <tr> <th>Key</th><th>Value/Remarks</th></tr> <tr> <td>docMimeType</td><td>"text/xml"</td></tr> <tr> <td>hash</td><td>The computed hash value of the unencrypted text.</td></tr> <tr> <td>key1</td><td>"" (Empty string)</td></tr> <tr> <td>key2</td><td>"" (Empty string)</td></tr> <tr> <td>iv</td><td>The initialization vector used in the AES encryption</td></tr> <tr> <td>doc</td><td>The encrypted text of the JSON object containing the detailed result of the PBEF inquiry.</td></tr> </table> <p>When decrypted, the JSON object contains the following key-value pairs</p> <table> <tr> <th>Key</th><th>Value/Remarks</th></tr> </table>	Key	Value/Remarks	docMimeType	"text/xml"	hash	The computed hash value of the unencrypted text.	key1	"" (Empty string)	key2	"" (Empty string)	iv	The initialization vector used in the AES encryption	doc	The encrypted text of the JSON object containing the detailed result of the PBEF inquiry.	Key	Value/Remarks
Key	Value/Remarks																
docMimeType	"text/xml"																
hash	The computed hash value of the unencrypted text.																
key1	"" (Empty string)																
key2	"" (Empty string)																
iv	The initialization vector used in the AES encryption																
doc	The encrypted text of the JSON object containing the detailed result of the PBEF inquiry.																
Key	Value/Remarks																

	isok	Eligibility response (YES NO)
	referenceno	Reference number for printing PBEF
	trackingno	Tracking number
	asof	Response Date

Sample encrypted JSON input payload:

```
{
  "result": {
    "docMimeType": "text/xml",
    "hash": "dc8f4d74d977dfe701c0c9bbca0678300540591fac928f71a3841317e7a99aec",
    "key1": "",
    "key2": "",
    "iv": "y1jPMxvQE2aJPVnqqn1pDQ==",
    "doc": "PMS1FWFZT+odAp0qf2zManmroSUr3lYgDFnhYeJqBkuhJNMJU5geEN=="
  }
}
```

Sample decrypted JSON input payload:

```
{
  "hospitalCode": "",
  "isForOPDHemodialysisClaim": "",
  "memberPIN": "",
  "memberBasicInformation": {
    "lastname": "",
    "firstname": "",
    "middlename": "",
    "maidenname": "",
    "suffix": "",
    "sex": "",
    "dateOfBirth": ""
  },
  "patientIs": "",
  "admissionDate": "",
  "patientPIN": "",
  "patientBasicInformation": {
    "lastname": "",
    "firstname": "",
    "middlename": "",
    "maidenname": "",
    "suffix": "",
    "sex": "",
    "dateOfBirth": ""
  },
  "membershipType": "",
  "pEN": "",
  "employerName": "",
  "isFinal": ""
}
```

Sample JSON output:

```
{
  "result": {
    "docMimeType": "text/xml",
    "hash": "dc8f4d74d977dfe701c0c9bbca06783005405",
    "key1": "",
    "key2": "",
    "iv": "y1jPMxvQE2aJPVnqqn",
  }
}
```

```
        "doc": " odAp0qf2zManmroSUr3lYgDFnhYeJqBkuhJNMJU5geEN=="
    }
    "success": true,
    "message": ""
}
```

Note: Decrypt the result which is an encrypted base64 PDF document. You may render it as shown below.

```
<embed src="data:application/pdf;base64,decrypted_base64_value_here"
type="application/pdf" />
```


19. Generate PBEF PDF Method

This method generates a PDF file of the PhilHealth Benefit Eligibility Form (PBEF) using the Reference Number obtained from the eligibility check.

Endpoint

`https://{pecws.domain}/PHIC/Claims3.0/generatePBEFPDF`

Method

POST

Header

Key	Value
token	PECWS authentication token

Parameter

Key	Value
accreno	Health facility accreditation number
referenceno	PBEF reference number from isClaimEligible method

Sample JSON Input Payload

```
{
  "accreno": "",
  "referenceno": ""
}
```

Output

This method returns a JSON object with the key-value pairs described below:

Key	Value/Remarks								
success	A value of 'true' indicates a successful operation								
message	If an error was encountered during the execution of this method, this will contain the error message								
result	An encrypted JSON object which contains the following key-value pairs <table><tr><th>Key</th><th>Value/Remarks</th></tr><tr><td>docMimeType</td><td>"text/xml"</td></tr><tr><td>hash</td><td>The computed hash value of the unencrypted text.</td></tr><tr><td>key1</td><td>"" (Empty string)</td></tr></table>	Key	Value/Remarks	docMimeType	"text/xml"	hash	The computed hash value of the unencrypted text.	key1	"" (Empty string)
Key	Value/Remarks								
docMimeType	"text/xml"								
hash	The computed hash value of the unencrypted text.								
key1	"" (Empty string)								

	key2	"" (Empty string)
	iv	The initialization vector used in the AES encryption
	doc	The encrypted text of the JSON object containing the detailed result of the PBEF inquiry.

When decrypted, it contains a base64 string representation of a PDF document. You may render the PDF using any base64 renderer.

Example:

```
<embed src="data:application/pdf;base64,decrypted_base64_value_here" type="application/pdf" />
```

Sample JSON output:

```
{
  "result": {
    "docMimeType": "text/xml",
    "hash": "dc8f4d74d977dfe701c0c9bbca06783005405",
    "key1": "",
    "key2": "",
    "iv": "y1jPMxvQE2aJPVnqqn",
    "doc": " odAp0qf2zManmroSUR3lYgDFnhYeJqBkuhJNMJU5geEN=="
  }
  "success": true,
  "message": ""
}
```

Annex A

Guidelines for the Data Encryption Using the Cipher Key of the Health Facility

Data will be encrypted using AES-256-CBC algorithm. PhilHealth issues a *cipher key* to the health facility for each certified software. The cipher key is to be used in generating the secret key for the AES encryption. System developers of service providers/health facilities may use their preferred programming language and tools for the data encryption/decryption.

```
{  
  "docMimeType": "{MIME type of the data}",  
  "hash": "{SHA-256 hash of the data before encryption}",  
  "key1": "",  
  "key2": "",  
  "iv": "{The initialization vector encoded as base-64 string}",  
  "doc": "{The encrypted data encoded as base-64 string}"  
}
```

Figure 1: Sample Format/Layout of the Resulting Encrypted Data

1. Get the hash total of the data to be decrypted using SHA-256 hash algorithm. After decrypting the data, a new hash total can be computed using the same algorithm, which can then be compared to the original hash total generated before decryption. If the two hash totals match, the data integrity is preserved. However, if the hash totals do not match, the data may have been tampered with or corrupted during decryption.
2. Set the value of the secret key to be used for the AES encryption.
 - a. Hash the given cipher key using SHA-256 algorithm.
 - b. Use the first 32 bytes of the resulting hash value as the passphrase for the AES encryption.
 - c. If the length of the resulting hash value is less than 32 bytes, pad the hash value with null character (with hexadecimal value of '0x00') to have a value having 32 bytes.
3. Prepare the initialization vector (IV) to be used
 - a. Generate an array of random 16 bytes (for a total of 128 bits).
 - b. Encode the array of bytes as base 64 string
4. Encrypt the data
 - a. Convert the data as array of bytes

- b. Encrypt the data using AES encryption using the values the array of bytes of the password, the IV and the data. Pad the data with null character (with hexadecimal value of '0x00') if it is not a multiple of 16 bytes.
 - c. Encode the resulting encrypted array of bytes as base-64 string.
- 5. Build the JSON string of the encrypted data using the following names and their corresponding values:
 - a. docMimeType: The value for this key will be the MIME type of the target data. For example, if the target data is a string in JSON format, the value for this key should be "application/json".
 - b. hash: Use the resulting value of Step #1 as the value of this key
 - c. key1: Set the value of this key as empty string.
 - d. key2: Set the value as empty string.
 - e. iv: Use the resulting value of Step #3.b as the value of this key
 - f. doc: Use the resulting value of Step #4.c as the value of this key

Annex B

Document Type Code and Description

Document Code	Description
CAB	Clinical Abstract
CAE	Certification of Approval/Agreement from the Employer
CF1	Claim Form 1
CF2	Claim Form 2
CF3	Claim Form 3
CF4	Claim Form 4
COE	Certificate of Eligibility
CSF	Claim Signature Form
CTR	Confirmatory Test Results by SACCL or RITM
DTR	Diagnostic Test Result
ITB	Itemized Billing (PDF format)
ITX	Itemized Billing (MS Excel Format)
MBC	Member's Birth Certificate
MDR	Proof of MDR with Payment Details
MEF	Member Empowerment Form
MMC	Member's Marriage Contract
MSR	Malarial Smear Results
NHC	Newborn Hearing Registry Card (Blue Card)
NHT	Newborn Hearing Screening Test Result
MWV	Waiver for Consent for Release of Confidential Patient Health Information
NTP	NTP Registry Card
OPR	Operative Record
ORS	Official Receipts
PAC	Pre-Authorization Clearance
PBC	Patient's Birth Certificate

PIC	Valid Philhealth Indigent ID
POR	PhilHealth Official Receipts
SOA	Statement of Account
STR	HIV Screening Test Result
TCC	TB-Diagnostic Committee Certification (-) Sputum
TYP	Three Years Payment of (2400 x 3 years of proof of payment)
MRF	PhilHealth Member Registration Form
ANR	Anesthesia Record
HDR	Hemodialysis Record
CF4	Claim Form 4
CF5	Electronic Claim Form 5 (for DRG Shadow Billing)
ESA	Electronic Statement of Account (eSOA)
NIR	Neuroimaging Result
RGR	Radiographic Result
BCR	Blood Culture Result
CBC	Complete Blood Count Result
OTH	Other documents

Annex C

Data Dictionary eClaimsUpload

Name	Length	Description	Valid Values
pUserName	String(20)	Provider user id	To be provided by PhilHealth
pUserPassword	String(20)	Provider user password	To be provided by PhilHealth
pHospitalCode	String(12)	Facility Accreditation Number	For now PMCC number should be used
pHospitalEmail	String(150)	Hospital Email Address where communication will be sent	Must not be blank
pHospitalTransmittalNo		Hospital Transmittal Number	Generated by the Hospital own batching system. This should be unique per hospital.
pTotalClaims	String(3)	Claims counter	Integer format
pClaimNumber	String(12)	Hospital Claim Number	Hospital Generated Claim Case #, this should be unique per hospital
pTrackingNumber	String(20)	The Claims Eligibility Tracking number assigned if undergone the Online Eligibility Checking	Formatted as: '####-####-####-####' Can be blank
pCataractPreAuth	String(20)	Cataract Pre-Authorization Application Number	
pPhilhealthClaimType	String(20)	Flag whether Claims Payment Mechanism	<ul style="list-style-type: none"> • 'ALL-CASE-RATE' • 'Z-BENEFIT'
pPatientType	String(1)	Patient Type	<ul style="list-style-type: none"> • 'I' – Inpatient • 'O' – Outpatient
pIsEmergency	String(1)	Flag if Emergency Case	<ul style="list-style-type: none"> • 'Y' – Yes • 'N' – No
pClaimSeriesLhio	String(15)	Philhealth Generated and Assigned Unique Number per Claim	Can be used by the hospital to reconcile their records with Philhealth This will be returned after the claim are uploaded to Philhealth
pMemberPIN	String(12)	PhilHealth Identification Number – a unique 12 digit number assigned to a member.	The last character in the PIN is a modulus 11 check digit.
pMemberLastName	String(60)	Member's Complete Surname	Any value consisting of : 'A' to 'Z', 'N'. Can include a space in between characters
pMemberFirstname	String(60)	Member's Complete First name	
pMemberMiddleName	String(60)	Member's Complete Middle name	
pMemberSuffix	String(5)	Member's Suffix name	'JR', 'SR', 'III', ...etc Suffixes can be blank
pMemberBirthDate	String(10)	Member's Birth Date	Date Format should be : 'MM-DD-YYYY'
pMemberShipType	String(2)	Philhealth membership type of the member	(Not limited to the following :) <ul style="list-style-type: none"> • 'S' – Employed Private • 'G' – Employer Government • 'I' – Indigent • 'NS' – Individually Paying • 'NO' – OFW • 'PS' – Non Paying Private • 'PG' – Non Paying Government • 'P' – Lifetime Member
pMailingAddress	String(150)	Mailing Address (address where the benefit payment notice will be sent)	This is where the notices will be mailed.
pZipCode	String(4)	Philippine Zip Code of the municipality	4 digit Philippine zip code value of the municipality/city
pMemberSex	String(1)	Member Sex	<ul style="list-style-type: none"> • 'M' – Male • 'F' – Female
pLandlineNo	String(20)	Members Landline Number	Can be blank
pMobileNo	String(20)	Members Cell Number	Can be blank
pEmailAddress	String(150)	Email Address	Can be blank
pPatientIs	String(1)	Flag whether patient is the member or if dependent the relationship of patient with the member.	<ul style="list-style-type: none"> • 'M' – patient is member(Self) • 'S' – patient is spouse • 'C' – patient is child • 'P' – patient is parent
pPatientPIN	String(12)	PhilHealth Identification Number – a unique 12 digit number assigned to a patient.	The last character in the PIN is a modulus 11 check digit.
pPatientLastName	String(60)	Patient's Complete Surname	

pPatientFirstName	String(60)	Patient's Complete First name	Same as for the member. These can be blank since these are disregarded if the value of pPatientIs is 'M'
pPatientMiddleName	String(60)	Patient's Complete Middle name	
pPatientSuffix	String(5)	Patient's Suffix name	'JR', 'SR', 'III', ...etc Suffixes can be blank
pPatientBirthDate	String(10)	Patient's Birth Date	Date Format should be : 'MM-DD-YYYY'
pPatientSex	String(1)	Member Sex	<ul style="list-style-type: none"> • 'M' – Male • 'F' – Female
pPEN	String(12)	Philhealth Employer Number – a unique 12 digit number assigned to an employer	These are disregarded if pMembershipType is not ('S' or 'G')
pEmployerName	String(100)	The Registered name of the employer	
pPatientReferred	String(1)	Referred Patient	<ul style="list-style-type: none"> • 'Y' – Yes • 'N' – No
pReferredIHCPAccreCode	String(12)	Referring Facility Accreditation Code	Required if the patient is referred by another IHCP
pAdmissionDate	String(10)	Admission Date	Date Format should be : 'MM-DD-YYYY'
pDischargeDate	String(10)	Discharge Date	
pAdmissionTime	String(10)	Admission Time	Time Format should be : 'HH:MM:SSAM/PM'
pDischargeTime	String(10)	Discharge Time	
pDisposition	String(1)	Patient's Disposition	<ul style="list-style-type: none"> • 'I' – Improved • 'R' – Recovered • 'H' – Home/Discharged Against Medical Advise • 'A' – Absconded • 'E' – Expired • 'T' – Transferred/Referred
pExpiredDate	String(10)	Date of Death of patient	Date Format should be : 'MM-DD-YYYY'
pExpiredTime	String(10)	Time of Death of patient	Time Format should be : 'HH:MM:SSAM/PM'
pReferralIHCPAccreCode	String(12)	Referral Facility Accreditation Code	Required when pDisposition = 'E'
pReferralReasons	String(150)	Reason/s for referral/transfer	
pAccommodationType	String(1)	Type of Accommodation	<ul style="list-style-type: none"> • 'P' – Private • 'N' – Non-Private (Charity/Service)
pAdmissionDiagnosis	String(500)	Admission Diagnosis	Can be multiple lines
pDischargeDiagnosis	String(500)	Discharge Diagnosis	
pICDCode	String(15)	ICD 10 Code of the illness	<i>Refer to ICD10 library</i>
pRelatedProcedure	String(150)	Related Procedure	Any value consisting of : 'A' to 'Z', 'N'. Can include a space in between characters
pRVSCode	String(6)	Relative Value Scale Code of the procedure/operation performed	<i>See RVS Library</i> Required for Operating Room, Surgeons and Anesthesiologist claims only
pProcedureDate	String(10)	Date of Procedure	Date Format should be : 'MM-DD-YYYY'
pLaterality	String(1)	Laterality	<ul style="list-style-type: none"> • 'L' – Left • 'R' – Right • 'B' – Both • 'N' – N/A
pSessionDate	String(10)	Inclusive Date of Session	Date Format should be : 'MM-DD-YYYY'
pCheckUpDate1	String(10)	1 st Checkup Date for MCP Package	Date Format should be : 'MM-DD-YYYY'
pCheckUpDate2	String(10)	2 nd Checkup Date for MCP Package	
pCheckUpDate3	String(10)	3 rd Checkup Date for MCP Package	
pCheckUpDate4	String(10)	4 th Checkup Date for MCP Package	
pTBType	String(1)	Type of TB-Dots claim	<ul style="list-style-type: none"> • 'I' – Intensive Phase • 'M' – Maintenance <p>NOTE For 'I', admission and discharge dates should be the first and last days of treatment in the intensive phase respectively. For 'M', admission and discharge dates should be the first and last days of treatment in the maintenance phase respectively. Required for TB-DOTS claims only</p>
pNTPCardNo	String(10)	NTP Card No	Required for TB-DOTS claims only

pDay0ARV	String(10)	Day 0 ARV (Anti Rabies Vaccine)	Date Format should be: 'MM-DD-YYYY' Required for Animal Bite Package claims only
pDay3ARV	String(10)	Day 3 ARV (Anti Rabies Vaccine)	
pDay7ARV	String(10)	Day 7 ARV (Anti Rabies Vaccine)	
pRIG	String(10)	RIG (Rabies Immunoglobulin)	
pABPOthers	String(10)	Other Date	
pABPSpecify	String(50)	Others (Specify)	Required for Animal Bite Package claims only
pEssentialNewbornCare	String(1)	Essential Newborn Care	Flag whether Yes (Y) or No (N) Required for Newborn Care Package
pNewbornHearingScreeningTest	String(1)	Newborn Hearing Screening Test	
pNewbornScreeningTest	String(1)	Newborn Screening Test	
pFilterCardNo	String(20)	Filter Card Number	
pDrying	String(1)	Immediate drying of newborn	
pSkinToSkin	String(1)	Early skin-to-skin contact	Required when the pEssentialNewbornCare = 'Y'
pCordClamping	String(1)	Timely cord clamping	
pProphylaxis	String(1)	Eye prophylaxis	
pWeighing	String(1)	Weighing of the newborn	
pVitaminK	String(1)	Vitamin K administration	
pBCG	String(1)	BCG vaccination	
pNonSeparation	String(1)	Non-separation of mother/baby for early breastfeeding initiation	
pHepatitisB	String(1)	Hepatitis B vaccination	
pLaboratoryNumber	String(20)	Laboratory Number	Required for Outpatient HIV/AIDS Treatment Package
pDoctorAccreCode	String(12)	Doctor's Accreditation Number	Formatted as: '#####-##'
pDoctorLastName	String(60)	Doctor's Complete Surname	Same as for the member.
pDoctorFirstName	String(60)	Doctor's Complete First name	
pDoctorMiddleName	String(60)	Doctor's Complete Middle name	
pDoctorSuffix	String(5)	Doctor's Suffix name	'JR', 'SR', 'III', ...etc Suffixes can be blank
pWithCoPay	String(1)	Flag whether the professional fee is no co-pay or with co-pay	<ul style="list-style-type: none"> • 'Y' – With co-pay • 'N' – Noco-pay
pDoctorCoPay	String(12)	Amount of co-pay	Formatted as: '#####.##' Required when pWithCoPay='Y'
pEnoughBenefits	String(1)	(Y) – if the PhilHealth benefit is enough to cover HCI and PF charges. No purchases of drugs/medicines, supplies, diagnostics, and co-pay for professional fees by the member/patient. (N) –if the benefit of the member/patient was completely consumed prior to co-pay OR the benefit of the member/patient is not completely consumed BUT with purchases/expenses for drugs/medicines, supplies, diagnostics and others.	Flag whether Yes(Y) or No(N)
pTotalHCIFees	String(12)	Total Health Care Institution Fees	Formatted as: '#####.##' Required when pEnoughBenefits = 'Y'
pTotalProfFees	String(12)	Total Professional Fees	
pGrandTotal	String(12)	Grand Total is equal to the HCI and Prof Fees	
pTotalActualCharges	String(12)	Total Actual Charges	Formatted as: '#####.##' Required when pEnoughBenefits = 'N'
pDiscount	String(12)	Amount after Application of Discount	
pPhilhealthBenefit	String(12)	PhilHealth Benefit	
pTotalAmount	String(12)	Amount after PhilHealth Deduction	
pMemberPatient	String(1)	Member/Patient	
pHMO	String(1)	HMO	Flag whether Yes(Y) or No(N) if applicable
pOthers	String(1)	Others (i.e., PCSO, Promissory note, etc.)	
pDrugsMedicinesSupplies	String(1)	(Y) –if there is a purchase/s for drugs/medicines and/or medical supplies bought by the patient/member within/outside the HCI during confinement (N) – None	

pDMSTotalAmount	String(12)	Total Amount for Drugs, Medicines and Supplies	Formatted as: '#####.##' Required when pDMSTotalAmount = 'Y'
pExaminations	String(1)	(Y) – if there is a diagnostic/laboratory examinations paid for by the patient/member done within/outside the HCI during confinement (N) – None	Flag whether Yes(Y) or No(N)
pExamTotalAmount	String(12)	Total Amount for Diagnostic	Formatted as: '#####.##' Required when pExamTotalAmount = 'Y'
pCaseRateCode	String(6)	Case Rate Codes for All Case Rates	See Case Rate Library
pZBenefitCode	String(7)	Z-Benefit Codes	<p>Standard Risk Acute Lymphocytic (lymphoblastic) Leukemia for Children</p> <ul style="list-style-type: none"> • Z0011 – 1st tranche • Z0012 – 2nd tranche • Z0013 – 3rd tranche <p>Early Stage Breast Cancer (Stage 0 to III-A)</p> <ul style="list-style-type: none"> • Z0021 –1st tranche • Z0022 –2nd tranche <p>Low to Intermediate Risk Prostate Cancer Requiring Prostatectomy</p> <ul style="list-style-type: none"> • Z003 – full payment <p>End Stage Renal Disease Eligible for Kidney Transplant (Low Risk)</p> <ul style="list-style-type: none"> • Z0041 –1st tranche • Z0042 –2nd tranche <p>Elective Surgery for Standard Risk Coronary Artery Bypass Graft (CABG)</p> <ul style="list-style-type: none"> • Z0051 – 1st tranche • Z0052 – 2nd tranche <p>Tetralogy of Fallot (TOF)</p> <ul style="list-style-type: none"> • Z0061 – 1st tranche • Z0062 – 2nd tranche <p>Ventricular Septal Defect(VSD)</p> <ul style="list-style-type: none"> • Z0071 – 1st tranche • Z0072 – 2nd tranche <p>Cervical Cancer Chemoradiation with Cobalt & Brachytherapy (Low Dose) or Primary Surgery for Stage IA1, IA2 – IIA1</p> <ul style="list-style-type: none"> • Z0081 – 1st tranche • Z0082 – 2nd tranche <p>Cervical Cancer Chemoradiation with Linear Accelerator & Brachytherapy (High Dose)</p> <ul style="list-style-type: none"> • Z0091 – 1st tranche • Z0092 – 2nd tranche
pChiefComplaint	String(200)	Chief Complaint or Reason for Admission	Text/Memo
pBriefHistory	String(2500)	Brief History of Present Illness (OB Score/OB History)	Text/Memo
pCourseWard	String(500)	Course in the Wards	Text/Memo
pPertinentFindings	String(500)	Pertinent Laboratory and Diagnostic Findings	Text/Memo
pPrenatalConsultation	String(10)	Initial Prenatal Consultation	Date Format should be : 'MM-DD-YYYY'
pMCPOrientation	String(1)	Orientation to MCP/Availment of Benefits	Flag whether Yes(Y) or No(N)
pExpectedDeliveryDate	String(10)	Expected data of delivery	Date Format should be : 'MM-DD-YYYY'
pVitalSigns	String(1)	Vital Signs are Normal	Flag whether Yes(Y) or No(N)
pPregnancyLowRisk	String(1)	Ascertain the Present. Pregnancy is low-risk	
pLMP	String(10)	Last Menstrual Period	Date Format should be :

			'MM-DD-YYYY'
pMenarcheAge	String(2)	Age of Menarche	Must be an integer
pObstetricG	String(10)	Obstetric History	Text
pObstetricP	String(10)		Text
pObstetric_T	String(10)		Text
pObstetric_P	String(10)		Text
pObstetric_A	String(10)		Text
pObstetric_L	String(10)		Text
pMultiplePregnancy	String(1)	Multiple Pregnancy	Obstetric Risk Factors Flag whether Yes(Y) or No(N)
pOvarianCyst	String(1)	Ovarian Cyst	
pMyomaUteri	String(1)	Myoma Uteri	
pPlacentaPrevia	String(1)	Placenta Previa	
pMiscarriages	String(1)	History of 3 miscarriages	
pStillBirth	String(1)	History of stillbirth	
pPreEclampsia	String(1)	History of pre-eclampsia	
pEclampsia	String(1)	History of eclampsia	
pPrematureContraction	String(1)	Premature contraction	
pHypertension	String(1)	Hypertension	Medical/Surgical risk factors Flag whether Yes(Y) or No(N)
pHeartDisease	String(1)	Heart Disease	
pDiabetes	String(1)	Diabetes	
pThyroidDisaster	String(1)	Thyroid Disorder	
pObesity	String(1)	Obesity	
pAsthma	String(1)	Moderate to severe asthma	
pEpilepsy	String(1)	Epilepsy	
pRenalDisease	String(1)	Renal Disease	
pBleedingDisorders	String(1)	Bleeding disorders	
pPreviousCS	String(1)	History of previous caesarian section	
pUrineMyomectomy	String(1)	History of uterine myomectomy	
pVisitDate	String(10)	Date of visit	Date Format should be : 'MM-DD-YYYY'
pAOGWeeks	String(3)	AOG in weeks	Text
pWeight	String(10)	Weight	Weight & Vital Signs
pCardiacRate	String(10)	Cardiac Rate	
pRespiratoryRate	String(10)	Respiratory Rate	
pBloodPressure	String(10)	Blood Pressure	
pTemperature	String(10)	Temperature	
pDeliveryDate	String(10)	Date of Delivery	Date Format should be : 'MM-DD-YYYY'
pDeliveryTime	String(10)	Time of Delivery	Time Format should be : 'HH:MM:SSAM/PM'
pObstetricIndex	String(50)	Obstetric Index	Text
pAOG LMP	String(50)	AOG by LMP	Text
pDeliveryManner	String(50)	Manner of Delivery	Text
pPresentation	String(50)	Presentation	Text
pFetalOutcome	String(50)	Fetal Outcome	Text
pSex	String(1)	Sex	<ul style="list-style-type: none"> • 'M' – Male • 'F' –Female
pBirthWeight	String(10)	Birth Weight (gram)	Must be an integer
pAPGARScore	String(10)	APGAR Score	
pPostpartum	String(10)	Scheduled Pospartum follow-up consultation 1 week after delivery	Date Format should be : 'MM-DD-YYYY'
pPerinealWoundCare	String(1)	Perineal wound care	Flag whether Yes(Y) or No(N)

pPerinealRemarks	String(100)	Remarks	Text/Memo
pMaternalComplications	String(1)	Signs of Maternal Postpartum Complications	Flag whether Yes(Y) or No(N)
pMaternalRemarks	String(100)	Remarks	Text/Memo
pBreastFeeding	String(1)	Breastfeeding and Nutrition	Flag whether Yes(Y) or No(N)
pBreastFeedingRemarks	String(100)	Remarks	Text/Memo
pFamilyPlanning	String(1)	Family Planning	Flag whether Yes(Y) or No(N)
pFamilyPlanningRemarks	String(100)	Remarks	Text/Memo
pPlanningService	String(1)	Provided family planning service to patient (as requested by patient)	Flag whether Yes(Y) or No(N)
pPlanningServiceRemarks	String(100)	Remarks	Text/Memo
pSurgicalSterilization	String(1)	Referred to partner physician for Voluntary Surgical Sterilization (as requested by patient)	Flag whether Yes(Y) or No(N)
pSterilizationRemarks	String(100)	Remarks	Text/Memo
pFollowupSchedule	String(1)	Schedule the next postpartum follow-up	Flag whether Yes(Y) or No(N)
pFollowupScheduleRemarks		Remarks	Text/Memo
pReferredReason	String(500)	Reason for referral from other HCI	Text/Memo
pIntensive	String(1)	Completed Intensive Phase	Flag whether Yes(Y) or No(N)
pMaintenance	String(1)	Completed Maintenance Phase	
pCriteria	String(200)	Clinical Criteria and/or Laboratories/Diagnostics Criteria	Text/Memo
pBP	String(20)	Blood Pressure	Text
pCR	String(20)		Text
pRR	String(20)	Respiratory Rate	Text
pTemp	String(20)	Temperature	Text
pHEENT	String(20)	Head, Ears, Eyes, Nose and Throat	Text
pChestLungs	String(20)	Chest/Lungs	Text
pCVS	String(20)		Text
pAbdomen	String(20)	Abdomen	Text
pGUIE	String(20)		Text
pSkinExtremities	String(20)	Skin Extremities	Text
pNeuroExam	String(20)	Neuro Examination	Text
pCourseDate	String(10)	Date in the Course in the Ward	Date Format should be : 'MM-DD-YYYY'
pFindings	String(200)	Pertinent PE/Lab Findings	Text/Memo
pAction	String(200)	Doctors Order/Action	Text/Memo
pPurchaseDate	String(10)	Date of Purchased	Date Format should be : 'MM-DD-YYYY'
pDrugCode	String(20)	Hospital Assigned Drug Code	Can be any format
pPNDFCode	String(20)	PNDF Code (Blank until PNDF lib is available)	Can be Blank for now
pGenericName	String(50)	Generic Name of Medicines/Drugs taken	Must not be blank
pBrandName	String(50)	Brand Name of Medicines/Drugs taken	Must not be blank
pPreparation	String(30)	Dose/ Cap/ Syrup/ Injectible/ Tab with ml/mg/gm content	Must not be blank
pQuantity	String(10)	Unit quantity of item	Integer format
pDiagnosticDate	String(10)	Date of Diagnostic	Date Format should be : 'MM-DD-YYYY'
pDiagnosticType	String(20)	Type of diagnostic/test done	<ul style="list-style-type: none"> • 'IMAGING' • 'LABORATORY' • 'SUPPLIES' • 'OTHERS'
pDiagnosticName	String(50)	Name of Imaging procedure for Imaging, Name of Laboratory procedure for Laboratory, Name of Supplies for Supplies or Others for Supplies and Others	Must not be blank
pCompanyName	String(100)	Company's Name	Text

pCompanyTIN	String(15)	Company's TIN	Formatted as: '###-###-###-###'
pBIRPermitNumber	String(20)	BIR Permit Number	
pReceiptNumber	String(20)	Official Receipt Number	
pReceiptDate	String(12)	Official Receipt Date	Date Format should be : 'MM-DD-YYYY'
pVATExemptSale	String(10)	VAT Exempt Sale	Formatted as: '#####.##'
pVAT	String(10)	VAT – 12%	
pTotal	String(10)	Total Amount in the Receipt	
pUnitPrice	String(10)	Unit price of item	
pDescription	String(100)	Item Description	Formatted as: '#####.##'
pAmount	String(10)	Total Amount of the Specific Item	Formatted as: '#####.##'
pDocumentType	String(3)	Document to support the claim	See Document Library
pDocumentURL	String(250)	URL of the document accessible via https.	The document must first be encrypted using philhealth public key before publishing online. Please see the Annex for the guidelines for encryption.
pHospitalTransmittalNo	String(20)	Hospital Generated Transmittal Number	
pTransmissionControlNumber	String(18)	Philhealth Generated Transmittal file control number.	Will be blank if the transmission is failed
pErrCode	String(3)	Claim file error Code	
pErrDescription	String(100)	Claim file error Description	
pReceiptTicketNumber	String(18)	Philhealth Generated Upload Confirmation Receipt ticket number	
pTransmissionDate	String(10)		Date Format should be : 'MM-DD-YYYY'
pTransmissionTime	String(10)		TIME Format should be : 'HH:MM:SSAM/PM'
pReceivedDate	String(10)	Date when the transmitted file received by PhilHealth	Date Format should be : 'MM-DD-YYYY'
pHasAttachedSOA	String(1)	Type of Accommodation	<ul style="list-style-type: none"> • 'Y' – With attached SOA • 'N' –Without attached SOA
pDateSigned	String(10)	Date signed	Date Format should be : 'MM-DD-YYYY'
pRelCode	String(1)	Code indicating the relation of the representative who signed the consent to access patient record	Valid values: <ul style="list-style-type: none"> • S = Spouse • C = Child • P = Parent • I = Sibling • O = Others (pRelDesc should have a value)
pRelDesc	String(50)	Specified relation of the representative who signed the consent on behalf of the patient to access patient records. This include relation not included in the list of defined values of relation of pRelCode	
pReasonCode	String(1)	Code indicating the reason why a representative signed the consent to access patient record	Valid values: <ul style="list-style-type: none"> • I = Patient is incapacitated • O = Patient is incapacitated (pReasonDesc should have a value)
pReasonDesc	String(50)	Text indicating reason why a reason why a representative signed the consent to access patient record	
pPreAuthDate	String(10)	Z-Benefit pre-authorization date	Date Format should be : 'MM-DD-YYYY'
pThumbmarkedBy	String(1)	Indicates whether the thumbmark is from the member/patient or from the representative of the member	Valid values: <ul style="list-style-type: none"> • P = of the patient/member • R = of a representative
pServiceProvider	String(50)	The acronym or short name of the company/institution that provided the system used in submitting the e-claims XML file	
pCertificateId	String(50)	The software validation certificate number issued upon passing the software compliance testing	

Annex D

Data Dictionary ValidateEsoa

Element	Name	Length	Description	Valid Values
eSOA	pHciPan	Varchar2 (9)	The accreditation number issued by PhilHealth to the health facility	
	pTransmittalId	Varchar2 (50)	A UNIQUE reference number assigned to the claim by the submitting health facility	
Summary of Fees	PChargesNetOfApplicableVat	Number (10,2)	Refed to the total charges net of applicable vat	0.00 to 99999999.99
	PSeniorCitizenDiscount	Number (10,2)	The amount of senior citizen (SC) discount, if applicable	0.00 to 99999999.99
	pPWDDIscount	Number (10,2)	The amount of Person with Disability (PWD) discount, if applicable	0.00 to 99999999.99
	pPCSO	Number (10,2)	The amount of Philippine Charity Sweepstakes Office (PCSO) discount, if applicable	0.00 to 99999999.99
	pDSWD	Number (10,2)	The amount of Department of Social Welfare and Development (DSWD) discount, if applicable	0.00 to 99999999.99
	pDOHMAP	Number (10,2)	The amount of Department of Health (DOHMAP) discount, if applicable	0.00 to 99999999.99
	pHMO	Number (10,2)	The amount of HMO discount, if applicable	0.00 to 99999999.99
Other Funding Source	pDescription	Varchar (10,2)	The description of other funding source, if applicable	
	pAmount	Number (10,2)	The amount of other funding source discount, if applicable	0.00 to 99999999.99
PhilHealth	PTotalCaseRateAmount	Number (8,2)	The amount of case rate package being claimed	0.00 to 999999.99
Balance	pAmount	Number (10,2)	The difference between the total amount of charges, net of applicable vat, for the room and board, drugs and medicine, laboratory and diagnostic, operating room and medical supplies charges, less applicable discounts and the case rate amount.	0.00 to 99999999.99
Professional Information	pPAN	Varchar (14)	The corresponding PhilHealth Accreditation Number (PAN) of the attending Physician	
	pFirstName	Varchar (60)	The first name of the attending physician	
	pMiddleName	varchar	The middle name of the attending physician	
	pLastName	Varchar (60)	The last name of the attending physician	

	pSuffixName	Varchar (60)	The suffix or extension name of the attending physician (ex. Jr., Sr. II, IIII and so on)	
Itemized Billing	pServiceDate	Date	the date of service for room and board, laboratory and diagnostic service and operating room fees; the date of issuance for drugs and medicine and medical supplies	mm-dd-yyyy
	pItemCode	Varchar (30)	Refer to the code of the item in the electronic Statement of Account (eSOA)	as listed in eSOA Item Library <i>For Drug and Medicine Category, use the Drug Code (30 characters) from the Medicine Library.</i> <i>For other item not included in the library, keep it blank</i>
	pItemName	Varchar (300)	Refer to the description of an item in the electronic Statement of Account (eSOA)	as listed in eSOA and Item Library <i>For other item not included in the library, provide the complete description.</i>
	pUnitOfMeasurement	Varchar (50)	Refer to the unit of measurement used for an item in the electronic Statement of Account (eSOA)	<i>For Drug and Medicine Category, keep it blank.</i>
	pUnitPrice	Number (8,2)	Refer to the unit amount of an item in the electronic Statement of Account (eSOA)	0.00 to 999999.99
	pQuantity	Number (4,0)	Refer to the total number or quantity of an item in the electronic Statement of Account (eSOA)	0 to 9999
	pTotalAmount	Number (11,2)	Refer to the total amount of of an item in the electronic Statement of Account (eSOA) i.e pTotalAmount = pQuantity x pUnitPrice)	0.00 to 999999999.99
	pCategory	Varchar (50)	Refer to eSOA category	as listed in eSOA and Category Library

Annex E

Data Dictionary validateCF5

Element	Name	Type/Length	Description	Valid Values
CF5	pHospitalCode	Varchar (6)	A unique code assigned by PhilHealth to the health facility, also known as the PMCC No.	6-digits alphanumeric code Format: 999999 or X99999
DRG Claim	PrimaryCode	Varchar (15)	The corresponding ICD10 code for primary diagnosis	Valid ICD10 codes
	NewBornAdm Weight	Numeric (2,1)	This is applicable to newborn patients (0-27 days old), refers to the admission weight of newborn in kilogram (kg)	Should be greater than 0.3 kg Formatted as: ##.# Up to one (1) decimal places only
	Remarks	Varchar (2000)	Logs for internal key per value validation	Internal logs
	ClaimNumber	Varchar (13)	A reference number assigned to the claim upon successful submission of claim via eClaims API	Format: 9999999999999
Secondary Diagnosis	SecondaryCode	Varchar (15)	The ICD-10 codes corresponding to the secondary conditions diagnosed, up to a total of 12 codes	Valid ICD10 codes
	Remarks	Varchar (2000)	Logs for internal key per value validation	Internal logs
Procedures	RvsCode	Varchar (6)	The RVS codes corresponding to the procedure performed, up to a total of 20 codes	Valid RVS code
	Laterality	Varchar (1)	Refer to the side of the body on which a procedure was performed, either left, right or both	L = Left R = Right B = Both N = None
	Ext1	Numeric (1)	Refer to the number of body sites	1-9
	Ext2	Numeric (1)	Refer to the number of times the procedure was done	1-9
	Remarks	Varchar (2000)	Logs for internal key per value validation	Internal logs

Annex F

eSOA Item Library

A. Item Library

pItemCode	pCategory	pItemName
1	MedicalSupplies	BATH SOAP
2	MedicalSupplies	ALCOHOL ISOPROPHYL 70%, 150ML.
3	MedicalSupplies	ALCOHOL ISOPROPHYL 70%, 500ML.
4	MedicalSupplies	BALANCE SALT SOLUTION 500ML
5	MedicalSupplies	DEPROTEINIZER SOLUTION
6	MedicalSupplies	ENZYMATIC ENZYME INSTRUMENT CLEANER"METRIZIME"
7	MedicalSupplies	ETHYL-PROPANOL DISINFECTANT
8	MedicalSupplies	FLOURIDE GEL (480ML)
9	MedicalSupplies	HYDROGEN PEROXIDE, 500ML
10	MedicalSupplies	ISOPROPANOL CHLORIDE 250ML SPRAY
11	MedicalSupplies	ISOPROPANOL CHLORIDE 50ML SPRAY
12	MedicalSupplies	LICOCAINE 10MG/DOSE 10% 50ML PUMP SPRAY (XYLOCAINE)
13	MedicalSupplies	LIDOCAINE, PUMP SPRAY 10% 100ML
14	MedicalSupplies	MACHINE SURFACE DISENFECTANT OR ITS EQUIVALENT
15	MedicalSupplies	PLAIN NSS 1 LITER PLASTIC BOTTLE (NIR-NS)
16	MedicalSupplies	PLASMA STERILIZING AGENT 140ML SHOULD BE COMPATIBLE
17	MedicalSupplies	POVIDONE IODINE GARGLE,60ML,1% SOL
18	MedicalSupplies	STERILE WATER WITH HUMIDIFIER 500-540ML
19	MedicalSupplies	STERILLIUM / HAND DISINFECTANT PROPAN-2OL, PROPAN-1OL
20	MedicalSupplies	TERPOLYMER BASED ALCOHOL FREE BARRIER FILM SKIN PROTECTANT 28ML SPRAY 3M
21	MedicalSupplies	TRANSEPTIC CLEANSING SOL. 250ML SPRAY
22	MedicalSupplies	WOUND IRRIGATION SOLUTION POLYHEXIDE + BETAINE 350ML "B BRUAN PRONTOSAN SOLUTION"
23	MedicalSupplies	ALCOHOL 70% (60 ML)
24	MedicalSupplies	CHLORHEXIDINE DIGLUCONATE 02-1%
25	MedicalSupplies	DISTILLED WATER 50 ML
26	MedicalSupplies	HYDROGEN PEROXIDE
27	MedicalSupplies	NAIL POLISH REMOVER (AT LEAST 10ML)
28	MedicalSupplies	POVIDONE IODINE- (ZETADONE) 30ML
29	MedicalSupplies	SPECIMEN BOTTLE
30	MedicalSupplies	THORACOSTOMY BOTTLE 2000ML (PLASTIC DISPOSABLE)
31	MedicalSupplies	3 SPRING RESERVOIR/HEMOVAC SILICONE ROUND DRAIN W/ TROCAR 400ML/BOX OF 10
32	MedicalSupplies	3M ATTEST STEAM BIOLOGICAL INDICATOR
33	MedicalSupplies	BLADE SURGICAL S-11
34	MedicalSupplies	GLOVES EXAMINATION, SMALL
35	MedicalSupplies	HANDPIECE OIL
36	MedicalSupplies	HEMOSTAT COLLAGEN 5 X 8 CM, LYOSTYPT
37	MedicalSupplies	HEMOSTAT, LAYERED OXIDIZED REGENERATED CELLULOSE, ABS. 1X2 IN, BY 10S BOX
38	MedicalSupplies	HYDROFIBER DRESSING WITH IONIC SILVER 10X10CM
39	MedicalSupplies	HYDROFIBER DRESSING WITH IONIC SILVER 15X15CM
40	MedicalSupplies	HYDROFIBER NON ADHESIVE FOAM DRESSING 10X10CM BY 10S BOX

41	MedicalSupplies	NASAL PACK 4.5CMX1.5CMX2CM(10PER PIECES/BOX)
42	MedicalSupplies	NASAL PACK 8CMX1.5CMX2CM(20PER PIECES/BOX)
43	MedicalSupplies	N-BUTYL-2-CYANOACRYLATE SKIN GLUE TWIST OFF CAP 0.5ML/TUBE BOX OF 5
44	MedicalSupplies	NEEDLE CORE BIOPSY (BARD MAX-CORE) G18
45	MedicalSupplies	NEEDLE DISP. G. 18 X 1 1/2
46	MedicalSupplies	NEEDLE DISPOSABLE G. 21 X 1
47	MedicalSupplies	APPLICATOR TIP (MICROFILL TIP 100'S)
48	MedicalSupplies	AQUASONIC ULTRASOUND GEL IN FOIL 20G
49	MedicalSupplies	ARTICULATING PAPER(STRAIGHT)
50	MedicalSupplies	BLADE SURGICAL S-10 100'S
51	MedicalSupplies	BLADE SURGICAL S-12, 100'S
52	MedicalSupplies	BLADE SURGICAL S-15
53	MedicalSupplies	BLADE SURGICAL S-22, 100S
54	MedicalSupplies	BONDING AGENT
55	MedicalSupplies	BOWIE AND DICK TEST STEAM/VACUUM STERILIZER INDICATOR
56	MedicalSupplies	CALIBRATION PACK
57	MedicalSupplies	CAMERA COVER 40PER PIECES/PACK
58	MedicalSupplies	CELLULOID STRIP
59	MedicalSupplies	CONDOM PLAIN
60	MedicalSupplies	DIAMOND BURS (REUSABLE) INVERTED CONE 33
61	MedicalSupplies	DIAMOND BURS (REUSABLE) INVERTED CONE 33 1/2
62	MedicalSupplies	DIAMOND BURS (REUSABLE) INVERTED CONE 33 1/4
63	MedicalSupplies	DIAMOND BURS. INVERTED CONE 1/8
64	MedicalSupplies	DIAMOND BURS. ROUND BURS 1
65	MedicalSupplies	DIAMOND BURS. ROUND BURS 1/2 (DENTAL)
66	MedicalSupplies	DIAMOND BURS. ROUND BURS 1/4
67	MedicalSupplies	DIAMOND BURS. ROUND BURS 1/8
68	MedicalSupplies	DIAMOND BURS. STRAIGHT BURS
69	MedicalSupplies	DIAMOND BURS. STRAIGHT BURS 1/2
70	MedicalSupplies	DIAMOND BURS. STRAIGHT BURS 1/4
71	MedicalSupplies	DIAMOND BURS. STRAIGHT BURS 1/8
72	MedicalSupplies	ECG PAPER (SCHILLER)
73	MedicalSupplies	ENDOCLIPS II ML 10MM AUTO SUTURE
74	MedicalSupplies	ETCHANT 3 GRAMS
75	MedicalSupplies	FLOURESCIN DYE 100MG/ML VIAL (12PER PIECES/BOX)
76	MedicalSupplies	FLOURIDE GEL 500ML
77	MedicalSupplies	FLOURIDE VANISH 10ML
78	MedicalSupplies	GLASS SLIDES, PLAIN (FROSTED)
79	MedicalSupplies	GLOVES EXAMINATION S, 100'S POWDER FREE
80	MedicalSupplies	GLOVES EXAMINATION, LARGE 100'S
81	MedicalSupplies	GLOVES EXAMINATION, MEDIUM 100'S
82	MedicalSupplies	GLOVES ORTHO 6.5 50'S STERILE
83	MedicalSupplies	GLOVES ORTHO 7.0 40'S STERILE
84	MedicalSupplies	GLOVES ORTHO 7.5 OR 7 1/2 40'S STERILE
85	MedicalSupplies	GLOVES ORTHO 8 50'S STERILE
86	MedicalSupplies	HYDROFIBER SURGICAL COVER DRESSING WITH SILVER 9X15CM BY 10'S
87	MedicalSupplies	HYDROFIBER SURGICAL COVER DRESSING WITH SILVER 9X25CM BY 10'S BOX

88	MedicalSupplies	HYDROFIBER W/ ADHESIVE SILICON BORDER DRESSING 12.5X12.5 BY 10'S BOX
89	MedicalSupplies	HYDROFIBER/LIPIDOCOLLOID ENHANCED DRESSING WITH SILVER AND STRENGTHENING FIBER
90	MedicalSupplies	LIDOCAINE EPINEPHRINE CARPULE 2%
91	MedicalSupplies	LUBRICANT DIFFUSER 4'S
92	MedicalSupplies	MASK, TIE-ON DISPOSABLE 50'S
93	MedicalSupplies	MATRIX BAND PACKAGE OF 10'S
94	MedicalSupplies	NEEDLE DISPOSABLE G. 26 X 1/2", 100'S
95	MedicalSupplies	NEEDLE DISPOSABLE G. 27 X 1/2
96	MedicalSupplies	NEEDLE DISPOSABLE G. 30 X 1/2
97	MedicalSupplies	NEEDLE, DENTAL G. 30 SHORT
98	MedicalSupplies	NEEDLE, DENTAL G.27, 100S
99	MedicalSupplies	NEEDLE, DISP. G.21 X1", 100'S PUNCTURED PROFF CAPS SHARPS
100	MedicalSupplies	NEEDLE, DISP. G19 X 1 1/2 PUNCTURED CAPS SHARP, 100'S
101	MedicalSupplies	NEEDLE, DISPOSABLE G. 23 X 1
102	MedicalSupplies	NEEDLE, DISPOSABLE G. 25 X 1, 100S PUNCTURED PROOF CAPS SHARPS
103	MedicalSupplies	NEEDLE, DISPOSABLE G.25 X 5/8, 100'S
104	MedicalSupplies	PITS & FISSURE SEALANT
105	MedicalSupplies	PLASMA BIOLOGICAL INDICATOR
106	MedicalSupplies	PLASTER ADHESIVE POROUS TAPE, WHITE, 1" X 10YARDS, HYPOALLERGENIC SILK-LIKE TAPE GOOD ADHESIVE PROPERTIES COMFORTABLE DURABLE 12'S
107	MedicalSupplies	POLISHING AND FINISHING STRIP 100"S
108	MedicalSupplies	POLYURETHANE FOAM WOUND DRESSING WITH 1% 10X10X0.5CM
109	MedicalSupplies	POLYURETHANE FOAM WOUND DRESSING WITH 3% 10X10X0.5CM
110	MedicalSupplies	POLYURETHANE FOAM WOUND DRESSING WITH 3% 20X20X0.5CM
111	MedicalSupplies	PROPHY PASTE 50 GRAMS
112	MedicalSupplies	PUMICE POWDER 200G
113	MedicalSupplies	PUMP COMPLETE TUBING ASSY
114	MedicalSupplies	SALIVA EJECTOR TIP 100'S
115	MedicalSupplies	STREAM INDICATOR STRIP FOR STEAM CLASS PMS
116	MedicalSupplies	SUTURE, RETENSION S-2.0, LR NEEDLE LONG
117	MedicalSupplies	SYRINGE DISP.5CC G. 23 X 1, LUER LOCK
118	MedicalSupplies	SYRINGE DISP. 3CC G23X1, LUER LOCK
119	MedicalSupplies	SYRINGE DISPOSABLE 10CC G23X1, LUER LOCK
120	MedicalSupplies	SYRINGE DISPOSABLE 20CC W/OUT NEEDLE, 50'S/BX
121	MedicalSupplies	SYRINGE DISPOSABLE 30CC W/OUT NEEDLE, 50'S/BX
122	MedicalSupplies	SYRINGE DISPOSABLE 50CC W/OUT NEEDLE, 25'S/BX
123	MedicalSupplies	SYRINGE FOR BLOOD EXTRACTION W/ NEEDLE 10ML
124	MedicalSupplies	SYRINGE FOR BLOOD EXTRACTION W/ NEEDLE 5ML
125	MedicalSupplies	SYRINGE INSULIN 0.5ML SHARP G. 31X6MM
126	MedicalSupplies	SYRINGE TUBERCULINE 1CC W/ NEEDLE G.25 X 5/8, 100'S
127	MedicalSupplies	TULLE GRAS DRESSING IMPREGNATED WITH 0.5 CHLORHEXIDINE ACETATE BP 10X10CM 50PER PIECES/BOX "BACTIGRAS"
128	MedicalSupplies	VAGINAL SPECULUM, DISPOSABLE (LARGE)
129	MedicalSupplies	VAGINAL SPECULUM, DISPOSABLE (MEDIUM)
130	MedicalSupplies	VAGINAL SPECULUM, DISPOSABLE (SMALL)
131	MedicalSupplies	DISINFECTANT SPRAY (CIF)
132	MedicalSupplies	CARESTREAM LASER IMAGING FILM
133	MedicalSupplies	O.R WRAPPING SHEET LIGHT BLUE 100CM X 100CM 160PER PIECES/CARTON

134	MedicalSupplies	O.R WRAPPING SHEET LIGHT BLUE 120CMX120CM
135	MedicalSupplies	O.R WRAPPING SHEET LIGHT BLUE 50CM X 50CM 640PER PIECES/CARTON
136	MedicalSupplies	O.R WRAPPING SHEET LIGHT BLUE 60CM X 60CM 480PER PIECES/CARTON
137	MedicalSupplies	O.R WRAPPING SHEET LIGHT BLUE 75CMX75CM
138	MedicalSupplies	O.R WRAPPING SHEET LIGHT BLUE 90CM X 90CM 160PER PIECES/CARTON
139	MedicalSupplies	ACID CONCENTRATE
140	MedicalSupplies	LOW CALCIUM ACID CONCENTRATES 5 LITER CONTAINER
141	MedicalSupplies	ANTIBACTERIAL BRAIDED SYNTHETIC ABSORBABLE SUTURE SIZE 0 ROUND NEEDLE 36-40MM BOX OF 36 "ETHICON"
142	MedicalSupplies	ANTIBACTERIAL BRAIDED SYNTHETIC ABSORBABLE SUTURE SIZE 1 ROUND NEEDLE 36-40MM BOX OF 36 "ETHICON"
143	MedicalSupplies	BARBED BIDIRECTIONAL SUTURE SIZE 1
144	MedicalSupplies	CATHETER RUBBER FR. 10
145	MedicalSupplies	CATHETER RUBBER FR. 12
146	MedicalSupplies	CATHETER RUBBER FR. 14
147	MedicalSupplies	CATHETER RUBBER FR. 16
148	MedicalSupplies	CATHETER RUBBER FR. 18
149	MedicalSupplies	CATHETER RUBBER FR. 20
150	MedicalSupplies	CATHETER RUBBER FR. 5
151	MedicalSupplies	CATHETER RUBBER FR. 8
152	MedicalSupplies	CHROMIC 3.0 W/ 1/2 ROUND NEEDLE
153	MedicalSupplies	DISPOSABLE GOWN AND DRAPES / ENT PACK
154	MedicalSupplies	DISPOSABLE GOWN AND DRAPES CRANIOTOMY PACK / NEURO PACK CONVERTORS
155	MedicalSupplies	DISPOSABLE GOWN AND DRAPES LAPAROTOMY PACK
156	MedicalSupplies	DISPOSABLE GOWN AND DRAPES ORTHO PACK
157	MedicalSupplies	DRESSING KIT
158	MedicalSupplies	MONOPOLAR NEEDLE ELECTRODES PRO 37 SAF DISPOSABLE FEMALE PLUG PEACH CAP
159	MedicalSupplies	ORAL CARE KIT
160	MedicalSupplies	POVIDONE IODINE PLEDGET (3 / PACK)
161	MedicalSupplies	SANITARY NAPKIN
162	MedicalSupplies	SELF-CARE PACKAGE FOR OB-GYNE PATIENTS
163	MedicalSupplies	SURGICAL GAUZE BALLS (CHERRY BALLS) NON STERILE 5PER PIECES/PACK
164	MedicalSupplies	SURGICAL GAUZE BALLS (CHERRY BALLS) STERILE 5PER PIECES/PACK
165	MedicalSupplies	SURGICAL PEANUT GAUZE STERILE (31/8"X3 1/2) 10PER PIECES/PER PIECEK
166	MedicalSupplies	SURGICAL PRE-CUT GAUZE STERILE 5PER PIECES/PACK X-RAY THREAD DETACTABLE (4X4X8 PLY)
167	MedicalSupplies	VITRECTOMY PACK
168	MedicalSupplies	ECG PAPER, FOR SCHILLER AT-1 (FOR AT-1, Z FOLDED)
169	MedicalSupplies	GLOVES I.E.
170	MedicalSupplies	GLOVES STERILE W/ EXCELLENT COMFORT AND FIT S-6.0
171	MedicalSupplies	GLOVES STERILE W/ EXCELLENT COMFORT AND FIT S-6.5
172	MedicalSupplies	GLOVES STERILE W/ EXCELLENT COMFORT AND FIT S-7.0
173	MedicalSupplies	GLOVES STERILE W/ EXCELLENT COMFORT AND FIT S-7.5
174	MedicalSupplies	GLOVES STERILE W/ EXCELLENT COMFORT AND FIT S-8.0
175	MedicalSupplies	GLOVES STERILE, S-6.0 NON LATEX
176	MedicalSupplies	GLOVES STERILE, S-6.5 NON LATEX
177	MedicalSupplies	GLOVES STERILE, S-7.0 NON LATEX
178	MedicalSupplies	GLOVES STERILE, S-7.5 NON LATEX

179	MedicalSupplies	GLOVES STERILE, S-8.0 NON LATEX
180	MedicalSupplies	GLOVES, NITRILE EXAMINATION, MEDIUM
181	MedicalSupplies	SHOE COVER, DISPOSABLE
182	MedicalSupplies	0.9% SODIUM CHLORIDE SOLUTION FOR IRRIGATION
183	MedicalSupplies	10-12MM BURRHOLE TYPE PUDENZ LOW PRESSURE SHUNT "MEDTRONIC"
184	MedicalSupplies	10-12MM TROCAR BLADED
185	MedicalSupplies	10MM TROCAR, BLADELESS
186	MedicalSupplies	10MM TROCAR, DILATING TIP
187	MedicalSupplies	12MM TROCAR, BLADELESS
188	MedicalSupplies	12MM TROCAR, DILATING TIP
189	MedicalSupplies	15MM TROCAR BLADELESS
190	MedicalSupplies	16MM BURRHOLE TYPE PUDENZ LOW PRESSURE SHUNT CONTOUR TYPE LOW PRESSURE SHUNT "MEDTRONIC"
191	MedicalSupplies	16MM BURRHOLE TYPE PUDENZ MEDIUM PRESSURE SHUNT "MEDTRONIC"
192	MedicalSupplies	2-3 LAYER POLYTETRAFLUOROETHYLENE (PTFE) VASCULAR GRAFT 6MMX50CM STANDARD
193	MedicalSupplies	2-3 LAYER POLYTETRAFLUOROETHYLENE (PTFE) VASCULAR GRAFT 6MMX80CM STANDARD
194	MedicalSupplies	2-3 LAYER POLYTETRAFLUOROETHYLENE (PTFE) VASCULAR GRAFT 8MMX50CM STANDARD
195	MedicalSupplies	2-3 LAYER POLYTETRAFLUOROETHYLENE (PTFE) VASCULAR GRAFT 8MMX80CM STANDARD
196	MedicalSupplies	3M N95 PARTICULATE RESPIRATORY MASK, REGULAR
197	MedicalSupplies	5MM TROCAR 5MM BLADED
198	MedicalSupplies	5MM TROCAR, BLADELESS
199	MedicalSupplies	5MM TROCAR, DILATING TIP
200	MedicalSupplies	9CM BALL 3.0MM DIAMOND BURR
201	MedicalSupplies	9CM BALL 4.0MM DIAMOND BURR
202	MedicalSupplies	9CM BALL 6.0MM CUTTING BURR "MEDTRONIC"
203	MedicalSupplies	9CM BALL CUTTING BURR 7.5MM "MEDTRONIC"
204	MedicalSupplies	ABBOCATH GA.16
205	MedicalSupplies	ABBOCATHGA.18
206	MedicalSupplies	ABBOCATHGA.20
207	MedicalSupplies	ABBOCATHGA.22
208	MedicalSupplies	ABDOMINAL BINDER
209	MedicalSupplies	ABSORBABLE GELATIN SPONGE HAEMOSTATIC 7CMX5CMX1CM
210	MedicalSupplies	CORD COVER ADHESIVE STRIPS
211	MedicalSupplies	CORE NEEDLE BIOPSY GAUGE 16
212	MedicalSupplies	CORE NEEDLE BIOPSY GAUGE 18
213	MedicalSupplies	COTTON BUDS (PAPER)X 108
214	MedicalSupplies	COTTON BUDS200 TIPS
215	MedicalSupplies	COUNTOUR CURVED 40MM
216	MedicalSupplies	COUNTOUR CURVED RELOAD 40MM
217	MedicalSupplies	COVER ALL
218	MedicalSupplies	COVERALL WITH HOOD, DISPOSABLE
219	MedicalSupplies	CR FILM 10 X 12 (D-HL)
220	MedicalSupplies	CR FILM 14 X 17 (D-HL)
221	MedicalSupplies	CR FILM 8 X 10 (D-HL)
222	MedicalSupplies	CSTD ADAPTER C100
223	MedicalSupplies	CSTD CONNECTOR LUER LOCK C45
224	MedicalSupplies	CSTD INJECTOR LUER LOCK N35

225	MedicalSupplies	CSTD PROTECTOR P50
226	MedicalSupplies	CT FILM ENVELOPE W/LOGO 14.5 X 17.5 INCH
227	MedicalSupplies	CURVED CUTTER STAPLER CARTRIDGE
228	MedicalSupplies	CURVED SHEAR SHAFT DIAMEER 5MM
229	MedicalSupplies	CUTTING LOOP ANGLED 24/26 IN COLOR YELLOW (SINGLE STEM)
230	MedicalSupplies	CVP MANOMETER (ANESTHESIA USE)
231	MedicalSupplies	CYTO BRUSH FOR PAP SMEAR
232	MedicalSupplies	DECONTAMINATION MAT 90 X 115CM W/BENZOILSOTHIASOLINE ANTIBACTERIAL AGENT
233	MedicalSupplies	DIALYZER LOPS 18
234	MedicalSupplies	DIALYZER, CELLULOSE HIGH EFFICIENCY AM-SD-750U
235	MedicalSupplies	DIALYZER, CELLULOSE MODERATE EFFICIENCY AM-SD-650U
236	MedicalSupplies	DIAMOND BURS (REUSABLE) ROUND 1
237	MedicalSupplies	DIAMOND BURS (REUSABLE) ROUND 1/2
238	MedicalSupplies	DIAMOND BURS ST. FISSURE S555
239	MedicalSupplies	DIAMOND BURS ST. FISSURE S556
240	MedicalSupplies	DIAPER LARGEADULT
241	MedicalSupplies	DIAPER NEWBORN
242	MedicalSupplies	DIAPER XLADULT
243	MedicalSupplies	DIAPER, ADULT (LARGE)
244	MedicalSupplies	DIAPER, ADULT (MEDIUM)
245	MedicalSupplies	DIGITAL, THERMOMETER AXILLARY
246	MedicalSupplies	DISINFECTANT HAND RUB (STERILLIUM)
247	MedicalSupplies	DISINFECTANT WIPES (BACIOL) PER SHEET
248	MedicalSupplies	DISPOSABLE DIGITAL ORAL THERMOMETER
249	MedicalSupplies	FRESHBATHROOM TISSUE PAPER X 8 X12
250	MedicalSupplies	FULL MASK LARGE (2 VENTED 1 NON VENTED)
251	MedicalSupplies	FULL MASK MEDIUM (6 VENTED 4 NON VENTED)
252	MedicalSupplies	FULL MASK SMALL (1 VENTED 1 NON VENTED)
253	MedicalSupplies	GAUZE 2X2 (2PCS/PACK)
254	MedicalSupplies	GAUZE 2X2 (5PCS/PACK)
255	MedicalSupplies	GAUZE 4X16 (2PCS/PACK)
256	MedicalSupplies	GAUZE 4X4 (5PCS/PACK)
257	MedicalSupplies	GAUZE 4X4X2
258	MedicalSupplies	GAUZE 4X8 (5PCS/PACK)
259	MedicalSupplies	GAUZE BANDAGE ROLLED 3" X 10
260	MedicalSupplies	GAUZE PAD 4X8 (2PCS/PACK)
261	MedicalSupplies	GAUZE PAD 4X8 (5PCS/PACK)
262	MedicalSupplies	GAUZE PAD 8X16 (2PCS/PACK)
263	MedicalSupplies	GAUZE STERILE 4X4
264	MedicalSupplies	GAUZE VASELINIZED 10X10
265	MedicalSupplies	GIGLI GUIDE SHIVA CHINA
266	MedicalSupplies	GIGLI WIRE 16" THICK COMPATIBLE FOR AMPUTATION (30X33 CM ALUMINUM HANDLE)
267	MedicalSupplies	GIGLI WIRE 20" THICK COMPATIBLE FOR AMPUTATION (30X33 CM ALUMINUM HANDLE)
268	MedicalSupplies	GIGLI WIRE DE MARTEL, 5X6 AESCULAP
269	MedicalSupplies	GLASS SLIDES
270	MedicalSupplies	GLASS TUBING 7MM BREAKABLE, AUTOCLAVABLE, REUSABLE
271	MedicalSupplies	GLOVES - SURGICAL 7

272	MedicalSupplies	GLOVES - SURGICAL 7 1/2
273	MedicalSupplies	GLOVES I.E. LARGE
274	MedicalSupplies	GLOVES- NITRRILE EXAMINATION,LARGE
275	MedicalSupplies	GLOVES, CHEMO, NON-LATEX 7.5
276	MedicalSupplies	GLUCOSE TEST STRIPS/TIP (INDIVIDUALLY WRAPPED) W/ 50 FREE AUTO CODED GLUCOMETER, 50 FREE EXTRA BATTERY, FREE LANCET EQUAL TO QTY. PURCHASED
277	MedicalSupplies	GOWN FOR PPE
278	MedicalSupplies	GOWN, DISPOSABLE STERILE
279	MedicalSupplies	GRAVITY FMS PACK, BASIC
280	MedicalSupplies	GRAVITY FMS PACK, WITH TIP
281	MedicalSupplies	GREEN FINE STONE
282	MedicalSupplies	GROUND PLATE ELECTRODES 48 FOR EMG-NCV
283	MedicalSupplies	GROUNDING/CABLE FOR CAUTERY PLATE
284	MedicalSupplies	HARMONIC ACE
285	MedicalSupplies	HARMONIC FOCUS PLUS 17CM
286	MedicalSupplies	HARMONIC FOCUS PLUS 9CM
287	MedicalSupplies	HEAT AND MOISTURE EXCHANGER
288	MedicalSupplies	HEAT MOISTURE EXCHANGER (HME)
289	MedicalSupplies	HEATED EXPIRATORY FILTER VENTILATOR ADULT
290	MedicalSupplies	HEATED EXPIRATORY FILTER VENTILATOR NEONATE
291	MedicalSupplies	HEMOCLIP, FOR EGD, ROTATABLE
292	MedicalSupplies	HEMOPERFUSION CARTRIDGE FOR SEPSIS CYTOKINE
293	MedicalSupplies	HEMORRHOID RUBBER BAND FOR LIGATION
294	MedicalSupplies	HEMORRHOIDAL CIRCULAR STAPLER 33MM
295	MedicalSupplies	HEMOSTAT, ABSORBABLE 1"*1" (SURGICEL)
296	MedicalSupplies	HEMOSTAT, ABSORBABLE 3"*4" (SURGICEL)
297	MedicalSupplies	HEMOSTAT, ABSORBABLE SPONGE (SPONGOSTAN)
298	MedicalSupplies	HEMOSTATIC OXIDIZED REGENERATED CELLULOSE 3 X 4
299	MedicalSupplies	HEPARIN
300	MedicalSupplies	HEPARIN LOCK
301	MedicalSupplies	HEPARIN LOCK (NEEDLE LESS)
302	MedicalSupplies	HEPLOCK
303	MedicalSupplies	HEPLOCK/IN STOPPERSET
304	MedicalSupplies	HIGH FLOW OXYGEN CANNULA
305	MedicalSupplies	HP BLUE HANDPIECE 100 USES FOR OPEN SHEARS HARMONIC FOCUS
306	MedicalSupplies	HP BOWIE & DICK TEST PACK
307	MedicalSupplies	HYDROCOLLOID DRESSING 4X4 CGF
308	MedicalSupplies	HYDROCOLLOID DRESSING 4X4 EXTRA THIN
309	MedicalSupplies	HYDROFIBER DRESSING WITH SILVER 4X4
310	MedicalSupplies	HYDROXYPROPYL METHYLCELLULOSE 2.0 %
311	MedicalSupplies	I.E GLOVESDISP.
312	MedicalSupplies	I.V. ADMINISTRATION SET FOR BLOOD
313	MedicalSupplies	I.V. ADMINISTRATION SET, ADULT (MACROSET)
314	MedicalSupplies	I.V. ADMINISTRATION SET, PEDIA (MICROSET)
315	MedicalSupplies	ICE BAG, 100% LEAK FREE RESEALABLE CAP
316	MedicalSupplies	ICP PROBE 3 XL
317	MedicalSupplies	ICP PROBE 3PN INTRAPHARENCHYMAL
318	MedicalSupplies	ICP TUNNELING SET PROBE 3

319	MedicalSupplies	IJ CATHETER FR. 12
320	MedicalSupplies	IMPLANTABLE VASCULAR ACCESS DEVICE PORT MEASUREMENT 12X13.2X35.27MM TITANIUM AND POLYSULPHONE MATERIAL
321	MedicalSupplies	IMPLANTABLE VASCULAR ACCESS DEVICE PORT MEASUREMENT 9.5X10.6X30X22MM TITANIUM AND POLYSULPHONE MATERIAL
322	MedicalSupplies	INFUSABLE PRESSURE/INFUSOR PUMP LARGE
323	MedicalSupplies	INFUSABLE PRESSURE/INFUSOR PUMP SMALL
324	MedicalSupplies	INFUSABLE PRESSURE/INFUSOR PUMP WITH INFLATE HOLD AIR DEFLATE CONTROL
325	MedicalSupplies	INFUSION SET/VENOFIXGAUGE 23
326	MedicalSupplies	INFUSOMAT SPACE LINE (DOSIFIX)NEEDLE
327	MedicalSupplies	INFUSOMAT SPACE LINE TYPE IV STANDAR
328	MedicalSupplies	INK CARTRIDGE
329	MedicalSupplies	INSULIN SYRINGE1ML
330	MedicalSupplies	INSYTE CATHETERGA.20
331	MedicalSupplies	INSYTE CATHETERGA.22
332	MedicalSupplies	INSYTE CATHETERGA.24
333	MedicalSupplies	INTER FOLDED TISSUE
334	MedicalSupplies	INTERNAL JUGULAR HEMODIALYSIS CATH. FR.13.5 X 16, DOUBLE LUMEN WITH LASER-CUT SIDE SLOTS, SYMMETRICAL & GREEN ATRIUMATIC TIP DESIGN, DEPTH MARKINGS, POLYURETHANE MATERIAL CURVED EXTENSIONS " MAHURKAR"
335	MedicalSupplies	INTERNAL JUGULAR HEMODIALYSIS CATHETER 11-12 FR. 16-25CM DOUBLE LUMEN
336	MedicalSupplies	INTERNAL JUGULAR HEMODIALYSIS CATHETER 12 FR. 16CM DOUBLE LUMEN
337	MedicalSupplies	INTERNAL JUGULAR HEMODIALYSIS CATHETER 12 FR. 20CM DOUBLE LUMEN
338	MedicalSupplies	INTERNAL JUGULAR HEMODIALYSIS CATHETER 12.5 FR. X 16, TRIPLE LUMEN
339	MedicalSupplies	INTERNAL JUGULAR HEMODIALYSIS CATHETER 12-13.5 FR. X 20CM TRIPLE LUMEN
340	MedicalSupplies	INTERNAL JUGULAR HEMODIALYSIS CATHETER 13.5 FR. 13CM DOUBLE LUMEN
341	MedicalSupplies	INTERNAL JUGULAR HEMODIALYSIS CATHETER, 11-12.5FR. X 24CM, DOUBLE LUMEN, LASER-CUT SIDE SLOTS, SYMMETRICAL & GREEN ATRIUMATIC TIP DESIGN, DEPTH MARKINGS, POLYURETHANE MATERIAL, STRAIGHT EXTENSIONS "MAHURKAR"
342	MedicalSupplies	INTERNAL JUGULAR HEMODIALYSIS CATHETER, 12-13.5FR. X 24CM, DOUBLE LUMEN, LASER-CUT SIDE SLOTS, SYMMETRICAL & GREEN ATRIUMATIC TIP DESIGN, DEPTH MARKINGS, POLYURETHANE MATERIAL, STRAIGHT EXTENSIONS "MAHURKAR"
343	MedicalSupplies	INTERNAL JUGULAR HEMODIALYSIS CATHETER, 12FR X 13CM, DOUBLE LUMEN, LASER-CUT SIDE SLOTS, SYMMETRICAL & GREEN ATRIUMATIC TIP DESIGN, DEPTH MARKINGS, POLYURETHANE MATERIAL, CURVED EXTENSIONS
344	MedicalSupplies	INTERNAL JUGULAR HEMODIALYSIS CATHETER, 13-13.5FR. X 24-25CM, DOUBLE LUMEN, LASER-CUT SIDE SLOTS, POLYURETHANE MATERIAL, STRAIGHT SILICON EXTENSIONS "MAHURKAR"
345	MedicalSupplies	INTRAOSSEOUS 45 MM ADULT NEEDLE WITH DRIVER
346	MedicalSupplies	INTRAOSSEOUS NEEDLE-ADULT
347	MedicalSupplies	ITO ES 320 ELECTRODE CABLE (BLUE)
348	MedicalSupplies	ITO ES 320 ELECTRODE CABLE (GRAY)
349	MedicalSupplies	IUD KITS (INTERVAL AND POST-PARTUM)
350	MedicalSupplies	IV CATHETER G. 16
351	MedicalSupplies	IV CATHETER G. 18
352	MedicalSupplies	IV CATHETER G. 20
353	MedicalSupplies	IV CATHETER G. 22
354	MedicalSupplies	IV CATHETER G. 24
355	MedicalSupplies	IV CATHETER G. 26
356	MedicalSupplies	IV EXTENSION TUBING 140-145CM

357	MedicalSupplies	IV EXTENSION TUBING 150CM
358	MedicalSupplies	IV EXTENSION TUBING 75CM
359	MedicalSupplies	IV EXTENSION TUBING 90-120CM
360	MedicalSupplies	IV INFUSION SET EXTENSION 100CM
361	MedicalSupplies	IV PUMP SET FOR PER PIECEA MACHINE "GEMSTAR"
362	MedicalSupplies	IV SET FOR INFUSION GRAVITY, W/LUER LOCK (IV INFUSION PUMP)
363	MedicalSupplies	IV START PACK
364	MedicalSupplies	IV STARTER KIT
365	MedicalSupplies	J GUIDE WIRE 0.38" (EXCHANGE WIRE)
366	MedicalSupplies	JACKSON REES 25CM TUBING W/ FLOW TROLL BLEED VALUE 150CM W/ 1L BAG
367	MedicalSupplies	JELLY, LUBRICATING 80GMS CLEAR, STERILE, WATER, SOLUBLE
368	MedicalSupplies	KELLY PAD W/ INFLATOR
369	MedicalSupplies	KMC BLOUSE LARGE
370	MedicalSupplies	KMC BLOUSE MEDIUM
371	MedicalSupplies	KMC BLOUSE XL
372	MedicalSupplies	KN-95 DONATION
373	MedicalSupplies	KNEE IMMOBILIZER, LARGE SIZE
374	MedicalSupplies	KNEE IMMOBILIZER, MEDIUM SIZE
375	MedicalSupplies	KNEE IMMOBILIZER, SMALL SIZE
376	MedicalSupplies	KNOTLESS TISSUE CONTROL DEVICE SIZE 0, CUTTING
377	MedicalSupplies	KNOTLESS TISSUE CONTROL DEVICE SIZE 1, CUTTING
378	MedicalSupplies	L SHAPED COAGULATING AND DESSECTING ELECTRODE WITH CONNECTOR PIN SIZE 5MM
379	MedicalSupplies	LACRIMAL TUBE, 7.5CM PROBE-23G, REGULAR TIP, SILICON
380	MedicalSupplies	LAPAROSCOPIC LINEAR CUTTING STAPLER RELOAD 45MM WITH TRI-STAPLE TECHNOLOGY
381	MedicalSupplies	LAPAROSCOPIC LINEAR CUTTING STAPLER RELOAD 45MM WITH TRI-STAPLE TECHNOLOGY, VASCULAR, MEDIUM
382	MedicalSupplies	LAPAROSCOPIC LINEAR CUTTING STAPLER RELOAD 45MM WITH TRI-STAPLE TECHNOLOGY, VASCULAR, MEDIUM, CURVED TIP
383	MedicalSupplies	LAPAROSCOPIC LINEAR CUTTING STAPLER RELOAD 60MM WITH TRI-STAPLE TECHNOLOGY
384	MedicalSupplies	LAPAROSCOPIC LINEAR CUTTING STAPLER STANDARD LENGTH COMPATIBLE WITH TRI-STAPLE TECHNOLOGY
385	MedicalSupplies	LAPAROSCOPIC LINEAR RELOAD 45MM
386	MedicalSupplies	LAPAROSCOPIC LINEAR RELOAD 60MM
387	MedicalSupplies	LAPAROSCOPIC LINEAR STAPLER 45MM
388	MedicalSupplies	LAPAROSCOPIC LINEAR STAPLER RELOAD 60MM "COVIDIEN"
389	MedicalSupplies	LARYNGEAL MASK AIRWAY SIZE 2
390	MedicalSupplies	LARYNGEAL MASK AIRWAY SIZE 3
391	MedicalSupplies	LARYNGEAL MASK AIRWAY SIZE 5
392	MedicalSupplies	LEGEND 8CM 2.3MM TAPERED
393	MedicalSupplies	LEGEND 9CM, 6MM
394	MedicalSupplies	LENS, ANTERIOR CHAMBER, INTRAOCULAR, POLYMETHYMETACRYLATE, (FDA APPROVED) RIGID, RANGE: 10.0D TO 26.0D (21.0D)
395	MedicalSupplies	LENS, POSTERIOR CHAMBER, INTRAOCULAR, HYDROPHILIC ACRYLIC SQUARE EDGE, 12.5MM LENGTH, OPTIC SIZE 6.0MM, EQUICONVEX, A CONSTANT 118.0. FDA APPROVED, FOLDABLE, RANGE: 1.0D TO 30.0D (20.0D)
396	MedicalSupplies	LIGATING CLIPS TITANIUM CLIPS, SMALL 30-36 CART./BOX 6CLIPS PER CART. WITH CLIP APPLIER.
397	MedicalSupplies	LIGATING CLIPS, MED-LARGE WITH LOCKING MECHANISM 6 CLIPS/14 CARTRIDGE/BOX W/ APPLIER "HOMOLOK LIGATING CLIPS, TELEFLEX USA"
398	MedicalSupplies	LINEAR CUTTER 100MM GIA AUTO SUTURE

399	MedicalSupplies	LINEAR CUTTER 80MM RELOAD (GREEN) MDR-01051
400	MedicalSupplies	LINEAR CUTTING RELOAD 100MM
401	MedicalSupplies	LINEAR CUTTING RELOAD 55MM
402	MedicalSupplies	LINEAR CUTTING RELOAD 60MM
403	MedicalSupplies	LINEAR CUTTING RELOAD 75MM
404	MedicalSupplies	LINEAR CUTTING STAPLER 100MM
405	MedicalSupplies	LINEAR CUTTING STAPLERS 55MM
406	MedicalSupplies	LINEAR CUTTING STAPLERS 60MM
407	MedicalSupplies	LINEAR CUTTING STAPLERS 75MM
408	MedicalSupplies	LINEAR CUTTING STAPLERS 80MM
409	MedicalSupplies	LINEAR RELOAD 45MM "COVIDIEN"
410	MedicalSupplies	LINEAR STAPLER 45MM TO 454.8
411	MedicalSupplies	LMA PROSEAL S. 1 LATEX FREE HYPOALLERGENIC DISP.
412	MedicalSupplies	LMA PROSEAL S. 1.5 LATEX FREE HYPOALLERGENIC DISP.
413	MedicalSupplies	LMA PROSEAL S. 2 1/2 LATEX FREE HYPOALLERGENIC DISP.
414	MedicalSupplies	LMA PROSEAL S. 2 LATEX FREE HYPOALLERGENIC DISP.
415	MedicalSupplies	LMA PROSEAL S. 3 LATEX FREE HYPOALLERGENIC DISP.
416	MedicalSupplies	LMA PROSEAL S. 4 LATEX FREE HYPOALLERGENIC DISP.
417	MedicalSupplies	LMA PROSEAL S. 5 LATEX FREE HYPOALLERGENIC DISP.
418	MedicalSupplies	LOW FLUX DIALYZER OR ITS EQUIVALENT
419	MedicalSupplies	LUBRICANT DIFFUSER FOR PNEUMATIC DRILL
420	MedicalSupplies	LUER LOCK ADAPTER
421	MedicalSupplies	LYNESTRENOL 500MCG
422	MedicalSupplies	LYSOL SPRAY (FOR BABYS ROOM SCENT)340G
423	MedicalSupplies	LYSOL SPRAY340G
424	MedicalSupplies	MACROSET
425	MedicalSupplies	MACROSET/ VENOSET/INTRAFIXADULT SET
426	MedicalSupplies	MALE BEDSIDE URINAL
427	MedicalSupplies	MALE CONDOM
428	MedicalSupplies	MANUAL VACUUM ASPIRATOR ENDOMETRIAL BIOPSY KIT
429	MedicalSupplies	MASK DISPOSABLE EARLOOP
430	MedicalSupplies	MASK, FLEX, SIZE 1
431	MedicalSupplies	MASK, FLEX, SIZE 2
432	MedicalSupplies	MASK, FLEX, SIZE 4
433	MedicalSupplies	MASK, FLEX, SIZE 5
434	MedicalSupplies	MDI CONNECTION PORT 22M-22F
435	MedicalSupplies	MEDICAL GRADE ISOLATION SUIT
436	MedicalSupplies	MEDICINE CAP 30ML (PLASTIC)
437	MedicalSupplies	MEDIPLAST BAND-AIDBANTAM
438	MedicalSupplies	MEDROXYPROGESTERONE ACETATE 50MG/3ML
439	MedicalSupplies	MESH 7.5CM X 15CM (LIGHT WEIGHT)
440	MedicalSupplies	MESH POLYPROPYLENE 6"X6"
441	MedicalSupplies	MG350 (FOR COVID19- SEPSIS ACUTE PATIENT)
442	MedicalSupplies	MICROSET
443	MedicalSupplies	MICROSETPEDIA SET
444	MedicalSupplies	MOUTH MIRROR HEAD #5
445	MedicalSupplies	MUCUS EXTRACTOR
446	MedicalSupplies	N POLYURETHANE FOAM WOUND DRESSING

447	MedicalSupplies	N-95, MASK
448	MedicalSupplies	NAIL POLISH REMOVER30ML
449	MedicalSupplies	NAIL POLISH REMOVER60ML
450	MedicalSupplies	NAME TAG- ADULT WHITE
451	MedicalSupplies	NANOCRYSTALLINE SILVER DRESSING
452	MedicalSupplies	NASAL AIRWAY, S-26
453	MedicalSupplies	NASAL AIRWAY, S-28
454	MedicalSupplies	NASAL AIRWAY, S-30
455	MedicalSupplies	NASAL CANNULA (MEDIUM) FOR HIGH FLOW
456	MedicalSupplies	NASAL CANNULA 2 PRONG ADULT,
457	MedicalSupplies	NASAL CANNULA 2 PRONG NEONATE
458	MedicalSupplies	NASAL CANNULA 2 PRONG PEDIA
459	MedicalSupplies	NASAL CANNULA 2 PRONG PRE TERM
460	MedicalSupplies	NASO GASTRIC TUBE FR 10 SILICON
461	MedicalSupplies	NASO GASTRIC TUBE FR 14 SILICON
462	MedicalSupplies	NASO GASTRIC TUBE FR 16 SILICON
463	MedicalSupplies	NASO GASTRIC TUBE FR 18 SILICON
464	MedicalSupplies	NASO GASTRIC TUBE, ADULT (FR. 14)
465	MedicalSupplies	NASOPHARYGEAL AIRWAY, LATEX (NASAL)
466	MedicalSupplies	NEBULIZER FILTER
467	MedicalSupplies	NEBULIZING KIT COMPLETE SET W/T-PIECE & MOUTHPIECE
468	MedicalSupplies	NEBULIZING KIT WITH MASK (ADULT)
469	MedicalSupplies	NEBULIZING KIT WITH MASK (PEDIA)
470	MedicalSupplies	NEBULIZING KIT WITH OXYGEN TUBING TEE ADAPTER
471	MedicalSupplies	NEEDLE CORE BIOPSY (BARD MAX-CORE) G14
472	MedicalSupplies	NEEDLE CORE BIOPSY G. 16/15CM LENGTH
473	MedicalSupplies	NEEDLE CORE BIOPSY G. 18 X 15CM
474	MedicalSupplies	NEEDLE CORE BIOPSY G. 18 X 20CM
475	MedicalSupplies	NEEDLE CORE BIOPSY G. 18/25CM LENGTH
476	MedicalSupplies	NEEDLE CORE BIOPSY G. 20/15CM LENGTH
477	MedicalSupplies	NEEDLE HYPO G 23
478	MedicalSupplies	NEEDLE, HYPO G. 25
479	MedicalSupplies	NEEDLE NERVE STIMULATOR
480	MedicalSupplies	NEEDLE, SPINAL G. 25
481	MedicalSupplies	NEEDLE SPINAL G.23
482	MedicalSupplies	NEEDLE SPINAL G-26
483	MedicalSupplies	NEEDLE, BIOPSY CORE G-14 PROGRAMMABLE (RADIOLOGY)
484	MedicalSupplies	NEEDLE, BUTTERLY, G 23
485	MedicalSupplies	NEEDLE, CHIBA G. 18 X 20CM
486	MedicalSupplies	NEEDLE, CHIBA G. 20 X 15CM
487	MedicalSupplies	NEEDLE, CHIBA G. 21 X 15CM
488	MedicalSupplies	NEEDLE, CHIBA G. 22 X 15CM
489	MedicalSupplies	NEEDLE, DISPOSABLE G. 18 X 3.5
490	MedicalSupplies	NEEDLE, MULTISAMPLE BD
491	MedicalSupplies	NEEDLE, SPINAL G. 25 X 3.5
492	MedicalSupplies	NEEDLE, SPINAL G. 26 X 3.5
493	MedicalSupplies	NEEDLE, SPINAL G.18 X 3.5 TOP
494	MedicalSupplies	NEEDLE, SPINAL G20/89MM 3.5" LENGTH

495	MedicalSupplies	NEEDLE, SPINAL G21/89MM 3.5" LENGTH
496	MedicalSupplies	NEEDLE, SPINAL G22/89MM 3.5" LENGTH
497	MedicalSupplies	NEEDLE, TEMNO BIOPSY 22 X 15CM
498	MedicalSupplies	NEEDLE, TOUGHY G. 18
499	MedicalSupplies	NEEDLES HYPO G 18
500	MedicalSupplies	NEEDLES HYPO G 30
501	MedicalSupplies	NEGATIVE WOUND PRESSURE THERAPHY MACHINE DISP. CANISTER 600CC
502	MedicalSupplies	NEPHROSTOMY KIT - TUBE FR. 10
503	MedicalSupplies	NEWBORN HEARING SCREENING REGISTRY CARD
504	MedicalSupplies	NKFS-02 LIQUID CONCENTRATION (ACID) 5L
505	MedicalSupplies	NON ADHESIVE HYDROCELLULAR FOAM DRESSING 10CMX10CM
506	MedicalSupplies	NON INVASIVE MASK MEDIUM (NON VENTED)
507	MedicalSupplies	NON INVASIVE MASK MEDIUM (VENTED)
508	MedicalSupplies	NON REBREATHER MASK ADULT
509	MedicalSupplies	PREGNANCY KIT
510	MedicalSupplies	PROBE COVER OR NON LUBRICATED CONDOM TRANSVAGINAL
511	MedicalSupplies	RADIAL INTRODUCER KIT FR 5 - TERUMO
512	MedicalSupplies	RADIAL INTRODUCER KIT FR 6 - TERUMO
513	MedicalSupplies	RAE TUBE ORAL WITH CUFF SIZE 4.0
514	MedicalSupplies	RAE TUBE ORAL WITH CUFF SIZE 4.5
515	MedicalSupplies	RAE TUBE ORAL WITH CUFF SIZE 6.5
516	MedicalSupplies	RECOMBINANT HUMAN GROWTH FACTOR SPRAY 10ML
517	MedicalSupplies	REINFORCED ADHESIVE SKIN CLOSURE, 1/2 "X4"
518	MedicalSupplies	RIGID LENS POSTERIOR INTRAOCULAR POWER 16.0
519	MedicalSupplies	RIGID LENS POSTERIOR INTRAOCULAR POWER 19.5
520	MedicalSupplies	RIGID LENS POSTERIOR INTRAOCULAR POWER 20.5
521	MedicalSupplies	RIGID LENS POSTERIOR INTRAOCULAR POWER 22.5
522	MedicalSupplies	RIGID LENS POSTERIOR INTRAOCULAR POWER 23.0
523	MedicalSupplies	RUBBER CONDOM
524	MedicalSupplies	SAFETY GOGGLES, POLYCARBONATE MATERIAL, HARD COAT LENSES
525	MedicalSupplies	SALIENT POWER INJECTOR KIT
526	MedicalSupplies	SAMPLE PUMP TUBE FOR 248 ABG MACHINE
527	MedicalSupplies	SAMPLE PUMP TUBING KIT
528	MedicalSupplies	SECURITY T-IUD TCU-380A
529	MedicalSupplies	SHARP/NEEDLE CONTAINER PLASTIC 1.4 QT
530	MedicalSupplies	SHAVER
531	MedicalSupplies	SHOE COVER
532	MedicalSupplies	SILICON BAND #242, 4.0MM WIDTH (OFFER:42)
533	MedicalSupplies	SILICON BUCKLE, #276, 9MM WIDTH
534	MedicalSupplies	SILICON OIL 1000 CENTISTOKES
535	MedicalSupplies	SILICON OIL 5000 CENTISTOKES
536	MedicalSupplies	SILICON OIL, ADATO, SYRINGE PACK, WITH TIPS G23 CX5720
537	MedicalSupplies	SILICONE GASTROSTOMY TUBE FR.14
538	MedicalSupplies	SILICONE GASTROSTOMY TUBE FR.16
539	MedicalSupplies	SILICONE OIL 1000 CENTIPOISE SPECIFIC GRAVITY 0.97 FRATIONATED 10ML
540	MedicalSupplies	SILICONE OIL 5000 CENTIPOISE SPECIFIC GRAVITY 0.97 FRATIONATED 10ML
541	MedicalSupplies	SILK 0, STRANDS

542	MedicalSupplies	SILK 1 STRANDS
543	MedicalSupplies	SILK 2-0 CUTTING
544	MedicalSupplies	SILK 2-0 CUTTING (WITH NEEDLE)
545	MedicalSupplies	SILK 2-0 ROUND
546	MedicalSupplies	SILK 2-0 ROUND NEEDLE
547	MedicalSupplies	SILK 2-0 STRANDS
548	MedicalSupplies	SILK 3-0 CUTTING NEEDLE
549	MedicalSupplies	SILK 3-0 ROUND NEEDLE
550	MedicalSupplies	SILK 3-0, STRANDS
551	MedicalSupplies	SILK 4-0 (W/O NEEDLE)
552	MedicalSupplies	SILK 4-0 ROUND NEEDLE
553	MedicalSupplies	SILK 4-0 W/ CUTTING NEEDLE
554	MedicalSupplies	SILK 5-0 CUTTING
555	MedicalSupplies	SILK 5-0 ROUND
556	MedicalSupplies	SILK 6-0 CUTTING (SINGLE ARM)
557	MedicalSupplies	SILK 6-0 DOUBLE ARMED
558	MedicalSupplies	SILVER HANDPIECE 95 USES FOR LAP SHEARS
559	MedicalSupplies	SINGLE CONVENTIONAL
560	MedicalSupplies	SKIN MARKER 4 IN 1 STERILE
561	MedicalSupplies	SKIN STAPLER
562	MedicalSupplies	SKIN STAPLER (6.9MMX 4.2)
563	MedicalSupplies	SKIN STAPLER DISPOSABLE W/STAPLE WIRE REMOVER
564	MedicalSupplies	SKIN STAPLER REMOVER, METAL
565	MedicalSupplies	SNARE: STANDARD OVAL SNARE:BRAIDED LOOP: 3 X 6 CM: FR 7 240CM
566	MedicalSupplies	SO2 CALIBRATOR (EQUIVALENT)
567	MedicalSupplies	SODIUM HYALURONATE 1.8% (CURAMED)
568	MedicalSupplies	SOLUSET
569	MedicalSupplies	SOLUSET (VOLUMETRIC INFUSION)
570	MedicalSupplies	SOLUSET SET120ML
571	MedicalSupplies	SOLUSET W/ CLEAR CALIBRATED CYLINDER, CLEAR TUBING, 100ML
572	MedicalSupplies	SOLUTION REF ELECTRODE M614
573	MedicalSupplies	SPINAL MANOMETER
574	MedicalSupplies	SPINAL NEEDLE G.21
575	MedicalSupplies	STERI STRIP
576	MedicalSupplies	STERILE OPHTHALMIC VISCOELASTIC GEL
577	MedicalSupplies	STERILE WATER 1LITER
578	MedicalSupplies	STERILIZING FLAT POUCH 150 X 350M
579	MedicalSupplies	STIMULATING NEEDLE SINGLE SHOT - 50MM
580	MedicalSupplies	STOCKINETTE BANDAGE SIZE 4 X 25M
581	MedicalSupplies	STOCKINETTE BANDAGE SIZE 6 X 25M
582	MedicalSupplies	STOCKINETTE BANDAGE SIZE 8 X 25M
583	MedicalSupplies	STONE BASKET FM, RUSCH STONE EXTRATOR, ROUND WIRE
584	MedicalSupplies	STYLE 240 (BAND)
585	MedicalSupplies	STYLE 270 (SLEEVE)
586	MedicalSupplies	STYLE 278 (TIRE)
587	MedicalSupplies	SUCTION AND IRRIGATOR SYSTEM 5MM
588	MedicalSupplies	SUCTION CONNECTING TUBE
589	MedicalSupplies	SUCTION SET. YANKAUER HANDLE

590	MedicalSupplies	SURGICAL BLADE NO. 15, UNSHARPENED BACK EDGE
591	MedicalSupplies	SURGICAL BLADE SIZE 10
592	MedicalSupplies	SURGICAL BLADE SIZE 11
593	MedicalSupplies	SURGICAL BLADE SIZE 11(ULTRA SHARP) OPHTHA USE
594	MedicalSupplies	SURGICAL FACEMASK WITH FACE SHIELD
595	MedicalSupplies	SURGICAL HANDBRUSH W/POVIDONE 7.5% DISPOSABLE
596	MedicalSupplies	SURGICAL HANDBRUSH, ORDINARY
597	MedicalSupplies	SURGICAL INCISION DRAPE POUCH INCISION FILM 30-40"
598	MedicalSupplies	SURGICAL INCISION DRAPE POUCH INCISION FILM, 22INCH X 17 INCH
599	MedicalSupplies	SURGICAL INCISION DRAPE POUCH INCISION FILM 13INCH X 13 INCH
600	MedicalSupplies	SURGICAL MASK
601	MedicalSupplies	SURGICAL NEEDLE / FERGUSON #10
602	MedicalSupplies	SURGICAL NEEDLE / FERGUSON #12
603	MedicalSupplies	SURGICAL NEEDLE / MAYO CUT GUT #3
604	MedicalSupplies	SURGICAL NEEDLE / MAYO CUT GUT #4
605	MedicalSupplies	SURGICAL NEEDLE / MAYO CUT GUT #5
606	MedicalSupplies	SURGICAL NEEDLE / MAYO CUT GUT #6
607	MedicalSupplies	SURGICAL NEEDLE / MURPHY INTESTIAL #3
608	MedicalSupplies	SURGICAL NEEDLE / ROUND POINT #2 FERGUSON
609	MedicalSupplies	SURGICAL NEEDLE / ROUND POINT #4 FERGUSON
610	MedicalSupplies	SURGICAL NEEDLE ANCHOR SN FERGUSON 1/2 CIRCLE #10
611	MedicalSupplies	SURGICAL NEEDLE ANCHOR SN FERGUSON 1/2 CIRCLE #12
612	MedicalSupplies	SURGICAL NEEDLE ANCHOR SN FERGUSON 1/2 CIRCLE #18
613	MedicalSupplies	SURGICAL NEEDLE ANCHOR SN FISTULA 1/2 CIRCLE (HEAVY GAUGE CUTTING NEEDLE)
614	MedicalSupplies	SURGICAL NEEDLE ANCHOR SN HERNIA 1/2 CIRCLE BLUNT POINT
615	MedicalSupplies	SURGICAL NEEDLE ANCHOR SN INTESTINAL #3/4 3/8 CIRCLE FINE TAPER POINT
616	MedicalSupplies	SURGICAL NEEDLE ANCHOR SN INTESTINAL 1/2 CIRCLE FINE TAPER POINT #4
617	MedicalSupplies	SURGICAL NEEDLE ANCHOR SN INTESTINAL 1/2 CIRCLE FINE TAPER POINT #5
618	MedicalSupplies	SURGICAL NEEDLE ANCHOR SN INTESTINAL 1/2 CIRCLE FINE TAPER POINT #6
619	MedicalSupplies	SURGICAL NEEDLE ANCHOR SN MAYO 1/2 CIRCLE TAPER POINT #5
620	MedicalSupplies	SURGICAL NEEDLE ANCHOR SN ROUND PINT 1/2 CIRCLE #2
621	MedicalSupplies	SURGICAL NEEDLE ANCHOR SN ROUND PINT 1/2 CIRCLE #3
622	MedicalSupplies	TRACHEOSTOMY TUBE WITHOUT FENESTRATION AND WITHOUT INNER CANNULA I.D.4.0MM;O.D. 5.9 MM;LENGTH: 34MM SIZE 4.0 (PEDIA CUFFLESS)
623	MedicalSupplies	TRACHEOSTOMY, PEDIATRIC, UNCUFFED, SIZE 4.0
624	MedicalSupplies	TRACHSEAL ADULT ENDOTRACHEAL CLOSED SUCTION SYSTEM, 72HOUR SIZE F.12
625	MedicalSupplies	TRACHSEAL ADULT ENDOTRACHEAL CLOSED SUCTION SYSTEM, 72HOUR SIZE F.14
626	MedicalSupplies	TRACHSEAL ADULT ENDOTRACHEAL CLOSED SUCTION SYSTEM, 72HOUR SIZE F.16
627	MedicalSupplies	TRANSPARENT DRESSING W/ ABSORBENT PAD, 15.5 X 8.5CM
628	MedicalSupplies	TRANSPARENT DRESSING W/ ABSORBENT PAD, 20X10CM
629	MedicalSupplies	TRANSPARENT DRESSING W/ ABSORBENT PAD, 25X10CM
630	MedicalSupplies	TRANSPARENT DRESSING W/ ABSORBENT PAD, 9.5 X 8.5CM
631	MedicalSupplies	TRANSPARENT DRESSING W/ NON-ADHERENT PAD 6CM X 10CM "3M"
632	MedicalSupplies	TRANSPARENT DRESSING W/ NON-ADHERENT PAD 9CM X 20CM
633	MedicalSupplies	TRANSPARENT DRESSING W/ NON-ADHERENT PAD 9CM X 25CM

634	MedicalSupplies	TRANSPARENT DRESSING W/ NON-ADHERENT PAD 9CM X 35CM "3M"
635	MedicalSupplies	TRANSPARENT DRESSING W/ABSORBENT PAD 6.5 X 5CM
636	MedicalSupplies	TRANSPARENT DRESSING W/NON-ADHERENT PAD 9CM X 15CM "3M"
637	MedicalSupplies	TRIANGULAR BANDAGE "WHITE STANDARD SIZE"
638	MedicalSupplies	TROCAR 10-15MM
639	MedicalSupplies	TRYPAN BLUE DYE, 1 MG/ML
640	MedicalSupplies	T-TUBE, DRAINAGE, LATEX F. 18
641	MedicalSupplies	T-TUBE, DRAINAGE, LATEX FR. 16
642	MedicalSupplies	TUBE AND CHAMBER KIT FOR NASAL HIGH FLOW
643	MedicalSupplies	TUBE CHAMBER KIT ADULT FOR NASAL HIGH FLOW SYSTEM
644	MedicalSupplies	TUBE CHAMBER KIT PEDIA FOR NASAL HIGH FLOW SYSTEM
645	MedicalSupplies	TUBE DOUBLE LUMEN (LEFT & RIGHT) S. 7.0
646	MedicalSupplies	TUBE DOUBLE LUMEN (LEFT & RIGHT) S. 7.5
647	MedicalSupplies	TUBE DOUBLE LUMEN (LEFT & RIGHT) S. 8.0
648	MedicalSupplies	TUBE ENDO S-2.5 UNCUFFED
649	MedicalSupplies	TUBE ENDO S-3.0 UNCUFFED
650	MedicalSupplies	TUBE ENDO S-3.5 UNCUFFED
651	MedicalSupplies	TUBE ENDO S-4.0 UNCUFFED
652	MedicalSupplies	TUBE ENDO S-4.5, UNCUFFED
653	MedicalSupplies	TUBE ENDO S-5.0, CUFFED
654	MedicalSupplies	TUBE ENDO S-5.5
655	MedicalSupplies	TUBE ENDO S-6.0, CUFFED
656	MedicalSupplies	TUBE ENDO S-6.5, CUFFED
657	MedicalSupplies	TUBE ENDO S-7.0 CUFFED
658	MedicalSupplies	TUBE ENDO S-7.5
659	MedicalSupplies	TUBE ENDO S-8.0 CUFFED
660	MedicalSupplies	TUBE ENDO S-8.5, CUFFED
661	MedicalSupplies	TUBE ENDOBRONCHIAL DOUBLE LUMEN S-35 LEFT
662	MedicalSupplies	TUBE ENDOBRONCHIAL DOUBLE LUMEN S-35 RIGHT
663	MedicalSupplies	TUBE ENDOBRONCHIAL DOUBLE LUMEN S-37 LEFT
664	MedicalSupplies	TUBE ENDOBRONCHIAL DOUBLE LUMEN S-37 RIGHT
665	MedicalSupplies	TUBE ENDOBRONCHIAL DOUBLE LUMEN S-39 LEFT
666	MedicalSupplies	TUBE ENDOBRONCHIAL DOUBLE LUMEN S-39 RIGHT
667	MedicalSupplies	TUBE FEEDING FR 8 X 100 CM STRAIGHT , DISP TRANSPARENT STERILE
668	MedicalSupplies	TUBE FEEDING FR 8 X 40 CM STRAIGHT , DISP TRANSPARENT STERILE
669	MedicalSupplies	TUBE FEEDING FR. 10 X 120 CM
670	MedicalSupplies	TUBE FEEDING FR. 12 X 120 CM
671	MedicalSupplies	TUBE FEEDING FR. 14 SILICONIZED
672	MedicalSupplies	TUBE FEEDING FR. 14 X 120 CM
673	MedicalSupplies	TUBE FEEDING FR. 16 SILICONIZED
674	MedicalSupplies	TUBE FEEDING FR. 16 X 125 CM
675	MedicalSupplies	TUBE FEEDING FR. 16 X 140 CM
676	MedicalSupplies	TUBE FEEDING FR. 5 X 100CM
677	MedicalSupplies	TUBE FEEDING FR. 5 X 40CM
678	MedicalSupplies	TUBE INTUBATING STYLET FR. 10 "RUSCH"
679	MedicalSupplies	TUBE INTUBATING STYLET FR. 12 "RUSCH"
680	MedicalSupplies	TUBE INTUBATING STYLET FR. 6 "RUSCH"
681	MedicalSupplies	TUBE INTUBATING STYLET, ADULT

682	MedicalSupplies	TUBE INTUBATING STYLET, CHILD
683	MedicalSupplies	TUBE INTUBATING STYLET, INFANT
684	MedicalSupplies	TUBE INTUBATING STYLET, NEONATE
685	MedicalSupplies	TUBE INTUBATING STYLET, PEDIA
686	MedicalSupplies	TUBE IV EXTENSION 50CM
687	MedicalSupplies	TUBE NGT FR 10 X 100CM
688	MedicalSupplies	TUBE NGT FR 12 X 100CM
689	MedicalSupplies	TUBE PREFORMED (NORTH & SOUTH) S-5.0
690	MedicalSupplies	TUBE PREFORMED (NORTH & SOUTH) S-5.5
691	MedicalSupplies	TUBE PREFORMED (NORTH & SOUTH) S-6.0
692	MedicalSupplies	TUBE PREFORMED (NORTH & SOUTH) S-6.5
693	MedicalSupplies	TUBE PREFORMED (NORTH & SOUTH) S-7.0
694	MedicalSupplies	TUBE PREFORMED (NORTH & SOUTH) S-7.5
695	MedicalSupplies	TUBE PREFORMED (NORTH & SOUTH) S-8.0
696	MedicalSupplies	TUBE PREFORMED (NORTH & SOUTH) S-8.5
697	MedicalSupplies	TUBE PREFORMED TRACHEAL W/O CUFF S 3.5
698	MedicalSupplies	TUBE PREFORMED TRACHEAL W/O CUFF S 4.0
699	MedicalSupplies	TUBE PREFORMED TRACHEAL W/O CUFF S 5.0
700	MedicalSupplies	TUBE RAE PREFORMED TRACHEAL W/ CUFF FR 4.5
701	MedicalSupplies	TUBE RAE PREFORMED TRACHEAL W/ CUFF FR 5.0
702	MedicalSupplies	TUBE RAE PREFORMED TRACHEAL W/ CUFF FR 6.0
703	MedicalSupplies	TUBE RAE PREFORMED TRACHEAL W/ CUFF FR 6.5
704	MedicalSupplies	TUBE RAE PREFORMED TRACHEAL W/ CUFF FR 7.0
705	MedicalSupplies	TUBE RAE PREFORMED TRACHEAL W/ CUFF FR 7.5
706	MedicalSupplies	TUBE RAE PREFORMED TRACHEAL W/ CUFF FR 8.5
707	MedicalSupplies	TUBE RAE PREFORMED TRACHEAL W/O CUFF FR 3.5
708	MedicalSupplies	TUBE RAE PREFORMED TRACHEAL W/O CUFF FR 4.0
709	MedicalSupplies	TUBE RAE PREFORMED TRACHEAL W/O CUFF FR 4.5
710	MedicalSupplies	TUBE RAE PREFORMED TRACHEAL W/O CUFF FR 5.0
711	MedicalSupplies	TUBE RAE PREFORMED TRACHEAL W/O CUFF FR 5.5
712	MedicalSupplies	TUBE RAE PREFORMED TRACHEAL W/O CUFF FR 6.0
713	MedicalSupplies	TUBE RAE PREFORMED TRACHEAL CUFF S-3.5
714	MedicalSupplies	TUBE RAE PREFORMED TRACHEAL CUFF S-4.0
715	MedicalSupplies	TUBE RAE/PREFORMED CUFFED S-5.5
716	MedicalSupplies	TUBE REINFORCED TRACHEAL / ARMoured S-3.5
717	MedicalSupplies	TUBE REINFORCED TRACHEAL / ARMoured S-4.0
718	MedicalSupplies	TUBE REINFORCED TRACHEAL / ARMoured S-4.5
719	MedicalSupplies	TUBE REINFORCED TRACHEAL / ARMoured S-5.0
720	MedicalSupplies	TUBE REINFORCED TRACHEAL / ARMoured S-5.5
721	MedicalSupplies	TUBE REINFORCED TRACHEAL / ARMoured S-6.0
722	MedicalSupplies	TUBE REINFORCED TRACHEAL / ARMoured S-6.5
723	MedicalSupplies	TUBE REINFORCED TRACHEAL / ARMoured S-7.0
724	MedicalSupplies	TUBE REINFORCED TRACHEAL / ARMoured S-7.5
725	MedicalSupplies	TUBE REINFORCED TRACHEAL / ARMoured S-8.0
726	MedicalSupplies	TUBE REINFORCED TRACHEAL / ARMoured S-8.5
727	MedicalSupplies	TUBE T. FR.16 LATEX FREE, STERILE
728	MedicalSupplies	TUBE T. FR.18 LATEX FREE, STERILE
729	MedicalSupplies	TUBE T-DRAINAGE FR.14

730	MedicalSupplies	TUBE TRACHEOSTOMY W/ SOFT SEAL CUFF S-5.0
731	MedicalSupplies	TUBE TRACHEOSTOMY W/ SOFT SEAL CUFF S-5.5
732	MedicalSupplies	TUBE TRACHEOSTOMY W/ SOFT SEAL CUFF S-6.0
733	MedicalSupplies	TUBE TRACHEOSTOMY W/ SOFT SEAL CUFF S-7.0
734	MedicalSupplies	TUBE TRACHEOSTOMY W/ SOFT SEAL CUFF S-7.5
735	MedicalSupplies	TUBE TRACHEOSTOMY W/ SOFT SEAL CUFF S-8.0
736	MedicalSupplies	TUBE TRACHEOSTOMY WITH FENESTRATION AND INNER CANNULA SIZE. 6 SOFT SWIVEL NECK FLANGE WITH LASER ETCHED MARKINGS
737	MedicalSupplies	TUBE TRACHEOSTOMY WITH FENESTRATION AND INNER CANNULA SIZE. 8
738	MedicalSupplies	TUBE, THORACOSTOMY G 28
739	MedicalSupplies	TUBE, THORACOSTOMY G 32
740	MedicalSupplies	TUBERCULIN SYRINGE 1ML
741	MedicalSupplies	TUBINGS FOR MECH. VENTILLATOR (ADULT)
742	MedicalSupplies	TUBINGS FOR MECH. VENTILLATOR (NEONATE)
743	MedicalSupplies	TUBINGS FOR MECH. VENTILLATOR (PEDIA)
744	MedicalSupplies	TUMOR STENT WITH HYDROGEL COATING FR. 6 X 26CM
745	MedicalSupplies	TUMOR STENT WITH HYDROGEL COATING FR. 7 X 24CM
746	MedicalSupplies	TUMOR STENT WITH HYDROGEL COATING FR. 7 X 26CM
747	MedicalSupplies	TYPANIC THERMOMETER
748	MedicalSupplies	ULTRASONIC OPEN SHEARS 17CM CURVED COMPATIBLE
749	MedicalSupplies	ULTRASONIC OPEN SHEARS 9CM CURVED COMPATIBLE
750	MedicalSupplies	ULTRASOUND FLUIDIC MGT SYSTEM
751	MedicalSupplies	ULTRASOUND TRANSMISSION GEL IN FOIL 20G/ 48FOILS PER BOX
752	MedicalSupplies	UMBILICAL CATHETER FR. 5
753	MedicalSupplies	UNDERPADS, ABSORBENT
754	MedicalSupplies	URETER CATHETER ER(STERILE) SIZE FR.4
755	MedicalSupplies	URETER CATHETER ER(STERILE) SIZE FR.5
756	MedicalSupplies	URETER CATHETER ER(STERILE) SIZE FR.6
757	MedicalSupplies	URETERAL CATHETER OPEN TIP FR. 6
758	MedicalSupplies	URINAL PLASTIC WHITE
759	MedicalSupplies	URINE BAG (ADULT)
760	MedicalSupplies	URINE BAG (WEE BAG / PEDIA URINE COLLECTOR)
761	MedicalSupplies	URINE BAG, CALIBRATION OF 1000-2000ML
762	MedicalSupplies	URINE COLLECTOR PEDIA
763	MedicalSupplies	VAGINAL SPECULUM, STAINLESS (LARGE)
764	MedicalSupplies	VAGINAL SPECULUM, STAINLESS (MEDIUM)
765	MedicalSupplies	VAGINAL SPECULUM, STAINLESS (SMALL)
766	MedicalSupplies	VAGINAL SPECULUM, STAINLESS (VAGINAL)
767	MedicalSupplies	VASCULAR PROSTHESIS IMPREGNATED WITH ABSORBABLE BOVINE GELATIN 14X7-10MM, 40CM BIFURCATED "UNIGRAFT"
768	MedicalSupplies	VASCULAR PROSTHESIS IMPREGNATED WITH ABSORBABLE BOVINE GELATIN 16X8-10MM, 40CM BIFURCATED "UNIGRAFT"
769	MedicalSupplies	VASCULAR PROSTHESIS IMPREGNATED WITH ABSORBABLE BOVINE GELATIN 18X9-10MM, 40CM BIFURCATED "UNIGRAFT"
770	MedicalSupplies	VASCULAR PROSTHESIS IMPREGNATED WITH ABSORBABLE BOVINE GELATIN 20X10MM, 40CM BIFURCATED "UNIGRAFT"
771	MedicalSupplies	VASCULAR PROSTHESIS IMPREGNATED WITH ABSORBABLE BOVINE GELATIN 22X10-11MM, 40CM BIFURCATED "UNIGRAFT"
772	MedicalSupplies	VASCUTEK GELWEAVE BIFURCATED GRAFT VASCULAR GRAFT 24MM X 12MM X 45MM
773	MedicalSupplies	VENTILATOR TUBING, ADULT
774	MedicalSupplies	VENTILATOR TUBING, NEONATE
775	MedicalSupplies	VENTILATOR TUBING, PEDIA

776	MedicalSupplies	VESSEL SEALER 18.8CM LENGTH JAW SURFACE
777	MedicalSupplies	VESSEL SEALER 18.8CM LENGTH RECTANGULAR SHAFT
778	MedicalSupplies	VESSEL SEALER 20CM LENGTH 10MM 359 DEGREES
779	MedicalSupplies	VESSEL SEALER 23CM LENGTH 5MM 360 DEGREES
780	MedicalSupplies	VESSEL SEALER 37CM LENGTH 5MM 180 DEGREES
781	MedicalSupplies	VESSEL SEALER 37CM LENGTH 5MM 350 DEGREES
782	MedicalSupplies	VESSEL SEALER/CORDLESS DISSECTOR 23-26CM LENGTH
783	MedicalSupplies	VITRECTOMY CUTTER G.23
784	MedicalSupplies	VITRECTOMY- POSTERIOR PROCEDURAL PACK G.23 BL5223
785	MedicalSupplies	VITRECTOMY-POSTERIOR PROCEDURAL PACK G.25 BL5225WVX
786	MedicalSupplies	VOLUMETRIC SET (SOLUSET)
787	MedicalSupplies	VP SHUNT, FLAT BOTTON SMALL, CONTOURED, MEDIUM PRESSURE
788	MedicalSupplies	WADDING SHEET 4" X 5 YARDS
789	MedicalSupplies	WADDING SHEET 5" X 5 YARDS
790	MedicalSupplies	WADDING SHEET 6" X 5 YARDS
791	MedicalSupplies	WOUND IRRIGATION SOLUTION (PRONTOSAN)
792	MedicalSupplies	WRIST SPLINT ADULT
793	MedicalSupplies	WRIST SPLINT MEDIUM
794	MedicalSupplies	WRIST SPLINT NEONATE
795	MedicalSupplies	WRIST SPLINT PEDIATRIC
796	MedicalSupplies	WRIST TAG, BLUE, PEDIA
797	MedicalSupplies	WRIST TAG, PINK, PEDIA
798	MedicalSupplies	X-RAY FILM ENVELOPE W/LOGO 14-5 X 17.5 INCH
799	MedicalSupplies	Y CONNECTOR W/ELBOW, ADULT
800	MedicalSupplies	Y CONNECTOR, PEDIA
801	MedicalSupplies	YARSAGIL PERM MINI CLIP STR 7MM
802	MedicalSupplies	YARSAGIL PERM STD-CLIP STR 6.6MM
803	MedicalSupplies	YASARGIL ANEURYSM CLIP 6.4MM
804	MedicalSupplies	YASARGIL ANEURYSM CLIP 6.5MM
805	MedicalSupplies	YASARGIL ANEURYSM CLIP 7MM
806	MedicalSupplies	YASARGIL ANEURYSM CLIP 8.3MM
807	MedicalSupplies	YASARGIL ANEURYSM CLIP 9MM
808	MedicalSupplies	MICROSTREAM CIRCUIT NASAL FILTER LINE
809	MedicalSupplies	AUTOCLAVE TAPE 1" STERILIZA INDICATOR FOR STEAM, GOOD ADHESIVE
810	MedicalSupplies	BANDAGE ELASTIC 6 X 5 YARDS
811	MedicalSupplies	BANDAGE, ELASTIC 2 X 5 YRDS
812	MedicalSupplies	BANDAGE, ELASTIC 3 X 5 YRDS
813	MedicalSupplies	BANDAGE, ELASTIC 4 X 5 YRDS PLIANT ADAPTABLE
814	MedicalSupplies	COTTON ABSORBENT 400 GRAMS, 100% COTTON
815	MedicalSupplies	COTTON ABSORBENT STERILE 12 GRAMS (100% COTTON)
816	MedicalSupplies	ECG PAPER FOR COMEN, 80 X 20 MM (3CHANNEL)
817	MedicalSupplies	GAUZE ABSORBENT 24X328 WITH X-RAY
818	MedicalSupplies	GAUZE ABSORBENT 28X24 MESH 36X100 YARDS
819	MedicalSupplies	GAUZE ROLL BANDAGE 4X10 YARD
820	MedicalSupplies	HEAVY DUTY PORTABLE BEDPAN FOR MEN AND WOMEN
821	MedicalSupplies	KINESEO TAPE (BLACK, BLUE & PINK)
822	MedicalSupplies	KNOTLESS TISSUE CONTROL DEVICE SIZE 2.0 30CM W/ 36MM CIRCLE TAPER NEEDLE 12'S
823	MedicalSupplies	NEEDLE FERGUSON EYE CUTTING EDGE #1

824	MedicalSupplies	NEEDLE FERGUSON EYE CUTTING EDGE#2
825	MedicalSupplies	NEEDLE FERGUSON EYE TAPER POINT #3
826	MedicalSupplies	NEEDLE FERGUSON EYE TAPER POINT #4
827	MedicalSupplies	NEEDLE FERGUSON TAPER POINT #10
828	MedicalSupplies	NEEDLE FERGUSON TAPER POINT #12
829	MedicalSupplies	NEEDLE FERGUSON TAPER POINT #8
830	MedicalSupplies	NEEDLE MAYO INTESTINAL TAPER POINT #3.0 ANCHOR
831	MedicalSupplies	NEEDLE REGULAR SURGEON CUTTING EDGE #1
832	MedicalSupplies	NEEDLE, SPINAL G. 27 X 3.5
833	MedicalSupplies	PARTIALLY ABSORBABLE POLYPROPYLENE MESH, 15CM X 15CM
834	MedicalSupplies	POLY 4-HYDROXYBUTYRATE MONOFILAMENT 1-0 HALF CIRCLE ROUND,
835	MedicalSupplies	POLY 4-HYDROXYBUTYRATE MONOFILAMENT EXTRA LONG TERM, SYNTHETIC ABSORBABLE SUTURE SIZE 1.0 W/HALF CIRCLE ROUND NEEDLE 40MM LOOP
836	MedicalSupplies	POLYBUTILATE COATED POLYESTER BRAIDED NON ABSORBABLE SUTURE SIZE 1 75CM WITH 30MM J NEEDLE, BOX OF 12 "ETHICON"
837	MedicalSupplies	POLYBUTILATE COATED POLYESTER BRAIDED NON ABSORBABLE SUTURE SIZE 5 4-75CM WITH 40MM TAPECUT-37 NEEDLE 12'S PER, BOX "ETHICON"
838	MedicalSupplies	POLYBUTILATE COATED POLYESTER BRAIDED NON ABSORBABLE SUTURE SIZE 5 4-75CM WITH 48MM 12'S PER, BOX "ETHICON"
839	MedicalSupplies	POLYBUTILATE COATED POLYESTER SUTURE SIZE 2, 75CM WITH 40MM TAPER CUT NEEDLE BOX OF 12 "ETHICON"
840	MedicalSupplies	POLYDIAAXANON 2.0 R 20-40MM 1/2 CIRCLE NEEDLE
841	MedicalSupplies	POLYDIAAXANONE WITH TRICLOSAN, SYNTHETIC, ABSORBABLE KNOTLESS TISSUE CONTROL DEVICE (SPIRAL DESIGN)WITH ADJUSTABLE FIXATION LOOP, SIZE 2-0, 30CM WITH CIRCLE TAPER NEEDLE, BOX OF 12
842	MedicalSupplies	POLYDIOXANON MONOFILAMENT SYNTHETIC ABS. SUTURE SIZE 2 W/ 26MM
843	MedicalSupplies	POLYIMIDE 4.0 NON ABSORBABLE W/ NEEDLE CUTTING (NYLON)
844	MedicalSupplies	POLYIMIDE 5.0 NON ABSORBABLE W/ NEEDLE CUTTING (NYLON)
845	MedicalSupplies	POLYIMIDE 6.0 NON ABSORBABLE W/ NEEDLE CUTTING (NYLON)
846	MedicalSupplies	POLYPROPHYLENE 2.0 SYNTHETIC NON ABSORBABLE MONOFILAMENT DOUBLE ARM ROUND NEEDLE
847	MedicalSupplies	POLYPROPHYLENE 3.0 SYNTHETIC NON ABSORBABLE MONOFILAMENT DOUBLE ARM ROUND NEEDLE
848	MedicalSupplies	POLYPROPHYLENE 4.0 SYNTHETIC NON ABSORBABLE MONOFILAMENT SUTURE DOUBLE ARM ROUND NEEDLE
849	MedicalSupplies	POLYPROPHYLENE 5.0 SYNTHETIC NON ABSORBABLE MONOFILAMENT SUTURE DOUBLE ARM ROUND NEEDLE
850	MedicalSupplies	POLYPROPHYLENE 6.0 SYNTHETIC 24-60CM, 9MM DOUBLE ARM NEEDLE
851	MedicalSupplies	PROCUREMENT OF LIQUID BASED PAP TEST KIT THINPREP LIQUID BASED CYTOLOGY
852	MedicalSupplies	SUTURE BRAIDED NATURAL SILK STRANDS 0, 15X60CM
853	MedicalSupplies	SUTURE MIDTERM SYNTHETIC 0, ABSORBABLE BRAIDED AND COATED MADE OF POLYGLYCOLIDE CO L LACTID 90CM ROUND NEEDLE
854	MedicalSupplies	SUTURE MIDTERM SYNTHETIC 1, ABSORBABLE BRAIDED AND COATED MADE OF POLYGLYCOLIDE CO L LACTID 90CM ROUND NEEDLE
855	MedicalSupplies	SUTURE MONOFILAMENT MID TERM 0, ABS. SUTURE MADE OF GLUCONATE ROUND NEEDLE
856	MedicalSupplies	SUTURE POLYGLECAPRONE SYNTHETIC MONOFILAMENT ABS. SUTURE 4-0 27-70CM 19MM
857	MedicalSupplies	SUTURE, CHROMIC 0 W/ NEEDLE ABSORBABLE 20-40MM
858	MedicalSupplies	SUTURE, CHROMIC 1.0 W/NEEDLE 40MM, CIRCLE ROUND
859	MedicalSupplies	SUTURE, CHROMIC 2.0 W/NEEDLE 37MM, CIRCLE ROUND
860	MedicalSupplies	SUTURE, CHROMIC 3.0 W/CUTTING NEEDLE
861	MedicalSupplies	SUTURE, CHROMIC 4.0 ROUND NEEDLE
862	MedicalSupplies	SUTURE, CHROMIC 5.0 W/ NEEDLE, ROUND, CUTGUT
863	MedicalSupplies	SUTURE, MONOFILAMENT NON ABSORBABLE POLYAMIDE S-10.0

864	MedicalSupplies	SUTURE, POLYDIOXANON MONOFILAMENT SYNTHETIC ABSORBABLE S-o 150CM W/ CIRCLE MAYO TAPER POINT NEEDLE
865	MedicalSupplies	SUTURE, POLYDIOXANON MONOFILAMENT SYNTHETIC ABSORBABLE S-o 70CM LONG, ROUND NEEDLE
866	MedicalSupplies	SUTURE, SILK 1.0, PRE-CUT
867	MedicalSupplies	SUTURE, SILK 2.0 W/25-26MM ROUND NEEDLE, NON ABSORBABLE, 75-80CM LONG
868	MedicalSupplies	SUTURE, SILK 2.0 W/OUT NEEDLE PRE-CUT BRAIDED NON ABSORBABLE 13-15 STRANDS BY 60-80CM LONG
869	MedicalSupplies	SUTURE, SILK 2.0 WITH 25-26MM, CUTTING NEEDLE NON ABSORBABLE 75-80CM LONG
870	MedicalSupplies	SUTURE, SILK 3.0 W/ 26MM CUTTING NEEDLE
871	MedicalSupplies	SUTURE, SILK 3.0 W/25-26MM ROUND BODIED NEEDLE NON ABSORBABLE
872	MedicalSupplies	SUTURE, SILK 3.0 W/OUT NEEDLE PRE-CUT BRAIDED NON ABSORBABLE 13-15 STRANDS 60-80CM LONG
873	MedicalSupplies	SUTURE, SILK 4.0 W/ CUTTING NEEDLE, NON-ABSORBABLE
874	MedicalSupplies	SUTURE, SILK 4.0 W/25-26MM ROUND BODIED NEEDLE NON ABSORBABLE
875	MedicalSupplies	SUTURE, SILK 4.0 W/OUT NEEDLE PRE-CUT BRAIDED NON ABSORBABLE 13-15 STRANDS BY 60-80CM LONG
876	MedicalSupplies	SUTURE, SILK 5.0 W/ 10MM ROUND NEEDLE
877	MedicalSupplies	SUTURE, SILK 5.0 W/ 3/8 CUTTING NEEDLE 19MM
878	MedicalSupplies	SUTURE, SILK 6.0 W/ 3/8 CIRCLE CUTTING NEEDLE
879	MedicalSupplies	SUTURE, SYNTHETIC BRAIDED ABS. 20MM 10.0, VICRYL (ETHILON)
880	MedicalSupplies	SUTURE, SYNTHETIC BRAIDED ABSORBABLE 6.0, SPATULATED NEEDLE, DOUBLE ARM 1/2
881	MedicalSupplies	SUTURE, SYNTHETIC BRAIDED ABSORBABLE S-o, ROUND NEEDLE 36-40MM
882	MedicalSupplies	SUTURE, SYNTHETIC BRAIDED ABSORBABLE S-1.0 ROUND NEEDLE
883	MedicalSupplies	SUTURE, SYNTHETIC BRAIDED ABSORBABLE S-2.0, ROUND NEEDLE 36-40MM
884	MedicalSupplies	SUTURE, SYNTHETIC BRAIDED ABSORBABLE S-3.0, ROUND NEEDLE
885	MedicalSupplies	SUTURE, SYNTHETIC BRAIDED ABSORBABLE S-4.0 CUTTING NEEDLE
886	MedicalSupplies	SUTURE, SYNTHETIC BRAIDED ABSORBABLE S-4.0 ROUND NEEDLE
887	MedicalSupplies	SUTURE, SYNTHETIC BRAIDED ABSORBABLE S-5.0 16-19 CUTTING
888	MedicalSupplies	SUTURE, SYNTHETIC BRAIDED ABSORBABLE SUTURE SIZE 2.0 70M W/ UR-5 NEEDLE
889	MedicalSupplies	SUTURE, SYNTHETIC MONOFILAMENT S-1
890	MedicalSupplies	SUTURE,POLYMAIDE 3.0 NON ABSORBABLE W/ NEEDLE CUTTING MONOFILAMENT 18"45CM,20-26MM (NYLON)
891	MedicalSupplies	SUTURE,SYNTHETIC ,MONOFILAMENT ABSORBABLE SUTURE S-2.0 W/ROUND NEEDLE
892	MedicalSupplies	SUTURE,SYNTHETIC BRAIDED ABSORBABLE SUTURE SIZE 2.0 70M W/ UR-6 NEEDLE
893	MedicalSupplies	SUTURE,SYNTHETIC BRAIDED ABSORBABLE SUTURES S-7.0 45CM WITH 5.5CM TG 160-8 12'S
894	MedicalSupplies	SUTURE,SYNTHETIC BRAIDED ABSORBABLE SUTURES S-8.0 45CM WITH 6.5CM TG 140-8 12'S
895	MedicalSupplies	SUTURE,SYNTHETIC MONOFILAMENT ABSORBABLE SUTURE S-4.0 CUTTING NEEDLE
896	MedicalSupplies	SUTURE. BRAIDED SYNTHETIC ABSORBABLE SUTURE S-6.0, 30CM, WITH 7.6-8MM S-29, S-29 NEEDLE "ETHICON"
897	MedicalSupplies	4 CHLORHEXIDINE SURGICAL HANDWASH 5LITERS/GAL "BIOSCRUB"
898	MedicalSupplies	70% ISOPROPYL ALCOHOL WITH MOISTURIZER
899	MedicalSupplies	ALKYLAMINE ALDEHYDE FREE " DECONEX" DISINFECTANT
900	MedicalSupplies	ANIOSYME SYNERGY 5 ENZYMATIC INSTRUMENT CLEANER 5L/G
901	MedicalSupplies	BLEACHING SOLUTION
902	MedicalSupplies	CARBON DIOXIDE ABSORBER 37LBS
903	MedicalSupplies	CARBON REMOVER
904	MedicalSupplies	DELIMER

905	MedicalSupplies	DISHWASHING LIQUID 1 GAL/CONT. 1:20 DILUTION
906	MedicalSupplies	DISINFECTANT (SPRAY) LYSOL 170G
907	MedicalSupplies	DISINFECTANT SOLUTION FOR CLEANING EQUIPMENT
908	MedicalSupplies	FABRIC CONDITIONER
909	MedicalSupplies	GLUTARALDEHYDE SOLUTION 2.45%
910	MedicalSupplies	HAND SANITIZER AT LEAST 70% ISOPROPHYL
911	MedicalSupplies	HAND SANITIZER GAL SIZE W/FREE STAINLESS HOLDER & PERSONALIZED BOTTLE W/PUMP & W/QMMC LOGO
912	MedicalSupplies	HIGH LEVEL DISINFECTANT BASED ON PARACETIC ACID
913	MedicalSupplies	LIQUID HAND SOAP ANTIBAC 1GAL/CONT. (FOR DIETARY USE)
914	MedicalSupplies	LIQUID SOAP, DISHWASHING
915	MedicalSupplies	POVIDONE IODINE ANTISEPTIC 10%
916	MedicalSupplies	POVIDONE IODINE ANTISEPTIC 7.5%
917	MedicalSupplies	ULTRASOUND GEL
918	MedicalSupplies	VINEGAR 4% NATURAL ACIDITY
919	MedicalSupplies	PLASTER SILK ACETATE - DURAPORE 1" PER INCHES
920	MedicalSupplies	PLASTER SILK ACETATE - DURAPORE 2" PER INCHES
921	MedicalSupplies	EEG PASTE, TEN 20 CONDUCTIVE 8OZ. JAR
922	MedicalSupplies	LIDOCAINE 5% OINTMENT 50 GRAMS
923	MedicalSupplies	DETERGENT POWDER ALL PURPOSE, 1KG
924	MedicalSupplies	200/200ML DISPOSABLE SYRINGE KIT
925	MedicalSupplies	24 KIT OF PHOX/B 1 PLUS PACK ON BOARD QC
926	MedicalSupplies	ADMISSION KIT FOR SUITE
927	MedicalSupplies	CATHETER NEPHROSTOMY, PIGTAIL FR.8
928	MedicalSupplies	DISPOSABLE SYRINGE KIT FOR (64 SLICE)
929	MedicalSupplies	HIGH PRESSURE SYRINGE INJECTOR KIT FOR (16 SLICE) MUST BE FIT FOR IMAXEON POWER INJECTOR)
930	MedicalSupplies	NEPHROSTOMY CATHETER KIT FR. 10
931	MedicalSupplies	NEPHROSTOMY CATHETER KIT FR. 12
932	MedicalSupplies	NON INVASIVE PRE-NATAL TEST
933	MedicalSupplies	QUALITY CONTROL PACK (QC)
934	MedicalSupplies	REFERENCE ELECTRODE
935	MedicalSupplies	TUBE AND CHAMBER KIT (900PT501)
936	MedicalSupplies	ALCOHOL 70% (500ML) PER ML
937	MedicalSupplies	DISINFECTANT BLEACH PER CC/ML
938	MedicalSupplies	LIQUID SOAP (ML)
939	MedicalSupplies	NON-IONIC SURFACTANT, ESEMTAN(100ML, SKIN CLEANSER)
940	MedicalSupplies	POVIDONE IODINE - PER ML/CC
941	MedicalSupplies	ADMISSION KIT
942	MedicalSupplies	ADMISSION KIT FOR CNS PT
943	MedicalSupplies	ADMISSION KIT FOR WARD
944	MedicalSupplies	APPLICATOR 3" PLASTIC STEM WITH COTTON TIP STERILE DOUBLE END
945	MedicalSupplies	BICARBONATE POWDER FOR HEMODIALYSIS
946	MedicalSupplies	CALIBRATOR PACK
947	MedicalSupplies	COTONIDS SURGICAL PATTIES 25C 75MM STERILE W/X-RAY DETECTABLE BLUE EYED STRING 100% ABSORBENT
948	MedicalSupplies	COTTON BALL 5 PCS/PACK
949	MedicalSupplies	COTTON BALLS, 3'S/PACK
950	MedicalSupplies	COTTON PLEDGET, 100'S
951	MedicalSupplies	DETERGENT SOAP

952	MedicalSupplies	DISINFECTION WIPES / TISSUE FOR MACHINE
953	MedicalSupplies	DISPOSABLE BIB, 100'S
954	MedicalSupplies	ECG ELECTRODE, ADULT 50'S
955	MedicalSupplies	ECG ELECTRODE, PEDIA
956	MedicalSupplies	ECG PAPER Z-FOLD 80MMX70MMX300 PAGES
957	MedicalSupplies	EEG DISK ELECTRODES, COLODION NE134A
958	MedicalSupplies	HEAD CAPS DISPOSABLE, GREEN
959	MedicalSupplies	MATERNITY PADS, 8'S (FOR PACKING)
960	MedicalSupplies	TONGUE DEPRESSOR, WOODEN 100'S
961	MedicalSupplies	VESSEL LOOPS MAXI 10 X 406MM RED MAXI SILICONE VESSEL LOOPS STERILE 10'S/BX
962	MedicalSupplies	AV FISTULA NEEDLE G. 16 2'S
963	MedicalSupplies	AV FISTULA NEEDLE G. 17 2'S
964	MedicalSupplies	ABSORBABLE GELATIN SPONGE HEMOSTATIC ANAL TAMPON 80MMX50MMX10MM 10'S
965	MedicalSupplies	ACTIVE IRRIGATION FMS BASIC (W/ 2 BAGS ALCON BSS)
966	MedicalSupplies	ADULT ISOTHERMAL OMNI FLEX CONNECTOR "AIRLIFE"
967	MedicalSupplies	ADVANCE WOUND DRESSING LARGE 260MMX 150MMX30MM
968	MedicalSupplies	ADVANCE WOUND DRESSING MEDIUM 160MMX 125MMX30MM
969	MedicalSupplies	ADVANCE WOUND DRESSING SMALL 100MMX75MMX30MM
970	MedicalSupplies	ADVANCED BIPOLAR SHEARS ARTICULATING CURVED
971	MedicalSupplies	AIR LIFE STERILE HUMIDIFIER500ML
972	MedicalSupplies	AIRLIFE RT ADULT FILTERED HYDROSCOPIC CONDENSER
973	MedicalSupplies	AIRWAY ORAL SIZE 60MM
974	MedicalSupplies	AIRWAY ORAL SIZE 70MM
975	MedicalSupplies	AIRWAY ORAL SIZE 80MM
976	MedicalSupplies	ALCOHOL 70% (100 ML)
977	MedicalSupplies	ALCOHOL 70% (500ML)/BOTTLE
978	MedicalSupplies	ALCOHOL 70%500ML
979	MedicalSupplies	ALDEHYDE AND PHENOL - FREE (DECONEX)
980	MedicalSupplies	AMPULE BREAKER
981	MedicalSupplies	AMPULE OPENER
982	MedicalSupplies	ANESTHESIA BREATHING CIRCUIT, ADULT 60"
983	MedicalSupplies	ANESTHESIA BREATHING CIRCUIT, PEDIA 72
984	MedicalSupplies	ANESTHESIA FACE MASK DISP. INFLATABLE SILICONE CHILD S-3
985	MedicalSupplies	ANESTHESIA FACE MASK DISP. INFLATABLE SILICONE INFANT S-2
986	MedicalSupplies	ANESTHESIA FACE MASK DISP., INFLATABLE SILICONE ADULT S-5.0
987	MedicalSupplies	ANESTHESIA FACE MASK, DISP INFLATABLE,SILICONE ADULT S-4
988	MedicalSupplies	ANESTHESIA FACE MASK, DISP INFLATABLE,SILICONE ADULT S-6
989	MedicalSupplies	ANESTHESIA FACE MASK, DISP INFLATABLE,SILICONE INFANT S-0
990	MedicalSupplies	ANESTHESIA FACE MASK, DISP INFLATABLE,SILICONE INFANT S-1
991	MedicalSupplies	ANESTHESIA MACHINE HOSE W/CYLINDER REGULATOR, MEDICAL OXYGEN HOSE, BRASS BLOCK,YOKE CYLINDER VALVE, FERRULES, OUTLET FIITNGS,THREADED NUT,HOSE NIPPLE
992	MedicalSupplies	ANTERIOR LENS CHAMBER POWER 21.0- RIGID
993	MedicalSupplies	ANTERIOR VITRECTOMY PACK- INFINITI
994	MedicalSupplies	ANTI FATIGUE MATT ORANGE SIZE 20INCH X 32 INCH,3/8 INCHES THICKNESS
995	MedicalSupplies	ANTIBACTERIAL BRAIDED SYNTHETIC ABSORBABLE SUTURE, SIZE 4-0, 30-75CM, 19MM PRECISION COSMETIC NEEDLE / BOX OF 36'S (CUTTING OR ROUND)
996	MedicalSupplies	APPLICATOR 6" WOODEN 2 PC/PACK, SINGLE END

997	MedicalSupplies	AQUACEL DRESSING SURGICAL COVER HYDROCOLLOID FIBER PAD 9X15
998	MedicalSupplies	AQUACEL HYDROFIBER 10 X 10CM
999	MedicalSupplies	AQUACEL HYDROFIBER 15 X 15CM
1000	MedicalSupplies	ARM SLING - MEDIUM
1001	MedicalSupplies	ARM SLING CLOTH , LARGE (PHILS.)
1002	MedicalSupplies	ARM SLING CLOTH- SMALL PHIL
1003	MedicalSupplies	ARM SLING W/BUE CLOTH & NYLON
1004	MedicalSupplies	ARMORED ENDOTRACHEAL TUBE SIZE 6.0
1005	MedicalSupplies	ARMORED ENDOTRACHEAL TUBE SIZE 7.0
1006	MedicalSupplies	ARMORED ENDOTRACHEAL TUBE SIZE 7.5
1007	MedicalSupplies	ASEPTO SYRINGE 60CC W/PLASTIC TIP COVER W/CONNECTOR
1008	MedicalSupplies	AUTOMATIC DOCUMENT SYSTEM LABEL FOR PLASMA
1009	MedicalSupplies	BABY WIPES80 SHEETS
1010	MedicalSupplies	BACILLUS CLAUSII 2 BILLION/5ML
1011	MedicalSupplies	BACTERIAL / VIRAL FILTER
1012	MedicalSupplies	BACTERIAL FILTER SHOULD BE COMPATIBLE W/ THE EXISTING SUCTION MACHINE (HOSPIVAC)
1013	MedicalSupplies	BAG COLOSTOMY WITH WAFER AND CLIP SIZE 56-57 MM
1014	MedicalSupplies	BAG VALVE MASK-ADULT
1015	MedicalSupplies	BAG VALVE MASK-NEONATE
1016	MedicalSupplies	BAG VALVE MASK-PEDIA
1017	MedicalSupplies	BAG, COLOSTOMY WITH WAFER AND CLIP SIZE 45MM
1018	MedicalSupplies	BAG, HOT WATER - BIG
1019	MedicalSupplies	BAG, HOT WATER - SMALL
1020	MedicalSupplies	BANDAGE, ELASTIC 2X4 YARDS
1021	MedicalSupplies	BANDAGE, ELASTIC 3"X 5YARDS
1022	MedicalSupplies	BANDAGE, ELASTIC 4" X 5YARDS
1023	MedicalSupplies	BANDAGE, ELASTIC 6" X 4 METERS
1024	MedicalSupplies	BIPAP MASK LARGE
1025	MedicalSupplies	BIPAP MASK MEDIUM
1026	MedicalSupplies	BIPAP MASK SMALL
1027	MedicalSupplies	BIPOLAR CORD
1028	MedicalSupplies	BIPOLAR FORCEP 0.7MM- 1.0MM TIP BAYONETTE TYPE, 2-PIN OR FLAT PLUG
1029	MedicalSupplies	BIPOLAR LOOP ELECTRODE
1030	MedicalSupplies	BIS MONITOR STRIPS (CONSUMABLE)
1031	MedicalSupplies	BLADE DERMATONE
1032	MedicalSupplies	BLADE SKIN GRAFT KNIFE 6.4" LONG X 7CM
1033	MedicalSupplies	BLADE SURGICAL NO. 11
1034	MedicalSupplies	BLADE SURGICAL NO. 15
1035	MedicalSupplies	BLOOD BAG TERUFLEX 450ML
1036	MedicalSupplies	BLOOD SETSANGOFIX
1037	MedicalSupplies	BLOOD TRANSFUSION SET
1038	MedicalSupplies	BLOOD TRANSFUSION SET-ADULT
1039	MedicalSupplies	BLOOD WARMER, L-INFUSION
1040	MedicalSupplies	BLOODLINES
1041	MedicalSupplies	BODY STRAP VELCO STRAP
1042	MedicalSupplies	BONE WAX
1043	MedicalSupplies	BONEWAX 2.5 GM (ETHICON)

1044	MedicalSupplies	BP ACCESSORIES W/ INFLATION BAG, ADULT
1045	MedicalSupplies	BP ACCESSORIES W/ INFLATION BAG, NEONATE
1046	MedicalSupplies	BP ACCESSORIES W/ INFLATION BAG, PEDIA
1047	MedicalSupplies	BP CUFFS, INFANT
1048	MedicalSupplies	BREATHING CIRCUIT SINGLE LIMB 72" WITH BREATHING BAG, ADULT
1049	MedicalSupplies	BUBBLE CPAP CANNULA SIZE 0
1050	MedicalSupplies	BUBBLE CPAP CANNULA SIZE 1
1051	MedicalSupplies	BUBBLE CPAP CANNULA SIZE 2
1052	MedicalSupplies	BUFFER REAGENT, FOR 248 ABG RAPID LAB ABG MACHINE (4'S)
1053	MedicalSupplies	CADAVER BAG WITH ZIPPER- ADULT
1054	MedicalSupplies	CADAVER BAG WITH ZIPPER-PEDIA
1055	MedicalSupplies	CAH SURGICAL CLIPPERS STANDARD PREP. BLADE
1056	MedicalSupplies	CAL LC 2G
1057	MedicalSupplies	CALCIUM HYDROXIDE (LIGHT CURE) 2.5G
1058	MedicalSupplies	CALIBRATED GLASS, PLASTIC, 210CC
1059	MedicalSupplies	CAP SURGICAL
1060	MedicalSupplies	CATARACT PACK
1061	MedicalSupplies	CATHETER CONDOM, S- LARGE
1062	MedicalSupplies	CATHETER CONDOM, S- MEDIUM
1063	MedicalSupplies	CATHETER EMBELECTOMY FR 4
1064	MedicalSupplies	CATHETER EMBELECTOMY FR 5 EMBOLOUS/THROMBUS EXTRACTOR STERILE HYPOALLERGENIC
1065	MedicalSupplies	CATHETER EMBELECTOMY FR 8 EMBOLOUS/THROMBUS EXTRACTOR STERILE HYPOALLERGENIC
1066	MedicalSupplies	CATHETER EMBELECTOMY FR. 3 EMBOLOUS/THROMBUS EXTRACTOR STERILE HYPOALLERGENIC
1067	MedicalSupplies	CATHETER EMBOLECTOMY FR.6
1068	MedicalSupplies	CATHETER EPIDURAL G.16 W/EPIDURAL NEEDLE KIT
1069	MedicalSupplies	CATHETER EPIDURAL W/ EPIDURAL NEEDLE KIT G.18 "PERIFIX"
1070	MedicalSupplies	CATHETER FOLEY 2 WAY FR. 12
1071	MedicalSupplies	CATHETER FOLEY 2 WAY FR. 12 W/SILICONE
1072	MedicalSupplies	CATHETER FOLEY 2 WAY FR. 14
1073	MedicalSupplies	CATHETER FOLEY 2 WAY FR. 14 W/SILICONE
1074	MedicalSupplies	CATHETER FOLEY 2 WAY FR. 16
1075	MedicalSupplies	CATHETER FOLEY 2 WAY FR. 16 W/SILICONE
1076	MedicalSupplies	CATHETER FOLEY 2 WAY FR. 18 W/SILICONE
1077	MedicalSupplies	CATHETER FOLEY 2 WAY FR. 20, 30-100CC CAP, BALLOON, LATEX
1078	MedicalSupplies	CATHETER FOLEY 2 WAY FR. 24, 30-100CC CAP, BALLOON, LATEX
1079	MedicalSupplies	CATHETER FOLEY 2-WAY FR 10, 15-30CC CAP BALLOON LATEX 100% SILICON
1080	MedicalSupplies	CATHETER FOLEY 2-WAY FR 8
1081	MedicalSupplies	CATHETER FOLEY 2-WAY FR.18 15-30CC CAP. BALLOON, LATEX
1082	MedicalSupplies	CATHETER FOLEY 3 WAY FR 22, 30-100CC
1083	MedicalSupplies	CATHETER FR. 4 8CM WITH ECG LEAD AND SAFESITE (PEDIATRIC DOUBLE LUMEN CENTRAL VENOUS)
1084	MedicalSupplies	CATHETER HEMODIALYSIS FR 11.5, CURVED
1085	MedicalSupplies	CATHETER HEMODIALYSIS FR 12, CURVED
1086	MedicalSupplies	CATHETER INTRA JUGULAR LUMER FR. 11-12X20CM
1087	MedicalSupplies	CATHETER INTRA JUGULAR TRIPLE LEMEN FR. 11.5X15CM PRE CURVED "MEDCOMP"
1088	MedicalSupplies	CATHETER INTRA JUGULAR VASCULAR DOUBLE LUMEN FR.12-14X20- 25CM (LEFT)
1089	MedicalSupplies	CATHETER INTRA JUGULAR VASCULAR DOUBLE LUMEN RIGHT SIDE FR.12-14X15-16CM

1090	MedicalSupplies	CATHETER INTRAJUGULAR DOUBLE LUMEN FR. 14 X 15CM W/BUTTERFLY RIGHT SIDED PRE CURVED
1091	MedicalSupplies	CATHETER INTRAJUGULAR FR. 11.5
1092	MedicalSupplies	CATHETER IV DRESSING 10 X 12CM
1093	MedicalSupplies	CATHETER IV DRESSING 4.5-5X4-4.5CM
1094	MedicalSupplies	CATHETER IV DRESSING 6 X 7CM
1095	MedicalSupplies	CATHETER IV DRESSING 7 X 9CM
1096	MedicalSupplies	CATHETER IV G. 16 WINGS
1097	MedicalSupplies	CATHETER IV G. 18 WINGS
1098	MedicalSupplies	CATHETER IV G. 20 W/ WINGS
1099	MedicalSupplies	CATHETER IV G. 22 W/ WINGS
1100	MedicalSupplies	CATHETER IV G. 24 W/WINGS
1101	MedicalSupplies	CATHETER IV G. 26 W/WINGS
1102	MedicalSupplies	CATHETER NEPHOSTOMY PIGTAIL FR. 12
1103	MedicalSupplies	CATHETER NEPHROSTOMY, PIGTAIL FR.10
1104	MedicalSupplies	CATHETER PERMANENT LEFT SIDED FR. 14 X23 CM CUFF TOM TIP SPLIT TIP TYPE STYLET SHEATHLESS TECHNIQUE RED AND BLUE EXTENSION/ATTACHMENTS "MEDCOMP"
1105	MedicalSupplies	CATHETER PERITOCRAT STYLET
1106	MedicalSupplies	CATHETER PERITONEAL 2 CUFF 47CM
1107	MedicalSupplies	CATHETER PERITONEAL DIALYSIS 48 HOLES, 2CUFF
1108	MedicalSupplies	CATHETER PERMANENT RIGHT SIDED FR. 14 X 19CM CUFF TOM TIP SPLIT TIP STYLET SHEATHLESS TECHNIQUE RED AND BLUE EXTENSIONS/ATTACHMENTS "MEDCOMP"
1109	MedicalSupplies	CATHETER PERMANENT RIGHT SIDED FR. 14.5 X 23CM CUFF TOM TIP SPLIT TIP STYLET SHEATHLESS TECHNIQUE RED AND BLUE EXTENSIONS/ATTACHMENTS
1110	MedicalSupplies	CATHETER PERMANENT RIGHT SIDED FR. 14.5X 19CM CUFF TOM TIP SPLIT TIP STYLET SHEATHLESS TECHNIQUE RED AND BLUE EXTENSIONS/ATTACHMENTS
1111	MedicalSupplies	CATHETER PERMANENT RIGHT SIDED FR. 14-15X 19CM CUFF TOM TIP SPLIT TIP STYLET SHEATHLESS TECHNIQUE RED AND BLUE EXTENSIONS/ATTACHMENTS "MEDCOMP"
1112	MedicalSupplies	CATHETER RUBBER STRAIGHT FR 12
1113	MedicalSupplies	CATHETER RUBBER STRAIGHT FR. 14
1114	MedicalSupplies	CATHETER RUBBER STRAIGHT FR. 16
1115	MedicalSupplies	CATHETER STRAIGHT SIZE 10
1116	MedicalSupplies	CATHETER SUCTION FR. 10
1117	MedicalSupplies	CATHETER SUCTION FR. 12
1118	MedicalSupplies	CATHETER SUCTION FR. 14
1119	MedicalSupplies	CATHETER SUCTION FR. 16
1120	MedicalSupplies	CATHETER SUCTION FR. 18
1121	MedicalSupplies	CATHETER SUCTION FR. 5
1122	MedicalSupplies	CATHETER SUCTION FR. 8
1123	MedicalSupplies	CATHETER SUCTION SIZE 5
1124	MedicalSupplies	CATHETER TENKHOFF FR. 15X62-63CM COILED SILICONE
1125	MedicalSupplies	CATHETER THORACIC FR. 24
1126	MedicalSupplies	CATHETER THORACIC FR. 28
1127	MedicalSupplies	CATHETER THORACIC FR. 32
1128	MedicalSupplies	CATHETER THORACIC FR. 36
1129	MedicalSupplies	CATHETER UMBILICAL FR.3.5
1130	MedicalSupplies	CATHETER UMBILICAL FR.4.0
1131	MedicalSupplies	CATHETER UMBILICAL FR.5.0
1132	MedicalSupplies	CATHETER URETERAL OPEN TIP FR. 4 X 70CM

1133	MedicalSupplies	CATHETER URETERAL OPEN TIP FR. 6X70CM SHIFT DIAMETER DISP. RADIO OPAQUE CATH STERILE "RUSCH"
1134	MedicalSupplies	CATHETER URETRAL OPEN TIP FR. 5 X 70CM SHIFT DIAMETER DISPOSABLE RADIO OPAQUE CATH, STERILE "RUSCH"
1135	MedicalSupplies	CATHETER VENOUS CATHETER DUO FR.7 X 12 30CM
1136	MedicalSupplies	CATHETER VENOUS CATHETER TRIO FR.7 X 12 30CM
1137	MedicalSupplies	CATHETER, CENTRAL VENOUS, 2 WAY LUMEN 4FT X 8CM
1138	MedicalSupplies	CATHETER, CENTRAL VENOUS, 3 WAY LUMEN FR 5 X 8 CM
1139	MedicalSupplies	CATHETER, INDWELLING FOLEY (2 WAY) FR. 18
1140	MedicalSupplies	CATHETER, INDWELLING FOLEY CATHETER (2 WAY) FR. 10
1141	MedicalSupplies	CATHETER, INDWELLING FOLEY CATHETER (2 WAY) FR. 12
1142	MedicalSupplies	CATHETER, INDWELLING FOLEY CATHETER (2 WAY) FR. 14
1143	MedicalSupplies	CATHETER, INDWELLING FOLEY CATHETER (2 WAY) FR. 16
1144	MedicalSupplies	CATHETER, INDWELLING FOLEY CATHETER (2 WAY) FR. 8
1145	MedicalSupplies	CATHETER, INDWELLING FOLEY CATHETER (2-WAY) FR. 20
1146	MedicalSupplies	CATHETER, INDWELLING FOLEY FR. 22 (3 WAY)
1147	MedicalSupplies	CATHETER, INDWELLING FOLEY FR. 24 (3 WAY)
1148	MedicalSupplies	CATHETER, INTRA JUGULAR DOUBLE LUMEN FR. 11-11.5X20CM BUTTERFLY LEFT SIDED "MEDCOMP"
1149	MedicalSupplies	CATHETER, INTRAJUGULAR VASCULAR SPLIT STEAM TRIO CATHETER RIGHT SIDED 11.5F X 15CM PRE CURVED "MEDCOMP"
1150	MedicalSupplies	CATHETER, TENCKHOFF SWAN NECK CURL-ADULT
1151	MedicalSupplies	CAUTERY PENCIL, W/ CORD/DETACHABLE TIPS
1152	MedicalSupplies	CAUTERY TIPS (POINTED TIP)
1153	MedicalSupplies	CAUTERY TIPS, FINE "VALLEYLAB"
1154	MedicalSupplies	CAUTERY TIPS, LONG "VALLEYLAB"
1155	MedicalSupplies	CAVIWIPES DISINFECTANT WIPES)PER SHEET
1156	MedicalSupplies	CENTRAL VENOUS CATH. PEDIATRIC DOUBLE LUMEN FR.4, 8CM WITH ECG LEAD
1157	MedicalSupplies	CENTRAL VENOUS CATH.SET PEDIATRIC G22 X 10CM
1158	MedicalSupplies	CENTRAL VENOUS CATHETER PEDIATRIC JUGULAR CATH. SET, SINGLE LUMEN, G20 X 12CM
1159	MedicalSupplies	CENTRAL VENOUS CATHETER SET PEDIATRIC G24X9CM "ARROW"
1160	MedicalSupplies	CERVICAL COLLAR (SOFT) LARGE
1161	MedicalSupplies	CERVICAL COLLAR (SOFT) MEDIUM
1162	MedicalSupplies	CERVICAL COLLAR, PHILADELPHIA, LARGE SIZE
1163	MedicalSupplies	CERVICAL COLLAR, PHILADELPHIA, MEDIUM SIZE
1164	MedicalSupplies	CERVICAL COLLAR, PHILADELPHIA, SMALL SIZE
1165	MedicalSupplies	CHEMICAL INDICATOR LABEL FOR PLASMA STERILIZER 40M/ROLL
1166	MedicalSupplies	CHEMICAL INDICATOR STRIPS FOR PLASMA STERILIZER 250 STRIPS/BOX
1167	MedicalSupplies	CIRCULAR 29MM STAPLER
1168	MedicalSupplies	CIRCULAR STAPLER 28MM, EEA TM
1169	MedicalSupplies	CIRCULAR STAPLER 31MM, EEA TM
1170	MedicalSupplies	CIRCULAR STAPLER 33MM
1171	MedicalSupplies	CLOSED CIRCUIT SUCTION CATHETER FR. 16
1172	MedicalSupplies	CLOSED SUCTION RESERVOIR KIT 20CM X 7-10MM 8CM TUBING
1173	MedicalSupplies	COAGULATING ELECTRODE BALL SHAPED DIAMETER
1174	MedicalSupplies	COAGULATING ELECTRODE BALL SHAPED DIAMETER 3MM 24/26 IN YELLOW (SINGLE STEM)
1175	MedicalSupplies	CO-AXIAL CORE BIOPSY GUN G-18- WITH G. 17 INTRODUCER NEEDLE
1176	MedicalSupplies	CO-AXIAL CORE BIOPSY GUN G-20- WITH G. 19 INTRODUCER NEEDLE
1177	MedicalSupplies	"COLOR INK CARTRIDGE (CODONICS VIRTUA) FOR MRI (MUST BE COMPATIBLE WTH THE EXISTING PRINTER CODONICS VIRTUA)"

1178	MedicalSupplies	COLOSTOMY SET 24MM, WAFER CLIP "ACTIVE LIFE DRAINAGE POUCH 19-64MM/ 10'S
1179	MedicalSupplies	COLOSTOMY SET 32MM POUCH, WAFER CLIP
1180	MedicalSupplies	COLOSTOMY SET 45MM-50MM POUCH, WAFER CLIP
1181	MedicalSupplies	COLOSTOMY SET 57MM-60MM POUCH, WAFER CLIP
1182	MedicalSupplies	COLOSTOMY SET 70MM POUCH, WAFER CLIP
1183	MedicalSupplies	COMBINED SPINAL EPIDURAL NEEDLES
1184	MedicalSupplies	COMPLETE T-PIECE SET WITH APN
1185	MedicalSupplies	COMPOSITE RESIN VISIBLE FOR LIGHT CURE, 4.5G
1186	MedicalSupplies	COMPRESSION STOCKING MEDIUM
1187	MedicalSupplies	CONDITIONER
1188	MedicalSupplies	CONDITIONER KIT
1189	MedicalSupplies	CONNECTOR TUBING
1190	MedicalSupplies	CONTOUR CURVE CUTTER RELOAD COMPATIBLE WITH NTLC75
1191	MedicalSupplies	CONTOUR CURVED 40MM STAPLER
1192	MedicalSupplies	CONTOUR CUVETTE SUTTER STAPLER 2.5MM/64MM
1193	MedicalSupplies	CONTOUR LINEAR 40MM RELOAD
1194	MedicalSupplies	CONTOUR TYPE LOW PRESSURE SHUNT "MEDTRONIC"
1195	MedicalSupplies	CORD CLAMP
1196	MedicalSupplies	CORD CLAMP, 100'S
1197	MedicalSupplies	DISPOSABLE FLOURIDE TRAYS 50'S
1198	MedicalSupplies	DISPOSABLE INCENTIVE SPIROMETER "AIRLIFE"
1199	MedicalSupplies	DISPOSABLE NEEDLEGA.19 X1
1200	MedicalSupplies	DISPOSABLE NEEDLEGA.23
1201	MedicalSupplies	DISPOSABLE NEEDLEGA.25X1
1202	MedicalSupplies	DISPOSABLE NEEDLEGA.26
1203	MedicalSupplies	DISPOSABLE NEEDLEGA.27 1/2
1204	MedicalSupplies	DISPOSABLE NEEDLES GA.21GA.21X1
1205	MedicalSupplies	DISPOSABLE PEDIATRIC VENTILATOR TUBING (W/ 2 WATER TRAP)
1206	MedicalSupplies	DISPOSABLE SURGICAL GOWN
1207	MedicalSupplies	DISPOSABLE SYRINGE 10ML /GAUGE 21 NEEDLE
1208	MedicalSupplies	DISPOSABLE SYRINGE 3ML
1209	MedicalSupplies	DISPOSABLE SYRINGE50ML
1210	MedicalSupplies	DISPOSABLE SYRINGE5ML
1211	MedicalSupplies	DOCTOR J ALCOHOL250ML
1212	MedicalSupplies	DOCTOR J ALCOHOL70 % SOLUTION (PUMP)/500ML
1213	MedicalSupplies	DOUBLE J STENT 6 X 22
1214	MedicalSupplies	DOUBLE J URETERAL STENT SIZE 6 FR 22-32
1215	MedicalSupplies	DOUBLE J-STENT SINGLE PACKAGE 6 X 24
1216	MedicalSupplies	DOUBLE J-STENT SINGLE PACKAGE 6 X 26
1217	MedicalSupplies	DRAIN, CLOSED SUCTION
1218	MedicalSupplies	DRAINAGE POUCH, 57MM
1219	MedicalSupplies	DRESSING - HYDROCOLLOID CONTROLLED 20X20
1220	MedicalSupplies	DRESSING PAD - HYDROCOLLOID 8 X 8 (DUODERM)
1221	MedicalSupplies	DRESSING PAD - TEGADERM 10X12
1222	MedicalSupplies	DRESSING PAD - TEGADERM 9X15
1223	MedicalSupplies	DRESSING PAD - TEGADERM 9X25
1224	MedicalSupplies	DRESSING PAD - TEGADERM IV W/ NOTCHED 8.5CM X 7CM
1225	MedicalSupplies	DRESSING PAD - TEGADERM/YAFOH 9X20

1226	MedicalSupplies	DRESSING PAD- HYDROCOLLOID 4X4, THIN
1227	MedicalSupplies	DRESSING PAD, HYDROCOLOID (LARGE) 8X8
1228	MedicalSupplies	DRESSING PAD-TEGADERM IV WITH NOTCHED 5CMX5,7CM
1229	MedicalSupplies	DRESSING TRANSPARENT FILM FRAME STYLE 10CM X 12CM
1230	MedicalSupplies	DRESSING TRANSPARENT FILM FRAME STYLE 6CM X 7CM
1231	MedicalSupplies	DUAL ANGEL BARB SUTURE 26MM ROUND NEEDLE 30CM 2.0 "V-LOC"
1232	MedicalSupplies	DUAL ANGEL BARB SUTURE 26MM ROUND NEEDLE 30CM 3.0 "V-LOC"
1233	MedicalSupplies	DUAL ANGEL BARB SUTURE 26MM ROUND NEEDLE 30CM 4.0 "V-LOC"
1234	MedicalSupplies	DUAL ANGEL BARB SUTURE 37MM ROUND NEEDLE 30CM 0 "V-LOC"
1235	MedicalSupplies	ECG ELECTRODES ADULT
1236	MedicalSupplies	EGIA 4 HANDLE STANDARD
1237	MedicalSupplies	EGIA 60MED THICK SULU
1238	MedicalSupplies	ELASTIC BANDAGE 2"
1239	MedicalSupplies	ELASTIC BANDAGE 2" X 5 YARDS
1240	MedicalSupplies	ELASTIC BANDAGE 3"
1241	MedicalSupplies	ELASTIC BANDAGE 3" X 5 YARDS
1242	MedicalSupplies	ELASTIC BANDAGE 4"
1243	MedicalSupplies	ELASTIC BANDAGE 4" X 5" YARDS
1244	MedicalSupplies	ELASTIC BANDAGE 6' X 5 YARDS
1245	MedicalSupplies	ELASTIC BANDAGE 6"
1246	MedicalSupplies	ELASTOMETRIC INFUSION PUMP, DISPOSABLE (EASY PUMP)
1247	MedicalSupplies	ELECTRO GROUNDING PAD W/CORD DISP., HYDROGEL ADHESIVE LATEX FREE COMP. FOR PEDIA
1248	MedicalSupplies	ELECTRO GROUNDING PAD, BI-POLAR W/CORD DISP. HYDROGEL, ADULT
1249	MedicalSupplies	ELECTRO SURGICAL GROUNDING PAD W/OUT CORD DISP
1250	MedicalSupplies	ELECTRODE CABLE FOR ELECTRICAL STIMULATOR
1251	MedicalSupplies	ELECTRODE CLEANING SOLUTION 6 X 8ML
1252	MedicalSupplies	ELECTRODE, FOR ABG MACHINE
1253	MedicalSupplies	ENDO GIA STAPLER UNIVERSAL
1254	MedicalSupplies	ENDOCATH GOLD SPECIMEN RETRIEVAL
1255	MedicalSupplies	ENDOCLIP MEDIUM/LARGE APPLIER
1256	MedicalSupplies	ENDOILLUMINATOR
1257	MedicalSupplies	ENDOTRACHEAL INTUBATION TRAY SIZE 7.0
1258	MedicalSupplies	ENDOTRACHEAL INTUBATION TRAY SIZE 7.5
1259	MedicalSupplies	ENDOTRACHEAL TUBE 2.5
1260	MedicalSupplies	ENDOTRACHEAL TUBE 3.0
1261	MedicalSupplies	ENDOTRACHEAL TUBE 3.5
1262	MedicalSupplies	ENDOTRACHEAL TUBE 4.0
1263	MedicalSupplies	ENDOTRACHEAL TUBE 4.5
1264	MedicalSupplies	ENDOTRACHEAL TUBE 5.0
1265	MedicalSupplies	ENDOTRACHEAL TUBE 5.5
1266	MedicalSupplies	ENDOTRACHEAL TUBE 6.0
1267	MedicalSupplies	ENDOTRACHEAL TUBE 6.5
1268	MedicalSupplies	ENDOTRACHEAL TUBE 7.0
1269	MedicalSupplies	ENDOTRACHEAL TUBE 7.5
1270	MedicalSupplies	ENDOTRACHEAL TUBE WITH SUBGLOTTIC SIZE 7.5
1271	MedicalSupplies	ENEMA CAN 1L CPA STAINLESS AUTOCLAVABLE, OLIO IRRIGATOR 130 X 190MM
1272	MedicalSupplies	ENNO MOIST HOTPACKS, CERVICAL

1273	MedicalSupplies	ENNO MOIST HOTPACKS, STANDARD
1274	MedicalSupplies	ENSEAL ARTICULATING SEALER 35CM
1275	MedicalSupplies	ENSEAL G2 ART CURVE SEALER 35CM
1276	MedicalSupplies	ENTERAL FEEDING BOTTLE
1277	MedicalSupplies	ENTERAL FEEDING TUBING
1278	MedicalSupplies	EPIDURAL G.16
1279	MedicalSupplies	EPIDURAL G.18
1280	MedicalSupplies	EPIDURAL SET G. 20, PEDIA CATHETHER G. 24
1281	MedicalSupplies	EPIDURAL SET, COMBINED SPINAL G. 18
1282	MedicalSupplies	EQ NEWBORN DIAPER 22 PER PIECES/PACK
1283	MedicalSupplies	EQ NEWBORNDIAPER/12
1284	MedicalSupplies	ET DOUBLE LUMEN FR 37
1285	MedicalSupplies	EXPANDED NEWBORN SCREENING
1286	MedicalSupplies	EXTRACTOR STONE BASKET METAL
1287	MedicalSupplies	F2 8CM TAPER 2.3MM "MEDTRONIC"
1288	MedicalSupplies	FACE MASK (EAR LOOP) DISPOSABLE
1289	MedicalSupplies	FACE SHIELD
1290	MedicalSupplies	FACESHIELD DISPOSABLE
1291	MedicalSupplies	FETAL MONITOR PAPER, FOR FM 20
1292	MedicalSupplies	FIBER GLASS 3"
1293	MedicalSupplies	FIBER GLASS 3X4 CAST
1294	MedicalSupplies	FIBER GLASS 4X4, CAST
1295	MedicalSupplies	FILTER LINE SET, MICROSTREAM
1296	MedicalSupplies	FILTER PAPER FOR NEWBORN SCREENING
1297	MedicalSupplies	FILTERED VENOSSET
1298	MedicalSupplies	FISTULA KIT
1299	MedicalSupplies	FIT TESTING KIT "3M"
1300	MedicalSupplies	FLOURIDE (TROPICAL)
1301	MedicalSupplies	FOAM DRESSING NON-ADHESIVE, 5MM 10X60MM
1302	MedicalSupplies	FOAM DRESSING NON-ADHESIVE, 7MM 4X70MM
1303	MedicalSupplies	FOLDABLE LENS POSTERIOR CHAMBER POWER 18.5
1304	MedicalSupplies	FOLDABLE LENS POSTERIOR CHAMBER POWER 19.00
1305	MedicalSupplies	FOLDABLE LENS POSTERIOR CHAMBER POWER 20.0
1306	MedicalSupplies	FOLDABLE LENS POSTERIOR CHAMBER POWER 20.5
1307	MedicalSupplies	FOLDABLE LENS POSTERIOR CHAMBER POWER 21.0
1308	MedicalSupplies	FOLDABLE LENS POSTERIOR CHAMBER POWER 21.5
1309	MedicalSupplies	FOLDABLE LENS POSTERIOR CHAMBER POWER 22.5
1310	MedicalSupplies	FOLDABLE LENS POSTERIOR CHAMBER POWER 23.0
1311	MedicalSupplies	FOLDABLE LENS POSTERIOR CHAMBER POWER 23.5
1312	MedicalSupplies	FOLDABLE LENS POSTERIOR CHAMBER POWER 26.5
1313	MedicalSupplies	FORCED AIR BLANKET "SMITH'S MEDICAL"
1314	MedicalSupplies	HYDROCELLULAR FOAM DRESSING, ADHESIVE 12.5 X 12.5CM 10'S ALLEVYN ADHESIVE
1315	MedicalSupplies	HYDROCELLULAR FOAM DRESSING, ADHESIVE 17.5 X 17.5CM 10'S ALLEVYN ADHESIVE
1316	MedicalSupplies	HYDROCELLULAR FOAM DRESSING, ADHESIVE 22.5 X 22.5CM 10'S ALLEVYN ADHESIVE
1317	MedicalSupplies	HYDROCELLULAR FOAM DRESSING, ADHESIVE 7.5 X 7.5CM 10'S ALLEVYN ADHESIVE
1318	MedicalSupplies	INCISED DRAPE, ANTIMICROBIAL 15 X 28, 10'S
1319	MedicalSupplies	INCISED DRAPE, ANTIMICROBIAL 15.5-35CM X 8.5-35CM

1320	MedicalSupplies	INCISED DRAPE, ANTIMICROBIAL 30 X 28
1321	MedicalSupplies	INCISED DRAPE, ANTIMICROBIAL 35-45 X 28-35CM, 10'S 3M
1322	MedicalSupplies	INCISED DRAPE, ANTIMICROBIAL 40-56 X 42-45CM, 10'S 3M
1323	MedicalSupplies	INCISED DRAPE, ANTIMICROBIAL 44-56CM X 45-55CM
1324	MedicalSupplies	INCISED DRAPE, ANTIMICROBIAL 50-56CM X 84-85CM
1325	MedicalSupplies	INCISED DRAPE, ANTIMICROBIAL 55-56 X 40-45CM, 10'S 3M
1326	MedicalSupplies	NON-ADHESIVE WOUND DRESSING 10CM X 20CM 10'S
1327	MedicalSupplies	NON-CONDUCTIVE RESPIRATORY BACTERIAL FILTER"AIRLIFE"
1328	MedicalSupplies	NON-CORING ACCESS PORT HUBER NEEDLE WITH FIXATION BASE, G20X15MM
1329	MedicalSupplies	NON-CORING ACCESS PORT HUBER NEEDLE WITH FIXATION BASE, G20X25(1MM)
1330	MedicalSupplies	NON-CORING ACCESS PORT HUBER NEEDLE WITH FIXATION BASE, G22 X25(1MM)
1331	MedicalSupplies	NON-CORING ACCESS PORT HUBER NEEDLE WITH FIXATION BASE, G22X15MM
1332	MedicalSupplies	NON-CORING ACCESS PORT HUBER NEEDLE WITH FIXATION BASE, G22X20MM
1333	MedicalSupplies	NYLON 3-0
1334	MedicalSupplies	NYLON 4-0
1335	MedicalSupplies	OGT FR. 5 - 100 CM
1336	MedicalSupplies	OGT FR. 5 - 20CM-40CM
1337	MedicalSupplies	OGT FR. 8 - 100 CM
1338	MedicalSupplies	OGT FR. 8 -20CM- 40CM
1339	MedicalSupplies	OPHTHALMIC CRESCENT KNIFE 2.5DEGREES
1340	MedicalSupplies	OPHTHALMIC KNIFE 15 DEGREES
1341	MedicalSupplies	OPHTHALMIC KNIFE 2.2 DEGREES
1342	MedicalSupplies	OPTIFLOW NASAL CANNULA (LARGE) OPT846
1343	MedicalSupplies	OPTIFLOW NASAL CANNULA (MEDIUM) OPT844
1344	MedicalSupplies	OPTIFLOW NASAL CANNULA (SMALL) OPT842
1345	MedicalSupplies	OPTIFLOW NASAL CANNULA HIGH FLOW INFANT
1346	MedicalSupplies	OPTIFLOW NASAL CANNULA HIGH FLOW PEDIA
1347	MedicalSupplies	ORAL AIRWAY S-0, PEDIA
1348	MedicalSupplies	ORAL AIRWAY S-00, PEDIA
1349	MedicalSupplies	ORAL AIRWAY S-1, ADULT
1350	MedicalSupplies	ORAL AIRWAY S-1, PEDIA
1351	MedicalSupplies	ORAL AIRWAY S-2, ADULT (GREEN)
1352	MedicalSupplies	ORAL AIRWAY S-2, PEDIA
1353	MedicalSupplies	ORAL AIRWAY S-3, ADULT (YELLOW)
1354	MedicalSupplies	ORAL AIRWAY S-3, PEDIA
1355	MedicalSupplies	ORAL AIRWAY S-30MM
1356	MedicalSupplies	ORAL AIRWAY S-4, ADULT (YELLOW)
1357	MedicalSupplies	ORAL AIRWAY S-40MM
1358	MedicalSupplies	ORAL AIRWAY S-70MM
1359	MedicalSupplies	ORAL AIRWAY S-80MM
1360	MedicalSupplies	ORAL AIRWAY S-90MM
1361	MedicalSupplies	ORBITAL IMPLANT WITH CONFORMER PMMA SIZE 18MM
1362	MedicalSupplies	OXYGEN CANNULA (ADULT)
1363	MedicalSupplies	OXYGEN CANNULA (NEONATAL)
1364	MedicalSupplies	OXYGEN CANNULA (PEDIA)
1365	MedicalSupplies	OXYGEN CONNECTING TUBE
1366	MedicalSupplies	OXYGEN FACE MASK ADULT

1367	MedicalSupplies	OXYGEN FACE MASK NEONATE
1368	MedicalSupplies	OXYGEN FACE MASK PEDIA
1369	MedicalSupplies	OXYGEN MASK W/TUBING, HYPOALLERGENIC AND TIE ON DISP, ADULT
1370	MedicalSupplies	OXYGEN MASK W/TUBING, HYPOALLERGENIC AND TIE ON DISP, PEDIA
1371	MedicalSupplies	OXYGEN MASK, WITH RESERVIOR BAG, ADULT
1372	MedicalSupplies	OXYGEN SENSOR, NELLCOR (PEDIA)
1373	MedicalSupplies	PAD
1374	MedicalSupplies	PAMPERS NEWBORN X 40
1375	MedicalSupplies	PAMPERS NEWBORN X 20`S
1376	MedicalSupplies	PAP SMEAR KIT
1377	MedicalSupplies	PARAFFIN WAX, W/ MINERAL OIL
1378	MedicalSupplies	PATIENT MOVER COVER 170X50 NYLON
1379	MedicalSupplies	PATIENT TRANSFER STATIC INSULATIVE RADIO LUCENT (PATSLIDE)
1380	MedicalSupplies	PD DISINFECTION MINI CAP WITH POVIDONE IODINE
1381	MedicalSupplies	PD DISINFECTION MINI CAP WITH POVIDONE IODINE SET
1382	MedicalSupplies	PD EXTENSION CATHETER
1383	MedicalSupplies	PEDIATRIC URINE COLLECTOR
1384	MedicalSupplies	PENROSE DRAIN 1/2
1385	MedicalSupplies	PENROSE DRAIN 1/4
1386	MedicalSupplies	PENROSE DRAIN 3/4
1387	MedicalSupplies	PENROSE DRAIN 3/8
1388	MedicalSupplies	PENROSE DRAIN S -1"
1389	MedicalSupplies	PER PIECE O ₂ ELECTRODE
1390	MedicalSupplies	PERCUTANEOUS NEPHROSTOMY PUNCTURE SET FR.12
1391	MedicalSupplies	PERCUTANEOUS NEPHROSTOMY PUNCTURE SET FR.14
1392	MedicalSupplies	PERCUTANEOUS NEPHROSTOMY SET
1393	MedicalSupplies	PERITONEAL DIALYSIS CATHETER FR. 15 X 47CM WITH 2 CUFFS , STRAIGHT " ARGYLE"
1394	MedicalSupplies	PERITONEAL DIALYSIS CATHETER FR. 15 X 31-32CM WITH 2 CUFFS , STRAIGHT " ARGYLE"
1395	MedicalSupplies	PERITONEAL DIALYSIS SOLUTION 2.3% 2L (1.75MMOL/L DIALYSATE CALCIUM)
1396	MedicalSupplies	PERMANENT HEMODIALYSIS CATHETER 14.5 FR. 19CM IMPLANT LENGTH & 36CM OVERALL
1397	MedicalSupplies	PERMANENT HEMODIALYSIS CATHETER 14.5 FR. 23CM IMPLANT LENGTH & 40CM OVERALL
1398	MedicalSupplies	PERMANENT HEMODIALYSIS CATHETER, 14.5-15.5FR X 23CM IMPLANT LENGTH & 40CM OVERALL LENGTH, WITH LASER-CUT SIDE SLOTS & SYMMETRICAL /STAGGERED TIP, CARBOTHANE MATERIAL, STRAIGHT EXTENSIONS WITH SHEATHLESS TECHNIQUE
1399	MedicalSupplies	PFT MOUTHPIECE WITH FILTER
1400	MedicalSupplies	PH ELECTRODE FOR 248 ABG MACHINE
1401	MedicalSupplies	PH ELECTRIC FILL SOLUTION (3UDC)
1402	MedicalSupplies	PHOTOTHERAPY EYE SHIELD
1403	MedicalSupplies	PLASTER - SILK ACETATE, (DURAPORE) 2
1404	MedicalSupplies	PLASTER OF PARIS 4"
1405	MedicalSupplies	PLASTER OF PARIS 4" X 5 YARDS
1406	MedicalSupplies	PLASTER OF PARIS 6" X 5 YARDS
1407	MedicalSupplies	PLASTER OF PARIS 6" X 5 YARDS (LONG BONE)
1408	MedicalSupplies	PLASTER SILK ACETATE - DURAPORE 1" PER 60 INCHES
1409	MedicalSupplies	PLASTER TRANSPARENT, (TRANSPORE) 1" PER 60 INCHES
1410	MedicalSupplies	PLATELET ADMINISTRATION SET

1411	MedicalSupplies	PLEUR EVAC CHEST DRAINAGE SYSTEM
1412	MedicalSupplies	PO2 ELECTRODE
1413	MedicalSupplies	POLYESTER 1,
1414	MedicalSupplies	POLYESTER 2, TAPERCUT
1415	MedicalSupplies	POLYESTER 5-0, DOUBLE SPATULATED
1416	MedicalSupplies	POLYGLACTIN 2-0 W/ "26MM" SH ROUND NEEDLE ORTHO/SURGE
1417	MedicalSupplies	POLYGLACTIN 2-0 W/ "36MM" CT-1ROUND NEEDLE OB
1418	MedicalSupplies	POLYGLACTIN 2-0 W/ UR-6 ROUND NEEDLE
1419	MedicalSupplies	POLYGLACTIN 2-0,CUTTING
1420	MedicalSupplies	POLYGLACTIN 3-0 ROUND NEEDLE
1421	MedicalSupplies	POLYGLACTIN 4-0 CUTTING
1422	MedicalSupplies	POLYGLACTIN 4-0 ROUND NEEDLE
1423	MedicalSupplies	POLYGLECAPRONE 3-0 CUTTING
1424	MedicalSupplies	POLYGLECAPRONE 4-0 CUTTING
1425	MedicalSupplies	POLYPROPELENE MESH NON-ABSORBABLE SYNTHETIC 8X15CM/BOX OF 6 SHEETS
1426	MedicalSupplies	POLYPROPYLENE 5.0 SYNTHETIC, NON ABSORBABLE MONOFILAMENT SUTURE 90CM
1427	MedicalSupplies	POOLE DRAIN ABDOMINAL SUCTION, SMOOTH TIP, YANKUER
1428	MedicalSupplies	PREGNA IUD - COPPER T, TW 38oA 1'S
1429	MedicalSupplies	SURGICAL NU KNIT 3" X 4" HEMOSTATIC OXIDIZED REGENERATED 12'S
1430	MedicalSupplies	SUTURE NYLON 10- 0
1431	MedicalSupplies	SUTURE NYLON 2-0
1432	MedicalSupplies	SUTURE NYLON 3-0
1433	MedicalSupplies	SUTURE NYLON 4-0
1434	MedicalSupplies	SUTURE NYLON 6-0
1435	MedicalSupplies	SUTURE NYLON 8-0 DOUBLE ARM
1436	MedicalSupplies	SUTURE NYLON 9-0 TAPER POINT SPATULATED NEEDLE
1437	MedicalSupplies	SUTURE POLYDIOXANONE 0, PDS PLUS 0
1438	MedicalSupplies	SUTURE POLYDIOXANONE 6-0 DOUBLE ARM, PDS II
1439	MedicalSupplies	SUTURE POLYGLACTIN 0, ROUND NEEDLE
1440	MedicalSupplies	SUTURE POLYGLACTIN 1-0 ROUND
1441	MedicalSupplies	SUTURE POLYGLACTIN 3-0 ROUND NEEDLE
1442	MedicalSupplies	SUTURE POLYGLACTIN 4-0 ROUND NEEDLE
1443	MedicalSupplies	SUTURE POLYGLACTIN 5-0 RB1- SURGERY
1444	MedicalSupplies	SUTURE POLYGLACTIN 6-0, SPATULATED- OPHTHA
1445	MedicalSupplies	SUTURE POLYGLACTIN 7-0, SPATULATED
1446	MedicalSupplies	SUTURE POLYPROPYLENE 3-0 DOUBLE ARM
1447	MedicalSupplies	SUTURE POLYPROPYLENE 4-0 RB1 , 17MM DOUBLE ARM
1448	MedicalSupplies	SUTURE POLYPROPYLENE 5-0 DOUBLE ARM
1449	MedicalSupplies	SUTURE POLYPROPYLENE 6-0 BV1
1450	MedicalSupplies	SUTURE POLYPROPYLENE 6-0 RB2 DOUBLE ARM
1451	MedicalSupplies	SUTURE, POLYGLACTIN (VICRYL) 2-0
1452	MedicalSupplies	SYRINGE ASEPTO
1453	MedicalSupplies	SYRINGE DISPOSABLE (20 CC)
1454	MedicalSupplies	SYRINGE DISPOSABLE (3.0 CC)
1455	MedicalSupplies	SYRINGE DISPOSABLE (30 CC)
1456	MedicalSupplies	SYRINGE DISPOSABLE (5.0 CC)
1457	MedicalSupplies	SYRINGE DISPOSABLE (50 CC)

1458	MedicalSupplies	SYRINGE DISPOSABLE (10 CC)
1459	MedicalSupplies	SYRINGE DISPOSABLE 1 CC
1460	MedicalSupplies	SYRINGE DISPOSABLE INSULIN G-30
1461	MedicalSupplies	SYRINGE INSULIN 1CC
1462	MedicalSupplies	SYRINGE TOMEY 80CC SHARP (SLIP TIP)
1463	MedicalSupplies	SYRINGE, TOOMEY
1464	MedicalSupplies	TEST STRIP (CONTROL)
1465	MedicalSupplies	TEST STRIP (GLUCOSE)
1466	MedicalSupplies	TEST STRIP2(CONTROL)
1467	MedicalSupplies	THERABAND RED, MEDIUM RESISTANT X 2 YRDS
1468	MedicalSupplies	THERMAL GUN SCANNER, INFRARED THERMOMETER
1469	MedicalSupplies	THERMAL PAPER FOR AUTOCLAVE
1470	MedicalSupplies	THORACIC TROCAR 10.5MM FOR INSTRUMENT UP TO 11MM, SINGLE USE TROCAR, NON-CONDUCTIVE SLEEVE "COVIDIEN"
1471	MedicalSupplies	THORACIC TROCAR 11.5MM FOR INSTRUMENT UP TO 12MM, SINGLE USE TROCAR, NON- CONDUCTIVE SLEEVE "COVIDIEN"
1472	MedicalSupplies	THORACIC TROCAR 15MM FOR INSTRUMENT UP TO 16.5MM, SINGLE USE TROCAR, NON- CONDUCTIVE SLEEVE "COVIDIEN"
1473	MedicalSupplies	THORACIC TROCAR 5.5MM FOR INSTRUMENT UP TO 6MM, SINGLE USE TROCAR NON-CONDUCTIVE SLEEVE "COVIDIEN"
1474	MedicalSupplies	THREE WAY STOPER PIECEOCK W/ 10 CM TUBING STERILE SINGLE
1475	MedicalSupplies	THREE WAY STOPER PIECEOCK W/ 5-7 CM TUBING STERILE SINGLE
1476	MedicalSupplies	TISSUE PAPER
1477	MedicalSupplies	TITANIUM LIGATING CLIPS, MEDIUM-LARGE
1478	MedicalSupplies	TITANIUM LIGATING CLIPS, POLYMER EXTRA LARGE 6 CLIPS/CART 14 CARTRIDGES
1479	MedicalSupplies	TITANIUM LIGATING CLIPS, POLYMER LARGE 6 CLIPS/CART, 14 CARTRIDGES
1480	MedicalSupplies	TONGUE DEPRESSOR
1481	MedicalSupplies	T-PIECE SET UP W/NEB
1482	MedicalSupplies	T-PIECE, VENTURI SETUP W/ FLOW METER ADAPTOR 500ML
1483	MedicalSupplies	TRACHE MASK "AIRLIFE"
1484	MedicalSupplies	TRACHEOSTOMY TUBE WITHOUT FENESTRATION AND WITHOUT INNER CANNULA I.D.3.0MM;O.D. 4.5 MM;LENGTH: 39MM SIZE 3.0 (PEDIA CUFFLESS)
1485	MedicalSupplies	TUBING SET ENDOMAT STERILE 10'S (LAPAROSCOPY)
1486	MedicalSupplies	TUBING SET HYSTEROMAT II STERILE 10'S
1487	MedicalSupplies	WOUND PROTECTOR 2.5-6CM 40'S SMALL
1488	MedicalSupplies	WOUND PROTECTOR 2-4CM (EXTRA SMALL)
1489	MedicalSupplies	WOUND PROTECTOR 5X9CM 40'S MEDIUM
1490	MedicalSupplies	PLASTER ADHESIVE POROUS TAPE 1" X 10 YARDS, 12'S/BX
1491	MedicalSupplies	PLASTER OF PARIS 4X4 YARDS
1492	MedicalSupplies	PLASTER OF PARIS 6 X 5YRDS
1493	MedicalSupplies	PLASTER, ADHESIVE PRE-CUT
1494	MedicalSupplies	POLYVINYL FILM "CLING WRAP 12X300M"
1495	MedicalSupplies	POUCHES, 350MMX80MMX100MM 14" X 3.5" X 330" (FOR GAS STERILIZATION)
1496	MedicalSupplies	RUBBER SHEET 25 YRDS REUSABLE, WASHABLE, STAINED FREE
1497	MedicalSupplies	RUBBER TUBING 1/4 X 1/6 X 3/8 50FT, SMALL
1498	MedicalSupplies	RUBBER TUBING 1/4 X 1/8 X 50FT, MEDIUM
1499	MedicalSupplies	SURGICAL PAPER TAPE 2X10 YARD (5CMX9.1INCH)
1500	MedicalSupplies	THERMAL PAPER
1501	MedicalSupplies	THERMAL PAPER FOR ECG (210MM X 140MMX143 SHEET)
1502	MedicalSupplies	THERMAL PAPER FOR STEAM STERILIZER

1503	MedicalSupplies	THERMAL PAPER II HD
1504	MedicalSupplies	THERMAL PAPER TYPE V (UPP-110HD 110MMX20M) COMPATIBLE WITH EXISTING PRINTER
1505	MedicalSupplies	TYVEC PAPER/POUCHES 100MM (W) X 80M
1506	MedicalSupplies	TYVEC PAPER/POUCHES 150MM (W) X 80M
1507	MedicalSupplies	TYVEC PAPER/POUCHES 200MM (W) X 80M
1508	MedicalSupplies	TYVEC PAPER/POUCHES 300MM (W) X 80M
1509	MedicalSupplies	TYVEK PAPER/POUCHES 250MM(W) X 70-80M "PMS"
1510	MedicalSupplies	TYVEK PAPER/POUCHES 75MM(W) X 70-80M "PMS"
1511	MedicalSupplies	WADDING SHEET, 4X6 YARDS 100% COTTON
1512	MedicalSupplies	WADDING SHEET, 5X6 YARDS 100% COTTON
1513	MedicalSupplies	ECG RECORDING PAPER SINGLE CHANNEL 50 X 30
1514	MedicalSupplies	LUBRICATING JELLY, 3 G/SACHET
1515	MedicalSupplies	STERILIZATION FEE (STEAM) 2X4X6
1516	MedicalSupplies	STERILIZATION FEE (STEAM) 4X8X12
1517	MedicalSupplies	BAG, COLOSTOMY WITH WAFER AND CLIP SIZE 32-40MM
1518	MedicalSupplies	BLOOD PRESSURE BULB, NIBP AIR RELEASE VALVE AND SMALL BULB DISPOSABLE BONE MARROW ASPIRATION SET SPECIAL PROC JAMSHID
1519	MedicalSupplies	11G X 10CM / BONE MARROW BIOSY
1520	MedicalSupplies	INCENTIVE SPIROMETER
1521	MedicalSupplies	JACKSON REESE SYSTEM, 1000ML
1522	MedicalSupplies	JACKSON REESE SYSTEM, 500ML
1523	MedicalSupplies	PPE FOR QMMC EMPLOYEES
1524	MedicalSupplies	PPE FOR QMMC EMPLOYEES FOR ANCILLARY
1525	MedicalSupplies	TRACHEOSTOMY KIT SIZE 6
1526	MedicalSupplies	TRACHEOSTOMY KIT SIZE 8
1527	MedicalSupplies	TRACHEOSTOMY TUBE, PEDIATRIC SIZE 4.5
1528	MedicalSupplies	TRI-PORT NEEDLELESS
1529	MedicalSupplies	WAX PAPER
1530	MedicalSupplies	PERACETIC ACID POTENCY STRIP (GENERIC)
1531	MedicalSupplies	PRESEPT- DISINFECTANT TABLET
1532	MedicalSupplies	POWER INJECTOR KIT (190ML SYRINGE W/QFT)
1533	MedicalSupplies	COMPOSITE RESIN VISIBLE FOR LIGHT CURE 3.5GMS.
1534	MedicalSupplies	EMAIL PREPARATORY BLUE 10ML. (ACID ETCH)
1535	MedicalSupplies	FLOWABLE COMPOSITE
1536	MedicalSupplies	INTERFACE 1.2ML TEMREX CAVITY LINER
1537	MedicalSupplies	POLIHAXANIDE AND BETAINE WOUND GEL 250GRAMS/TUBE "B BRAUN PTONTSAN SOLUTION"
1538	MedicalSupplies	THREE LINE LABELER WITH 12 ALPHANUMERIC CHARACTER
1539	MedicalSupplies	HEPARIN SODIUM
1540	LaboratoryAndDiagnostic	ABG: ABG
1541	LaboratoryAndDiagnostic	ABG: VBG
1542	LaboratoryAndDiagnostic	ABG: COVID PANEL: ABG
1543	LaboratoryAndDiagnostic	HEMATOLOGY: PTT
1544	LaboratoryAndDiagnostic	HEMATOLOGY: BLEEDING TIME (BT)
1545	LaboratoryAndDiagnostic	HEMATOLOGY: CBC
1546	LaboratoryAndDiagnostic	HEMATOLOGY: CBC WITH VCS (VOLUME CONDUCTIVITY SCATTER)
1547	LaboratoryAndDiagnostic	HEMATOLOGY: CLOT RETRACTION TIME (CRT)
1548	LaboratoryAndDiagnostic	HEMATOLOGY: CLOTTING TIME (CT)
1549	LaboratoryAndDiagnostic	HEMATOLOGY: D-DIMER

1550	LaboratoryAndDiagnostic	HEMATOLOGY: DIFFERENTIAL COUNT
1551	LaboratoryAndDiagnostic	HEMATOLOGY: ESR
1552	LaboratoryAndDiagnostic	HEMATOLOGY: FACTOR VIII
1553	LaboratoryAndDiagnostic	HEMATOLOGY: FACTOR IX
1554	LaboratoryAndDiagnostic	HEMATOLOGY: FILARIAL SMEAR
1555	LaboratoryAndDiagnostic	HEMATOLOGY: HCT
1556	LaboratoryAndDiagnostic	HEMATOLOGY: CSF-HCT
1557	LaboratoryAndDiagnostic	HEMATOLOGY: PLEURAL FLUID-HCT
1558	LaboratoryAndDiagnostic	HEMATOLOGY: HGB
1559	LaboratoryAndDiagnostic	HEMATOLOGY: LE CELL PREPARATION
1560	LaboratoryAndDiagnostic	HEMATOLOGY: LUPUS SCREEN
1561	LaboratoryAndDiagnostic	HEMATOLOGY: DRVV SCREEN
1562	LaboratoryAndDiagnostic	HEMATOLOGY: LUPUS CONFIRM
1563	LaboratoryAndDiagnostic	HEMATOLOGY: DRVV CONFIRM
1564	LaboratoryAndDiagnostic	HEMATOLOGY: MALARIAL SMEAR
1565	LaboratoryAndDiagnostic	HEMATOLOGY: OSMOTIC FRAGILITY TEST (OFT)
1566	LaboratoryAndDiagnostic	HEMATOLOGY: PLATELET COUNT
1567	LaboratoryAndDiagnostic	HEMATOLOGY: PROTEIN C
1568	LaboratoryAndDiagnostic	HEMATOLOGY: PROTEIN S
1569	LaboratoryAndDiagnostic	HEMATOLOGY: PROTIME
1570	LaboratoryAndDiagnostic	HEMATOLOGY: RBC/WBC MORPHOLOGY (PBS)
1571	LaboratoryAndDiagnostic	HEMATOLOGY: PBS, RETICULOCYTE COUNT
1572	LaboratoryAndDiagnostic	HEMATOLOGY: RETICULOCYTE COUNT
1573	LaboratoryAndDiagnostic	HEMATOLOGY: SILICA CLOTTING TIME (SCT)
1574	LaboratoryAndDiagnostic	HEMATOLOGY: WHITE BLOOD CELL COUNT
1575	LaboratoryAndDiagnostic	HEMATOLOGY: FIBRINOGEN
1576	LaboratoryAndDiagnostic	HEMATOLOGY: BAT PANEL: CBC WITH PC, PT/PTT
1577	LaboratoryAndDiagnostic	HEMATOLOGY: COVID PANEL: CBC WITH PC, ESR, PT/PTT, D-DIMER
1578	LaboratoryAndDiagnostic	CHEMISTRY: SERUM ALBUMIN
1579	LaboratoryAndDiagnostic	CHEMISTRY: ALKALINE PHOSPHATASE
1580	LaboratoryAndDiagnostic	CHEMISTRY: ALP
1581	LaboratoryAndDiagnostic	CHEMISTRY: ALT
1582	LaboratoryAndDiagnostic	CHEMISTRY: AMMONIA
1583	LaboratoryAndDiagnostic	CHEMISTRY: AMYLASE
1584	LaboratoryAndDiagnostic	CHEMISTRY: AST
1585	LaboratoryAndDiagnostic	CHEMISTRY: DIRECT BILIRUBIN (DB)
1586	LaboratoryAndDiagnostic	CHEMISTRY: TOTAL BILIRUBIN (TB)
1587	LaboratoryAndDiagnostic	CHEMISTRY: INDIRECT BILIRUBIN (IB)
1588	LaboratoryAndDiagnostic	CHEMISTRY: TB-DB-IB
1589	LaboratoryAndDiagnostic	CHEMISTRY: BUN
1590	LaboratoryAndDiagnostic	CHEMISTRY: SERUM CALCIUM
1591	LaboratoryAndDiagnostic	CHEMISTRY: CHLORIDE
1592	LaboratoryAndDiagnostic	CHEMISTRY: CK TOTAL
1593	LaboratoryAndDiagnostic	CHEMISTRY: CKMB
1594	LaboratoryAndDiagnostic	CHEMISTRY: CREATININE
1595	LaboratoryAndDiagnostic	CHEMISTRY: DRUG ASSAY-ACETAMINOPHEN
1596	LaboratoryAndDiagnostic	CHEMISTRY: DRUG ASSAY-DIGOXIN
1597	LaboratoryAndDiagnostic	CHEMISTRY: DRUG ASSAY-PHENOBARBITAL

1598	LaboratoryAndDiagnostic	CHEMISTRY: DRUG ASSAY-LITHIUM
1599	LaboratoryAndDiagnostic	CHEMISTRY: DRUG ASSAY-SALICYLATES
1600	LaboratoryAndDiagnostic	CHEMISTRY: DRUG ASSAY-VALPROIC ACID
1601	LaboratoryAndDiagnostic	CHEMISTRY: DRUG ASSAY-VANCOMYCIN
1602	LaboratoryAndDiagnostic	CHEMISTRY: SERUM ALCOHOL
1603	LaboratoryAndDiagnostic	CHEMISTRY: NA, K, CL
1604	LaboratoryAndDiagnostic	CHEMISTRY: NA, K, CL - URINE
1605	LaboratoryAndDiagnostic	CHEMISTRY: GGT
1606	LaboratoryAndDiagnostic	CHEMISTRY: HBA1C
1607	LaboratoryAndDiagnostic	CHEMISTRY: IONIZED CA
1608	LaboratoryAndDiagnostic	CHEMISTRY: LDH
1609	LaboratoryAndDiagnostic	CHEMISTRY: LIPASE
1610	LaboratoryAndDiagnostic	CHEMISTRY: HDL
1611	LaboratoryAndDiagnostic	CHEMISTRY: LIPID PROFILE
1612	LaboratoryAndDiagnostic	CHEMISTRY: TRIGLYCERIDES
1613	LaboratoryAndDiagnostic	CHEMISTRY: TOTAL CHOLESTEROL
1614	LaboratoryAndDiagnostic	CHEMISTRY: LIVER FUNCTION TEST PANEL (AST, ALT, ALP, LDH, TPAG, TBIL, DBIL, IBIL)
1615	LaboratoryAndDiagnostic	CHEMISTRY: MAGNESIUM
1616	LaboratoryAndDiagnostic	CHEMISTRY: NEONATAL BILIRUBIN - TB, DB, IB (0-15 DAYS OLD)
1617	LaboratoryAndDiagnostic	CHEMISTRY: PHOSPHOROUS
1618	LaboratoryAndDiagnostic	CHEMISTRY: PLASMA K
1619	LaboratoryAndDiagnostic	CHEMISTRY: POTASSIUM
1620	LaboratoryAndDiagnostic	CHEMISTRY: RAPID HEP-HBSAG (QUALITATIVE)
1621	LaboratoryAndDiagnostic	CHEMISTRY: RAPID HBSAG (QUALITATIVE)
1622	LaboratoryAndDiagnostic	CHEMISTRY: RAPID HEP-HCV (QUALITATIVE)
1623	LaboratoryAndDiagnostic	CHEMISTRY: RAPID HCV (QUALITATIVE)
1624	LaboratoryAndDiagnostic	CHEMISTRY: SODIUM
1625	LaboratoryAndDiagnostic	CHEMISTRY: FBS / RBS, 2 HRS PPBS
1626	LaboratoryAndDiagnostic	CHEMISTRY: OGCT (50 GMS)
1627	LaboratoryAndDiagnostic	CHEMISTRY: OGTT (100 GMS)
1628	LaboratoryAndDiagnostic	CHEMISTRY: OGTT (75 GMS)
1629	LaboratoryAndDiagnostic	CHEMISTRY: TCO ₂
1630	LaboratoryAndDiagnostic	CHEMISTRY: TOTAL PROTEIN
1631	LaboratoryAndDiagnostic	CHEMISTRY: TPAG/GLOBULIN
1632	LaboratoryAndDiagnostic	CHEMISTRY: URIC ACID
1633	LaboratoryAndDiagnostic	CHEMISTRY: CSF & OTHER BODY FLUIDS - ALBUMIN
1634	LaboratoryAndDiagnostic	CHEMISTRY: CSF & OTHER BODY FLUIDS - CHOLESTEROL
1635	LaboratoryAndDiagnostic	CHEMISTRY: CSF & OTHER BODY FLUIDS - LDH
1636	LaboratoryAndDiagnostic	CHEMISTRY: CSF & OTHER BODY FLUIDS - GLUCOSE
1637	LaboratoryAndDiagnostic	CHEMISTRY: CSF & OTHER BODY FLUIDS - TOTAL PROTEIN
1638	LaboratoryAndDiagnostic	CHEMISTRY: CSF & OTHER BODY FLUIDS - TOTAL PROTEIN & GLUCOSE (QUANTITATIVE)
1639	LaboratoryAndDiagnostic	CHEMISTRY: CSF & OTHER BODY FLUIDS - TRIGLYCERIDES
1640	LaboratoryAndDiagnostic	CHEMISTRY: DIRECT LDL
1641	LaboratoryAndDiagnostic	CHEMISTRY: ELECTROPHORESIS - HEMOGLOBIN
1642	LaboratoryAndDiagnostic	CHEMISTRY: ELECTROPHORESIS - SERUM PROTEIN
1643	LaboratoryAndDiagnostic	CHEMISTRY: ELECTROPHORESIS - URINE PROTEIN
1644	LaboratoryAndDiagnostic	CHEMISTRY: IMMUNOTYPING
1645	LaboratoryAndDiagnostic	CHEMISTRY: MA1B

1646	LaboratoryAndDiagnostic	CHEMISTRY: MICRAL
1647	LaboratoryAndDiagnostic	CHEMISTRY: MICROALBUMIN TEST
1648	LaboratoryAndDiagnostic	CHEMISTRY: PLASMA LACTATE
1649	LaboratoryAndDiagnostic	CHEMISTRY: PRE-ALBUMIN
1650	LaboratoryAndDiagnostic	CHEMISTRY: URINE - AMYLASE
1651	LaboratoryAndDiagnostic	CHEMISTRY: URINE - CALCIUM
1652	LaboratoryAndDiagnostic	CHEMISTRY: URINE - CHLORIDE
1653	LaboratoryAndDiagnostic	CHEMISTRY: URINE - CREATININE
1654	LaboratoryAndDiagnostic	CHEMISTRY: URINE - MG (FOR 24 HRS URINE COLLECTION)
1655	LaboratoryAndDiagnostic	CHEMISTRY: URINE - PHOSPHORUS
1656	LaboratoryAndDiagnostic	CHEMISTRY: URINE - K
1657	LaboratoryAndDiagnostic	CHEMISTRY: URINE - NA
1658	LaboratoryAndDiagnostic	CHEMISTRY: URINE - GLUCOSE
1659	LaboratoryAndDiagnostic	CHEMISTRY: URINE - TOTAL PROTEIN
1660	LaboratoryAndDiagnostic	CHEMISTRY: URINE - UREAN NITROGEN
1661	LaboratoryAndDiagnostic	CHEMISTRY: URINE - URIC ACID
1662	LaboratoryAndDiagnostic	CHEMISTRY: SERUM IRON
1663	LaboratoryAndDiagnostic	CHEMISTRY: URINE - NA (FOR 24 HRS URINE COLLECTION)
1664	LaboratoryAndDiagnostic	CHEMISTRY: URINE - K (FOR 24 HRS URINE COLLECTION)
1665	LaboratoryAndDiagnostic	CHEMISTRY: URINE - CHLORIDE (FOR 24 HRS URINE COLLECTION)
1666	LaboratoryAndDiagnostic	CHEMISTRY: URINE - ALBUMIN (FOR 24 HRS URINE COLLECTION)
1667	LaboratoryAndDiagnostic	CHEMISTRY: URINE - CALCIUM (FOR 24 HRS URINE COLLECTION)
1668	LaboratoryAndDiagnostic	CHEMISTRY: URINE - CREATININE (FOR 24 HRS URINE COLLECTION)
1669	LaboratoryAndDiagnostic	CHEMISTRY: URINE - PROTEIN (FOR 24 HRS URINE COLLECTION)
1670	LaboratoryAndDiagnostic	CHEMISTRY: FBS
1671	LaboratoryAndDiagnostic	CHEMISTRY: RBS
1672	LaboratoryAndDiagnostic	CHEMISTRY: G6PD
1673	LaboratoryAndDiagnostic	CHEMISTRY: FERRITIN
1674	LaboratoryAndDiagnostic	CHEMISTRY: TROPONIN I
1675	LaboratoryAndDiagnostic	CHEMISTRY: COVID TOTAL ANTIBODY TEST
1676	LaboratoryAndDiagnostic	CHEMISTRY: COVID IGG ANTIBODY TEST
1677	LaboratoryAndDiagnostic	CHEMISTRY: PERITONEAL FLUID, DIALYSATE, JP DRAIN (FROM ABDOMINAL CAVITY ONLY) - CREATININE
1678	LaboratoryAndDiagnostic	CHEMISTRY: NEPHROCHECK (URINARY TIMP-2 AND IGFBP-7)
1679	LaboratoryAndDiagnostic	CHEMISTRY: HSCRP
1680	LaboratoryAndDiagnostic	CHEMISTRY: TOTAL IRON BINDING CAPACITY (TBIC)
1681	LaboratoryAndDiagnostic	CHEMISTRY: DIRECT TBIC
1682	LaboratoryAndDiagnostic	CHEMISTRY: NT PRO BNP
1683	LaboratoryAndDiagnostic	CHEMISTRY: PRO-BNP
1684	LaboratoryAndDiagnostic	CHEMISTRY: BAT PANEL: NA, K, CA, CL, ALB, RBS, MG, BUN, CREA, ALT, AST
1685	LaboratoryAndDiagnostic	CHEMISTRY: COVID PANEL: HSCRP, LDH, RBS, BUN, CREATININE, NA, K, CL, CA, MG, ALB, TG, AST, ALT, TB, DB, IB, ALP, FERRITIN
1686	LaboratoryAndDiagnostic	MICROSCOPY: CSF-DIFF COUNT
1687	LaboratoryAndDiagnostic	MICROSCOPY: CSF-PH
1688	LaboratoryAndDiagnostic	MICROSCOPY: CSF-QUALITATIVE
1689	LaboratoryAndDiagnostic	MICROSCOPY: CSF-RBC COUNT
1690	LaboratoryAndDiagnostic	MICROSCOPY: CSF-RBC MORPHOLOGY
1691	LaboratoryAndDiagnostic	MICROSCOPY: CSF-SPECIFIC GRAVITY
1692	LaboratoryAndDiagnostic	MICROSCOPY: CSF-WBC COUNT

1693	LaboratoryAndDiagnostic	MICROSCOPY: CSF-BODY FLUID WET MOUNT EXAMINATION FOR PARASITES
1694	LaboratoryAndDiagnostic	MICROSCOPY: DIRECT FECAL SMEAR WITH FECAL CONCENTRATION TECHNIQUE
1695	LaboratoryAndDiagnostic	MICROSCOPY: DIRECT FECAL SMEAR WITH FECAL CONCENTRATION TECHNIQUE AND OCCULT BLOOD (GUAIAC METHOD)
1696	LaboratoryAndDiagnostic	MICROSCOPY: FERNING TEST
1697	LaboratoryAndDiagnostic	MICROSCOPY: SEMENALYSIS
1698	LaboratoryAndDiagnostic	MICROSCOPY: POLARIZED MICROSCOPY IN SYNOVIAL FLUID
1699	LaboratoryAndDiagnostic	MICROSCOPY: STOOL-FIT
1700	LaboratoryAndDiagnostic	MICROSCOPY: STOOL EXAMINATION FOR ENTAMOEBA HISTOLYTICA ANTIGEN
1701	LaboratoryAndDiagnostic	MICROSCOPY: STOOL EXAMINATION FOR ROTAVIRUS ANTIGEN
1702	LaboratoryAndDiagnostic	MICROSCOPY: STOOL-FECT (STOOL CONCENTRATION)
1703	LaboratoryAndDiagnostic	MICROSCOPY: STOOL-OCCULT BLOOD
1704	LaboratoryAndDiagnostic	MICROSCOPY: URINE-URINALYSIS
1705	LaboratoryAndDiagnostic	MICROSCOPY: URINE-MYOGLOBIN
1706	LaboratoryAndDiagnostic	MICROSCOPY: URINE-PREGNANCY TEST
1707	LaboratoryAndDiagnostic	MICROSCOPY: URINE-RBC MORPHOLOGY
1708	LaboratoryAndDiagnostic	MICROSCOPY: VAGINAL WASHING ANALYSIS FOR SPERM CELLS
1709	LaboratoryAndDiagnostic	MICROSCOPY: CLOSTRIDIUM DIFFICILE ANTIGEN AND TOXIN
1710	LaboratoryAndDiagnostic	MICROSCOPY: PLEURAL FLUID (QUALITATIVE)
1711	LaboratoryAndDiagnostic	MICROSCOPY: SYNOVIAL FLUID (QUALITATIVE)
1712	LaboratoryAndDiagnostic	MICROSCOPY: ASCITIC FLUID (QUALITATIVE)
1713	LaboratoryAndDiagnostic	MICROSCOPY: OTHER BODY FLUIDS (QUALITATIVE)
1714	LaboratoryAndDiagnostic	MICROSCOPY: URINE PH
1715	LaboratoryAndDiagnostic	MICROSCOPY: URINE SPECIFIC GRAVITY
1716	LaboratoryAndDiagnostic	MICROSCOPY: URINE KETONES
1717	LaboratoryAndDiagnostic	MICROSCOPY: BAT PANEL: URINALYSIS
1718	LaboratoryAndDiagnostic	PATHOLOGY: AFP
1719	LaboratoryAndDiagnostic	PATHOLOGY: AMACR
1720	LaboratoryAndDiagnostic	PATHOLOGY: BCL6
1721	LaboratoryAndDiagnostic	PATHOLOGY: CA 19-9
1722	LaboratoryAndDiagnostic	PATHOLOGY: CALDESMON
1723	LaboratoryAndDiagnostic	PATHOLOGY: CD79A
1724	LaboratoryAndDiagnostic	PATHOLOGY: CD21
1725	LaboratoryAndDiagnostic	PATHOLOGY: CD1A
1726	LaboratoryAndDiagnostic	PATHOLOGY: CD4
1727	LaboratoryAndDiagnostic	PATHOLOGY: CD8
1728	LaboratoryAndDiagnostic	PATHOLOGY: CDX2
1729	LaboratoryAndDiagnostic	PATHOLOGY: CELL BLOCKS / CYTOLOGY (CB)
1730	LaboratoryAndDiagnostic	PATHOLOGY: DAK-PAX5
1731	LaboratoryAndDiagnostic	PATHOLOGY: GATA3
1732	LaboratoryAndDiagnostic	PATHOLOGY: HBME1
1733	LaboratoryAndDiagnostic	PATHOLOGY: HER2 FISH
1734	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - ALK
1735	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - BCL2
1736	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CALRETININ
1737	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CD3
1738	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CD5
1739	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CD10

1740	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CD15
1741	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CD20
1742	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CD30
1743	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CD34
1744	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CD45
1745	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - LCA
1746	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CD56
1747	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - NCAM
1748	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CD68
1749	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CD99
1750	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CD117
1751	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CD138
1752	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CEA
1753	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CK (PANCYTOKERATIN)
1754	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - PANCK
1755	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CK7
1756	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CK19
1757	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CK20
1758	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CK56
1759	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CKLMW
1760	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CYCLIN D1
1761	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - DESMIN
1762	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - E-CADHERIN
1763	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - EMA
1764	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - ER & PR ASSAY
1765	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - ER OR PR ONLY
1766	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - ER/PR WITH HER2NEU
1767	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - GFAP
1768	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - HEPPAR1
1769	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - HER2NEU
1770	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - HERCEP 2
1771	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - HMB 45
1772	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - IDH1
1773	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - KI67
1774	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - MAMMAGLOBIN
1775	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - MELANIN A
1776	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - MUSCLE ACTIN (MSA)
1777	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - NSE
1778	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - P40
1779	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - P53
1780	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - P57
1781	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - P63
1782	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - PAX5
1783	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - PLAP
1784	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - S-100
1785	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - SMA
1786	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - SYNAPTOPHYSIN
1787	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - TDT

1788	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - TOXOPLASMA
1789	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - TTF-1
1790	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - VIMENTIN
1791	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - WT 1
1792	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - THYROGLOBULIN
1793	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - MUM1
1794	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - MYELOPEROXIDASE
1795	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - MYOD1
1796	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - NEUROFILAMENT
1797	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - NKX3
1798	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - P16
1799	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - PAX-8
1800	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - PD-LHHC
1801	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CAM 5.2 (CK LMW)
1802	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CD2
1803	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CD7
1804	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CD21
1805	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CD31
1806	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CHROMOGRANIN
1807	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - CK 8/18
1808	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - C-MYC
1809	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - DOG1
1810	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - ERG
1811	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - GATA3
1812	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - GLYPICAN 3
1813	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - GRANZYME B
1814	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - HBME1
1815	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - INHIBIN
1816	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - INI1
1817	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - KAPPA LIGHT CHAIN
1818	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - LAMBDA LIGHT CHAIN
1819	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - MDM2
1820	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - MUC4
1821	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - HMLH1
1822	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - MURL PROTEIN HOMOLOG 1
1823	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - MUTS PROTEIN HOMOLOG 2
1824	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - HMSH2
1825	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - MUTS PROTEIN HOMOLOG 6
1826	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - NAPSIN A
1827	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - PSMA (PROSTATE-SPECIFIC)
1828	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - SATB2
1829	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - STAT6
1830	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - SOX10
1831	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - TFE1
1832	LaboratoryAndDiagnostic	PATHOLOGY: IMMUNOSTAINS - TFE3
1833	LaboratoryAndDiagnostic	PAP SMEAR, CONVENTIONAL:
1834	LaboratoryAndDiagnostic	LIQUID-BASED CYTOLOGY FOR NON-GYNECOLOGICAL SAMPLES:
1835	LaboratoryAndDiagnostic	LIQUID-BASED CYTOLOGY FOR GYNECOLOGICAL SAMPLES: SALL4

1836	LaboratoryAndDiagnostic	LIQUID-BASED CYTOLOGY FOR GYNECOLOGICAL SAMPLES: SLIDE REVIEW (SR)
1837	LaboratoryAndDiagnostic	LIQUID-BASED CYTOLOGY FOR GYNECOLOGICAL SAMPLES: SPECIAL STAINS - GIEMSA
1838	LaboratoryAndDiagnostic	LIQUID-BASED CYTOLOGY FOR GYNECOLOGICAL SAMPLES: SPECIAL STAINS - PERIODIC ACID SCHIFF
1839	LaboratoryAndDiagnostic	LIQUID-BASED CYTOLOGY FOR GYNECOLOGICAL SAMPLES: SPECIAL STAINS - RETICULIN
1840	LaboratoryAndDiagnostic	LIQUID-BASED CYTOLOGY FOR GYNECOLOGICAL SAMPLES: SPECIAL STAINS - WRIGHTS
1841	LaboratoryAndDiagnostic	LIQUID-BASED CYTOLOGY FOR GYNECOLOGICAL SAMPLES: SPECIAL STAINS - CONGO RED
1842	LaboratoryAndDiagnostic	LIQUID-BASED CYTOLOGY FOR GYNECOLOGICAL SAMPLES: SPECIAL STAINS - ELASTIC STAIN
1843	LaboratoryAndDiagnostic	LIQUID-BASED CYTOLOGY FOR GYNECOLOGICAL SAMPLES: SPECIAL STAINS - FITE FARACO
1844	LaboratoryAndDiagnostic	LIQUID-BASED CYTOLOGY FOR GYNECOLOGICAL SAMPLES: SPECIAL STAINS - MASSON TRICHOME
1845	LaboratoryAndDiagnostic	LIQUID-BASED CYTOLOGY FOR GYNECOLOGICAL SAMPLES: SPECIAL STAINS - PAPANICOLAU STAIN
1846	LaboratoryAndDiagnostic	LIQUID-BASED CYTOLOGY FOR GYNECOLOGICAL SAMPLES: SPECIAL STAINS - DIFF QUICK STAIN
1847	LaboratoryAndDiagnostic	LIQUID-BASED CYTOLOGY FOR GYNECOLOGICAL SAMPLES: TISSUE BIOPSY - EXTRA LARGE
1848	LaboratoryAndDiagnostic	LIQUID-BASED CYTOLOGY FOR GYNECOLOGICAL SAMPLES: TISSUE BIOPSY - FROZEN SECTION
1849	LaboratoryAndDiagnostic	LIQUID-BASED CYTOLOGY FOR GYNECOLOGICAL SAMPLES: TISSUE BIOPSY - LARGE (>2 CM)
1850	LaboratoryAndDiagnostic	LIQUID-BASED CYTOLOGY FOR GYNECOLOGICAL SAMPLES: TISSUE BIOPSY - MEDIUM (1 CM <= 2 CM)
1851	LaboratoryAndDiagnostic	LIQUID-BASED CYTOLOGY FOR GYNECOLOGICAL SAMPLES: TISSUE BIOPSY - SMALL (<1 CM)
1852	LaboratoryAndDiagnostic	MICROBIOLOGY: AFB (AURAMINE)
1853	LaboratoryAndDiagnostic	MICROBIOLOGY: AFB (MODIFIED)
1854	LaboratoryAndDiagnostic	MICROBIOLOGY: AFB CS
1855	LaboratoryAndDiagnostic	MICROBIOLOGY: ANAEROBIC CULTURE
1856	LaboratoryAndDiagnostic	MICROBIOLOGY: BLOOD C/S
1857	LaboratoryAndDiagnostic	MICROBIOLOGY: CSF GS/CS
1858	LaboratoryAndDiagnostic	MICROBIOLOGY: EXUDATES GS/CS
1859	LaboratoryAndDiagnostic	MICROBIOLOGY: FUNGAL C/S
1860	LaboratoryAndDiagnostic	MICROBIOLOGY: GRAM STAIN
1861	LaboratoryAndDiagnostic	MICROBIOLOGY: INDIA INK
1862	LaboratoryAndDiagnostic	MICROBIOLOGY: KOH MOUNT
1863	LaboratoryAndDiagnostic	MICROBIOLOGY: MAKI/ENVIRONMENTAL
1864	LaboratoryAndDiagnostic	MICROBIOLOGY: MRSA SCREENING
1865	LaboratoryAndDiagnostic	MICROBIOLOGY: MTB/RIF RT-PCR
1866	LaboratoryAndDiagnostic	MICROBIOLOGY: MTB/RIF RT-PCR + AFB CS (PACKAGE)
1867	LaboratoryAndDiagnostic	MICROBIOLOGY: RT-PCR C. DIFFICILE
1868	LaboratoryAndDiagnostic	MICROBIOLOGY: RT-PCR C. TRACHOMATIS
1869	LaboratoryAndDiagnostic	MICROBIOLOGY: RESPIRATORY GS/CS
1870	LaboratoryAndDiagnostic	MICROBIOLOGY: STOOL GS/CS
1871	LaboratoryAndDiagnostic	MICROBIOLOGY: TMG SMEAR
1872	LaboratoryAndDiagnostic	MICROBIOLOGY: TPN SOLUTION
1873	LaboratoryAndDiagnostic	MICROBIOLOGY: URINE GS/CS
1874	LaboratoryAndDiagnostic	MICROBIOLOGY: URINE GS
1875	LaboratoryAndDiagnostic	MICROBIOLOGY: URINE CS
1876	LaboratoryAndDiagnostic	MICROBIOLOGY: ETA GS/CS
1877	LaboratoryAndDiagnostic	MICROBIOLOGY: ETA AFB
1878	LaboratoryAndDiagnostic	MICROBIOLOGY: SPUTUM GS/CS

1879	LaboratoryAndDiagnostic	MICROBIOLOGY: SPUTUM AFB
1880	LaboratoryAndDiagnostic	MICROBIOLOGY: URINE AFB
1881	LaboratoryAndDiagnostic	MICROBIOLOGY: STOOL AFB
1882	LaboratoryAndDiagnostic	MICROBIOLOGY: WOUND DISCHARGE GS/CS
1883	LaboratoryAndDiagnostic	MICROBIOLOGY: WOUND DISCHARGE AFB
1884	LaboratoryAndDiagnostic	MICROBIOLOGY: PERITONEAL FLUID GS/CS
1885	LaboratoryAndDiagnostic	MICROBIOLOGY: PERITONEAL FLUID AFB
1886	LaboratoryAndDiagnostic	MICROBIOLOGY: PLEURAL FLUID GS/CS
1887	LaboratoryAndDiagnostic	MICROBIOLOGY: PLEURAL FLUID AFB
1888	LaboratoryAndDiagnostic	MICROBIOLOGY: SYNOVIAL FLUID GS/CS
1889	LaboratoryAndDiagnostic	MICROBIOLOGY: SYNOVIAL FLUID AFB
1890	LaboratoryAndDiagnostic	MICROBIOLOGY: MILKBANK CS
1891	LaboratoryAndDiagnostic	MICROBIOLOGY: BODY FLUID CS
1892	LaboratoryAndDiagnostic	MICROBIOLOGY: VP SHUNT TIP GS/CS
1893	LaboratoryAndDiagnostic	MICROBIOLOGY: TISSUE GS/CS
1894	LaboratoryAndDiagnostic	MICROBIOLOGY: TISSUE AFB
1895	LaboratoryAndDiagnostic	MICROBIOLOGY: COVID PANEL: BLOOD CS
1896	LaboratoryAndDiagnostic	BLOOD TYPING: BLOOD TYPING
1897	RoomAndBoard	ROOM AND BOARD
1898	OperatingRoomFees	OPERATING ROOM

B. Category Library

cat_id	cat_desc	billtype_id
RoomAndBoard	ROOM AND BOARD	I
DrugsAndMedicine	DRUGS AND MEDICINE	I
LaboratoryAndDiagnostic	LABORATORY AND DIAGNOSTIC	I
OperatingRoomFees	OPERATING ROOM FEES	I
MedicalSupplies	MEDICAL SUPPLIES	I
Others	OTHERS	I
PAN: (FN MN LN Suffix)	PAN: (FN MN LN SUFFIX)	P
PhilHealth	PHILHEALTH	
Balance	BALANCE	

C. Sub-Category Library

subcat_id	subcat_desc
1	ACTUAL CHARGES
2	SENIOR CITIZEN DISCOUNT
3	PWD DISCOUNT
4	PCSO DISCOUNT
5	DSWD DISCOUNT
6	DOH MEDICAL ASSISTANCE PROGRAM
7	HMO DISCOUNT

D. Bill Type Library

billtype_id	billtype_desc
I	HEALTH CARE INSTITUTION
P	HEALTH CARE PROFESSIONAL

E. Medicine Library (Drug and Medicine Category)

Drug Code	Drug Description
0196100000000000SOL3200195AMPUL	VITAMINS INTRAVENOUS, FAT-SOLUBLE SOLUTION 10 mL AMPULE
0196200000000000SOL3200195AMPUL	VITAMINS INTRAVENOUS, TRACE ELEMENTS SOLUTION 10 mL AMPULE
0196400000000361SOL3200000AMPUL	VITAMINS INTRAVENOUS, WATER-SOLUBLE 2 mL SOLUTION AMPULE
0196400000000516SOL3200000VIALX	VITAMINS INTRAVENOUS, WATER-SOLUBLE 1 mL SOLUTION VIAL
09SOD00000000000SOL3200011AMPUL	0.9% SODIUM CHLORIDE SOLUTION 20 mL AMPULE
09SOD00000000000SOL3200020BOTTL	0.9% SODIUM CHLORIDE SOLUTION 250 mL BOTTLE
09SOD00000000000SOL3200044BOTTL	0.9% SODIUM CHLORIDE SOLUTION 500 mL BOTTLE
09SOD00000000000SOL3200055BOTTL	0.9% SODIUM CHLORIDE SOLUTION 50 mL BOTTLE
09SOD00000000000SOL3200060AMPUL	0.9% SODIUM CHLORIDE SOLUTION 5 mL AMPULE
09SOD00000000000SOL3200178BOTTL	0.9% SODIUM CHLORIDE SOLUTION 100 mL BOTTLE
09SOD00000000000SOL3200195AMPUL	0.9% SODIUM CHLORIDE SOLUTION 10 mL AMPULE
09SOD00000000000SOL3200247BOTTL	0.9% SODIUM CHLORIDE SOLUTION 1 L BOTTLE
09SOD00000000000SOL3200361AMPUL	0.9% SODIUM CHLORIDE SOLUTION 2 mL AMPULE
09SOD00000000000SOL3201124AMPUL	0.9% SODIUM CHLORIDE SOLUTION 2.5 mL AMPULE
09SOD00000000000SOL3201125BOTTL	0.9% SODIUM CHLORIDE SOLUTION 200 mL BOTTLE
10DEX00000000000SOL3200020BOTTL	10% DEXTROSE IN WATER SOLUTION 250 mL BOTTLE
10DEX00000000000SOL3200044BOTTL	10% DEXTROSE IN WATER SOLUTION 500 mL BOTTLE
10DEX00000000000SOL3200247BOTTL	10% DEXTROSE IN WATER SOLUTION 1 L BOTTLE
10DEX00000000000SOL3200388AMPUL	10% DEXTROSE IN WATER SOLUTION 3 mL AMPULE
5D3SC00000000000SOL3200020BOTTL	5% DEXTROSE IN 0.3% SODIUM CHLORIDE SOLUTION 250 mL BOTTLE
5D3SC00000000000SOL3200044BOTTL	5% DEXTROSE IN 0.3% SODIUM CHLORIDE SOLUTION 500 mL BOTTLE
5D3SC00000000000SOL3200247BOTTL	5% DEXTROSE IN 0.3% SODIUM CHLORIDE SOLUTION 1 L BOTTLE
5D45S00000000000SOL3200247BAGXX	5 % DEXTROSE IN SODIUM CHLORIDE SOLUTION 1 L BAG
5D9SC00000000000SOL3200020BOTTL	5% DEXTROSE IN 0.9% SODIUM CHLORIDE SOLUTION 250 mL BOTTLE
5D9SC00000000000SOL3200044BOTTL	5% DEXTROSE IN 0.9% SODIUM CHLORIDE SOLUTION 500 mL BOTTLE
5D9SC00000000000SOL3200247BOTTL	5% DEXTROSE IN 0.9% SODIUM CHLORIDE SOLUTION 1 L BOTTLE
5DEXW00000000000SOL3200020BOTTL	5% DEXTROSE IN WATER SOLUTION 250 mL BOTTLE
5DEXW00000000000SOL3200044BOTTL	5% DEXTROSE IN WATER SOLUTION 500 mL BOTTLE
5DEXW00000000000SOL3200247BOTTL	5% DEXTROSE IN WATER SOLUTION 1 L BOTTLE
5DLRS00000000000SOL3200020BOTTL	5% DEXTROSE IN LACTATED RINGER'S SOLUTION 250 mL BOTTLE
5DLRS00000000000SOL3200044BOTTL	5% DEXTROSE IN LACTATED RINGER'S SOLUTION 500 mL BOTTLE
5DLRS00000000000SOL3200247BOTTL	5% DEXTROSE IN LACTATED RINGER'S SOLUTION 1 L BOTTLE
ABA010000000028TAB490000000000	ABACAIR 300 mg TABLET
AC0010000000000POW2701266USP02	ACTIVATED CHARCOAL POWDER 1 kg USP
ACARB00000000050TAB490000000000	ACARBOSE 50 mg TABLET
ACARB00000000139TAB490000000000	ACARBOSE 100 mg TABLET
ACET10000000000SOL3200044BAGXX	ACETATED RINGER'S SOLUTION SOLUTION 500 mL BAG
ACET10000000000SOL3200044BOTTL	ACETATED RINGER'S SOLUTION SOLUTION 500 mL BOTTLE
ACET10000000000SOL3200247BAGXX	ACETATED RINGER'S SOLUTION SOLUTION 1 L BAG
ACET10000000000SOL3200247BOTTL	ACETATED RINGER'S SOLUTION SOLUTION 1 L BOTTLE
ACET20000000015TAB490000000000	ACETAZOLAMIDE 250 mg TABLET
ACEY100000000040000000000SAC01	ACETYLCYSTEINE 200 mg SACHET
ACEY100000000061TAB200000000000	ACETYLCYSTEINE 600 mg TABLET EFFERVESCENT

ACEY100000001390000000000SAC01	ACETYLCYSTEINE 100 mg SACHET
ACEY10000000162SUS1401374BOTTL	ACETYLCYSTEINE 100 mg/5 mL SUSPENSION 150 mL BOTTLE
ACEY100000001123SOL3200195AMPUL	ACETYLCYSTEINE 200 mg/mL SOLUTION 10 mL AMPULE
ACEY100000001123SOL3201375BOTTL	ACETYLCYSTEINE 200 mg/mL SOLUTION 25 mL BOTTLE
ACICL0000000004TAB490000000000	ACICLOVIR 200 mg TABLET
ACICL0000000006SUS1400201BOTTL	ACICLOVIR 200 mg/5 mL SUSPENSION 120 mL BOTTLE
ACICL0000000006SUS1400469BOTTL	ACICLOVIR 200 mg/5 mL SUSPENSION 60 mL BOTTLE
ACICL00000000034TAB490000000000	ACICLOVIR 400 mg TABLET
ACICL000000000364EYO0101376TUB01	ACICLOVIR 3% EYE OINTMENT 4.5 g TUBE
ACICL00000000495TAB490000000000	ACICLOVIR 800 mg TABLET
ACICL00000000985SOL3200195VIALX	ACICLOVIR 25 mg/mL SOLUTION 10 mL VIAL
ADENO00000000387SOL3200361VIALX	ADENOSINE 3 mg/mL SOLUTION 2 mL VIAL
AHMHY00000000320SUS1400020BOTTL	ALUMINUM HYDROXIDE + MAGNESIUM HYDROXIDE 225 mg + 200 mg/5 mL SUSPENSION 250 mL BOTTLE
AHMHY00000000320SUS1400201BOTTL	ALUMINUM HYDROXIDE + MAGNESIUM HYDROXIDE 225 mg + 200 mg/5 mL SUSPENSION 120 mL BOTTLE
AHMHY00000000320SUS1400469BOTTL	ALUMINUM HYDROXIDE + MAGNESIUM HYDROXIDE 225 mg + 200 mg/5 mL SUSPENSION 60 mL BOTTLE
AHMHY00000000320SUS1401377BOTTL	ALUMINUM HYDROXIDE + MAGNESIUM HYDROXIDE 225 mg + 200 mg/5 mL SUSPENSION 180 mL BOTTLE
AHMHY00000000320SUS1401378BOTTL	ALUMINUM HYDROXIDE + MAGNESIUM HYDROXIDE 225 mg + 200 mg/5 mL SUSPENSION 355 mL BOTTLE
AKAPU00000000512LOT1100469BOTTL	AKAPULKO 50% LOTION 60 mL BOTTLE
ALBED0000000006SUS1400195BOTTL	ALBENDAZOLE 200 mg/5 mL SUSPENSION 10 mL BOTTLE
ALBED0000000006SUS1400231BOTTL	ALBENDAZOLE 200 mg/5 mL SUSPENSION 15 mL BOTTLE
ALBED0000000006SUS1400379BOTTL	ALBENDAZOLE 200 mg/5 mL SUSPENSION 30 mL BOTTLE
ALBED0000000006SUS1400469BOTTL	ALBENDAZOLE 200 mg/5 mL SUSPENSION 60 mL BOTTLE
ALBED00000000034TAB490000000000	ALBENDAZOLE 400 mg TABLET
ALBUM00000000289SOL3200055BOTTL	ALBUMIN, HUMAN 20% SOLUTION 50 mL BOTTLE
ALBUM00000000289SOL3200178BOTTL	ALBUMIN, HUMAN 20% SOLUTION 100 mL BOTTLE
ALBUM00000000588SOL3200055BAGXX	ALBUMIN, HUMAN 25% SOLUTION 50 mL BAG
ALBUM00000000588SOL3200055BOTTL	ALBUMIN, HUMAN 25% SOLUTION 50 mL BOTTLE
ALBUM00000000588SOL3200178BOTTL	ALBUMIN, HUMAN 25% SOLUTION 100 mL BOTTLE
ALCOE00000000000SOL3200516AMPUL	ALCOHOL, ETHYL SOLUTION 1 mL AMPULE
ALCOE00000000479SOL3200020BOTTL	ALCOHOL, ETHYL 70% SOLUTION 250 mL BOTTLE
ALCOE00000000479SOL3200044BOTTL	ALCOHOL, ETHYL 70% SOLUTION 500 mL BOTTLE
ALCOE00000000479SOL3200178BOTTL	ALCOHOL, ETHYL 70% SOLUTION 100 mL BOTTLE
ALCOE00000000479SOL3200469BOTTL	ALCOHOL, ETHYL 70% SOLUTION 60 mL BOTTLE
ALCOE00000000479SOL3201374BOTTL	ALCOHOL, ETHYL 70% SOLUTION 150 mL BOTTLE
ALE000005300518TAB490000000000	ALENDRONATE + CHOLECALCIFEROL (VIT. D3) (as SODIUM SALT) 70 mg + 2800 IU TABLET
ALEN20005300482TAB490000000000	ALENDRONATE SODIUM (as SODIUM SALT) 70 mg TABLET
ALEND0005300069TAB490000000000	ALENDRONATE (as SODIUM SALT) 10 mg TABLET
ALFUZ0000200003TAB490000000000	ALFUZOSIN (as HYDROCHLORIDE) 2.5 mg TABLET
ALFUZ0000200056TAB490000000000	ALFUZOSIN (as HYDROCHLORIDE) 5 mg TABLET
ALFUZ0000200069TAB490000000000	ALFUZOSIN (as HYDROCHLORIDE) 10 mg TABLET
ALL-I0000000072SOL3200000BOTTL	ALL-IN-ONE ADMIXTURES ("3-IN-1" OR "DUAL ENERGY" SOLUTIONS) 3 in 1 SOLUTION BOTTLE
ALL-I0000001350SOL3200000BOTTL	ALL-IN-ONE ADMIXTURES ("3-IN-1" OR "DUAL ENERGY" SOLUTIONS) 3 in 1 1400 Kcal SOLUTION BOTTLE
ALL-I0000001351SOL3200000BOTTL	ALL-IN-ONE ADMIXTURES ("3-IN-1" OR "DUAL ENERGY" SOLUTIONS) 3 in 1 1000 Kcal SOLUTION BOTTLE

ALL-I0000001354SOL3200000BOTTL	ALL-IN-ONE ADMIXTURES ("3-IN-1" OR "DUAL ENERGY" SOLUTIONS) 3 in 1 1900 Kcal SOLUTION BOTTLE
ALL-I0000001355SOL3200000BOTTL	ALL-IN-ONE ADMIXTURES ("3-IN-1" OR "DUAL ENERGY" SOLUTIONS) 3 in 1 1300 Kcal SOLUTION BOTTLE
ALL-I0000001356SOL3200000BOTTL	ALL-IN-ONE ADMIXTURES ("3-IN-1" OR "DUAL ENERGY" SOLUTIONS) 3 in 1 1500 Kcal SOLUTION BOTTLE
ALLOP0000000028TAB490000000000	ALLOPURINOL 300 mg TABLET
ALLOP0000000139TAB490000000000	ALLOPURINOL 100 mg TABLET
ALP0000000100194SUS1400195BOTTL	ALPHA-TOCOPHEROL (VITAMIN E) (FOR PREMATURES) (as ACETATE) 10 mg/mL SUSPENSION 10 mL BOTTLE
ALPRA0000000013TAB490000000000	ALPRAZOLAM 250 mcg TABLET
ALPRA0000000045TAB490000000000	ALPRAZOLAM 500 mcg TABLET
ALPRA00000000252TAB490000000000	ALPRAZOLAM 1 mg TABLET
ALUM200000000521GELXX00322BOTTL	ALUMINUM HYDROXIDE 600 mg/5 mL GEL 240 mL BOTTLE
AMACY0000000133SOL3200044GB001	AMINO ACIDS, CRYSTALLINE STANDARD 10% SOLUTION 500 mL GLASS BOTTLE
AMACY0000000133SOL3200044PLASB	AMINO ACIDS, CRYSTALLINE STANDARD 10% SOLUTION 500 mL PLASTIC BOTTLE
AMACY0000000133SOL3200178GB001	AMINO ACIDS, CRYSTALLINE STANDARD 10% SOLUTION 100 mL GLASS BOTTLE
AMACY0000000133SOL3200178PLASB	AMINO ACIDS, CRYSTALLINE STANDARD 10% SOLUTION 100 mL PLASTIC BOTTLE
AMACY0000000133SOL3200247GB001	AMINO ACIDS, CRYSTALLINE STANDARD 10% SOLUTION 1 L GLASS BOTTLE
AMACY0000000133SOL3200247PLASB	AMINO ACIDS, CRYSTALLINE STANDARD 10% SOLUTION 1 L PLASTIC BOTTLE
AMACY0000000366SOL3200044GB001	AMINO ACIDS, CRYSTALLINE STANDARD 3.50% SOLUTION 500 mL GLASS BOTTLE
AMACY0000000366SOL3200044PLASB	AMINO ACIDS, CRYSTALLINE STANDARD 3.50% SOLUTION 500 mL PLASTIC BOTTLE
AMACY0000000462SOL3200020GB001	AMINO ACIDS, CRYSTALLINE STANDARD 6% SOLUTION 250 mL GLASS BOTTLE
AMACY0000000462SOL3200020PLASB	AMINO ACIDS, CRYSTALLINE STANDARD 6% SOLUTION 250 mL PLASTIC BOTTLE
AMACY0000000462SOL3200178GB001	AMINO ACIDS, CRYSTALLINE STANDARD 6% SOLUTION 100 mL GLASS BOTTLE
AMACY0000000462SOL3200178PLASB	AMINO ACIDS, CRYSTALLINE STANDARD 6% SOLUTION 100 mL PLASTIC BOTTLE
AMACY0000000475SOL3200044GB001	AMINO ACIDS, CRYSTALLINE STANDARD 7% SOLUTION 500 mL GLASS BOTTLE
AMACY0000000475SOL3200044PLASB	AMINO ACIDS, CRYSTALLINE STANDARD 7% SOLUTION 500 mL PLASTIC BOTTLE
AMACY0000000475SOL3200178GB001	AMINO ACIDS, CRYSTALLINE STANDARD 7% SOLUTION 100 mL GLASS BOTTLE
AMACY0000000475SOL3200178PLASB	AMINO ACIDS, CRYSTALLINE STANDARD 7% SOLUTION 100 mL PLASTIC BOTTLE
AMACY0000000475SOL3200247GB001	AMINO ACIDS, CRYSTALLINE STANDARD 7% SOLUTION 1 L GLASS BOTTLE
AMACY0000000475SOL3200247PLASB	AMINO ACIDS, CRYSTALLINE STANDARD 7% SOLUTION 1 L PLASTIC BOTTLE
AMACY0000000493SOL3200044GB001	AMINO ACIDS, CRYSTALLINE STANDARD 8.50% SOLUTION 500 mL GLASS BOTTLE
AMACY0000000493SOL3200044PLASB	AMINO ACIDS, CRYSTALLINE STANDARD 8.50% SOLUTION 500 mL PLASTIC BOTTLE
AMACY0000000493SOL3200178GB001	AMINO ACIDS, CRYSTALLINE STANDARD 8.50% SOLUTION 100 mL GLASS BOTTLE
AMACY0000000493SOL3200178PLASB	AMINO ACIDS, CRYSTALLINE STANDARD 8.50% SOLUTION 100 mL PLASTIC BOTTLE
AMACY0000000493SOL3200247GB001	AMINO ACIDS, CRYSTALLINE STANDARD 8.50% SOLUTION 1 L GLASS BOTTLE
AMACY0000000493SOL3200247PLASB	AMINO ACIDS, CRYSTALLINE STANDARD 8.50% SOLUTION 1 L PLASTIC BOTTLE
AMACY0000000500SOL3200011AMPUL	AMINO ACIDS, CRYSTALLINE STANDARD 9.12% SOLUTION 20 mL AMPULE

AMACY0000000500SOL3201125GB001	AMINO ACIDS, CRYSTALLINE STANDARD 9.12% SOLUTION 200 mL GLASS BOTTLE
AMACY0000000500SOL3201125PLASB	AMINO ACIDS, CRYSTALLINE STANDARD 9.12% SOLUTION 200 mL PLASTIC BOTTLE
AMACY0000000541SOL3200044GB001	AMINO ACIDS, CRYSTALLINE STANDARD 11.40% SOLUTION 500 mL GLASS BOTTLE
AMACY0000000541SOL3200044PLASB	AMINO ACIDS, CRYSTALLINE STANDARD 11.40% SOLUTION 500 mL PLASTIC BOTTLE
AMACY0000000541SOL3200178GB001	AMINO ACIDS, CRYSTALLINE STANDARD 11.40% SOLUTION 100 mL GLASS BOTTLE
AMACY0000000541SOL3200178PLASB	AMINO ACIDS, CRYSTALLINE STANDARD 11.40% SOLUTION 100 mL PLASTIC BOTTLE
AMACY0000000541SOL3200247GB001	AMINO ACIDS, CRYSTALLINE STANDARD 11.40% SOLUTION 1 L GLASS BOTTLE
AMACY0000000541SOL3200247PLASB	AMINO ACIDS, CRYSTALLINE STANDARD 11.40% SOLUTION 1 L PLASTIC BOTTLE
AMACY0000000611SOL3200044GB001	AMINO ACIDS, CRYSTALLINE STANDARD 8% SOLUTION 500 mL GLASS BOTTLE
AMACY0000000611SOL3200044PLASB	AMINO ACIDS, CRYSTALLINE STANDARD 8% SOLUTION 500 mL PLASTIC BOTTLE
AMACY0000001315SOL3200020GB001	AMINO ACIDS, CRYSTALLINE STANDARD 5% SOLUTION 250 mL GLASS BOTTLE
AMACY0000001315SOL3200020PLASB	AMINO ACIDS, CRYSTALLINE STANDARD 5% SOLUTION 250 mL PLASTIC BOTTLE
AMACY0000001315SOL3200044GB001	AMINO ACIDS, CRYSTALLINE STANDARD 5% SOLUTION 500 mL GLASS BOTTLE
AMACY0000001315SOL3200044PLASB	AMINO ACIDS, CRYSTALLINE STANDARD 5% SOLUTION 500 mL PLASTIC BOTTLE
AMACY0000001315SOL3200178GB001	AMINO ACIDS, CRYSTALLINE STANDARD 5% SOLUTION 100 mL GLASS BOTTLE
AMACY0000001315SOL3200178PLASB	AMINO ACIDS, CRYSTALLINE STANDARD 5% SOLUTION 100 mL PLASTIC BOTTLE
AMI000003600526SOL3200011AMPUL	AMIDOTRIZOATE (DIATRIZOATE) (as MEGLUMINE AND/OR SODIUM SALT) 150 mg - 650 mg iodine/mL SOLUTION 2
AMI000003600526SOL3200011BOTTL	AMIDOTRIZOATE (DIATRIZOATE) (as MEGLUMINE AND/OR SODIUM SALT) 150 mg - 650 mg iodine/mL SOLUTION 2
AMI000003600526SOL3200055AMPUL	AMIDOTRIZOATE (DIATRIZOATE) (as MEGLUMINE AND/OR SODIUM SALT) 150 mg - 650 mg iodine/mL SOLUTION 5
AMI000003600526SOL3200055BOTTL	AMIDOTRIZOATE (DIATRIZOATE) (as MEGLUMINE AND/OR SODIUM SALT) 150 mg - 650 mg iodine/mL SOLUTION 5
AMI000003600526SOL3200178AMPUL	AMIDOTRIZOATE (DIATRIZOATE) (as MEGLUMINE AND/OR SODIUM SALT) 150 mg - 650 mg iodine/mL SOLUTION 1
AMI000003600526SOL3200178BOTTL	AMIDOTRIZOATE (DIATRIZOATE) (as MEGLUMINE AND/OR SODIUM SALT) 150 mg - 650 mg iodine/mL SOLUTION 1
AMI000003600526SOL3200379AMPUL	AMIDOTRIZOATE (DIATRIZOATE) (as MEGLUMINE AND/OR SODIUM SALT) 150 mg - 650 mg iodine/mL SOLUTION 3
AMI000003600526SOL3200379BOTTL	AMIDOTRIZOATE (DIATRIZOATE) (as MEGLUMINE AND/OR SODIUM SALT) 150 mg - 650 mg iodine/mL SOLUTION 3
AMI000003600526SOL3201125AMPUL	AMIDOTRIZOATE (DIATRIZOATE) (as MEGLUMINE AND/OR SODIUM SALT) 150 mg - 650 mg iodine/mL SOLUTION 2
AMI000003600526SOL3201125BOTTL	AMIDOTRIZOATE (DIATRIZOATE) (as MEGLUMINE AND/OR SODIUM SALT) 150 mg - 650 mg iodine/mL SOLUTION 2
AMI000003600526SOL3201375AMPUL	AMIDOTRIZOATE (DIATRIZOATE) (as MEGLUMINE AND/OR SODIUM SALT) 150 mg - 650 mg iodine/mL SOLUTION 2
AMI000003600526SOL3201375BOTTL	AMIDOTRIZOATE (DIATRIZOATE) (as MEGLUMINE AND/OR SODIUM SALT) 150 mg - 650 mg iodine/mL SOLUTION 2
AMIKA0000000000POW2700047VIALX	AMIKACIN POWDER 500 mg VIAL
AMIKA0006000000POW2700139VIALX	AMIKACIN (as SULFATE) POWDER 100 mg VIAL
AMIKA0006000000SOL3200241VIALX	AMIKACIN (as SULFATE) SOLUTION 1 g VIAL
AMIKA0006000015POW2700000VIALX	AMIKACIN (as SULFATE) 250 mg POWDER VIAL
AMIKA0006000142SOL3200361AMPUL	AMIKACIN (as SULFATE) 100 mg/ml SOLUTION 2 mL AMPULE
AMIKA0006000142SOL3200361VIALX	AMIKACIN (as SULFATE) 100 mg/ml SOLUTION 2 mL VIAL
AMIKA0006000214SOL3200361AMPUL	AMIKACIN (as SULFATE) 125 mg/mL SOLUTION 2 mL AMPULE
AMIKA0006000214SOL3200361VIALX	AMIKACIN (as SULFATE) 125 mg/mL SOLUTION 2 mL VIAL

AMIKA0006000851SOL3200361AMPUL	AMIKACIN (as SULFATE) 50 mg/mL SOLUTION 2 mL AMPULE
AMIKA0006000851SOL3200361VIALX	AMIKACIN (as SULFATE) 50 mg/mL SOLUTION 2 mL VIAL
AMIKA0006000972SOL3200361AMPUL	AMIKACIN (as SULFATE) 250 mg/mL SOLUTION 2 mL AMPULE
AMIKA0006000972SOL3200361VIALX	AMIKACIN (as SULFATE) 250 mg/mL SOLUTION 2 mL VIAL
AMIN10000000133SOL3200044BOTTLE	AMINO ACID SOLUTIONS FOR HEPATIC FAILURE 10% SOLUTION 500 mL BOTTLE
AMIN10000000462SOL3200044BOTTLE	AMINO ACID SOLUTIONS FOR HEPATIC FAILURE 6% SOLUTION 500 mL BOTTLE
AMIN10000000475SOL3200044BOTTLE	AMINO ACID SOLUTIONS FOR HEPATIC FAILURE 7% SOLUTION 500 mL BOTTLE
AMIN10000000611SOL3200044BOTTLE	AMINO ACID SOLUTIONS FOR HEPATIC FAILURE 8% SOLUTION 500 mL BOTTLE
AMIN100000001315SOL3200044BOTTLE	AMINO ACID SOLUTIONS FOR HEPATIC FAILURE 5% SOLUTION 500 mL BOTTLE
AMIN100000001320SOL3200044BOTTLE	AMINO ACID SOLUTIONS FOR HEPATIC FAILURE 9% SOLUTION 500 mL BOTTLE
AMIN200000000000000000055BOTTLE	AMINO ACID SOLUTIONS FOR IMMUNONUTRITION/IMMUNOENHANCEMENT 50 mL BOTTLE
AMIN200000000000000000178BOTTLE	AMINO ACID SOLUTIONS FOR IMMUNONUTRITION/IMMUNOENHANCEMENT 100 mL BOTTLE
AMIN30000000133SOL3200178BOTTLE	AMINO ACID SOLUTIONS FOR INFANTS 10% SOLUTION 100 mL BOTTLE
AMIN30000000462SOL3200178BOTTLE	AMINO ACID SOLUTIONS FOR INFANTS 6% SOLUTION 100 mL BOTTLE
AMIN30000000475SOL3200178BOTTLE	AMINO ACID SOLUTIONS FOR INFANTS 7% SOLUTION 100 mL BOTTLE
AMIN30000000611SOL3200178BOTTLE	AMINO ACID SOLUTIONS FOR INFANTS 8% SOLUTION 100 mL BOTTLE
AMIN300000001320SOL3200178BOTTLE	AMINO ACID SOLUTIONS FOR INFANTS 9% SOLUTION 100 mL BOTTLE
AMIN400000003660000000044BOTTLE	AMINO ACID SOLUTIONS FOR RENAL CONDITIONS 3.50% 500 mL BOTTLE
AMIN400000004620000000044BOTTLE	AMINO ACID SOLUTIONS FOR RENAL CONDITIONS 6% 500 mL BOTTLE
AMIN400000004750000000044BOTTLE	AMINO ACID SOLUTIONS FOR RENAL CONDITIONS 7% 500 mL BOTTLE
AMIN400000001315000000044BOTTLE	AMINO ACID SOLUTIONS FOR RENAL CONDITIONS 5% 500 mL BOTTLE
AMINO00000000985SOL3200195AMPUL	AMINOPHYLLINE 25 mg/mL SOLUTION 10 mL AMPULE
AMIOD0000000004TAB490000000000	AMIODARONE 200 mg TABLET
AMIOD0000200851SOL3200388AMPUL	AMIODARONE (as HYDROCHLORIDE) 50 mg/mL SOLUTION 3 mL AMPULE
AMLO20000001437TAB490000000000	AMLODIPINE AS BESILATE + ATORVASTATIN CALCIUM 5 mg + 10 mg TABLET
AMLO20000001438TAB490000000000	AMLODIPINE AS BESILATE + ATORVASTATIN CALCIUM 5 mg + 20 mg TABLET
AMLO20000001439TAB490000000000	AMLODIPINE AS BESILATE + ATORVASTATIN CALCIUM 5 mg + 40 mg TABLET
AMLO20000001440TAB490000000000	AMLODIPINE AS BESILATE + ATORVASTATIN CALCIUM 5 mg + 80 mg TABLET
AMLO20000001441TAB490000000000	AMLODIPINE AS BESILATE + ATORVASTATIN CALCIUM 10 mg + 10 mg TABLET
AMLO20000001442TAB490000000000	AMLODIPINE AS BESILATE + ATORVASTATIN CALCIUM 10 mg + 20 mg TABLET
AMLO20000001443TAB490000000000	AMLODIPINE AS BESILATE + ATORVASTATIN CALCIUM 10 mg + 40 mg TABLET
AMLO20000001444TAB490000000000	AMLODIPINE AS BESILATE + ATORVASTATIN CALCIUM 10 mg + 80 mg TABLET
AMLO30000000003TAB490000000000	AMLODIPINE 2.5 mg TABLET
AMLO30000000056TAB490000000000	AMLODIPINE 5 mg TABLET
AMLO30000000069TAB490000000000	AMLODIPINE 10 mg TABLET
AMOX50005700015CAPSU0000000000	AMOXICILLIN (as TRIHYDRATE) 250 mg CAPSULE
AMOX50005700047CAPSU0000000000	AMOXICILLIN (as TRIHYDRATE) 500 mg CAPSULE

AMOX50005700142SUS1400195DRO01	AMOXICILLIN (as TRIHYDRATE) 100 mg/ml SUSPENSION 10 mL DROPS
AMOX50005700142SUS1400231DRO01	AMOXICILLIN (as TRIHYDRATE) 100 mg/ml SUSPENSION 15 mL DROPS
AMOX50005700209SUS1400379BOTTL	AMOXICILLIN (as TRIHYDRATE) 125 mg/5 mL SUSPENSION 30 mL BOTTLE
AMOX50005700209SUS1400469BOTTL	AMOXICILLIN (as TRIHYDRATE) 125 mg/5 mL SUSPENSION 60 mL BOTTLE
AMOX50005701265SUS1400379BOTTL	AMOXICILLIN (as TRIHYDRATE) 250 mg/5 mL SUSPENSION 30 mL BOTTLE
AMOX50005701265SUS1400469BOTTL	AMOXICILLIN (as TRIHYDRATE) 250 mg/5 mL SUSPENSION 60 mL BOTTLE
AMPH10000000050SOL3200000VIALX	AMPHOTERICIN B (LIPID COMPLEX) 50 mg SOLUTION VIAL
AMPH10000000139SOL3200000VIALX	AMPHOTERICIN B (LIPID COMPLEX) 100 mg SOLUTION VIAL
AMPH20000000050SOL3200000VIALX	AMPHOTERICIN B (NON-LIPID COMPLEX) 50 mg SOLUTION VIAL
AMPI100000000241SOL3200000VIALX	AMPICILLIN 1 g SOLUTION VIAL
AMPI10005300015SOL3200241VIALX	AMPICILLIN (as SODIUM SALT) 250 mg SOLUTION 1 g VIAL
AMPI10005300047SOL3200241VIALX	AMPICILLIN (as SODIUM SALT) 500 mg SOLUTION 1 g VIAL
AMPI10005300205SOL3200241VIALX	AMPICILLIN (as SODIUM SALT) 125 mg SOLUTION 1 g VIAL
AMPI200000000120SOL3200000VIALX	AMPICILLIN + SULBACTAM 1.5 g SOLUTION VIAL
AMPI200000001207SOL3200000VIALX	AMPICILLIN + SULBACTAM 2 g + 1 g SOLUTION VIAL
AMPI20005300334SOL3200000VIALX	AMPICILLIN + SULBACTAM (as SODIUM SALT) 250 mg + 125 mg SOLUTION VIAL
AMPI20005300431SOL3200000VIALX	AMPICILLIN + SULBACTAM (as SODIUM SALT) 500mg + 250 mg SOLUTION VIAL
ANITG00000000548SOL3200000VIALX	ANTITHYMOCYTE IMMUNOGLOBULIN ATG (RABBIT) 25 mg/5 mL SOLUTION VIAL
ANTID00000001018SOL3200122AMPUL	ANTI-D IMMUNOGLOBULIN (HUMAN ANTI-D IMMUNOGLOBULIN) 200 mcg/mL SOLUTION 1.5 mL AMPULE
ANTIL00000000162SOL3200000VIALX	ANTILYMPHOCYTE IMMUNOGLOBULIN (ALG) (EQUINE) 100 mg/5 mL SOLUTION VIAL
ANTIR00000000293SOL3200060VIALX	ANTI-RABIES SERUM (EQUINE) 200 IU/mL SOLUTION 5 mL VIAL
ANTIR00000000394SOL3200060VIALX	ANTI-RABIES SERUM (EQUINE) 400 IU/ml SOLUTION 5 mL VIAL
ANTIT00000000216SOL3200122AMPUL	ANTI-TETANUS SERUM (ATS) (EQUINE) 1500 IU/mL SOLUTION 1.5 mL AMPULE
ANTIT00000000216SOL3200122VIALX	ANTI-TETANUS SERUM (ATS) (EQUINE) 1500 IU/mL SOLUTION 1.5 mL VIAL
ANTIT00000000216SOL3200516AMPUL	ANTI-TETANUS SERUM (ATS) (EQUINE) 1500 IU/mL SOLUTION 1 mL AMPULE
ANTIT00000000216SOL3200516VIALX	ANTI-TETANUS SERUM (ATS) (EQUINE) 1500 IU/mL SOLUTION 1 mL VIAL
ANTIT00000000549SOL3201124VIALX	ANTI-TETANUS SERUM (ATS) (EQUINE) 4000 IU/mL SOLUTION 2.5 mL VIAL
ARTEM00000000550TAB490000000000	ARTEMETHER + LUMEFANTRIN 20 mg + 120 mg TABLET
ASCO10000000015TAB490000000000	ASCORBIC ACID (VITAMIN C) 250 mg TABLET
ASCO10000000047TAB490000000000	ASCORBIC ACID (VITAMIN C) 500 mg TABLET
ASCO100000000139TAB490000000000	ASCORBIC ACID (VITAMIN C) 100 mg TABLET
ASCO100000000142DROPS00231BOTTL	ASCORBIC ACID (VITAMIN C) 100 mg/ml DROPS 15 mL BOTTLE
ASCO100000000142DROPS00379BOTTL	ASCORBIC ACID (VITAMIN C) 100 mg/ml DROPS 30 mL BOTTLE
ASCO100000000142DROPS00469BOTTL	ASCORBIC ACID (VITAMIN C) 100 mg/ml DROPS 60 mL BOTTLE
ASCO100000000162SYRUP00020BOTTL	ASCORBIC ACID (VITAMIN C) 100 mg/5 mL SYRUP 250 mL BOTTLE
ASCO100000000162SYRUP00201BOTTL	ASCORBIC ACID (VITAMIN C) 100 mg/5 mL SYRUP 120 mL BOTTLE
ASCO100000000162SYRUP00469BOTTL	ASCORBIC ACID (VITAMIN C) 100 mg/5 mL SYRUP 60 mL BOTTLE
ASCO100000000972SOL3200361AMPUL	ASCORBIC ACID (VITAMIN C) 250 mg/mL SOLUTION 2 mL AMPULE
ASPO100000000135SOL3200000VIALX	ASPARAGINASE 10,000 IU SOLUTION VIAL
ASPIR00000000028TAB490000000000	ASPIRIN 300 mg TABLET
ASPIR00000000066TAB490000000000	ASPIRIN 80 mg TABLET

ASPIR0000000139TAB490000000000	ASPIRIN 100 mg TABLET
ASPIR00000000551TAB490000000000	ASPIRIN 325 mg TABLET
ATEN10000000050TAB490000000000	ATENOLOL 50 mg TABLET
ATEN100000000139TAB490000000000	ATENOLOL 100 mg TABLET
ATOR10000000007TAB490000000000	ATORVASTATIN 20 mg TABLET
ATOR10000000036TAB490000000000	ATORVASTATIN 40 mg TABLET
ATOR10000000066TAB490000000000	ATORVASTATIN 80 mg TABLET
ATOR10000000069TAB490000000000	ATORVASTATIN 10 mg TABLET
ATRA10000900194SOL3200060AMPUL	ATRACURIUM (as BESILATE) 10 mg/mL SOLUTION 5 mL AMPULE
ATRA10000900194SOL3201124AMPUL	ATRACURIUM (as BESILATE) 10 mg/mL SOLUTION 2.5 mL AMPULE
ATROE0006001062EYD0100060BOTTL	ATROPINE (as SULFATE) 1% EYE DROPS 5 mL BOTTLE
ATROE0006001062EYD0100195BOTTL	ATROPINE (as SULFATE) 1% EYE DROPS 10 mL BOTTLE
ATROP0006000261SOL3200516AMPUL	ATROPINE (as SULFATE) 1 mg/mL SOLUTION 1 mL AMPULE
ATROP0006000552TAB490000000000	ATROPINE (as SULFATE) 600 mcg TABLET
ATROP0006000553SOL3200516AMPUL	ATROPINE (as SULFATE) 600 mcg/mL SOLUTION 1 mL AMPULE
ATROP0006000668SOL3200516AMPUL	ATROPINE (as SULFATE) 500 mcg/mL SOLUTION 1 mL AMPULE
AZATH0000000050SOL3200000VIALX	AZATHIOPRINE 50 mg SOLUTION VIAL
AZATH0000000050TAB490000000000	AZATHIOPRINE 50 mg TABLET
AZITH0001700006POW2400231BOTTL	AZITHROMYCIN (as DIHYDRATE) 200 mg/5 mL POWDER FOR SUSPENSION 15 mL BOTTLE
AZITH0001700006POW2400379BOTTL	AZITHROMYCIN (as DIHYDRATE) 200 mg/5 mL POWDER FOR SUSPENSION 30 mL BOTTLE
AZITH0001700006SUS1401379BOTTL	AZITHROMYCIN (as DIHYDRATE) 200 mg/5 mL SUSPENSION 22.5 mL BOTTLE
AZITH0001700015CAPSU0000000000	AZITHROMYCIN (as DIHYDRATE) 250 mg CAPSULE
AZITH0001700015TAB490000000000	AZITHROMYCIN (as DIHYDRATE) 250 mg TABLET
AZITH0001700025GRA0100000BOTTL	AZITHROMYCIN (as DIHYDRATE) 2 g GRANULES BOTTLE
AZITH0001700047POW2700000VIALX	AZITHROMYCIN (as DIHYDRATE) 500 mg POWDER VIAL
AZITH0001700047TAB490000000000	AZITHROMYCIN (as DIHYDRATE) 500 mg TABLET
BACI20006000554OINTX01380TUB01	BACITRACIN + NEOMYCIN + POLYMYXIN B (as SULFATE) 200 units + 3 mg + 4000 units/g OINTMENT 10 g TUBE
BACI20006000556OINTX01245TUB01	BACITRACIN + NEOMYCIN + POLYMYXIN B (as SULFATE) 400 units + 5 mg + 5000 units/g OINTMENT 5 g TUBE
BACI20006000556OINTX01381JAR01	BACITRACIN + NEOMYCIN + POLYMYXIN B (as SULFATE) 400 units + 5 mg + 5000 units/g OINTMENT 500 g JAR
BACLO0000000069TAB490000000000	BACLOFEN 10 mg TABLET
BALAN0000000000SOL3200044BOTTL	BALANCED MULTIPLE REPLACEMENT SOLUTION WITH PH 7.4 SOLUTION 500 mL BOTTLE
BALAN0000000000SOL3200247BOTTL	BALANCED MULTIPLE REPLACEMENT SOLUTION WITH PH 7.4 SOLUTION 1 L BOTTLE
BAMR500000001315SOL3200044BOTTL	BALANCED MULTIPLE REPLACEMENT SOLUTION WITH 5% DEXTROSE 5% SOLUTION 500 mL BOTTLE
BAMR500000001315SOL3200247BOTTL	BALANCED MULTIPLE REPLACEMENT SOLUTION WITH 5% DEXTROSE 5% SOLUTION 1 L BOTTLE
BARIU00000000557POW2700000PLB01	BARIUM SULFATE 340 g POWDER POUCHED LDPE BAG
BARIU00000000558POW2700000PLB01	BARIUM SULFATE 454 g POWDER POUCHED LDPE BAG
BASIL0000000007SOL3200000VIALX	BASILIXIMAB 20 mg SOLUTION VIAL
BCGVA00000000561SOL3200122VIALX	BCG VACCINE 100 mcg/0.1 mL SOLUTION 1.5 mL VIAL
BCGVA00000000561SOL3200361VIALX	BCG VACCINE 100 mcg/0.1 mL SOLUTION 2 mL VIAL
BCGVA00000000561SOL3200516VIALX	BCG VACCINE 100 mcg/0.1 mL SOLUTION 1 mL VIAL
BCGVA00000000562SOL3200000VIALX	BCG VACCINE 500 mcg/mL + 1 mL diluent SOLUTION VIAL
BECL00001900563DPIO100000DIS01	BECLOMETASONE (as DIPROPIONATE) 100 mcg/dose DRY POWDER DISPENSER

BECL00001900564DPI0100000DIS01	BECLOMETASONE (as DIPROPIONATE) 200 mcg/dose DRY POWDER DISPENSER
BECL00001900565DPI0100000INH01	BECLOMETASONE (as DIPROPIONATE) 400 mcg/dose DRY POWDER INHALER
BECL0000190056600000000000MDI01	BECLOMETASONE (as DIPROPIONATE) 50 mcg/dose x 100 doses METERED DOSE INHALER
BECL0000190056700000000000BAM01	BECLOMETASONE (as DIPROPIONATE) 50 mcg/dose x 200 doses BREATH ACTUATED MDI
BECL0000190056700000000000MDI01	BECLOMETASONE (as DIPROPIONATE) 50 mcg/dose x 200 doses METERED DOSE INHALER
BECL0000190056800000000000MDI01	BECLOMETASONE (as DIPROPIONATE) 250 mcg/dose x 200 doses METERED DOSE INHALER
BECL0000190056900000000000BAM01	BECLOMETASONE (as DIPROPIONATE) 100 mcg/dose x 200 doses BREATH ACTUATED MDI
BEN100000000133GELXX01380TUB01	BENZOYL PEROXIDE 10% GEL 10 g TUBE
BEN100000000133GELXX01382BOTTL	BENZOYL PEROXIDE 10% GEL 40 g BOTTLE
BEN100000000133GELXX01383TUB01	BENZOYL PEROXIDE 10% GEL 60 g TUBE
BEN100000000279GELXX01380TUB01	BENZOYL PEROXIDE 2.50% GEL 10 g TUBE
BEN100000000279GELXX01382BOTTL	BENZOYL PEROXIDE 2.50% GEL 40 g BOTTLE
BEN100000000658CREA301382TUB01	BENZOYL PEROXIDE 4% CREAM 40 g TUBE
BEN100000001315GELXX00227BOTTL	BENZOYL PEROXIDE 5% GEL 15 g BOTTLE
BEN100000001315GELXX01380BOTTL	BENZOYL PEROXIDE 5% GEL 10 g BOTTLE
BEN100000001315GELXX01380TUB01	BENZOYL PEROXIDE 5% GEL 10 g TUBE
BEN100000001315GELXX01382TUB01	BENZOYL PEROXIDE 5% GEL 40 g TUBE
BEN100000001315LOTI100201BOTTL	BENZOYL PEROXIDE 5% LOTION 120 mL BOTTLE
BEN100000001315LOTI100379BOTTL	BENZOYL PEROXIDE 5% LOTION 30 mL BOTTLE
BEN100000001315SOAPX01385SOA01	BENZOYL PEROXIDE 5% SOAP 75 g SOAP
BEN100000001315WAS0101064BOTTL	BENZOYL PEROXIDE 5% WASH 100 g BOTTLE
BEN100000001315WAS0101064TUB01	BENZOYL PEROXIDE 5% WASH 100 g TUBE
BEN100000001315WAS0101384BOTTL	BENZOYL PEROXIDE 5% WASH 50 g BOTTLE
BEN100000001315WAS0101384TUB01	BENZOYL PEROXIDE 5% WASH 50 g TUBE
BENZ70000000570CREA300227TUB01	BENZOIC ACID + SALICYLIC ACID 6% + 3% CREAM 15 g TUBE
BENZ70000000570CREA300367TUB01	BENZOIC ACID + SALICYLIC ACID 6% + 3% CREAM 30 g TUBE
BENZ70000000570OINTX00227TUB01	BENZOIC ACID + SALICYLIC ACID 6% + 3% OINTMENT 15 g TUBE
BENZ70000000570OINTX00367TUB01	BENZOIC ACID + SALICYLIC ACID 6% + 3% OINTMENT 30 g TUBE
BENZY0000000588LOTI100201BOTTL	BENZYL BENZOATE 25% LOTION 120 mL BOTTLE
BENZY0000000588LOTI100469BOTTL	BENZYL BENZOATE 25% LOTION 60 mL BOTTLE
BERAC0000000985SUS1400516VIALX	BERACTANT 25 mg/mL SUSPENSION 1 mL VIAL
BERAC0000000985SUS1401386VIALX	BERACTANT 25 mg/mL SUSPENSION 4 mL VIAL
BETAH0005800063TAB4900000000000	BETAHISTINE (as TYDROCHLORIDE) 6 mg TABLET
BETAH0005800068TAB4900000000000	BETAHISTINE (as TYDROCHLORIDE) 8 mg TABLET
BETAH0005800323TAB4900000000000	BETAHISTINE (as TYDROCHLORIDE) 24 mg TABLET
BETAM0000000078OINTX00227TUB01	BETAMETHASONE 0.05% OINTMENT 15 g TUBE
BETAM0001900045TAB4900000000000	BETAMETHASONE (as DIPROPIONATE) 500 mcg TABLET
BETAM0001900078CREA301245TUB01	BETAMETHASONE (as DIPROPIONATE) 0.05% CREAM 5 g TUBE
BETAM0001900078CREA301380TUB01	BETAMETHASONE (as DIPROPIONATE) 0.05% CREAM 10 g TUBE
BETAM0001900078LOTI100379BOTTL	BETAMETHASONE (as DIPROPIONATE) 0.05% LOTION 30 mL BOTTLE
BETAM0001900082OINTX01245TUB01	BETAMETHASONE (as DIPROPIONATE) 0.10% OINTMENT 5 g TUBE
BETAM0001900590SOL3200361VIALX	BETAMETHASONE (as DIPROPIONATE) 5 mg (as dipropionate) + 2 mg (as sodium phosphate) per mL SOLUTION
BETAM0001900590SOL3200516AMPUL	BETAMETHASONE (as DIPROPIONATE) 5 mg (as dipropionate) + 2 mg (as sodium phosphate) per mL SOLUTION

BETAM0005900082CREA301245TUB01	BETAMETHASONE (as VALERATE) 0.10% CREAM 5 g TUBE
BETAX0000200085EDS0100060BOTTL	BETAXOLOL (as HYDROCHLORIDE) 0.25% EYE DROPS SUSPENSION 5 mL BOTTLE
BETAX0000200649EDS0100060BOTTL	BETAXOLOL (as HYDROCHLORIDE) 0.50% EYE DROPS SUSPENSION 5 mL BOTTLE
BETAX0000200649EDS0100195BOTTL	BETAXOLOL (as HYDROCHLORIDE) 0.50% EYE DROPS SUSPENSION 10 mL BOTTLE
BICAL0000000050TAB490000000000	BICALUTAMIDE 50 mg TABLET
BICAL00000000215TAB490000000000	BICALUTAMIDE 150 mg TABLET
BIPER0000200026TAB490000000000	BIPERIDEN (as HYDROCHLORIDE) 2 mg TABLET
BIPER0003200459SOL3200516AMPUL	BIPERIDEN (as LACTATE) 5 mg/mL SOLUTION 1 mL AMPULE
BIPHA00000000593SOL3200195VIALX	BIPHASIC ISOPHANE HUMAN INSULIN 70/30 (RECOMBINANT DNA) 70% isophane suspension + 30% soluble insuli
BIPHA00000000593SOL3200388DS001	BIPHASIC ISOPHANE HUMAN INSULIN 70/30 (RECOMBINANT DNA) 70% isophane suspension + 30% soluble insuli
BIPHA00000000593SOL3200388GC001	BIPHASIC ISOPHANE HUMAN INSULIN 70/30 (RECOMBINANT DNA) 70% isophane suspension + 30% soluble insuli
BIPHA00000000761SOL3200388GC001	BIPHASIC ISOPHANE HUMAN INSULIN 70/30 (RECOMBINANT DNA) 100 IU/mL SOLUTION 3 mL GLASS CARTRIDGE
BIPHA00000001339SOL3200060VIALX	BIPHASIC ISOPHANE HUMAN INSULIN 70/30 (RECOMBINANT DNA) 100 IU/mL (Equiv. to 3.571 mg/mL) SOLUTION 5
BISAC0000000056SUPP400000000000	BISACODYL 5 mg SUPPOSITORY
BISAC0000000056TAB490000000000	BISACODYL 5 mg TABLET
BISAC0000000069SUPP400000000000	BISACODYL 10 mg SUPPOSITORY
BISOP0002300003TAB490000000000	BISOPROLOL (as FUMARATE) 2.5 mg TABLET
BISOP0002300056TAB490000000000	BISOPROLOL (as FUMARATE) 5 mg TABLET
BLEOM0006000000POW2700228AMPUL	BLEOMYCIN (as SULFATE) POWDER 15 mg AMPULE
BLEOM0006000000POW2700228VIALX	BLEOMYCIN (as SULFATE) POWDER 15 mg VIAL
BMSDA00000001315SOL3200247BOTTL	BALANCED MULTIPLE MAINTENANCE SOLUTION WITH 5% DEXTROSE FOR ADULT 5% SOLUTION 1 L BOTTLE
BMSDP00000001315SOL3200020BOTTL	BALANCED MULTIPLE MAINTENANCE SOLUTION WITH 5% DEXTROSE FOR PEDIA 5% SOLUTION 250 mL BOTTLE
BMSDP00000001315SOL3200044BOTTL	BALANCED MULTIPLE MAINTENANCE SOLUTION WITH 5% DEXTROSE FOR PEDIA 5% SOLUTION 500 mL BOTTLE
BRIMO00000000084OPS0100060BOTTL	BRIMONIDINE TARTRATE 0.15% OPHTHALMIC SOLUTION 5 mL BOTTLE
BRINZ00000001062OPS0200060BOTTL	BRINZOLAMIDE 1% OPHTHALMIC SUSPENSION 5 mL BOTTLE
BROMA00000000121TAB490000000000	BROMAZEPAM 1.5 mg TABLET
BROMO0003700003TAB490000000000	BROMOCRIPTINE (as MESILATE) 2.5 mg TABLET
BUDE10002300605DPI0101388DIS01	BUDESONIDE + FORMOTEROL (as FUMARATE) 320 mcg + 9 mcg DRY POWDER 60 Doses DISPENSER
BUDE10002400602DPI0101387DIS01	BUDESONIDE + FORMOTEROL (as FUMARATE DIHYDRATE) 80 mcg + 4.5 mcg DRY POWDER 120 Doses DISPENSER
BUDE10002400602DPI0101388DIS01	BUDESONIDE + FORMOTEROL (as FUMARATE DIHYDRATE) 80 mcg + 4.5 mcg DRY POWDER 60 Doses DISPENSER
BUDE10002400604DPI0101387DIS01	BUDESONIDE + FORMOTEROL (as FUMARATE DIHYDRATE) 160 mcg + 4.5 mcg DRY POWDER 120 Doses DISPENSER
BUDE10002400604DPI0101388DIS01	BUDESONIDE + FORMOTEROL (as FUMARATE DIHYDRATE) 160 mcg + 4.5 mcg DRY POWDER 60 Doses DISPENSER
BUDES00000000427SOL3200361NEB01	BUDESONIDE 500mcg/ml SOLUTION 2 mL NEBULE
BUDES00000000563DPI0100000INH01	BUDESONIDE 100 mcg/dose DRY POWDER INHALER
BUDES00000000596DPI0100000INH01	BUDESONIDE 200 mcg/dose x 300 doses DRY POWDER INHALER
BUDES00000000597POW2700000INH01	BUDESONIDE 400 mcg/dose x 50 doses POWDER INHALER
BUDES00000000598NAS010000000000	BUDESONIDE 100 mcg/dose x 50 metered doses NASAL AQUEOUS SOLUTION
BUDES00000000599NAS010000000000	BUDESONIDE 100 mcg/dose x 200 metered doses NASAL AQUEOUS SOLUTION
BUDES00000000600DPI0100000INH01	BUDESONIDE 200 mcg/dose x 100 doses DRY POWDER INHALER

BUDES00000000600MDIXX00000INH01	BUDESONIDE 200 mcg/dose x 100 doses METERED DOSE INHALER INHALER
BUDES00000000771SOL3200361NEB01	BUDESONIDE 250 mcg/mL SOLUTION 2 mL NEBULE
BUDES00000001208DPI0100000INH01	BUDESONIDE 200 mcg/dose x 200 doses DRY POWDER INHALER
BUMET00000000252TAB4900000000000	BUMETANIDE 1 mg TABLET
BUMET00000000668SOL3201386AMPUL	BUMETANIDE 500 mcg/mL SOLUTION 4 mL AMPULE
BUPIV0000200649SOL3200011AMPUL	BUPIVACAINE (as HYDROCHLORIDE) 0.50% SOLUTION 20 mL AMPULE
BUPIV0000200649SOL3200011VIALX	BUPIVACAINE (as HYDROCHLORIDE) 0.50% SOLUTION 20 mL VIAL
BUPIV0000200649SOL3200055AMPUL	BUPIVACAINE (as HYDROCHLORIDE) 0.50% SOLUTION 50 mL AMPULE
BUPIV0000200649SOL3200055VIALX	BUPIVACAINE (as HYDROCHLORIDE) 0.50% SOLUTION 50 mL VIAL
BUPIV0000200649SOL3200060AMPUL	BUPIVACAINE (as HYDROCHLORIDE) 0.50% SOLUTION 5 mL AMPULE
BUPIV0000200649SOL3200060VIALX	BUPIVACAINE (as HYDROCHLORIDE) 0.50% SOLUTION 5 mL VIAL
BUPIV0000200649SOL3200195AMPUL	BUPIVACAINE (as HYDROCHLORIDE) 0.50% SOLUTION 10 mL AMPULE
BUPIV0000200649SOL3200195VIALX	BUPIVACAINE (as HYDROCHLORIDE) 0.50% SOLUTION 10 mL VIAL
BUPIV0000200649SOL3201389AMPUL	BUPIVACAINE (as HYDROCHLORIDE) 0.50% SOLUTION 4 mL (spinal) with 8% Dextrose AMPULE
BUTAM0001500050TAB4900000000000	BUTAMIRATE (as CITRATE) 50 mg TABLET
BUTAM0001500478SYRUP00201BOTTL	BUTAMIRATE (as CITRATE) 7.5 mg/5 mL SYRUP 120 mL BOTTLE
BUTAM0001500478SYRUP00469BOTTL	BUTAMIRATE (as CITRATE) 7.5 mg/5 mL SYRUP 60 mL BOTTLE
BUTOR0005600357SOL3200361VIALX	BUTORPHANOL (as TARTRATE) 2 mg/ml SOLUTION 2 mL VIAL
BUTOR0005600357SOL3200516VIALX	BUTORPHANOL (as TARTRATE) 2 mg/ml SOLUTION 1 mL VIAL
CAB010002800616OINTX00367TUB01	CALCIPOTRIOL + BETAMETHASONE (as HYDRATE + DIPROPIONATE) 50 mcg + 500 mcg/g OINTMENT 30 g TUBE
CABO200000000000POW2700040VIALX	CARBOPLATIN POWDER 450 mg VIAL
CABO200000000000POW2700215VIALX	CARBOPLATIN POWDER 150 mg VIAL
CABO200000000194SOL3200055VIALX	CARBOPLATIN 10 mg/mL SOLUTION 50 mL VIAL
CABO200000000194SOL3200060VIALX	CARBOPLATIN 10 mg/mL SOLUTION 5 mL VIAL
CABO200000000194SOL3200231VIALX	CARBOPLATIN 10 mg/mL SOLUTION 15 mL VIAL
CABO200000000194SOL3201390VIALX	CARBOPLATIN 10 mg/mL SOLUTION 45 mL VIAL
CALA100000000611LOT1100201BOTTL	CALAMINE, PLAIN 8% LOTION 120 mL BOTTLE
CALA100000000611LOT1100469BOTTL	CALAMINE, PLAIN 8% LOTION 60 mL BOTTLE
CALC10000000087CAPSU00000000000	CALCITRIOL 0.25 mcg CAPSULE
CALC100000000617CAPSU00000000000	CALCITRIOL 0.5 mcg CAPSULE
CALC20003200047TAB4900000000000	CALCIUM (as LACTATE) 500 mg TABLET
CALC300000000047TAB4900000000000	CALCIUM CARBONATE 500 mg TABLET
CALC300000000061TAB4900000000000	CALCIUM CARBONATE 600 mg TABLET
CALC400000000196TAB4900000000000	CALCIUM CARBONATE + CHOLECALCIFEROL (VIT. D3) 1.25 g (equiv. to 500 mg elemental calcium) + 250 IU T
CALC400000000485TAB4900000000000	CALCIUM CARBONATE + CHOLECALCIFEROL (VIT. D3) 750 mg (equiv. to 300 mg elemental calcium) + 150 IU T
CALC400000001267TAB4900000000000	CALCIUM CARBONATE + CHOLECALCIFEROL (VIT. D3) Equiv. to 500 mg elemental Ca + 200 IU Vit D3 TABLET
CALC400000001275TAB4900000000000	CALCIUM CARBONATE + CHOLECALCIFEROL (VIT. D3) 600 mg elementa Ca + 400 IU Vit. D3 TABLET
CALC400000001278TAB4900000000000	CALCIUM CARBONATE + CHOLECALCIFEROL (VIT. D3) 600 mg elemental Ca + 200 IU Vit. D3 TABLET
CALC400000001306TAB4900000000000	CALCIUM CARBONATE + CHOLECALCIFEROL (VIT. D3) 500 mg elemental Ca + 400 IU Vit. D3 TABLET
CALC60000000050POW2700000VIALX	CALCIUM FOLINATE (LEUCOVORIN CALCIUM) 50 mg POWDER VIAL
CALC600000000139POW2700000VIALX	CALCIUM FOLINATE (LEUCOVORIN CALCIUM) 100 mg POWDER VIAL
CALC600000000194SOL3200060AMPUL	CALCIUM FOLINATE (LEUCOVORIN CALCIUM) 10 mg/mL SOLUTION 5 mL AMPULE

CALC60000000194SOL3200060VIALX	CALCIUM FOLINATE (LEUCOVORIN CALCIUM) 10 mg/mL SOLUTION 5 mL VIAL
CALC60000000194SOL3200195AMPUL	CALCIUM FOLINATE (LEUCOVORIN CALCIUM) 10 mg/mL SOLUTION 10 mL AMPULE
CALC60000000194SOL3200195VIALX	CALCIUM FOLINATE (LEUCOVORIN CALCIUM) 10 mg/mL SOLUTION 10 mL VIAL
CALC60000000194SOL3200388AMPUL	CALCIUM FOLINATE (LEUCOVORIN CALCIUM) 10 mg/mL SOLUTION 3 mL AMPULE
CALC60000000194SOL3200388VIALX	CALCIUM FOLINATE (LEUCOVORIN CALCIUM) 10 mg/mL SOLUTION 3 mL VIAL
CALC60000000387SOL3200195AMPUL	CALCIUM FOLINATE (LEUCOVORIN CALCIUM) 3 mg/mL SOLUTION 10 mL AMPULE
CALC60000000387SOL3200516AMPUL	CALCIUM FOLINATE (LEUCOVORIN CALCIUM) 3 mg/mL SOLUTION 1 mL AMPULE
CALC60000000622SOL3200361AMPUL	CALCIUM FOLINATE (LEUCOVORIN CALCIUM) 7.5 mg/mL SOLUTION 2 mL AMPULE
CALC60000600228CAPSU0000000000	CALCIUM FOLINATE (LEUCOVORIN CALCIUM) (as ANHYDROUS) 15 mg CAPSULE
CALC60000600228TAB490000000000	CALCIUM FOLINATE (LEUCOVORIN CALCIUM) (as ANHYDROUS) 15 mg TABLET
CALC60000600619TAB490000000000	CALCIUM FOLINATE (LEUCOVORIN CALCIUM) (as ANHYDROUS) 25 mg (equiv. to 25 mg folinic acid) TABLET
CALC70000000133SOL3200011BOTTL	CALCIUM GLUCONATE 10% SOLUTION 20 mL BOTTLE
CALC70000000133SOL3200195AMPUL	CALCIUM GLUCONATE 10% SOLUTION 10 mL AMPULE
CALC70000000133SOL3200195VIALX	CALCIUM GLUCONATE 10% SOLUTION 10 mL VIAL
CALC70000000133SOL3201375BOTTL	CALCIUM GLUCONATE 10% SOLUTION 25 mL BOTTLE
CALCT0000000614CREA300367BOTTL	CALCIPOTRIOL 50 mcg/g CREAM 30 g BOTTLE
CALCT0000000614CREA300367TUB01	CALCIPOTRIOL 50 mcg/g CREAM 30 g TUBE
CALCT0000000614OINTX00367BOTTL	CALCIPOTRIOL 50 mcg/g OINTMENT 30 g BOTTLE
CALCT0000000614OINTX00367TUB01	CALCIPOTRIOL 50 mcg/g OINTMENT 30 g TUBE
CALCT0000000615SOL3200379BOTTL	CALCIPOTRIOL 50 mcg/mL SOLUTION 30 mL BOTTLE
CAN010001400068TAB490000000000	CANDESARTAN (as CILEXETIL) 8 mg TABLET
CAN010001400235TAB490000000000	CANDESARTAN (as CILEXETIL) 16 mg TABLET
CAPEC0000000047TAB490000000000	CAPECITABINE 500 mg TABLET
CAPEC00000000215TAB490000000000	CAPECITABINE 150 mg TABLET
CAPTO0000000021TAB490000000000	CAPTOPRIL 25 mg TABLET
CAPTO0000000050TAB490000000000	CAPTOPRIL 50 mg TABLET
CARB10000000004TAB490000000000	CARBAMAZEPINE 200 mg TABLET
CARB10000000034TAB490000000000	CARBAMAZEPINE 400 mg TABLET
CARB10000000162SYRUP00201BOTTL	CARBAMAZEPINE 100 mg/5 mL SYRUP 120 mL BOTTLE
CARBA0000000075SOL3200122VIALX	CARBACHOL 0.01% SOLUTION 1.5 mL VIAL
CARBI0000000007TAB490000000000	CARBIMAZOLE 20 mg TABLET
CARBI0000000056TAB490000000000	CARBIMAZOLE 5 mg TABLET
CARMU0000000630POW2700000VIALX	CARMUSTINE 100 mg + 3 mL diluent POWDER VIAL
CARVE0000000021TAB490000000000	CARVEDILOL 25 mg TABLET
CARVE0000000466TAB490000000000	CARVEDILOL 6.25 mg TABLET
CASTO00000000000000000469BOTTL	CASTOR OIL 60 mL BOTTLE
CASTO00000000000USP0100000BOTTL	CASTOR OIL USP GRADE BOTTLE
CEFAD0003800015CAPSU0000000000	CEFADROXIL (as MONOHYDRATE) 250 mg CAPSULE
CEFAD0003800047CAPSU0000000000	CEFADROXIL (as MONOHYDRATE) 500 mg CAPSULE
CEFAD0003800209SYRUP00469BOTTL	CEFADROXIL (as MONOHYDRATE) 125 mg/5 mL SYRUP 60 mL BOTTLE
CEFAD0003801265SUS1400469BOTTL	CEFADROXIL (as MONOHYDRATE) 250 mg/5 mL SUSPENSION 60 mL BOTTLE
CEFAL0003800015CAPSU0000000000	CEFALEXIN (as MONOHYDRATE) 250 mg CAPSULE

CEFAL0003800047CAPSU0000000000	CEFALEXIN (as MONOHYDRATE) 500 mg CAPSULE
CEFAL0003800142DROPS00195BOTTL	CEFALEXIN (as MONOHYDRATE) 100 mg/ml DROPS 10 mL BOTTLE
CEFAL0003800142DROPS00231BOTTL	CEFALEXIN (as MONOHYDRATE) 100 mg/ml DROPS 15 mL BOTTLE
CEFAL000380020900000000469BOTTL	CEFALEXIN (as MONOHYDRATE) 125 mg/5 mL 60 mL BOTTLE
CEFAL0003800209SYRUP00379BOTTL	CEFALEXIN (as MONOHYDRATE) 125 mg/5 mL SYRUP 30 mL BOTTLE
CEFAL0003800241TAB4900000000000	CEFALEXIN (as MONOHYDRATE) 1 g TABLET
CEFAL000380126500000000379BOTTL	CEFALEXIN (as MONOHYDRATE) 250 mg/5 mL 30 mL BOTTLE
CEFAL000380126500000000469BOTTL	CEFALEXIN (as MONOHYDRATE) 250 mg/5 mL 60 mL BOTTLE
CEFAZ0005300047SOL3200000VIALX	CEFAZOLIN (as SODIUM SALT) 500 mg SOLUTION VIAL
CEFAZ0005300241SOL3200000VIALX	CEFAZOLIN (as SODIUM SALT) 1 g SOLUTION VIAL
CEFEP0000200025SOL3200000VIALX	CEFEPIME (as HYDROCHLORIDE) 2 g SOLUTION VIAL
CEFEP0000200047SOL3200000VIALX	CEFEPIME (as HYDROCHLORIDE) 500 mg SOLUTION VIAL
CEFEP0000200241SOL3200000VIALX	CEFEPIME (as HYDROCHLORIDE) 1 g SOLUTION VIAL
CEFIX0000000004CAPSU0000000000	CEFIXIME 200 mg CAPSULE
CEFIX0000000034CAPSU0000000000	CEFIXIME 400 mg CAPSULE
CEFIX0000000139CAPSU0000000000	CEFIXIME 100 mg CAPSULE
CEFIX0000000162SUS1400379BOTTL	CEFIXIME 100 mg/5 mL SUSPENSION 30 mL BOTTLE
CEFIX0000000162SUS1400469BOTTL	CEFIXIME 100 mg/5 mL SUSPENSION 60 mL BOTTLE
CEFIX00000001199SUS1400195BOTTL	CEFIXIME 20 mg/mL SUSPENSION 10 mL BOTTLE
CEFOT0005300636SOL3200000VIALX	CEFOTAXIME (as SODIUM SALT) 250 mg + 2 mL diluent SOLUTION VIAL
CEFOT0005300637SOL3200000VIALX	CEFOTAXIME (as SODIUM SALT) 500 mg + 2 mL diluent SOLUTION VIAL
CEFOT0005300638SOL3200000VIALX	CEFOTAXIME (as SODIUM SALT) 1 g + 4 mL diluent SOLUTION VIAL
CEFOX0005300241SOL3200000VIALX	CEFOXITIN (as SODIUM SALT) 1 g SOLUTION VIAL
CEFOX0005301321SOL3200000DICO1	CEFOXITIN (as SODIUM SALT) 1 g with 4% dextrose SOLUTION DUPLEX III CONTAINER
CEFTA0004600015SOL3200000VIALX	CEFTAZIDIME (as PENTAHYDRATE) 250 mg SOLUTION VIAL
CEFTA0004600025SOL3200000VIALX	CEFTAZIDIME (as PENTAHYDRATE) 2 g SOLUTION VIAL
CEFTA0004600047SOL3200000VIALX	CEFTAZIDIME (as PENTAHYDRATE) 500 mg SOLUTION VIAL
CEFTA0004600241SOL3200000VIALX	CEFTAZIDIME (as PENTAHYDRATE) 1 g SOLUTION VIAL
CEFTR00000001322SOL3600000DICO1	CEFTRIAXONE 1 g with 3.74% dextrose SOLUTION OF LIDOCAINE DUPLEX III CONTAINER
CEFTR0002000639SOL3600000VIALX	CEFTRIAXONE (as DISODIUM/SODIUM SALT) 250 mg + 2 mL 1% solution of lidocaine SOLUTION OF LIDOCAINE
CEFTR0002000640SOL3600000VIALX	CEFTRIAXONE (as DISODIUM/SODIUM SALT) 500 mg + 2 mL 1% solution of lidocaine SOLUTION OF LIDOCAINE
CEFTR0002000641SOL3200000VIALX	CEFTRIAXONE (as DISODIUM/SODIUM SALT) 250 mg + 5 mL diluent SOLUTION VIAL
CEFTR0002000642SOL3200000VIALX	CEFTRIAXONE (as DISODIUM/SODIUM SALT) 500 mg + 5 mL diluent SOLUTION VIAL
CEFTR0002000643SOL3600000VIALX	CEFTRIAXONE (as DISODIUM/SODIUM SALT) 1 g + 3.5 mL 1% solution of lidocaine SOLUTION OF LIDOCAINE
CEFTR0002000644SOL3200000VIALX	CEFTRIAXONE (as DISODIUM/SODIUM SALT) 1 g + 10 mL diluent SOLUTION VIAL
CEFUR0000000015SOL3200000VIALX	CEFUROXIME 250 mg SOLUTION VIAL
CEFUR0000000064SOL3200000VIALX	CEFUROXIME 750 mg SOLUTION VIAL
CEFUR0000001265SUS1400469BOTTL	CEFUROXIME 250 mg/5 mL SUSPENSION 60 mL BOTTLE
CEFUR0000700015TAB4900000000000	CEFUROXIME (as AXETIL) 250 mg TABLET
CEFUR0000700047TAB4900000000000	CEFUROXIME (as AXETIL) 500 mg TABLET
CEFUR0000700205SAC0100000SAC01	CEFUROXIME (as AXETIL) 125 mg SACHET SACHET
CEFUR0000700209SUS1400055BOTTL	CEFUROXIME (as AXETIL) 125 mg/5 mL SUSPENSION 50 mL BOTTLE
CEFUR0000700209SUS1401391BOTTL	CEFUROXIME (as AXETIL) 125 mg/5 mL SUSPENSION 70 mL BOTTLE

CEFUR0005300120SOL3200000VIALX	CEFUROXIME (as SODIUM SALT) 1.5 g SOLUTION VIAL
CELEC0000000004CAPSU0000000000	CELECOXIB 200 mg CAPSULE
CELEC0000000034CAPSU0000000000	CELECOXIB 400 mg CAPSULE
CELEC0000000139CAPSU0000000000	CELECOXIB 100 mg CAPSULE
CETIR0001800069TAB490000000000	CETIRIZINE (as DIHYDROCHLORIDE) 10 mg TABLET
CETIR0001800194DROPS00195BOTTL	CETIRIZINE (as DIHYDROCHLORIDE) 10 mg/mL DROPS 10 mL BOTTLE
CETIR0001800261SOL3200469BOTTL	CETIRIZINE (as DIHYDROCHLORIDE) 1 mg/mL SOLUTION 60 mL BOTTLE
CETIR0001800454SYRUP00379BOTTL	CETIRIZINE (as DIHYDROCHLORIDE) 5 mg/5 mL SYRUP 30 mL BOTTLE
CHICK0000000281SOL3200000VIALX	RABIES CHICK EMBRYO CELL (PURIFIED INACTIVATED) 2.5 IU/mL, 1 dose + 1 mL diluent SOLUTION VIAL
CHLO10003400041TAB490000000000	CHLORPHENAMINE (as MALEATE) 4 mg TABLET
CHLO10003400194SOL3200195VIALX	CHLORPHENAMINE (as MALEATE) 10 mg/mL SOLUTION 10 mL VIAL
CHLO10003400194SOL3200516AMPUL	CHLORPHENAMINE (as MALEATE) 10 mg/mL SOLUTION 1 mL AMPULE
CHLO10003400194SOL3200516VIALX	CHLORPHENAMINE (as MALEATE) 10 mg/mL SOLUTION 1 mL VIAL
CHLO10003400349SYRUP00469BOTTL	CHLORPHENAMINE (as MALEATE) 2 mg/5 mL SYRUP 60 mL BOTTLE
CHLO10003400663SYRUP00469BOTTL	CHLORPHENAMINE (as MALEATE) 2.5 mg/5 mL SYRUP 60 mL BOTTLE
CHLO300000001062EYO0100025TUB01	CHLORAMPHENICOL 1% EYE OINTMENT 2 g TUBE
CHLO300000001062EYO0100682TUB01	CHLORAMPHENICOL 1% EYE OINTMENT 4 g TUBE
CHLO300000001062EYO0101393TUB01	CHLORAMPHENICOL 1% EYE OINTMENT 3.5 g TUBE
CHLO40000000649EDS0200060BOTTL	CHLORAMPHENICOL 0.50% EYE DROPS SOLUTION 5 mL BOTTLE
CHLO40000000649EDS0200195BOTTL	CHLORAMPHENICOL 0.50% EYE DROPS SOLUTION 10 mL BOTTLE
CHLO40000000649EDS0201392BOTTL	CHLORAMPHENICOL 0.50% EYE DROPS SOLUTION 7.5 mL BOTTLE
CHLO50000000649SOL3700060BOTTL	CHLORAMPHENICOL EAR DROPS 0.50% EAR DROPS SOLUTION 5 mL BOTTLE
CHLO10000000645SYRUP000000BOTTL	CHLORAL HYDRATE 500 mg/5 mL SYRUP BOTTLE
CHLO20000000026TAB490000000000	CHLORAMBUCIL 2 mg TABLET
CHLO30000000015CAPSU0000000000	CHLORAMPHENICOL 250 mg CAPSULE
CHLO30000000047CAPSU0000000000	CHLORAMPHENICOL 500 mg CAPSULE
CHLO30004400209SUS1400379BOTTL	CHLORAMPHENICOL (as PALMITATE) 125 mg/5 mL SUSPENSION 30 mL BOTTLE
CHLO30004400209SUS1400469BOTTL	CHLORAMPHENICOL (as PALMITATE) 125 mg/5 mL SUSPENSION 60 mL BOTTLE
CHLO30005400241SOL3200000VIALX	CHLORAMPHENICOL (as SODIUM SUCCINATE) 1 g SOLUTION VIAL
CHLO50002600655SOL3200044BOTTL	CHLORHEXIDINE (as GLUCONATE) 0.12% SOLUTION 500 mL BOTTLE
CHLO50002600655SOL3200055BOTTL	CHLORHEXIDINE (as GLUCONATE) 0.12% SOLUTION 50 mL BOTTLE
CHLO50002600655SOL3200201BOTTL	CHLORHEXIDINE (as GLUCONATE) 0.12% SOLUTION 120 mL BOTTLE
CHLO50002600655SOL3201347BOTTL	CHLORHEXIDINE (as GLUCONATE) 0.12% SOLUTION 5 L BOTTLE
CHLO50002600655SOL3201394BOTTL	CHLORHEXIDINE (as GLUCONATE) 0.12% SOLUTION 380 mL BOTTLE
CHLO50002600658SOL3200044BOTTL	CHLORHEXIDINE (as GLUCONATE) 4% SOLUTION 500 mL BOTTLE
CHLO50002600658SOL3200055BOTTL	CHLORHEXIDINE (as GLUCONATE) 4% SOLUTION 50 mL BOTTLE
CHLO50002600658SOL3200201BOTTL	CHLORHEXIDINE (as GLUCONATE) 4% SOLUTION 120 mL BOTTLE
CHLO50002600658SOL3201347BOTTL	CHLORHEXIDINE (as GLUCONATE) 4% SOLUTION 5 L BOTTLE
CHLO50002600658SOL3201394BOTTL	CHLORHEXIDINE (as GLUCONATE) 4% SOLUTION 380 mL BOTTLE
CHLO50002600658SOL3201395BOTTL	CHLORHEXIDINE (as GLUCONATE) 4% SOLUTION 4 L BOTTLE
CHLO70004800015TAB490000000000	CHLOROQUINE (as PHOSPHATE OR DIPHOSPHATE) 250 mg TABLET
CHLO70004800851SOL3200011VIALX	CHLOROQUINE (as PHOSPHATE OR DIPHOSPHATE) 50 mg/mL SOLUTION 20 mL VIAL
CHLOR0000200004TAB490000000000	CHLORPROMAZINE (as HYDROCHLORIDE) 200 mg TABLET

CHLOR0000200050TAB490000000000	CHLORPROMAZINE (as HYDROCHLORIDE) 50 mg TABLET
CHLOR0000200139TAB490000000000	CHLORPROMAZINE (as HYDROCHLORIDE) 100 mg TABLET
CHLOR0000200985SOL3200516AMPUL	CHLORPROMAZINE (as HYDROCHLORIDE) 25 mg/mL SOLUTION 1 mL AMPULE
CHLP20000000015TAB490000000000	CHLORPROPAMIDE 250 mg TABLET
CICLO0000000021CAPSU0000000000	CICLOSPORIN 25 mg CAPSULE
CICLO0000000050CAPSU0000000000	CICLOSPORIN 50 mg CAPSULE
CICLO0000000078EMUL501396BOTTL	CICLOSPORIN 0.05% EMULSION 0.4 mL BOTTLE
CICLO0000000139CAPSU0000000000	CICLOSPORIN 100 mg CAPSULE
CICLO0000000142SOL3200055BOTTL	CICLOSPORIN 100 mg/ml SOLUTION 50 mL BOTTLE
CICLO00000000851SOL3200060AMPUL	CICLOSPORIN 50 mg/mL SOLUTION 5 mL AMPULE
CILOS0000000050TAB490000000000	CILOSTAZOL 50 mg TABLET
CILOS0000000139TAB490000000000	CILOSTAZOL 100 mg TABLET
CINNA0000000021CAPSU0000000000	CINNARIZINE 25 mg CAPSULE
CINNA0000000050CAPSU0000000000	CINNARIZINE 50 mg CAPSULE
CINNA0000000065CAPSU0000000000	CINNARIZINE 75 mg CAPSULE
CIPRO0000000047TAB520000000000	CIPROFLOXACIN 500 mg TABLET MODIFIED RELEASE
CIPRO0000000241TAB490000000000	CIPROFLOXACIN 1 g TABLET
CIPRO00000001456EYD010000000000	CIPROFLOXACIN 0.3% x 5 mL EYE DROPS
CIPRO00000200015TAB490000000000	CIPROFLOXACIN (as HYDROCHLORIDE) 250 mg TABLET
CIPRO00000200047TAB490000000000	CIPROFLOXACIN (as HYDROCHLORIDE) 500 mg TABLET
CIPRO00003200357SOL3200055VIALX	CIPROFLOXACIN (as LACTATE) 2 mg/ml SOLUTION 50 mL VIAL
CIPRO00003200357SOL3200178VIALX	CIPROFLOXACIN (as LACTATE) 2 mg/ml SOLUTION 100 mL VIAL
CIPRO00003200357SOL3201125VIALX	CIPROFLOXACIN (as LACTATE) 2 mg/ml SOLUTION 200 mL VIAL
CISPL0000000050POW2700000VIALX	CISPLATIN 50 mg POWDER VIAL
CISPL0000000069POW2700000VIALX	CISPLATIN 10 mg POWDER VIAL
CISPL0000000139SOL3200000VIALX	CISPLATIN 100 mg SOLUTION VIAL
CISPL0000000261SOL3200011VIALX	CISPLATIN 1 mg/mL SOLUTION 20 mL VIAL
CISPL0000000261SOL3200055VIALX	CISPLATIN 1 mg/mL SOLUTION 50 mL VIAL
CISPL0000000261SOL3200195VIALX	CISPLATIN 1 mg/mL SOLUTION 10 mL VIAL
CISPL0000000668SOL3200011VIALX	CISPLATIN 500 mcg/mL SOLUTION 20 mL VIAL
CISPL0000000668SOL3200055VIALX	CISPLATIN 500 mcg/mL SOLUTION 50 mL VIAL
CISPL0000000668SOL3200178VIALX	CISPLATIN 500 mcg/mL SOLUTION 100 mL VIAL
CITIC00000000470000000000VIALX	CITICOLINE 500 mg VIAL
CITIC00000001445SOL3200000VIALX	CITICOLINE 1000 mg SOLUTION VIAL
CLARI0000000015TAB490000000000	CLARITHROMYCIN 250 mg TABLET
CLARI0000000047TAB490000000000	CLARITHROMYCIN 500 mg TABLET
CLARI0000000209SUS1400055BOTTL	CLARITHROMYCIN 125 mg/5 mL SUSPENSION 50 mL BOTTLE
CLARI0000000209SUS1401375BOTTL	CLARITHROMYCIN 125 mg/5 mL SUSPENSION 25 mL BOTTLE
CLARI0000000209SUS1401391BOTTL	CLARITHROMYCIN 125 mg/5 mL SUSPENSION 70 mL BOTTLE
CLARI00000001265SUS1401391BOTTL	CLARITHROMYCIN 250 mg/5 mL SUSPENSION 70 mL BOTTLE
CLEVU0000000031TAB490000000000	CLEVUDINE 30 mg TABLET
CLIND0000200028CAPSU0000000000	CLINDAMYCIN (as HYDROCHLORIDE) 300 mg CAPSULE
CLIND0000200215CAPSU0000000000	CLINDAMYCIN (as HYDROCHLORIDE) 150 mg CAPSULE
CLIND0004500491SUS1400469BOTTL	CLINDAMYCIN (as PALMITATE HYDROCHLORIDE) 75 mg/5 mL SUSPENSION 60 mL BOTTLE
CLIND0004700224SOL3200361AMPUL	CLINDAMYCIN (as PHOSPHATE) 150 mg/mL SOLUTION 2 mL AMPULE
CLIND0004700224SOL3200361VIALX	CLINDAMYCIN (as PHOSPHATE) 150 mg/mL SOLUTION 2 mL VIAL

CLIND0004700224SOL3201386AMPUL	CLINDAMYCIN (as PHOSPHATE) 150 mg/mL SOLUTION 4 mL AMPULE
CLOBE0005000078CREA300227TUB01	CLOBETASOL (as PROPIONATE) 0.05% CREAM 15 g TUBE
CLOBE0005000078CREA301245TUB01	CLOBETASOL (as PROPIONATE) 0.05% CREAM 5 g TUBE
CLOBE0005000078CREA301380TUB01	CLOBETASOL (as PROPIONATE) 0.05% CREAM 10 g TUBE
CLOBE0005000078OINTX00227TUB01	CLOBETASOL (as PROPIONATE) 0.05% OINTMENT 15 g TUBE
CLOBE0005000078OINTX01245TUB01	CLOBETASOL (as PROPIONATE) 0.05% OINTMENT 5 g TUBE
CLOBE0005000078OINTX01380TUB01	CLOBETASOL (as PROPIONATE) 0.05% OINTMENT 10 g TUBE
CLOBE0005000078SHAMX01397BOTTL	CLOBETASOL (as PROPIONATE) 0.05% SHAMPOO 27 mL BOTTLE
CLOFA0000000050CAPSU0000000000	CLOFAZIMINE 50 mg CAPSULE
CLOFA0000000139CAPSU0000000000	CLOFAZIMINE 100 mg CAPSULE
CLOMC0001500050TAB490000000000	CLOMIFENE (as CITRATE) 50 mg TABLET
CLONA0000000026TAB490000000000	CLONAZEPAM 2 mg TABLET
CLONI0000200218TAB490000000000	CLONIDINE (as HYDROCHLORIDE) 150 mcg TABLET
CLONI0000200220SOL3200516AMPUL	CLONIDINE (as HYDROCHLORIDE) 150 mcg/mL SOLUTION 1 mL AMPULE
CLONI0000200486TAB490000000000	CLONIDINE (as HYDROCHLORIDE) 75 mcg TABLET
CLOPI0000000065TAB490000000000	CLOPIDOGREL 75 mg TABLET
CLOTZ0000000000LOT1100195BOTTL	CLOTRIMAZOLE LOTION 10 mL BOTTLE
CLOTZ0000000000SHAMX00195SAC01	CLOTRIMAZOLE SHAMPOO 10 mL SACHET
CLOTZ0000000000SHAMX00469BOTTL	CLOTRIMAZOLE SHAMPOO 60 mL BOTTLE
CLOTZ0000000000SHAMX00842SAC01	CLOTRIMAZOLE SHAMPOO 6 mL SACHET
CLOTZ0000000028OVULE0000000000	CLOTRIMAZOLE 300 mg OVULE
CLOTZ0000000047OVULE0000000000	CLOTRIMAZOLE 500 mg OVULE
CLOTZ00000000139OVULE0000000000	CLOTRIMAZOLE 100 mg OVULE
CLOTZ0000000017SHAMX00178BOTTL	CLOTRIMAZOLE 100 mL SHAMPOO 100 mL BOTTLE
CLOTZ00000000194SOL3200231BOTTL	CLOTRIMAZOLE 10 mg/mL SOLUTION 15 mL BOTTLE
CLOTZ00000000194SOL3201375BOTTL	CLOTRIMAZOLE 10 mg/mL SOLUTION 25 mL BOTTLE
CLOTZ000000010620000001334TUB01	CLOTRIMAZOLE 1% 400 g TUBE
CLOTZ00000001062CREA300227TUB01	CLOTRIMAZOLE 1% CREAM 15 g TUBE
CLOTZ00000001062CREA301393TUB01	CLOTRIMAZOLE 1% CREAM 3.5 g TUBE
CLOTZ00000001314CREA300227TUB01	CLOTRIMAZOLE 2% CREAM 15 g TUBE
CLOTZ00000001314CREA301393TUB01	CLOTRIMAZOLE 2% CREAM 3.5 g TUBE
CLOTZ00000001314CREA301398TUB01	CLOTRIMAZOLE 2% CREAM 450 g TUBE
CLOXA0005300015CAPSU0000000000	CLOXACILLIN (as SODIUM SALT) 250 mg CAPSULE
CLOXA0005300047CAPSU0000000000	CLOXACILLIN (as SODIUM SALT) 500 mg CAPSULE
CLOXA0005300209SUS1400469BOTTL	CLOXACILLIN (as SODIUM SALT) 125 mg/5 mL SUSPENSION 60 mL BOTTLE
CLOXA0005301265SUS1400469BOTTL	CLOXACILLIN (as SODIUM SALT) 250 mg/5 mL SUSPENSION 60 mL BOTTLE
CLOZA0000000021TAB490000000000	CLOZAPINE 25 mg TABLET
CLOZA00000000139TAB490000000000	CLOZAPINE 100 mg TABLET
CML010005100649EYD0101386BOTTL	CARBOXYMETHYLCELLULOSE (as SODIUM) 0.50% EYE DROPS 4 mL BOTTLE
CO-AM0000000061POW2700000VIALX	CO-AMOXICLAV (AMOXICILLIN + POTASSIUM CLAVULANATE) 600 mg POWDER VIAL
CO-AM00000000334TAB490000000000	CO-AMOXICLAV (AMOXICILLIN + POTASSIUM CLAVULANATE) 250 mg + 125 mg TABLET
CO-AM00000000399SUS1401463BOTTL	CO-AMOXICLAV (AMOXICILLIN + POTASSIUM CLAVULANATE) 400 mg + 57 mg/5 ml SUSPENSION 35 mL BOTTLE
CO-AM00000001451POW2700000VIALX	CO-AMOXICLAV (AMOXICILLIN + POTASSIUM CLAVULANATE) 1.2 g POWDER VIAL

CO-AM0000001452SUS1400469BOTTL	CO-AMOXICLAV (AMOXICILLIN + POTASSIUM CLAVULANATE) 125 mg + 31.25 mg/5 mL SUSPENSION 60 mL BOTTLE
CO-AM0005700206SUS1400379BOTTL	CO-AMOXICLAV (AMOXICILLIN + POTASSIUM CLAVULANATE) (as TRIHYDRATE) 125 mg + 31 mg/5 mL SUSPENSION 3
CO-AM0005700206SUS1400469BOTTL	CO-AMOXICLAV (AMOXICILLIN + POTASSIUM CLAVULANATE) (as TRIHYDRATE) 125 mg + 31 mg/5 mL SUSPENSION 6
CO-AM0005700241TAB4900000000000	CO-AMOXICLAV (AMOXICILLIN + POTASSIUM CLAVULANATE) (as TRIHYDRATE) 1 g TABLET
CO-AM00057003010000001391BOTTL	CO-AMOXICLAV (AMOXICILLIN + POTASSIUM CLAVULANATE) (as TRIHYDRATE) 200 mg + 28.5 mg/5 mL 70 mL BOT
CO-AM0005700336SUS1400178BOTTL	CO-AMOXICLAV (AMOXICILLIN + POTASSIUM CLAVULANATE) (as TRIHYDRATE) 250 mg + 62.5 mg/5 mL SUSPENSION
CO-AM0005700336SUS1400469BOTTL	CO-AMOXICLAV (AMOXICILLIN + POTASSIUM CLAVULANATE) (as TRIHYDRATE) 250 mg + 62.5 mg/5 mL SUSPENSION
CO-AM0005700399SUS1400379BOTTL	CO-AMOXICLAV (AMOXICILLIN + POTASSIUM CLAVULANATE) (as TRIHYDRATE) 400 mg + 57 mg/5 ml SUSPENSION 3
CO-AM0005700399SUS1401391BOTTL	CO-AMOXICLAV (AMOXICILLIN + POTASSIUM CLAVULANATE) (as TRIHYDRATE) 400 mg + 57 mg/5 ml SUSPENSION 7
CO-AM0005701241TAB4900000000000	CO-AMOXICLAV (AMOXICILLIN + POTASSIUM CLAVULANATE) (as TRIHYDRATE) 500 mg + 125 mg TABLET
COAL10000000082SHAMX01374BOTTL	COAL TAR 0.10% SHAMPOO 150 mL BOTTLE
COAL100000000279SHAMX01399BOTTL	COAL TAR 2.50% SHAMPOO 130 ml BOTTLE
COAL100000000649SHAMX01399BOTTL	COAL TAR 0.50% SHAMPOO 130 ml BOTTLE
COAL100000001062SHAMX01374BOTTL	COAL TAR 1% SHAMPOO 150 mL BOTTLE
COAL100000001315SHAMX01399BOTTL	COAL TAR 5% SHAMPOO 130 ml BOTTLE
COBRA00000000679SOL3200516AMPUL	COBRA ANTIVENIN 800 MU/4.8 mL SOLUTION 1 mL AMPULE
COBRA00000000680SOL3200060AMPUL	COBRA ANTIVENIN 800 IU/mL SOLUTION 5 mL AMPULE
CODEI00000000184SUS1400469BOTTL	CODEINE 10 mg/5 mL SUSPENSION 60 mL BOTTLE
CODEI0004700031CAPSU00000000000	CODEINE (as PHOSPHATE) 30 mg CAPSULE
COLCH0000000045TAB4900000000000	COLCHICINE 500 mcg TABLET
COLES00000000682POW2700000SAC01	COLESTYRAMINE 4 g POWDER SACHET
COMAM0000000044SOL3200000BOTTL	COMBINED AMINO ACIDS SORBITOL 500 mL SOLUTION BOTTLE
COMAM00000000178SOL3200000BOTTL	COMBINED AMINO ACIDS SORBITOL 100 mL SOLUTION BOTTLE
COMBI0000000044SOL3200000BOTTL	COMBINED GLUCOSE-AMINO ACID SOLUTIONS 500 mL SOLUTION BOTTLE
COMBI00000000178SOL3200000BOTTL	COMBINED GLUCOSE-AMINO ACID SOLUTIONS 100 mL SOLUTION BOTTLE
COMBI000000001125SOL3200000BOTTL	COMBINED GLUCOSE-AMINO ACID SOLUTIONS 200 mL SOLUTION BOTTLE
COMBI000000001309SOL3200000BOTTL	COMBINED GLUCOSE-AMINO ACID SOLUTIONS 300 mL SOLUTION BOTTLE
COMBI000000001310SOL3200000BOTTL	COMBINED GLUCOSE-AMINO ACID SOLUTIONS 400 mL SOLUTION BOTTLE
CONJ100000000684TAB4900000000000	CONJUGATED EQUINE ESTROGEN + MEDROXYPROGESTERONE ACETATE 625 mcg + 2.5 mg TABLET
CONJ100000000685TAB4900000000000	CONJUGATED EQUINE ESTROGEN + MEDROXYPROGESTERONE ACETATE 625 mcg + 5 mg TABLET
CONJ200000000470TAB4900000000000	CONJUGATED ESTROGENS 625 mcg TABLET
CONJ200000000683TAB4900000000000	CONJUGATED ESTROGENS 1.25 mg TABLET
CONJ200000000686POW2700000VIALX	CONJUGATED ESTROGENS 25 mg powder + 5 mL diluent POWDER VIAL
CONJ200000000687TAB4900000000000	CONJUGATED ESTROGENS 300 mcg TABLET
CONJU00000000470TAB4900000000000	CONJUGATED EQUINE ESTROGEN 625 mcg TABLET
CONJU00000000683TAB4900000000000	CONJUGATED EQUINE ESTROGEN 1.25 mg TABLET
CONJU00000000687TAB4900000000000	CONJUGATED EQUINE ESTROGEN 300 mcg TABLET
COTR100000000035CAPSU00000000000	COTRIMOXAZOLE (SULFAMETHOXAZOLE + TRIMETHOPRIM) 400 mg + 80 mg CAPSULE
COTR100000000035SOL3200060AMPUL	COTRIMOXAZOLE (SULFAMETHOXAZOLE + TRIMETHOPRIM) 400 mg + 80 mg SOLUTION 5 mL AMPULE

COTR10000000035TAB490000000000	COTRIMOXAZOLE (SULFAMETHOXAZOLE + TRIMETHOPRIM) 400 mg + 80 mg TABLET
COTR100000000400SUS1400379BOTTL	COTRIMOXAZOLE (SULFAMETHOXAZOLE + TRIMETHOPRIM) 400 mg + 80 mg/5 mL SUSPENSION 30 mL BOTTLE
COTR100000000400SUS1400469BOTTL	COTRIMOXAZOLE (SULFAMETHOXAZOLE + TRIMETHOPRIM) 400 mg + 80 mg/5 mL SUSPENSION 60 mL BOTTLE
COTR100000001220SUS1400178BOTTL	COTRIMOXAZOLE (SULFAMETHOXAZOLE + TRIMETHOPRIM) 200 mg + 40 mg/5 mL SUSPENSION 100 mL BOTTLE
COTR100000001220SUS1400379BOTTL	COTRIMOXAZOLE (SULFAMETHOXAZOLE + TRIMETHOPRIM) 200 mg + 40 mg/5 mL SUSPENSION 30 mL BOTTLE
COTR100000001220SUS1400469BOTTL	COTRIMOXAZOLE (SULFAMETHOXAZOLE + TRIMETHOPRIM) 200 mg + 40 mg/5 mL SUSPENSION 60 mL BOTTLE
COTR100000001220SUS1401391BOTTL	COTRIMOXAZOLE (SULFAMETHOXAZOLE + TRIMETHOPRIM) 200 mg + 40 mg/5 mL SUSPENSION 70 mL BOTTLE
COTR100000001357TAB490000000000	COTRIMOXAZOLE (SULFAMETHOXAZOLE + TRIMETHOPRIM) 800 mg + 160 mg TABLET
CROTA00000000133CREA301380TUB01	CROTAMITON 10% CREAM 10 g TUBE
CROTA00000000133LOT1100201BOTTL	CROTAMITON 10% LOTION 120 mL BOTTLE
CROTA00000000133LOT1100469BOTTL	CROTAMITON 10% LOTION 60 mL BOTTLE
CYCLO00000000000POW2700004VIALX	CYCLOPHOSPHAMIDE POWDER 200 mg VIAL
CYCLO00000000000POW2700047VIALX	CYCLOPHOSPHAMIDE POWDER 500 mg VIAL
CYCLO00000000000POW2700139VIALX	CYCLOPHOSPHAMIDE POWDER 100 mg VIAL
CYCLO00000000000SOL3200241VIALX	CYCLOPHOSPHAMIDE SOLUTION 1 g VIAL
CYCLO00000600050TAB490000000000	CYCLOPHOSPHAMIDE (as ANHYDROUS) 50 mg TABLET
CYPRO00000100050TAB490000000000	CYPROTERONE (as ACETATE) 50 mg TABLET
CYPRO00000100142SOL3200388AMPUL	CYPROTERONE (as ACETATE) 100 mg/ml SOLUTION 3 mL AMPULE
CYTAR00000000000POW2700139VIALX	CYTARABINE POWDER 100 mg VIAL
CYTAR00000000142SOL320006000000	CYTARABINE 100 mg/ml SOLUTION 5 mL
CYTAR00000000142SOL3200195VIALX	CYTARABINE 100 mg/ml SOLUTION 10 mL VIAL
CYTAR00000000142SOL3200516VIALX	CYTARABINE 100 mg/ml SOLUTION 1 mL VIAL
CYTAR000000001199SOL3200060AMPUL	CYTARABINE 20 mg/mL SOLUTION 5 mL AMPULE
D504500000000000SOL3200044BOTTL	5% DEXTROSE IN 0.45% SODIUM CHLORIDE SOLUTION 500 mL BOTTLE
DACAR00000000000POW2700004VIALX	DACARBAZINE POWDER 200 mg VIAL
DACAR00000000000POW2700139VIALX	DACARBAZINE POWDER 100 mg VIAL
DACTI00000000000POW2700045VIALX	DACTINOMYCIN (ACTINOMYCIN D) POWDER 500 mcg VIAL
DALTE0005100694SOL3200000PFS01	DALTEPARIN (as SODIUM) 2500 IU/0.2 mL SOLUTION PRE-FILLED SYRINGE
DALTE0005100695SOL3200000PFS01	DALTEPARIN (as SODIUM) 5000 IU/0.2 mL SOLUTION PRE-FILLED SYRINGE
DALTE0005100696SOL3200516AMPUL	DALTEPARIN (as SODIUM) 10000 IU/mL SOLUTION 1 mL AMPULE
DANAZ00000000004CAPSU0000000000	DANAZOL 200 mg CAPSULE
DANAZ00000000139CAPSU0000000000	DANAZOL 100 mg CAPSULE
DANTR0005300021CAPSU0000000000	DANTROLENE (as SODIUM SALT) 25 mg CAPSULE
DANTR0005300050CAPSU0000000000	DANTROLENE (as SODIUM SALT) 50 mg CAPSULE
DANTR0005300697SOL3200000VIALX	DANTROLENE (as SODIUM SALT) 20 mg (with 3 g mannitol/vial) SOLUTION VIAL
DAPSO00000000139TAB490000000000	DAPSONE 100 mg TABLET
DAUNO00000000357SOL3200195VIALX	DAUNORUBICIN 2 mg/ml SOLUTION 10 mL VIAL
DAUNO00000000357SOL3201375VIALX	DAUNORUBICIN 2 mg/ml SOLUTION 25 mL VIAL
DEFER10000000047TAB490000000000	DEFERIPRONE 500 mg TABLET
DEFER0003700025SOL3200000VIALX	DEFEROXAMINE (as MESILATE) 2 g SOLUTION VIAL
DEFER0003700047SOL3200000VIALX	DEFEROXAMINE (as MESILATE) 500 mg SOLUTION VIAL
DESF100000000322SOL3200000BOTTL	DESFLURANE 240 mL SOLUTION BOTTLE

DESMO0000100148TAB490000000000	DESMOPRESSIN (as ACETATE) 100 mcg TABLET
DESMO0000100508TAB490000000000	DESMOPRESSIN (as ACETATE) 200 mcg TABLET
DESMO0000100699SOL3200516AMPUL	DESMOPRESSIN (as ACETATE) 15 mcg/mL SOLUTION 1 mL AMPULE
DEX4000000000701SOL3200044BOTTL	DEXTRAN, LOW MOLECULAR WEIGHT (DEXTRAN 40) 10% dextran 40 in 0.9% NaCl SOLUTION 500 mL BOTTLE
DEX4000000000702SOL3200044BOTTL	DEXTRAN, LOW MOLECULAR WEIGHT (DEXTRAN 40) 10% dextran 40 in 5% dextrose SOLUTION 500 mL BOTTLE
DEX7000000000464SOL3200044BOTTL	DEXTRAN, HIGH MOLECULAR WEIGHT (DEXTRAN 70) 6%, Dextran 70 in 0.9% NaCl SOLUTION 500 mL BOTTLE
DEX7000000000465SOL3200044BOTTL	DEXTRAN, HIGH MOLECULAR WEIGHT (DEXTRAN 70) 6%, Dextran 70 in 5% Dextrose SOLUTION 500 mL BOTTLE
DEXAM0000000033TAB490000000000	DEXAMETHASONE 3 mg TABLET
DEXAM0000000041TAB490000000000	DEXAMETHASONE 4 mg TABLET
DEXAM0000000045TAB490000000000	DEXAMETHASONE 500 mcg TABLET
DEXAM0000000082EDS0100060BOTTL	DEXAMETHASONE 0.10% EYE DROPS SUSPENSION 5 mL BOTTLE
DEXAM00000000700TAB490000000000	DEXAMETHASONE 750 mcg TABLET
DEXAM0005200414SOL3200361AMPUL	DEXAMETHASONE (as SODIUM PHOSPHATE) 4 mg/mL SOLUTION 2 mL AMPULE
DEXAM0005200414SOL3200361VIALX	DEXAMETHASONE (as SODIUM PHOSPHATE) 4 mg/mL SOLUTION 2 mL VIAL
DEXAM0005200414SOL3200516AMPUL	DEXAMETHASONE (as SODIUM PHOSPHATE) 4 mg/mL SOLUTION 1 mL AMPULE
DEXAM0005200414SOL3200516VIALX	DEXAMETHASONE (as SODIUM PHOSPHATE) 4 mg/mL SOLUTION 1 mL VIAL
DEXAM0005200459SOL3200516AMPUL	DEXAMETHASONE (as SODIUM PHOSPHATE) 5 mg/mL SOLUTION 1 mL AMPULE
DEXTR0000200069TAB490000000000	DEXTROMETHORPHAN (as HYDROCHLORIDE) 10 mg TABLET
DEXTR0000200184SYRUP00469BOTTL	DEXTROMETHORPHAN (as HYDROCHLORIDE) 10 mg/5 mL SYRUP 60 mL BOTTLE
DEXTR0000200454SYRUP00201BOTTL	DEXTROMETHORPHAN (as HYDROCHLORIDE) 5 mg/5 mL SYRUP 120 mL BOTTLE
DEXTR0000200454SYRUP00379BOTTL	DEXTROMETHORPHAN (as HYDROCHLORIDE) 5 mg/5 mL SYRUP 30 mL BOTTLE
DEXTR0000200454SYRUP00469BOTTL	DEXTROMETHORPHAN (as HYDROCHLORIDE) 5 mg/5 mL SYRUP 60 mL BOTTLE
DEXTR0000200706SYRUP00201BOTTL	DEXTROMETHORPHAN (as HYDROCHLORIDE) 15 mg/5 mL SYRUP 120 mL BOTTLE
DEXTR0000200706SYRUP00379BOTTL	DEXTROMETHORPHAN (as HYDROCHLORIDE) 15 mg/5 mL SYRUP 30 mL BOTTLE
DEXTR0000200706SYRUP00469BOTTL	DEXTROMETHORPHAN (as HYDROCHLORIDE) 15 mg/5 mL SYRUP 60 mL BOTTLE
DIAZE0000000026TAB490000000000	DIAZEPAM 2 mg TABLET
DIAZE0000000056TAB490000000000	DIAZEPAM 5 mg TABLET
DIAZE0000000069TAB490000000000	DIAZEPAM 10 mg TABLET
DIAZE00000000459SOL3200361AMPUL	DIAZEPAM 5 mg/mL SOLUTION 2 mL AMPULE
DIAZE00000000708RTU0000000TUB01	DIAZEPAM 5 mg/2.5 mL RECTAL TUBE
DIAZE00000000709RTU0000000TUB01	DIAZEPAM 10 mg/2.5 mL RECTAL TUBE
DICLO0005100021CAPSU0000000000	DICLOFENAC (as SODIUM) 25 mg CAPSULE
DICLO0005100021TAB490000000000	DICLOFENAC (as SODIUM) 25 mg TABLET
DICLO0005100065CAPSU0000000000	DICLOFENAC (as SODIUM) 75 mg CAPSULE
DICLO0005100065TAB490000000000	DICLOFENAC (as SODIUM) 75 mg TABLET
DICLO0005100139CAPSU0000000000	DICLOFENAC (as SODIUM) 100 mg CAPSULE
DICLO0005100139TAB490000000000	DICLOFENAC (as SODIUM) 100 mg TABLET
DICLO0005300050CAPSU0000000000	DICLOFENAC (as SODIUM SALT) 50 mg CAPSULE
DICLO0005300050TAB490000000000	DICLOFENAC (as SODIUM SALT) 50 mg TABLET
DICLO0005300713SOL3200361AMPUL	DICLOFENAC (as SODIUM SALT) 37.5 mg/ml SOLUTION 2 mL AMPULE

DICLO0005300713SOL3200361PFS01	DICLOFENAC (as SODIUM SALT) 37.5 mg/ml SOLUTION 2 mL PRE-FILLED SYRINGE
DICLO0005300985SOL3200195VIALX	DICLOFENAC (as SODIUM SALT) 25 mg/mL SOLUTION 10 mL VIAL
DICLO0005300985SOL3200361AMPUL	DICLOFENAC (as SODIUM SALT) 25 mg/mL SOLUTION 2 mL AMPULE
DICLO0005300985SOL3200388AMPUL	DICLOFENAC (as SODIUM SALT) 25 mg/mL SOLUTION 3 mL AMPULE
DICY10000200069TAB4900000000000	DICYCLOVERINE (DICYCLOMINE) (as HYDROCHLORIDE) 10 mg TABLET
DICY10000200184SYRUP00201BOTTL	DICYCLOVERINE (DICYCLOMINE) (as HYDROCHLORIDE) 10 mg/5 mL SYRUP 120 mL BOTTLE
DICY10000200184SYRUP00379BOTTL	DICYCLOVERINE (DICYCLOMINE) (as HYDROCHLORIDE) 10 mg/5 mL SYRUP 30 mL BOTTLE
DICY10000200184SYRUP00469BOTTL	DICYCLOVERINE (DICYCLOMINE) (as HYDROCHLORIDE) 10 mg/5 mL SYRUP 60 mL BOTTLE
DIDAN0000000004TAB4900000000000	DIDANOSINE 200 mg TABLET
DIDAN0000000015CAPSU00000000000	DIDANOSINE 250 mg CAPSULE
DIETH0000000050TAB4900000000000	DIETHYLCARBAMAZINE 50 mg TABLET
DIETH00000000139TAB4900000000000	DIETHYLCARBAMAZINE 100 mg TABLET
DIGOX0000000013TAB4900000000000	DIGOXIN 250 mcg TABLET
DIGOX00000000615SOL3200469ELI01	DIGOXIN 50 mcg/mL SOLUTION 60 mL ELIXIR
DIGOX00000000771SOL3200361AMPUL	DIGOXIN 250 mcg/mL SOLUTION 2 mL AMPULE
DILOX0002500047TAB4900000000000	DILOXANIDE (as FUROATE) 500 mg TABLET
DILOX0002500209SYRUP00379BOTTL	DILOXANIDE (as FUROATE) 125 mg/5 mL SYRUP 30 mL BOTTLE
DILOX0002500209SYRUP00469BOTTL	DILOXANIDE (as FUROATE) 125 mg/5 mL SYRUP 60 mL BOTTLE
DILTIO0000200031TAB4900000000000	DILTIAZEM (as HYDROCHLORIDE) 30 mg TABLET
DILTIO0000200062CAPSU00000000000	DILTIAZEM (as HYDROCHLORIDE) 60 mg CAPSULE
DILTIO0000200062TAB4900000000000	DILTIAZEM (as HYDROCHLORIDE) 60 mg TABLET
DILTIO0000200199CAPSU00000000000	DILTIAZEM (as HYDROCHLORIDE) 120 mg CAPSULE
DILTIO0000200199TAB4900000000000	DILTIAZEM (as HYDROCHLORIDE) 120 mg TABLET
DILTIO0000200236CAPSU00000000000	DILTIAZEM (as HYDROCHLORIDE) 180 mg CAPSULE
DILTIO0000200236TAB4900000000000	DILTIAZEM (as HYDROCHLORIDE) 180 mg TABLET
DILTIO0000200321CAPSU00000000000	DILTIAZEM (as HYDROCHLORIDE) 240 mg CAPSULE
DILTIO0000200503CAPSU00000000000	DILTIAZEM (as HYDROCHLORIDE) 90 mg CAPSULE
DILTIO0000200503TAB4900000000000	DILTIAZEM (as HYDROCHLORIDE) 90 mg TABLET
DIME10000000142SOL3200195VIALX	DIMERCAPTOPROPANE SULFONATE (DMPS) 100 mg/ml SOLUTION 10 mL VIAL
DIME10000000142SOL3200516AMPUL	DIMERCAPTOPROPANE SULFONATE (DMPS) 100 mg/ml SOLUTION 1 mL AMPULE
DIMEG00000000719SOL3200011VIALX	DIMEGLUMINE GADOPENTETATE 469 mg/mL SOLUTION 20 mL VIAL
DIMEG00000000719SOL3200060VIALX	DIMEGLUMINE GADOPENTETATE 469 mg/mL SOLUTION 5 mL VIAL
DIMEG00000000719SOL3200195VIALX	DIMEGLUMINE GADOPENTETATE 469 mg/mL SOLUTION 10 mL VIAL
DIMEG00000000719SOL3200231VIALX	DIMEGLUMINE GADOPENTETATE 469 mg/mL SOLUTION 15 mL VIAL
DIMER00000000851SOL3200361AMPUL	DIMERCAPROL 50 mg/mL SOLUTION 2 mL AMPULE
DIPH100000000732SOL3200000AMPUL	DIPHThERIA-TETANUS TOXOIDS (DT) 30 IU diphtheria toxoid + 40 IU tetanus toxoid per 0.5 mL SOLUTION
DIPH200000000733SOL3200000AMPUL	DIPHThERIA-TETANUS TOXOIDS (TD) 2 IU diphtheria toxoid + 20 IU tetanus toxoid SOLUTION AMPULE
DIPH300000000094SOL3200000PFS01	DIPHThERIA-TETANUS TOXOIDS AND ACELLULAR PERTUSSIS VACCINE 0.5 mL SOLUTION PRE-FILLED SYRINGE
DIPH300000000094SOL3200000VIALX	DIPHThERIA-TETANUS TOXOIDS AND ACELLULAR PERTUSSIS VACCINE 0.5 mL SOLUTION VIAL
DIPH400000000094SOL3200000AMPUL	DIPHThERIA-TETANUS TOXOIDS AND PERTUSSIS VACCINE (DTP) 0.5 mL SOLUTION AMPULE
DIPH400000000094SOL3200000PFS01	DIPHThERIA-TETANUS TOXOIDS AND PERTUSSIS VACCINE (DTP) 0.5 mL SOLUTION PRE-FILLED SYRINGE
DIPH400000000094SOL3200060VIALX	DIPHThERIA-TETANUS TOXOIDS AND PERTUSSIS VACCINE (DTP) 0.5 mL SOLUTION 5 mL VIAL

DIPH40000000094SOL3200195VIALX	DIPHTHERIA-TETANUS TOXOIDS AND PERTUSSIS VACCINE (DTP) 0.5 mL SOLUTION 10 mL VIAL
DIPH40000000094SOL3201392VIALX	DIPHTHERIA-TETANUS TOXOIDS AND PERTUSSIS VACCINE (DTP) 0.5 mL SOLUTION 7.5 mL VIAL
DIPH500000000725SYRUP000379BOTTLE	DIPHENHYDRAMINE 12.5 mg/5 mL SYRUP 30 mL BOTTLE
DIPH500000200021CAPSU0000000000	DIPHENHYDRAMINE (as HYDROCHLORIDE) 25 mg CAPSULE
DIPH500000200050CAPSU0000000000	DIPHENHYDRAMINE (as HYDROCHLORIDE) 50 mg CAPSULE
DIPH500000200197SYRUP000469BOTTLE	DIPHENHYDRAMINE (as HYDROCHLORIDE) 12.5 mg/ 5 mL SYRUP 60 mL BOTTLE
DIPH500000200725SYRUP000201BOTTLE	DIPHENHYDRAMINE (as HYDROCHLORIDE) 12.5 mg/5 mL SYRUP 120 mL BOTTLE
DIPH500000200851SOL3200195VIALX	DIPHENHYDRAMINE (as HYDROCHLORIDE) 50 mg/mL SOLUTION 10 mL VIAL
DIPH500000200851SOL3200516AMPUL	DIPHENHYDRAMINE (as HYDROCHLORIDE) 50 mg/mL SOLUTION 1 mL AMPULE
DIPYR00000000021TAB490000000000	DIPYRIDAMOLE 25 mg TABLET
DIPYR00000000050TAB490000000000	DIPYRIDAMOLE 50 mg TABLET
DIPYR00000000065TAB490000000000	DIPYRIDAMOLE 75 mg TABLET
DITHR00000000737OINTX01384TUB01	DITHRANOL 0.1% - 2% OINTMENT 50 g TUBE
DOBUT00000000357SOL3201400DPB01	DOBUTAMINE 2 mg/ml SOLUTION 250 mL D5W PRE-MIXED BOTTLE
DOBUT00000000414SOL3201400DPB01	DOBUTAMINE 4 mg/mL SOLUTION 250 mL D5W PRE-MIXED BOTTLE
DOBUT00000000851SOL3200060AMPUL	DOBUTAMINE 50 mg/mL SOLUTION 5 mL AMPULE
DOBUT00000200198SOL3200011VIALX	DOBUTAMINE (as HYDROCHLORIDE) 12.5 mg/mL SOLUTION 20 mL VIAL
DOBUT00000200261SUS1400020DPB01	DOBUTAMINE (as HYDROCHLORIDE) 1 mg/mL SUSPENSION 250 mL D5W PRE-MIXED BOTTLE
DOCET00000000039SOL3200361VIALX	DOCETAXEL 40 mg/mL SOLUTION 2 mL VIAL
DOCET00000000194SOL3200361VIALX	DOCETAXEL 10 mg/mL SOLUTION 2 mL VIAL
DOCET00000000194SOL3201401VIALX	DOCETAXEL 10 mg/mL SOLUTION 8 mL VIAL
DOCET00000000313SOL3200094VIALX	DOCETAXEL 20 mg/0.5 mL SOLUTION 0.5 mL VIAL
DOCET00000001199SOL3201402VIALX	DOCETAXEL 20 mg/mL SOLUTION 7 mL VIAL
DOCET00000001200SOL3201402VIALX	DOCETAXEL 80 mg/4 mL SOLUTION 7 mL VIAL
DOMPE00000000069TAB490000000000	DOMPERIDONE 10 mg TABLET
DOMPE00000000261SUS1400379BOTTLE	DOMPERIDONE 1 mg/mL SUSPENSION 30 mL BOTTLE
DOMPE00000000261SUS1400469BOTTLE	DOMPERIDONE 1 mg/mL SUSPENSION 60 mL BOTTLE
DONEP00000000056TAB490000000000	DONEPEZIL HYDROCHLORIDE 5 mg TABLET
DONEP00000000069TAB490000000000	DONEPEZIL HYDROCHLORIDE 10 mg TABLET
DOPAM00000000123SUS1400020DPB01	DOPAMINE 1.6 mg/mL SUSPENSION 250 mL D5W PRE-MIXED BOTTLE
DOPAM00000000494SOL3200020DPB01	DOPAMINE 800 mcg/mL SOLUTION 250 mL D5W PRE-MIXED BOTTLE
DOPAM00000200039SOL3200060AMPUL	DOPAMINE (as HYDROCHLORIDE) 40 mg/mL SOLUTION 5 mL AMPULE
DOPAM00000200039SOL3200060VIALX	DOPAMINE (as HYDROCHLORIDE) 40 mg/mL SOLUTION 5 mL VIAL
DOPAM00000200741SOL3200060VIALX	DOPAMINE (as HYDROCHLORIDE) 80 mg/mL SOLUTION 5 mL VIAL
DOPAM00000200742SOL3200020DPB01	DOPAMINE (as HYDROCHLORIDE) 3.2 mg/mL SOLUTION 250 mL D5W PRE-MIXED BOTTLE
DOPAM00000201199SOL3200060AMPUL	DOPAMINE (as HYDROCHLORIDE) 20 mg/mL SOLUTION 5 mL AMPULE
DOPAM00000201199SOL3200195AMPUL	DOPAMINE (as HYDROCHLORIDE) 20 mg/mL SOLUTION 10 mL AMPULE
DORZO00000201314EDS0100060VIALX	DORZOLAMIDE (as HYDROCHLORIDE) 2% EYE DROPS SUSPENSION 5 mL VIAL
DOXOR00000200000POW2700050VIALX	DOXORUBICIN (as HYDROCHLORIDE) POWDER 50 mg VIAL
DOXOR00000200000POW2700069VIALX	DOXORUBICIN (as HYDROCHLORIDE) POWDER 10 mg VIAL
DOXOR00000200357SOL3200060VIALX	DOXORUBICIN (as HYDROCHLORIDE) 2 mg/ml SOLUTION 5 mL VIAL
DOXOR00000200357SOL3201375VIALX	DOXORUBICIN (as HYDROCHLORIDE) 2 mg/ml SOLUTION 25 mL VIAL

DOXYC0002700050CAPSU0000000000	DOXYCYCLINE (as HYCLATE) 50 mg CAPSULE
DOXYC0002700139CAPSU0000000000	DOXYCYCLINE (as HYCLATE) 100 mg CAPSULE
DPHo100000000094SOL3200000VIALX	DIPHThERIA, TETANUS, PERTUSSIS, HEPATITIS B RECOMBINANT AND HAEMOPHILUS INFLUENZA TYPE B (HIB) COMBI
DTIB00000000000SOL3200094PFS01	DTAP + HIB SOLUTION 0.5 mL PRE-FILLED SYRINGE
DTPHB00000000000SOL3200094VIALX	DTP + HEPATITIS B VACCINE (RECOMBINANT) SOLUTION 0.5 mL VIAL
DTPHI00000000744SOL3200000VIALX	DTP + HIB 0.5 mL DTP diluent + 10 mcg Hib SOLUTION VIAL
DTPI000000000094SOL3200000PFS01	DTP + INACTIVATED POLIO VACCINE 0.5 mL SOLUTION PRE-FILLED SYRINGE
DTPI000000000745SOL3200000VIALX	DTP + INACTIVATED POLIO VACCINE 0.5 mL monodose SOLUTION VIAL
DTPIH00000000094SOL3200000PFS01	DTP + IPV + HIB 0.5 mL SOLUTION PRE-FILLED SYRINGE
DUTAS00000000047CAPSU0000000000	DUTASTERIDE 500 mg CAPSULE
DWP010000000000SOL3201347CRRT1	DOUSOL WITHOUT POTASSIUM SOLUTION 5 L CRRT
DYDRO00000000069TAB490000000000	DYDROGESTERONE 10 mg TABLET
DZAC100000000000SOL3201346DRU01	DIALYSATE (ACETATE BASED) SOLUTION 5 gallon (approx. 20 L) DRUM
DZAC100000000000SOL3201347DRU01	DIALYSATE (ACETATE BASED) SOLUTION 5 L DRUM
DZAC100000001348SOL3201348DRU01	DIALYSATE (ACETATE BASED) 10 L SOLUTION 10 L DRUM
DZBC100000000000SOL3201346DRU01	DIALYSATE (BICARBONATE BASED) SOLUTION 5 gallon (approx. 20 L) DRUM
DZBC100000000000SOL3201347DRU01	DIALYSATE (BICARBONATE BASED) SOLUTION 5 L DRUM
DZBC100000000000SOL3201348DRU01	DIALYSATE (BICARBONATE BASED) SOLUTION 10 L DRUM
DZLC100000000000SOL3201346DRU01	DIALYSATE (LOW CALCIUM) SOLUTION 5 gallon (approx. 20 L) DRUM
DZLC100000000000SOL3201347DRU01	DIALYSATE (LOW CALCIUM) SOLUTION 5 L DRUM
DZLC100000000000SOL3201348DRU01	DIALYSATE (LOW CALCIUM) SOLUTION 10 L DRUM
ECONA00000000000LOT100195BOTTL	ECONAZOLE LOTION 10 mL BOTTLE
ECONA00000000000SHAMX00178BOTTL	ECONAZOLE SHAMPOO 100 mL BOTTLE
ECONA00000000000SHAMX00195SAC01	ECONAZOLE SHAMPOO 10 mL SACHET
ECONA00000000000SHAMX00469BOTTL	ECONAZOLE SHAMPOO 60 mL BOTTLE
ECONA00000000000SHAMX00842SAC01	ECONAZOLE SHAMPOO 6 mL SACHET
ECONA00000000028OVULE0000000000	ECONAZOLE 300 mg OVULE
ECONA00000000047OVULE0000000000	ECONAZOLE 500 mg OVULE
ECONA00000000139OVULE0000000000	ECONAZOLE 100 mg OVULE
ECONA00000000194SOL3200231BOTTL	ECONAZOLE 10 mg/mL SOLUTION 15 mL BOTTLE
ECONA00000000194SOL3201375BOTTL	ECONAZOLE 10 mg/mL SOLUTION 25 mL BOTTLE
ECONA00000001062CREA300227TUB01	ECONAZOLE 1% CREAM 15 g TUBE
ECONA00000001062CREA301393TUB01	ECONAZOLE 1% CREAM 3.5 g TUBE
ECONA00000001062CREA301398TUB01	ECONAZOLE 1% CREAM 450 g TUBE
ECONA00000001314CREA301393TUB01	ECONAZOLE 2% CREAM 3.5 g TUBE
ECONA00000001314CREA301403TUB01	ECONAZOLE 2% CREAM 450 g TUBE
EDROP0001300194SOL3200516AMPUL	EDROPHONIUM (as CHLORIDE) 10 mg/mL SOLUTION 1 mL AMPULE
EFA0000000000004CAPSU0000000000	EFAVIRENZ 200 mg CAPSULE
EFA0000000000004TAB490000000000	EFAVIRENZ 200 mg TABLET
EFA0000000000050CAPSU0000000000	EFAVIRENZ 50 mg CAPSULE
EFA0000000000050TAB490000000000	EFAVIRENZ 50 mg TABLET
EFA0000000000061TAB490000000000	EFAVIRENZ 600 mg TABLET
EFA0000000000139CAPSU0000000000	EFAVIRENZ 100 mg CAPSULE
EFA0000000000139TAB490000000000	EFAVIRENZ 100 mg TABLET
ENAL100000000311TAB490000000000	ENALAPRIL + HYDROCHLOROTHIAZIDE 20 mg + 12.5 mg TABLET

ENALA0003400007TAB4900000000000	ENALAPRIL (as MALEATE) 20 mg TABLET
ENALA0003400056TAB4900000000000	ENALAPRIL (as MALEATE) 5 mg TABLET
ENALA0003400069TAB4900000000000	ENALAPRIL (as MALEATE) 10 mg TABLET
ENDS20000000000000000001334CAN01	ENTERAL NUTRITION - DISEASE-SPECIFIC (DIABETES) 400 g CAN
ENDS30000000000000000001353CAN01	ENTERAL NUTRITION - DISEASE SPECIFIC (CANCER-INDUCED WEIGHT LOSS) 380 g CAN
ENOXA0005300142SOL3200516VIALX	ENOXAPARIN (as SODIUM SALT) 100 mg/ml SOLUTION 1 mL VIAL
ENOXA0005300142SOL3201396PFS01	ENOXAPARIN (as SODIUM SALT) 100 mg/ml SOLUTION 0.4 mL PRE-FILLED SYRINGE
ENOXA0005300142SOL3201404PFS01	ENOXAPARIN (as SODIUM SALT) 100 mg/ml SOLUTION 0.2 mL PRE-FILLED SYRINGE
ENOXA0005300142SOL3201405PFS01	ENOXAPARIN (as SODIUM SALT) 100 mg/ml SOLUTION 0.6 mL PRE-FILLED SYRINGE
ENOXA0005300142SOL3201406PFS01	ENOXAPARIN (as SODIUM SALT) 100 mg/ml SOLUTION 0.8 mL PRE-FILLED SYRINGE
ENPo10000000000000000001334CAN01	ENTERAL NUTRITION - POLYMERIC 400 g CAN
ENTE30000000000000000001334CAN01	ENTERAL NUTRITION - FIBER-CONTAINING 400 g CAN
ENTE50000000000000000001064CAN01	ENTERAL NUTRITION - PEDIATRIC POLYMERIC 100 g CAN
EPHED0006000851SOL3200516AMPUL	EPHEDRINE (as SULFATE) 50 mg/mL SOLUTION 1 mL AMPULE
EPINE0000001235SOL3200000AI001	EPINEPHRINE 0.3 mg SOLUTION AUTO-INJECTOR
EPINE0000200261SOL3200516AMPUL	EPINEPHRINE (as HYDROCHLORIDE) 1 mg/mL SOLUTION 1 mL AMPULE
EPIRU0000200050SOL3200000VIALX	EPIRUBICIN (as HYDROCHLORIDE) 50 mg SOLUTION VIAL
EPIRU0000200069SOL3200000VIALX	EPIRUBICIN (as HYDROCHLORIDE) 10 mg SOLUTION VIAL
EPOE10000000295SOL3200000PSN01	EPOETIN BETA (RECOMBINANT ERYTHROPOIETIN) 2000 IU/ 0.3 mL SOLUTION PRE-FILLED SYRINGE WITH NEEDLE
EPOE10000000395SOL3200000PSN01	EPOETIN BETA (RECOMBINANT ERYTHROPOIETIN) 4000 IU/0.3 mL SOLUTION PRE-FILLED SYRINGE WITH NEEDLE
EPOE10000000424SOL3200000PSN01	EPOETIN BETA (RECOMBINANT ERYTHROPOIETIN) 5000IU/0.3ml SOLUTION PRE-FILLED SYRINGE WITH NEEDLE
EPOE10000000749SOL3200000PFS01	EPOETIN BETA (RECOMBINANT ERYTHROPOIETIN) 10000 IU/0.6 mL SOLUTION PRE-FILLED SYRINGE
EPOE10000000750SOL3200516PFS01	EPOETIN BETA (RECOMBINANT ERYTHROPOIETIN) 20000 IU/mL SOLUTION 1 mL PRE-FILLED SYRINGE
EPOE10000000751SOL3200516PFS01	EPOETIN BETA (RECOMBINANT ERYTHROPOIETIN) 30000 IU/0.6 mL SOLUTION 1 mL PRE-FILLED SYRINGE
EPOET0000000296SOL3200000PFS01	EPOETIN ALFA (RECOMBINANT HUMAN ERYTHROPOIETIN) 2000 IU/0.5 mL SOLUTION PRE-FILLED SYRINGE
EPOET0000000296SOL3200000PFS02	EPOETIN ALFA (RECOMBINANT HUMAN ERYTHROPOIETIN) 2000 IU/0.5 mL SOLUTION PRE-FILLED SYRINGE (HAS FRE
EPOET0000000297SOL3200516PFG01	EPOETIN ALFA (RECOMBINANT HUMAN ERYTHROPOIETIN) 2000 IU/mL SOLUTION 1 mL PRE-FILLED GLASS SYRINGE
EPOET0000000297SOL3200516VIALX	EPOETIN ALFA (RECOMBINANT HUMAN ERYTHROPOIETIN) 2000 IU/mL SOLUTION 1 mL VIAL
EPOET0000000396SOL3200000PFS01	EPOETIN ALFA (RECOMBINANT HUMAN ERYTHROPOIETIN) 4000 IU/0.4 mL SOLUTION PRE-FILLED SYRINGE
EPOET0000000396SOL3200000PFS02	EPOETIN ALFA (RECOMBINANT HUMAN ERYTHROPOIETIN) 4000 IU/0.4 mL SOLUTION PRE-FILLED SYRINGE (HAS FRE
EPOET0000000549SOL3200516PFG01	EPOETIN ALFA (RECOMBINANT HUMAN ERYTHROPOIETIN) 4000 IU/mL SOLUTION 1 mL PRE-FILLED GLASS SYRINGE
EPOET0000000549SOL3200516VIALX	EPOETIN ALFA (RECOMBINANT HUMAN ERYTHROPOIETIN) 4000 IU/mL SOLUTION 1 mL VIAL
EPOET0000000747SOL3200000PFS01	EPOETIN ALFA (RECOMBINANT HUMAN ERYTHROPOIETIN) 10,000 IU/mL SOLUTION PRE-FILLED SYRINGE
EPOET0000000747SOL3200000PFS02	EPOETIN ALFA (RECOMBINANT HUMAN ERYTHROPOIETIN) 10,000 IU/mL SOLUTION PRE-FILLED SYRINGE (HAS FREE)
EPOET0000000748SOL3200000PFS02	EPOETIN ALFA (RECOMBINANT HUMAN ERYTHROPOIETIN) 40,000 IU/mL SOLUTION PRE-FILLED SYRINGE (HAS FREE)
EPR010003700061TAB4900000000000	EPROSARTAN (as MESILATE) 600 mg TABLET
ERGOC0000000752CAPSU00000000000	ERGOCALCIFEROL (CALCIFEROL, VITAMIN D2) 1.25 mg (50,000 IU) CAPSULE

ERGOC0000000752TAB490000000000	ERGOCALCIFEROL (CALCIFEROL, VITAMIN D2) 1.25 mg (50,000 IU) TABLET
ERGOC0000000753SOL3200469BOTTL	ERGOCALCIFEROL (CALCIFEROL, VITAMIN D2) 250 mcg/mL (10,000 IU/mL) SOLUTION 60 mL BOTTLE
ERGOT0005600252TAB490000000000	ERGOTAMINE (as TARTRATE) 1 mg TABLET
ERLOT0000000215TAB490000000000	ERLOTINIB 150 mg TABLET
ERT020005100241SOL3200000VIALX	ERTAPENEM (as SODIUM) 1 g SOLUTION VIAL
ERYTH0000000015CAPSU0000000000	ERYTHROMYCIN 250 mg CAPSULE
ERYTH0000000015TAB490000000000	ERYTHROMYCIN 250 mg TABLET
ERYTH0000000649EYO0101245TUB01	ERYTHROMYCIN 0.50% EYE OINTMENT 5 g TUBE
ERYTH0000000649EYO0101393TUB01	ERYTHROMYCIN 0.50% EYE OINTMENT 3.5 g TUBE
ERYTH0002200006SUS1400469BOTTL	ERYTHROMYCIN (as ETHYL SUCCINATE) 200 mg/5 mL SUSPENSION 60 mL BOTTLE
ERYTH0002200039SUS1400379BOTTL	ERYTHROMYCIN (as ETHYL SUCCINATE) 40 mg/mL SUSPENSION 30 mL BOTTLE
ERYTH0002200757SUS1400055BOTTL	ERYTHROMYCIN (as ETHYL SUCCINATE) 400 mg/5 mL SUSPENSION 50 mL BOTTLE
ERYTH0002200757SUS1400469BOTTL	ERYTHROMYCIN (as ETHYL SUCCINATE) 400 mg/5 mL SUSPENSION 60 mL BOTTLE
ERYTH0005500047TAB490000000000	ERYTHROMYCIN (as STEARATE) 500 mg TABLET
ESC010004300069TAB490000000000	ESCITALOPRAM (as OXALATE) 10 mg TABLET
ESMOL0000200142SOL3200195VIALX	ESMOLOL (as HYDROCHLORIDE) 100 mg/ml SOLUTION 10 mL VIAL
ESMOL0000200194SOL3200195VIALX	ESMOLOL (as HYDROCHLORIDE) 10 mg/mL SOLUTION 10 mL VIAL
ESOM10000000007CAPSU0000000000	OMEPRAZOLE 20 mg CAPSULE
ESOM10000000036CAPSU0000000000	OMEPRAZOLE 40 mg CAPSULE
ESOM10000000069CAPSU0000000000	OMEPRAZOLE 10 mg CAPSULE
ESOM10000001025LPS0100000AMPUL	OMEPRAZOLE 40 mg powder vial + 10 mL solvent LYOPHILIZED POWDER + SOLVENT AMPULE
ETHAM0000200004TAB490000000000	ETHAMBUTOL (as HYDROCHLORIDE) 200 mg TABLET
ETHAM0000200034TAB490000000000	ETHAMBUTOL (as HYDROCHLORIDE) 400 mg TABLET
ETHDE0000000375TAB490000000000	ETHINYLESTRADIOL + DESOGESTREL 30 mcg + 150 mcg TABLET
ETHI20000000758TAB490000000000	ETHINYLESTRADIOL + NORETHISTERONE 35 mcg + 400 mcg TABLET
ETHIN0000000375TAB490000000000	ETHINYLESTRADIOL + LEVONORGESTREL 30 mcg + 150 mcg TABLET
ETHNG0000000759TAB490000000000	ETHINYLESTRADIOL + NORGESTREL 30 mcg + 300 mcg TABLET
ETON0000000134000000000000IMP01	ETONOGESTREL 68 mg IMPLANT
ETOPO0000000021CAPSU0000000000	ETOPOSIDE 25 mg CAPSULE
ETOPO0000000050CAPSU0000000000	ETOPOSIDE 50 mg CAPSULE
ETOPO00000000139CAPSU0000000000	ETOPOSIDE 100 mg CAPSULE
ETOPO00000000139SOL3200000VIALX	ETOPOSIDE 100 mg SOLUTION VIAL
ETOPO000000001620000000000AMPUL	ETOPOSIDE 100 mg/5 mL AMPULE
ETOPO00000001199SOL3200060AMPUL	ETOPOSIDE 20 mg/mL SOLUTION 5 mL AMPULE
ETOPO00000001199SOL3200060VIALX	ETOPOSIDE 20 mg/mL SOLUTION 5 mL VIAL
ETOPO00000001199SOL3200195VIALX	ETOPOSIDE 20 mg/mL SOLUTION 10 mL VIAL
ETOPO00000001199SOL3201124VIALX	ETOPOSIDE 20 mg/mL SOLUTION 2.5 mL VIAL
EVER00000000013TAB490000000000	EVEROLIMUS 250 mcg TABLET
EVER00000000045TAB490000000000	EVEROLIMUS 500 mcg TABLET
EVER00000000700TAB490000000000	EVEROLIMUS 750 mcg TABLET
EZET100000001441TAB490000000000	EZETIMIBE + SIMVASTATIN 10 mg + 10 mg TABLET
EZET100000001442TAB490000000000	EZETIMIBE + SIMVASTATIN 10 mg + 20 mg TABLET
EZET100000001443TAB490000000000	EZETIMIBE + SIMVASTATIN 10 mg + 40 mg TABLET
EZETI00000000069TAB490000000000	EZETIMIBE 10 mg TABLET

FACT10000000764SOL3200000VIALX	FACTOR VIII CONCENTRATE 100 IU/g + diluent SOLUTION VIAL
FACTO0000000761SOL3200060VIALX	FACTOR IX COMPLEX (COAGULATION FACTORS II, VII, IX, X) CONCENTRATE 100 IU/mL SOLUTION 5 mL VIAL
FACTO0000000761SOL3200195VIALX	FACTOR IX COMPLEX (COAGULATION FACTORS II, VII, IX, X) CONCENTRATE 100 IU/mL SOLUTION 10 mL VIAL
FACTO0000000763SOL3200000VIALX	FACTOR IX COMPLEX (COAGULATION FACTORS II, VII, IX, X) CONCENTRATE 500 IU + diluent SOLUTION VIAL
FAMCI0000000015TAB490000000000	FAMCICLOVIR 250 mg TABLET
FAMCI00000000205TAB490000000000	FAMCICLOVIR 125 mg TABLET
FAMOT0000000007LYP0100000VIALX	FAMOTIDINE 20 mg LYOPHILIZED POWDER VIAL
FAMOT0000000007TAB490000000000	FAMOTIDINE 20 mg TABLET
FAMOT0000000036TAB490000000000	FAMOTIDINE 40 mg TABLET
FAMOT0000000069TAB490000000000	FAMOTIDINE 10 mg TABLET
FAMOT0000000194SOL3200361AMPUL	FAMOTIDINE 10 mg/mL SOLUTION 2 mL AMPULE
FAMOT0000000194SOL3200361VIALX	FAMOTIDINE 10 mg/mL SOLUTION 2 mL VIAL
FELOD0000000003TAB490000000000	FELODIPINE 2.5 mg TABLET
FELOD0000000056TAB490000000000	FELODIPINE 5 mg TABLET
FELOD0000000069TAB490000000000	FELODIPINE 10 mg TABLET
FEME10000001454TAB490000000000	FELODIPINE + METOPROLOL 5 mg + 50 mg TABLET
FENOF0000000004CAPSU0000000000	FENOFIBRATE 200 mg CAPSULE
FENOF0000000139CAPSU0000000000	FENOFIBRATE 100 mg CAPSULE
FENOF0000000233TAB490000000000	FENOFIBRATE 160 mg TABLET
FENOF00000001216TAB490000000000	FENOFIBRATE 145 mg TABLET
FENTA0001500443SOL3200195AMPUL	FENTANYL (as CITRATE) 50mcg/mL SOLUTION 10 mL AMPULE
FENTA0001500443SOL3200361AMPUL	FENTANYL (as CITRATE) 50mcg/mL SOLUTION 2 mL AMPULE
FENTA0001500615SOL3200011AMPUL	FENTANYL (as CITRATE) 50 mcg/mL SOLUTION 20 mL AMPULE
FERFU00000001280TAB490000000000	FERROUS FUMARATE 60 mg elemental iron TABLET
FERR200000001280TAB490000000000	FERROUS SULFATE 60 mg elemental iron TABLET
FERR200000001282DROPS00231BOTTL	FERROUS SULFATE 15 mg elemental iron/0.6 mL DROPS 15 mL BOTTLE
FERR200000001282DROPS00379BOTTL	FERROUS SULFATE 15 mg elemental iron/0.6 mL DROPS 30 mL BOTTLE
FERR200000001283SYRUP00469BOTTL	FERROUS SULFATE 30 mg elemental iron/5 mL SYRUP 60 mL BOTTLE
FERR300000001193TAB490000000000	FERROUS SULFATE + FOLIC ACID 60 mg elemental iron + 400 mcg folic acid TABLET
FERRO00000001189TAB490000000000	FERROUS SALT equiv. to 60 mg elemental iron TABLET
FERRO00000001190DROPS00231BOTTL	FERROUS SALT (equiv. to 15 mg elemental iron/0.6 mL) DROPS 15 mL BOTTLE
FERRO00000001190DROPS00379BOTTL	FERROUS SALT (equiv. to 15 mg elemental iron/0.6 mL) DROPS 30 mL BOTTLE
FERRO00000001192SYRUP00469BOTTL	FERROUS SALT (equiv. to 30 mg elemental iron/5 mL) SYRUP 60 mL BOTTLE
FILGR00000000219SOL3200000VIALX	FILGRASTIM 150 mcg/0.6 mL SOLUTION VIAL
FILGR00000000370SOL3200000VIALX	FILGRASTIM 300 mcg/1.2 mL SOLUTION VIAL
FILGR00000000371SOL3200000VIALX	FILGRASTIM 300 mcg/mL SOLUTION VIAL
FINAS0000000056TAB490000000000	FINASTERIDE 5 mg TABLET
FLU010001600985SOL3200195VIALX	FLUPHENAZINE (as DECANOATE) 25 mg/mL SOLUTION 10 mL VIAL
FLU010001600985SOL3200516AMPUL	FLUPHENAZINE (as DECANOATE) 25 mg/mL SOLUTION 1 mL AMPULE
FLUC00000000047TAB490000000000	FLUCYTOSINE 500 mg TABLET
FLUCO0000000004CAPSU0000000000	FLUCONAZOLE 200 mg CAPSULE
FLUCO0000000050CAPSU0000000000	FLUCONAZOLE 50 mg CAPSULE
FLUCO00000000215CAPSU0000000000	FLUCONAZOLE 150 mg CAPSULE
FLUCO00000000357SOL3200178VIALX	FLUCONAZOLE 2 mg/ml SOLUTION 100 mL VIAL

FLUMA00000000767SOL3200060AMPUL	FLUMAZENIL 100 mcg/mL SOLUTION 5 mL AMPULE
FLUMA00000000767SOL3200195AMPUL	FLUMAZENIL 100 mcg/mL SOLUTION 10 mL AMPULE
FLUNA00000200056CAPSU00000000000	FLUNARIZINE (as HYDROCHLORIDE) 5 mg CAPSULE
FLUO100000000000SOL3200060BOTTL	FLUOCINOLONE ACETONIDE + POLYMXIN B SULFATE NEOMYCIN SULFATE SOLUTION 5 mL BOTTLE
FLUO200000000078CREA300227TUB01	FLUOCINONIDE 0.05% CREAM 15 g TUBE
FLUO200000000078CREA301245TUB01	FLUOCINONIDE 0.05% CREAM 5 g TUBE
FLUO200000000078CREA301380TUB01	FLUOCINONIDE 0.05% CREAM 10 g TUBE
FLUO200000000078OINTX00227TUB01	FLUOCINONIDE 0.05% OINTMENT 15 g TUBE
FLUO200000000078OINTX01245TUB01	FLUOCINONIDE 0.05% OINTMENT 5 g TUBE
FLUO200000000078OINTX01380TUB01	FLUOCINONIDE 0.05% OINTMENT 10 g TUBE
FLUO30005300133SOL3200060AMPUL	FLUORESCEIN (as SODIUM SALT) 10% SOLUTION 5 mL AMPULE
FLUO30005300142SOL3200060AMPUL	FLUORESCEIN (as SODIUM SALT) 100 mg/ml SOLUTION 5 mL AMPULE
FLUO30005300252STRI100000000000	FLUORESCEIN (as SODIUM SALT) 1 mg STRIPS
FLUO30005300770STRI100000000000	FLUORESCEIN (as SODIUM SALT) 9 mg STRIPS
FLUO40005300013TAB490000000000	FLUORIDE (as SODIUM SALT) 250 mcg TABLET
FLUO40005300045TAB490000000000	FLUORIDE (as SODIUM SALT) 500 mcg TABLET
FLUO40005300771DROPS000000BOTTL	FLUORIDE (as SODIUM SALT) 250 mcg/mL DROPS BOTTLE
FLUO40005300772DROPS000000BOTTL	FLUORIDE (as SODIUM SALT) 500 mcg/15 mL DROPS BOTTLE
FLUO500000000851SOL3200011AMPUL	FLUOROURACIL 50 mg/mL SOLUTION 20 mL AMPULE
FLUO500000000851SOL3200011VIALX	FLUOROURACIL 50 mg/mL SOLUTION 20 mL VIAL
FLUO500000000851SOL3200060AMPUL	FLUOROURACIL 50 mg/mL SOLUTION 5 mL AMPULE
FLUO500000000851SOL3200060VIALX	FLUOROURACIL 50 mg/mL SOLUTION 5 mL VIAL
FLUO500000000851SOL3200178AMPUL	FLUOROURACIL 50 mg/mL SOLUTION 100 mL AMPULE
FLUO500000000851SOL3200178VIALX	FLUOROURACIL 50 mg/mL SOLUTION 100 mL VIAL
FLUO500000000851SOL3200195AMPUL	FLUOROURACIL 50 mg/mL SOLUTION 10 mL AMPULE
FLUO500000000851SOL3200195VIALX	FLUOROURACIL 50 mg/mL SOLUTION 10 mL VIAL
FLUO500000000985SOL3200195VIALX	FLUOROURACIL 25 mg/mL SOLUTION 10 mL VIAL
FLUOX00000000007CAPSU00000000000	FLUOXETINE 20 mg CAPSULE
FLUOX00000000007DISPE00000000000	FLUOXETINE 20 mg DISPERSIBLE TABLET
FLUP30001601199SOL3200195VIALX	FLUPENTIXOL (as DECANOATE) 20 mg/mL SOLUTION 10 mL VIAL
FLUP30001601199SOL3200516AMPUL	FLUPENTIXOL (as DECANOATE) 20 mg/mL SOLUTION 1 mL AMPULE
FLUP30001800013TAB490000000000	FLUPENTIXOL (as DIHYDROCHLORIDE) 250 mcg TABLET
FLUP30001800045TAB490000000000	FLUPENTIXOL (as DIHYDROCHLORIDE) 500 mcg TABLET
FLUP30001800056TAB490000000000	FLUPENTIXOL (as DIHYDROCHLORIDE) 5 mg TABLET
FLUP30001800252TAB490000000000	FLUPENTIXOL (as DIHYDROCHLORIDE) 1 mg TABLET
FLURA0003900228CAPSU00000000000	FLURAZEPAM (as MONOHYDROCHLORIDE) 15 mg CAPSULE
FLUT100000000137DPI01000000INH01	FLUTICASONE + SALMETEROL 100 mcg + 50 mcg x 28 doses DRY POWDER INHALER
FLUT100000000138DPI01000000INH01	FLUTICASONE + SALMETEROL 100 mcg + 50 mcg x 60 doses DRY POWDER INHALER
FLUT100000000202MDIXX000000INH01	FLUTICASONE + SALMETEROL 125 mcg + 25 mcg x 120 doses METERED DOSE INHALER INHALER
FLUT100000000327MDIXX000000INH01	FLUTICASONE + SALMETEROL 250 mcg + 25 mcg x 120 doses METERED DOSE INHALER INHALER
FLUT100000000328DPI01000000INH01	FLUTICASONE + SALMETEROL 250 mcg + 50 mcg x 28 doses DRY POWDER INHALER
FLUT100000000329DPI01000000INH01	FLUTICASONE + SALMETEROL 250 mcg + 50 mcg x 60 doses DRY POWDER INHALER
FLUT100000000421DPI01000000INH01	FLUTICASONE + SALMETEROL 500 mcg + 50 mcg x 28 doses DRY POWDER INHALER

FLUT10000000422DPI0100000INH01	FLUTICASONE + SALMETEROL 500 mcg + 50 mcg x 60 doses DRY POWDER INHALER
FLUT10000000438MDIXX00000INH01	FLUTICASONE + SALMETEROL 50mcg + 25mcg x 120 doses METERED DOSE INHALER INHALER
FLUT10000000783DPI0100000INH01	FLUTICASONE + SALMETEROL 50 mcg + 500 mcg x 60 doses DRY POWDER INHALER
FLUT80005000074OINTX01245TUB01	FLUTICASONE OINTMENT (as PROPIONATE) 0.01% OINTMENT 5 g TUBE
FLUTA0000000015TAB490000000000	FLUTAMIDE 250 mg TABLET
FLUTI00000000204MDIXX00000INH01	FLUTICASONE 125 mcg/dose X 120 doses METERED DOSE INHALER INHALER
FLUTI00000000771RSOLU0036100000	FLUTICASONE 250 mcg/mL RESPIRATORY SOLUTION 2 mL
FLUTI00000000774DPI0100000INH01	FLUTICASONE 50 mcg/dose DRY POWDER INHALER
FLUTI00000000775DPI0100000INH01	FLUTICASONE 250 mcg/dose DRY POWDER INHALER
FLUTI00000000776MDIXX00000INH01	FLUTICASONE 50 mcg/dose x 60 doses METERED DOSE INHALER INHALER
FLUTI00000000777MDIXX00000INH01	FLUTICASONE 50 mcg/dose x 120 doses METERED DOSE INHALER INHALER
FLUTI00000000778MDIXX00000INH01	FLUTICASONE 125 mcg/dose x 60 doses METERED DOSE INHALER INHALER
FLUTI00000000782NAS010000000000	FLUTICASONE 0.05% dose x 120 doses NASAL AQUEOUS SOLUTION
FLUTI0005000078CREA301245TUB01	FLUTICASONE (as PROPIONATE) 0.05% CREAM 5 g TUBE
FOLIC0000000056CAPSU0000000000	FOLIC ACID 5 mg CAPSULE
FOLIC0000000056TAB490000000000	FOLIC ACID 5 mg TABLET
FOLIC00000000252CAPSU0000000000	FOLIC ACID 1 mg CAPSULE
FOLIC00000000252TAB490000000000	FOLIC ACID 1 mg TABLET
FOLIC00000000454SYRUP00201BOTTL	FOLIC ACID 5 mg/5 mL SYRUP 120 mL BOTTLE
FOLIC00000000454SYRUP00469BOTTL	FOLIC ACID 5 mg/5 mL SYRUP 60 mL BOTTLE
FOLIC00000000980SOL3200379BOTTL	FOLIC ACID 2.5 mg/ml SOLUTION 30 mL BOTTLE
FOLIC0005300261SOL3200516AMPUL	FOLIC ACID (as SODIUM SALT) 1 mg/mL SOLUTION 1 mL AMPULE
FOM010000000784SOL3200122VIALX	FOMEPIZOLE 1 g/mL SOLUTION 1.5 mL VIAL
FOND10005300284SOL3200000BOTTL	FONDAPARINUX (as SODIUM SALT) 2.5 mg/0.5 mL SOLUTION BOTTLE
FUROS0000000007TAB490000000000	FUROSEMIDE 20 mg TABLET
FUROS0000000036TAB490000000000	FUROSEMIDE 40 mg TABLET
FUROS0000000066TAB490000000000	FUROSEMIDE 80 mg TABLET
FUROS00000000194SOL3200361AMPUL	FUROSEMIDE 10 mg/mL SOLUTION 2 mL AMPULE
FUROS00000000194SOL3201375AMPUL	FUROSEMIDE 10 mg/mL SOLUTION 25 mL AMPULE
FUSI100000001314CREA301245TUB01	FUSIDATE SODIUM/FUSIDIC ACID 2% CREAM 5 g TUBE
FUSI100000001314OINTX00227TUB01	FUSIDATE SODIUM/FUSIDIC ACID 2% OINTMENT 15 g TUBE
FUSI100000001314OINTX01245TUB01	FUSIDATE SODIUM/FUSIDIC ACID 2% OINTMENT 5 g TUBE
FUSI10006001062DROPS01245TUB01	FUSIDATE SODIUM/FUSIDIC ACID (as SULFATE) 1% DROPS 5 g TUBE
FUSS200000001301OINTX00120SAC01	FUSIDATE SODIUM/FUSIDIC ACID OINTMENT 2%, 1.5 g (medicated surgical dressing) OINTMENT 1.5 g SACHET
GABAP0000000028CAPSU0000000000	GABAPENTIN 300 mg CAPSULE
GABAP0000000034CAPSU0000000000	GABAPENTIN 400 mg CAPSULE
GABAP00000000139CAPSU0000000000	GABAPENTIN 100 mg CAPSULE
GADOD00000000788SOL3200011VIALX	GADODIAMIDE 287 mg/mL SOLUTION 20 mL VIAL
GADOD00000000788SOL3200060VIALX	GADODIAMIDE 287 mg/mL SOLUTION 5 mL VIAL
GADOD00000000788SOL3200195VIALX	GADODIAMIDE 287 mg/mL SOLUTION 10 mL VIAL
GADOD00000000788SOL3200231VIALX	GADODIAMIDE 287 mg/mL SOLUTION 15 mL VIAL
GAL010000000041TAB490000000000	GALANTAMINE 4 mg TABLET
GAL010000000068CAPSU0000000000	GALANTAMINE 8 mg CAPSULE

GAL010000000068TAB490000000000	GALANTAMINE 8 mg TABLET
GAL0100000000235CAPSU0000000000	GALANTAMINE 16 mg CAPSULE
GAL0100000000323CAPSU0000000000	GALANTAMINE 24 mg CAPSULE
GANC100000000084GELXX01245TUB01	GANCICLOVIR 0.15% GEL 5 g TUBE
GANC10005100047SOL3200000VIALX	GANCICLOVIR (as SODIUM) 500 mg SOLUTION VIAL
GEMC10000200004SOL3800000VIALX	GEMCITABINE (as HYDROCHLORIDE) 200 mg POWDER SOLUTION VIAL
GEMC10000200241SOL3800000VIALX	GEMCITABINE (as HYDROCHLORIDE) 1 g POWDER SOLUTION VIAL
GENTA0000000039SOL3200122AMPUL	GENTAMICIN 40 mg/mL SOLUTION 1.5 mL AMPULE
GENTA0000000039SOL3200122VIALX	GENTAMICIN 40 mg/mL SOLUTION 1.5 mL VIAL
GENTA0000000039SOL3200361AMPUL	GENTAMICIN 40 mg/mL SOLUTION 2 mL AMPULE
GENTA0000000039SOL3200361VIALX	GENTAMICIN 40 mg/mL SOLUTION 2 mL VIAL
GENTA0000000039SOL3200516AMPUL	GENTAMICIN 40 mg/mL SOLUTION 1 mL AMPULE
GENTA0000000039SOL3200516VIALX	GENTAMICIN 40 mg/mL SOLUTION 1 mL VIAL
GENTA0000000194SOL3200361VIALX	GENTAMICIN 10 mg/mL SOLUTION 2 mL VIAL
GENTA0000000194SOL3200516AMPUL	GENTAMICIN 10 mg/mL SOLUTION 1 mL AMPULE
GENTA0006000089EDSo200060BOTTL	GENTAMICIN (as SULFATE) 0.30% EYE DROPS SOLUTION 5 mL BOTTLE
GENTA0006000089EYO0101245TUB01	GENTAMICIN (as SULFATE) 0.30% EYE OINTMENT 5 g TUBE
GENTA0006000089EYO0101393TUB01	GENTAMICIN (as SULFATE) 0.30% EYE OINTMENT 3.5 g TUBE
GENTA0006000089EYO0101407TUB01	GENTAMICIN (as SULFATE) 0.30% EYE OINTMENT 3 g TUBE
GENTA0006000649EDSo200060BOTTL	GENTAMICIN (as SULFATE) 0.50% EYE DROPS SOLUTION 5 mL BOTTLE
GFA1000000001188TAB490000000000	GAS FORMING AGENT 470 mg NaHCO ₃ + 70 mg Glutamic Acid + 420 mg Tartaric Acid + 25 mg Silicon Resin T
GLIB0000000003TAB490000000000	GLIBENCLAMIDE 2.5 mg TABLET
GLIB00000000056TAB490000000000	GLIBENCLAMIDE 5 mg TABLET
GLICL00000000031TAB490000000000	GLICLAZIDE 30 mg TABLET
GLICL00000000062TAB490000000000	GLICLAZIDE 60 mg TABLET
GLICL00000000066TAB490000000000	GLICLAZIDE 80 mg TABLET
GLIPI0000000003TAB490000000000	GLIPIZIDE 2.5 mg TABLET
GLIPI00000000056TAB490000000000	GLIPIZIDE 5 mg TABLET
GLU0100000001243SOL3200201GAL01	GLUTARALDEHYDE (GLUTARAL) 2% (with alkaline activating solution) SOLUTION 120 mL GALLON
GLU0100000001243SOL3200247GAL01	GLUTARALDEHYDE (GLUTARAL) 2% (with alkaline activating solution) SOLUTION 1 L GALLON
GLUC100000001446CAPSU0000000000	GLUCOMETAMINE + GLUCODIAMINE + NICOTINAMIDE ASCORBATE 150 mg + 30 mg + 20 mg CAPSULE
GLUC100000001447SYRUP00201BOTTL	GLUCOMETAMINE + GLUCODIAMINE + NICOTINAMIDE ASCORBATE 187.5 mg + 50 mg + 25 mg/5 mL SYRUP 120 mL BOT
GLUC100000001447SYRUP00469BOTTL	GLUCOMETAMINE + GLUCODIAMINE + NICOTINAMIDE ASCORBATE 187.5 mg + 50 mg + 25 mg/5 mL SYRUP 60 mL BOTT
GLUCA0000200252LPS010000000000	GLUCAGON (as HYDROCHLORIDE) 1 mg LYOPHILIZED POWDER + SOLVENT
GLUCO00000000512SOL3200011AMPUL	GLUCOSE (DEXTROSE) 50% SOLUTION 20 mL AMPULE
GLUCO00000000512SOL3200055VIALX	GLUCOSE (DEXTROSE) 50% SOLUTION 50 mL VIAL
GLUCO00000000512SOL3200195AMPUL	GLUCOSE (DEXTROSE) 50% SOLUTION 10 mL AMPULE
GLUCO00000000802SOL3200011VIALX	GLUCOSE (DEXTROSE) 50% (85 Kcal) SOLUTION 20 mL VIAL
GLUCO00000000802SOL3200055VIALX	GLUCOSE (DEXTROSE) 50% (85 Kcal) SOLUTION 50 mL VIAL
GLUCO00000000802SOL3200195VIALX	GLUCOSE (DEXTROSE) 50% (85 Kcal) SOLUTION 10 mL VIAL
GLYC10000000025SUPP400000000000	GLYCEROL 2 g SUPPOSITORY
GLYC100000000805SUPP400000000000	GLYCEROL 2.5 g SUPPOSITORY
GOSER0000100809SOL3900000PFS01	GOSERELIN (as ACETATE) 3.6 mg DEPOT SOLUTION PRE-FILLED SYRINGE

GOSER00000100810SOL3900000PFS01	GOSERELIN (as ACETATE) 10.8 mg DEPOT SOLUTION PRE-FILLED SYRINGE
GRISE0000000047TAB490000000000	GRISEOFULVIN 500 mg TABLET
GRISE00000000205TAB490000000000	GRISEOFULVIN 125 mg TABLET
GTT010004200056PATC400000000000	GLYCERYL TRINITRATE (as NITROGLYCERINE) 5 mg PATCH
GTT010004200069PATC400000000000	GLYCERYL TRINITRATE (as NITROGLYCERINE) 10 mg PATCH
GTT010004200261SOL3200055GVIAL	GLYCERYL TRINITRATE (as NITROGLYCERINE) 1 mg/mL SOLUTION 50 mL GLASS VIAL
GTT010004200261SOL3200060AMPUL	GLYCERYL TRINITRATE (as NITROGLYCERINE) 1 mg/mL SOLUTION 5 mL AMPULE
GTT010004200261SOL3200195AMPUL	GLYCERYL TRINITRATE (as NITROGLYCERINE) 1 mg/mL SOLUTION 10 mL AMPULE
GTT010004200261SOL3201375AMPUL	GLYCERYL TRINITRATE (as NITROGLYCERINE) 1 mg/mL SOLUTION 25 mL AMPULE
GTT010004200808TAB490000000000	GLYCERYL TRINITRATE (as NITROGLYCERINE) 400 mcg TABLET
HACP100000000000000000134700000	HEMODIALYSIS ACID CONCENTRATE WITH POTASSIUM 5 L
HACP200000000000000000134700000	HEMODIALYSIS ACID CONCENTRATE WITHOUT POTASSIUM 5 L
HAEM100000000000000000134700000	HEMODIALYSIS BICARBONATE CONCENTRATE 5 L
HAEM100000000000000000134800000	HEMODIALYSIS BICARBONATE CONCENTRATE 10 L
HAEMO00000000000000000134700000	HEMODIALYSIS ACID CONCENTRATE (DIALYSATE ACETATE BASED) 5 L
HAEMO00000000000000000134800000	HEMODIALYSIS ACID CONCENTRATE (DIALYSATE ACETATE BASED) 10 L
HALOP0000000007TAB490000000000	HALOPERIDOL 20 mg TABLET
HALOP0000000026TAB490000000000	HALOPERIDOL 2 mg TABLET
HALOP0000000045TAB490000000000	HALOPERIDOL 500 mcg TABLET
HALOP0000000056TAB490000000000	HALOPERIDOL 5 mg TABLET
HALOP0000000069TAB490000000000	HALOPERIDOL 10 mg TABLET
HALOP0000000121TAB490000000000	HALOPERIDOL 1.5 mg TABLET
HALOP00000004590000000516AMPUL	HALOPERIDOL 5 mg/mL 1 mL AMPULE
HALOP00000008510000000516AMPUL	HALOPERIDOL 50 mg/mL 1 mL AMPULE
HALOT00000000200000000000BOTTL	HALOTHANE 250 mL BOTTLE
HEMO10000000813SOL3200000VIALX	HEMOPHILUS INFLUENZAE TYPE B (HIB) CONJUGATED VACCINE 10 mcg/0.5 mL, 1 dose + 0.5 mL diluent with te
HEMO10000000814SOL3200000VIALX	HEMOPHILUS INFLUENZAE TYPE B (HIB) CONJUGATED VACCINE 10 mcg/0.5 mL + 0.9% Sodium Chloride with Diph
HEMO10000000815SOL3200000VIALX	HEMOPHILUS INFLUENZAE TYPE B (HIB) CONJUGATED VACCINE 0.5 mL with Meningococcal Protein SOLUTION VI
HEMOD0000000811SOL3201347DRU01	HEMODIALYSIS SOLUTION 5 gallon (concentrate) SOLUTION 5 L DRUM
HEMOD0000000811SOL3201348DRU01	HEMODIALYSIS SOLUTION 5 gallon (concentrate) SOLUTION 10 L DRUM
HEP010005100145SOL3200060VIALX	HEPARIN (as SODIUM) 1000 IU/mL SOLUTION 5 mL VIAL
HEP010005100145SOL3200379VIALX	HEPARIN (as SODIUM) 1000 IU/mL SOLUTION 30 mL VIAL
HEP010005101285SOL3200060VIALX	HEPARIN (as SODIUM) 5000 IU/mL SOLUTION 5 mL VIAL
HEPA20000000816SOL3200000MV001	HEPATITIS A INACTIVATED VACCINE viral antigen not less than 720 EU in 0.5 mL (junior) SOLUTION MONO
HEPA20000000817SOL3200000MV001	HEPATITIS A INACTIVATED VACCINE viral antigen not less than 1440 EU in 1 mL (adult) SOLUTION MONODO
HEPA20000000818SOL3200000PFS03	HEPATITIS A INACTIVATED VACCINE 80 units/0.5 mL (pediatric) SOLUTION PRE-FILLED SYRINGE (SINGLE DOS
HEPA20000000819SOL3200000PFS03	HEPATITIS A INACTIVATED VACCINE 160 units/0.5 mL (adult) SOLUTION PRE-FILLED SYRINGE (SINGLE DOSE)
HEPA30000000094SOL3200000VIALX	HEPATITIS B IMMUNOGLOBULIN (HUMAN) 0.5 mL SOLUTION VIAL
HEPA30000000361SOL3200000VIALX	HEPATITIS B IMMUNOGLOBULIN (HUMAN) 2 mL SOLUTION VIAL
HEPA30000000516SOL3200000VIALX	HEPATITIS B IMMUNOGLOBULIN (HUMAN) 1 mL SOLUTION VIAL
HEPA50000000822SOL3200000MV001	HEPATITIS B VACCINE (RECOMBINANT DNA) 10 mcg/0.5 mL (pediatric) SOLUTION MONODOSE VIAL

HEPA50000000824SOL3200060VIALX	HEPATITIS B VACCINE (RECOMBINANT DNA) 20 mcg/mL SOLUTION 5 mL VIAL
HEPA50000000824SOL3200195VIALX	HEPATITIS B VACCINE (RECOMBINANT DNA) 20 mcg/mL SOLUTION 10 mL VIAL
HEPA50000000824SOL3200516VIALX	HEPATITIS B VACCINE (RECOMBINANT DNA) 20 mcg/mL SOLUTION 1 mL VIAL
HUMA10000001361LPS0200000VIALX	HUMAN GROWTH HORMONE (BIOSYNTHETIC) 4 mg + 2 mL LYOPHILIZED POWDER + DILUENT VIAL
HUMA10000001362LPS0200000VIALX	HUMAN GROWTH HORMONE (BIOSYNTHETIC) 5 mg + 5 mL LYOPHILIZED POWDER + DILUENT VIAL
HUMA20000000831LPS0100000AMPUL	HUMAN MENOPAUSAL GONADOTROPHIN (HMG, MENOTROPHIN) 75 IU FSH + 75 IU LH + 1 mL LYOPHILIZED POWDER + S
HUMAN0000001358LPS0100000AMPUL	HUMAN CHORIONIC GONADOTROPHIN (HCG) 1000 IU + 1 mL LYOPHILIZED POWDER + SOLVENT AMPULE
HUMAN0000001358LPS0100000VIALX	HUMAN CHORIONIC GONADOTROPHIN (HCG) 1000 IU + 1 mL LYOPHILIZED POWDER + SOLVENT VIAL
HUMAN0000001359LPS0100000AMPUL	HUMAN CHORIONIC GONADOTROPHIN (HCG) 500 IU + 1 mL LYOPHILIZED POWDER + SOLVENT AMPULE
HUMAN0000001360LPS0100000AMPUL	HUMAN CHORIONIC GONADOTROPHIN (HCG) 5000 IU + 1 mL LYOPHILIZED POWDER + SOLVENT AMPULE
HUMAN0000001360LPS0100000VIALX	HUMAN CHORIONIC GONADOTROPHIN (HCG) 5000 IU + 1 mL LYOPHILIZED POWDER + SOLVENT VIAL
HYDR10000000279LOT1101375BOTTL	HYDROCORTISONE 2.50% LOTION 25 mL BOTTLE
HYDR10000001062CREA301245TUB01	HYDROCORTISONE 1% CREAM 5 g TUBE
HYDR10000001062CREA301380TUB01	HYDROCORTISONE 1% CREAM 10 g TUBE
HYDR10000001062CREA301381JAR01	HYDROCORTISONE 1% CREAM 500 g JAR
HYDR10000001062LOT1101375BOTTL	HYDROCORTISONE 1% LOTION 25 mL BOTTLE
HYDR10000001062OINTX01245TUB01	HYDROCORTISONE 1% OINTMENT 5 g TUBE
HYDR10000001062OINTX01381JAR01	HYDROCORTISONE 1% OINTMENT 500 g JAR
HYDR10005400007TAB490000000000	HYDROCORTISONE (as SODIUM SUCCINATE) 20 mg TABLET
HYDR10005400015POW2700000VIALX	HYDROCORTISONE (as SODIUM SUCCINATE) 250 mg POWDER VIAL
HYDR10005400047POW2700000VIALX	HYDROCORTISONE (as SODIUM SUCCINATE) 500 mg POWDER VIAL
HYDR10005400056TAB490000000000	HYDROCORTISONE (as SODIUM SUCCINATE) 5 mg TABLET
HYDR10005400139POW2700000VIALX	HYDROCORTISONE (as SODIUM SUCCINATE) 100 mg POWDER VIAL
HYDR10005400214SOL3200361VIALX	HYDROCORTISONE (as SODIUM SUCCINATE) 125 mg/mL SOLUTION 2 mL VIAL
HYDR10005400214SOL3201386VIALX	HYDROCORTISONE (as SODIUM SUCCINATE) 125 mg/mL SOLUTION 4 mL VIAL
HYDR10005400851SOL3200361VIALX	HYDROCORTISONE (as SODIUM SUCCINATE) 50 mg/mL SOLUTION 2 mL VIAL
HYDR10005401062OINTX01380TUB01	HYDROCORTISONE (as SODIUM SUCCINATE) 1% OINTMENT 10 g TUBE
HYDR40000000364SOL3200201BOTTL	HYDROGEN PEROXIDE 3% SOLUTION 120 mL BOTTLE
HYDR40000000364SOL3200469BOTTL	HYDROGEN PEROXIDE 3% SOLUTION 60 mL BOTTLE
HYDR60000000013TAB490000000000	HYDROXOCOBALAMIN (VITAMIN B12) 250 mcg TABLET
HYDR600000000148TAB490000000000	HYDROXOCOBALAMIN (VITAMIN B12) 100 mcg TABLET
HYDR600000000261SOL3200516AMPUL	HYDROXOCOBALAMIN (VITAMIN B12) 1 mg/mL SOLUTION 1 mL AMPULE
HYDR600000000261SOL3200516VIALX	HYDROXOCOBALAMIN (VITAMIN B12) 1 mg/mL SOLUTION 1 mL VIAL
HYDR70000000004TAB490000000000	HYDROXYCHLOROQUINE 200 mg TABLET
HYDR800000000133SOL3200020BOTTL	HYDROXYETHYL STARCH 10% SOLUTION 250 mL BOTTLE
HYDR800000000133SOL3200044BOTTL	HYDROXYETHYL STARCH 10% SOLUTION 500 mL BOTTLE
HYDR800000000462SOL3200020BOTTL	HYDROXYETHYL STARCH 6% SOLUTION 250 mL BOTTLE
HYDR800000000462SOL3200044BOTTL	HYDROXYETHYL STARCH 6% SOLUTION 500 mL BOTTLE
HYDR90000000047CAPSU0000000000	HYDROXYUREA 500 mg CAPSULE
HYDRA0000200021TAB490000000000	HYDRALAZINE (as HYDROCHLORIDE) 25 mg TABLET
HYDRA0000200050TAB490000000000	HYDRALAZINE (as HYDROCHLORIDE) 50 mg TABLET

HYDRA0000200069TAB490000000000	HYDRALAZINE (as HYDROCHLORIDE) 10 mg TABLET
HYDRA0000201199SOL3200516AMPUL	HYDRALAZINE (as HYDROCHLORIDE) 20 mg/mL SOLUTION 1 mL AMPULE
HYDRO0000000021TAB490000000000	HYDROCHLOROTHIAZIDE 25 mg TABLET
HYDRO0000000050TAB490000000000	HYDROCHLOROTHIAZIDE 50 mg TABLET
HYDRO00000000974TAB490000000000	HYDROCHLOROTHIAZIDE 12.5 mg TABLET
HYDRX00000000851SOL3200516VIALX	HYDROXYZINE 50 mg/mL SOLUTION 1 mL VIAL
HYDRX0001800021TAB490000000000	HYDROXYZINE (as DIHYDROCHLORIDE) 25 mg TABLET
HYDRX0001800069TAB490000000000	HYDROXYZINE (as DIHYDROCHLORIDE) 10 mg TABLET
HYDRX0001800357SYRUP00178BOTTL	HYDROXYZINE (as DIHYDROCHLORIDE) 2 mg/ml SYRUP 100 mL BOTTLE
HYDRX0001800357SYRUP00469BOTTL	HYDROXYZINE (as DIHYDROCHLORIDE) 2 mg/ml SYRUP 60 mL BOTTLE
HYOSC0002900668SOL3200516AMPUL	HYOSCINE (as HYDROBROMIDE) 500 mcg/mL SOLUTION 1 mL AMPULE
HYOSC0002901018SOL3200516AMPUL	HYOSCINE (as HYDROBROMIDE) 200 mcg/mL SOLUTION 1 mL AMPULE
HYOSC0004100454SYRUP00469BOTTL	HYOSCINE (as N-BUTYL BROMIDE) 5 mg/5 mL SYRUP 60 mL BOTTLE
HYOSC0004101199SOL3200516AMPUL	HYOSCINE (as N-BUTYL BROMIDE) 20 mg/mL SOLUTION 1 mL AMPULE
HYP010000000364SOL3200000BOTTL	HYPERTONIC SALINE 3% SOLUTION BOTTLE
HYPRO00000000089OPS0100195BOTTL	HYPROMELLOSE 0.30% OPHTHALMIC SOLUTION 10 mL BOTTLE
HYPRO00000000194EDS0200195BOTTL	HYPROMELLOSE 10 mg/mL EYE DROPS SOLUTION 10 mL BOTTLE
HYPRO00000000194EDS0200231BOTTL	HYPROMELLOSE 10 mg/mL EYE DROPS SOLUTION 15 mL BOTTLE
HYPRO00000000459EDS0200195BOTTL	HYPROMELLOSE 5 mg/mL EYE DROPS SOLUTION 10 mL BOTTLE
HYPRO00000000459EDS0200231BOTTL	HYPROMELLOSE 5 mg/mL EYE DROPS SOLUTION 15 mL BOTTLE
HYS010004100069TAB490000000000	HYOSCINE (as N-BUTYL BROMIDE) 10 mg TABLET
IBUPR00000000004TAB490000000000	IBUPROFEN 200 mg TABLET
IBUPR00000000006SUS1400469BOTTL	IBUPROFEN 200 mg/5 mL SUSPENSION 60 mL BOTTLE
IBUPR00000000034TAB490000000000	IBUPROFEN 400 mg TABLET
IBUPR00000000162SYRUP00469BOTTL	IBUPROFEN 100 mg/5 mL SYRUP 60 mL BOTTLE
IBUPR00000000495TAB520000000000	IBUPROFEN 800 mg TABLET MODIFIED RELEASE
IDARU00000000056POW2700000VIALX	IDARUBICIN 5 mg POWDER VIAL
IFOSF00000000025POW2700000VIALX	IFOSFAMIDE 2 g POWDER VIAL
IFOSF000000000241POW2700000VIALX	IFOSFAMIDE 1 g POWDER VIAL
IIHo100000000761SOL3200060VIALX	ISOPHANE INSULIN HUMAN 100 IU/mL SOLUTION 5 mL VIAL
IIHo100000000761SOL3200195VIALX	ISOPHANE INSULIN HUMAN 100 IU/mL SOLUTION 10 mL VIAL
IIHo100000000761SOL3200388PFS01	ISOPHANE INSULIN HUMAN 100 IU/mL SOLUTION 3 mL PRE-FILLED SYRINGE
IIHo100000001339SOL3200060VIALX	ISOPHANE INSULIN HUMAN 100 IU/mL (Equiv. to 3.571 mg/mL) SOLUTION 5 mL VIAL
IMB690000000034TAB490000000000	IMATINIB 400 mg TABLET
IMB6900000000139TAB490000000000	IMATINIB 100 mg TABLET
IMIDA00000000000SHAMX00178BOTTL	IMIDAZOLES (TOPICAL) SHAMPOO 100 mL BOTTLE
IMIDA00000000000SHAMX00195SAC01	IMIDAZOLES (TOPICAL) SHAMPOO 10 mL SACHET
IMIDA00000000000SHAMX00469BOTTL	IMIDAZOLES (TOPICAL) SHAMPOO 60 mL BOTTLE
IMIDA00000000000SHAMX00842SAC01	IMIDAZOLES (TOPICAL) SHAMPOO 6 mL SACHET
IMIDA00000000028OVULE0000000000	IMIDAZOLES (TOPICAL) 300 mg OVULE
IMIDA00000000047OVULE0000000000	IMIDAZOLES (TOPICAL) 500 mg OVULE
IMIDA00000000139OVULE0000000000	IMIDAZOLES (TOPICAL) 100 mg OVULE
IMIDA00000000194SOL3200231BOTTL	IMIDAZOLES (TOPICAL) 10 mg/mL SOLUTION 15 mL BOTTLE
IMIDA00000000194SOL3201375BOTTL	IMIDAZOLES (TOPICAL) 10 mg/mL SOLUTION 25 mL BOTTLE

IMIDA0000000227CREA300000TUB01	IMIDAZOLES (TOPICAL) 15 g CREAM TUBE
IMIDA0000000841CREA300227TUB01	IMIDAZOLES (TOPICAL) 1%-2% CREAM 15 g TUBE
IMIDA0000000841CREA301393TUB01	IMIDAZOLES (TOPICAL) 1%-2% CREAM 3.5 g TUBE
IMIDA0000000841CREA301398TUB01	IMIDAZOLES (TOPICAL) 1%-2% CREAM 450 g TUBE
IMIGL00000001332SOL3200011VIALX	IMIGLUCERASE 400 U Powder SOLUTION 20 mL VIAL
IMIPR00000200021TAB490000000000	IMIPRAMINE (as HYDROCHLORIDE) 25 mg TABLET
IMMU100000000177SOL3200000VIALX	IMMUNOGLOBULIN NORMAL, HUMAN (IGIV) 100mg/mL/1g SOLUTION VIAL
IMMU100000000331SOL3200000VIALX	IMMUNOGLOBULIN NORMAL, HUMAN (IGIV) 2500 mg SOLUTION VIAL
IMMU100000000851SOL3200011VIALX	IMMUNOGLOBULIN NORMAL, HUMAN (IGIV) 50 mg/mL SOLUTION 20 mL VIAL
IMMU100000000851SOL3200055VIALX	IMMUNOGLOBULIN NORMAL, HUMAN (IGIV) 50 mg/mL SOLUTION 50 mL VIAL
IMMU100000000851SOL3200178VIALX	IMMUNOGLOBULIN NORMAL, HUMAN (IGIV) 50 mg/mL SOLUTION 100 mL VIAL
IMMU100000000851SOL3200195VIALX	IMMUNOGLOBULIN NORMAL, HUMAN (IGIV) 50 mg/mL SOLUTION 10 mL VIAL
IMMU100000000851SOL3201125VIALX	IMMUNOGLOBULIN NORMAL, HUMAN (IGIV) 50 mg/mL SOLUTION 200 mL VIAL
IMMU100000001160FPR0100000BOTTL	IMMUNOGLOBULIN NORMAL, HUMAN (IGIV) 1 g + 50 mL FREEZE-DRIED POWDER + RECONSTITUTION FLUID BOTTLE
IMMU100000001363FPR0100000BOTTL	IMMUNOGLOBULIN NORMAL, HUMAN (IGIV) 1 g + 33 mL FREEZE-DRIED POWDER + RECONSTITUTION FLUID BOTTLE
IMMU100000001364FDP0100000BOTTL	IMMUNOGLOBULIN NORMAL, HUMAN (IGIV) 10 g + 260 mL FREEZE-DRIED POWDER + DILUENT BOTTLE
IMMU100000001365FPR0100000BOTTL	IMMUNOGLOBULIN NORMAL, HUMAN (IGIV) 2.5 g + 100 mL FREEZE-DRIED POWDER + RECONSTITUTION FLUID BOTT
IMMU100000001366FDP0100000BOTTL	IMMUNOGLOBULIN NORMAL, HUMAN (IGIV) 5 g + 100 mL FREEZE-DRIED POWDER + DILUENT BOTTLE
IMMUN00000000847SOL3200060VIALX	IMMUNOGLOBULIN NORMAL, HUMAN (IGIM) 160 mg/mL SOLUTION 5 mL VIAL
IMMUN00000000847SOL3200195VIALX	IMMUNOGLOBULIN NORMAL, HUMAN (IGIM) 160 mg/mL SOLUTION 10 mL VIAL
IMMUN00000000847SOL3200361VIALX	IMMUNOGLOBULIN NORMAL, HUMAN (IGIM) 160 mg/mL SOLUTION 2 mL VIAL
INDAP00000000858TAB490000000000	INDAPAMIDE 1.5 mg MR TABLET
INDIN00000000004CAPSU0000000000	INDINAVIR 200 mg CAPSULE
INDIN00000000034CAPSU0000000000	INDINAVIR 400 mg CAPSULE
INDOM00000000021CAPSU0000000000	INDOMETACIN 25 mg CAPSULE
INDOM00000001018SOL3200060AMPUL	INDOMETACIN 200 mcg/mL SOLUTION 5 mL AMPULE
INPOL00000000094SOL3200060VIALX	INACTIVATED POLIOMYELITIS VACCINE (TYPE 1, 2 & 3) 0.5 mL SOLUTION 5 mL VIAL
INPOL00000000745SOL3200000VIALX	INACTIVATED POLIOMYELITIS VACCINE (TYPE 1, 2 & 3) 0.5 mL monodose SOLUTION VIAL
INSUL00000000761SOL3200195VIALX	INSULIN ZINC SUSPENSION (HUMAN) 100 IU/mL SOLUTION 10 mL VIAL
INTE100000000868SOL3200000VIALX	INTERFERON ALFA 2B (HUMAN) 3 million IU/mL + diluent SOLUTION VIAL
INTE100000000869SOL3200000VIALX	INTERFERON ALFA 2B (HUMAN) 5 million IU/mL + diluent SOLUTION VIAL
INTE100000000870SOL3200000VIALX	INTERFERON ALFA 2B (HUMAN) 10 million IU/mL + diluent SOLUTION VIAL
INTER00000000861TAB490000000000	INTERFERON ALFA 2A (HUMAN) 200 IU TABLET
INTER00000000862SOL3200000PFS01	INTERFERON ALFA 2A (HUMAN) 3 million IU SOLUTION PRE-FILLED SYRINGE
INTER00000000862SOL3200516VIALX	INTERFERON ALFA 2A (HUMAN) 3 million IU SOLUTION 1 mL VIAL
INTER00000000864SOL3200516VIALX	INTERFERON ALFA 2A (HUMAN) 4.5 million IU SOLUTION 1 mL VIAL
INTER00000000865SOL3200516VIALX	INTERFERON ALFA 2A (HUMAN) 9 million IU SOLUTION 1 mL VIAL
INTER00000000866SOL3200000PFS01	INTERFERON ALFA 2A (HUMAN) 4.5 million IU/0.5 mL SOLUTION PRE-FILLED SYRINGE

INTER00000000867SOL3200000PFS01	INTERFERON ALFA 2A (HUMAN) 6 million IU/0.5 mL SOLUTION PRE-FILLED SYRINGE
INTRA00000000000SOL3200020BOTTL	INTRAOCULAR IRRIGATING SOLUTION (BALANCED SALT SOLUTION) SOLUTION 250 mL BOTTLE
INTRA00000000000SOL3200044BOTTL	INTRAOCULAR IRRIGATING SOLUTION (BALANCED SALT SOLUTION) SOLUTION 500 mL BOTTLE
INTRA00000000000SOL3200231BOTTL	INTRAOCULAR IRRIGATING SOLUTION (BALANCED SALT SOLUTION) SOLUTION 15 mL BOTTLE
IODAM00000000874SOL3200055VIALX	IODAMIDE 495 mg/mL equiv. to 300 mg/mL iodine SOLUTION 50 mL VIAL
IODAM00000000874SOL3200178VIALX	IODAMIDE 495 mg/mL equiv. to 300 mg/mL iodine SOLUTION 100 mL VIAL
IODAM00000000874SOL3200379AMPUL	IODAMIDE 495 mg/mL equiv. to 300 mg/mL iodine SOLUTION 30 mL AMPULE
IODAM00000000875SOL3200055BOTTL	IODAMIDE 496 mg/mL equiv. to 300 mg/mL iodine SOLUTION 50 mL BOTTLE
IODAM00000000875SOL3200178BOTTL	IODAMIDE 496 mg/mL equiv. to 300 mg/mL iodine SOLUTION 100 mL BOTTLE
IODAM00000000878SOL3200055BOTTL	IODAMIDE 627.9 mg/mL equiv. to 380 mg/mL iodine SOLUTION 50 mL BOTTLE
IODAM00000000878SOL3200055VIALX	IODAMIDE 627.9 mg/mL equiv. to 380 mg/mL iodine SOLUTION 50 mL VIAL
IODAM00000000878SOL3200178BOTTL	IODAMIDE 627.9 mg/mL equiv. to 380 mg/mL iodine SOLUTION 100 mL BOTTLE
IODAM00000000878SOL3200178VIALX	IODAMIDE 627.9 mg/mL equiv. to 380 mg/mL iodine SOLUTION 100 mL VIAL
IODIZ00000000047CAP3100000000000	IODIZED OIL FLUID 500 mg CAPSULE SOFT GELATIN
IOHEX00000000880SOL3200055VIALX	IOHEXOL 180 mg iodine/mL SOLUTION 50 mL VIAL
IOHEX00000000881SOL3200055VIALX	IOHEXOL 240 mg iodine/mL SOLUTION 50 mL VIAL
IOHEX00000000882SOL3200055VIALX	IOHEXOL 300 mg iodine/mL SOLUTION 50 mL VIAL
IOHEX00000000882SOL3200178VIALX	IOHEXOL 300 mg iodine/mL SOLUTION 100 mL VIAL
IOHEX00000000884SOL3200055VIALX	IOHEXOL 350 mg iodine/mL SOLUTION 50 mL VIAL
IOHEX00000000884SOL3200178VIALX	IOHEXOL 350 mg iodine/mL SOLUTION 100 mL VIAL
IOHEX00000000884SOL3201125VIALX	IOHEXOL 350 mg iodine/mL SOLUTION 200 mL VIAL
IOPAM00000000887SOL3200195AMPUL	IOPAMIDOL 408 mg/mL equiv. to 200 mg iodine SOLUTION 10 mL AMPULE
IOPAM00000000890SOL3200055VIALX	IOPAMIDOL 612 mg/mL equiv. to 300 mg iodine SOLUTION 50 mL VIAL
IOPAM00000000890SOL3200178VIALX	IOPAMIDOL 612 mg/mL equiv. to 300 mg iodine SOLUTION 100 mL VIAL
IOPAM00000000890SOL3200195AMPUL	IOPAMIDOL 612 mg/mL equiv. to 300 mg iodine SOLUTION 10 mL AMPULE
IOPAM00000000890SOL3200379VIALX	IOPAMIDOL 612 mg/mL equiv. to 300 mg iodine SOLUTION 30 mL VIAL
IOPAM00000000893SOL3200055VIALX	IOPAMIDOL 755 mg/mL equiv. to 370 mg iodine SOLUTION 50 mL VIAL
IOPAM00000000893SOL3200178BOTTL	IOPAMIDOL 755 mg/mL equiv. to 370 mg iodine SOLUTION 100 mL BOTTLE
IOPAM00000000893SOL3200379AMPUL	IOPAMIDOL 755 mg/mL equiv. to 370 mg iodine SOLUTION 30 mL AMPULE
IOPRO00000000895SOL3200055VIALX	IOPROMIDE 240 mg/mL equiv. to 499 mg iodine SOLUTION 50 mL VIAL
IOPRO00000000900SOL3200044VIALX	IOPROMIDE 300 mg/mL equiv. to 623 mg iodine SOLUTION 500 mL VIAL
IOPRO00000000900SOL3200055VIALX	IOPROMIDE 300 mg/mL equiv. to 623 mg iodine SOLUTION 50 mL VIAL
IOPRO00000000900SOL3200178VIALX	IOPROMIDE 300 mg/mL equiv. to 623 mg iodine SOLUTION 100 mL VIAL
IOPRO00000000900SOL3200379VIALX	IOPROMIDE 300 mg/mL equiv. to 623 mg iodine SOLUTION 30 mL VIAL
IOPRO00000000900SOL3201125VIALX	IOPROMIDE 300 mg/mL equiv. to 623 mg iodine SOLUTION 200 mL VIAL
IOPRO00000000902SOL3200055VIALX	IOPROMIDE 370 mg/mL equiv. to 769 mg iodine SOLUTION 50 mL VIAL
IOPRO00000000902SOL3200178VIALX	IOPROMIDE 370 mg/mL equiv. to 769 mg iodine SOLUTION 100 mL VIAL
IOTHA0003500905SOL3200055VIALX	IOTHALAMATE (as MEGLUMINE) 600 mg/mL SOLUTION 50 mL VIAL
IOTHA0003500905SOL3200178VIALX	IOTHALAMATE (as MEGLUMINE) 600 mg/mL SOLUTION 100 mL VIAL
IOTHA0003500905SOL3200379VIALX	IOTHALAMATE (as MEGLUMINE) 600 mg/mL SOLUTION 30 mL VIAL

IOVER00000000908SOL3200055VIALX	IOVERSOL 636 mg/mL equiv. to 300 mg/mL iodine SOLUTION 50 mL VIAL
IOVER00000000908SOL3200178VIALX	IOVERSOL 636 mg/mL equiv. to 300 mg/mL iodine SOLUTION 100 mL VIAL
IOVER00000000908SOL3200379VIALX	IOVERSOL 636 mg/mL equiv. to 300 mg/mL iodine SOLUTION 30 mL VIAL
IOVER00000000910SOL3200055VIALX	IOVERSOL 678 mg/mL equiv. to 320 mg/mL iodine SOLUTION 50 mL VIAL
IOVER00000000910SOL3200178VIALX	IOVERSOL 678 mg/mL equiv. to 320 mg/mL iodine SOLUTION 100 mL VIAL
IOVER00000000912SOL3200055VIALX	IOVERSOL 741 mg/mL equiv. to 350 mg/mL iodine SOLUTION 50 mL VIAL
IOVER00000000912SOL3200178VIALX	IOVERSOL 741 mg/mL equiv. to 350 mg/mL iodine SOLUTION 100 mL VIAL
IOXIT00000000915SOL3200011VIALX	IOXITHALAMIC ACID 9.66 g sodium ioxithalamate + 65.09 g meglumine ioxithalamate (equiv. to 35 g iod
IOXIT00000000915SOL3200055VIALX	IOXITHALAMIC ACID 9.66 g sodium ioxithalamate + 65.09 g meglumine ioxithalamate (equiv. to 35 g iod
IOXIT00000000915SOL3200178VIALX	IOXITHALAMIC ACID 9.66 g sodium ioxithalamate + 65.09 g meglumine ioxithalamate (equiv. to 35 g iod
IOXIT00000000917SOL3200379BOTTL	IOXITHALAMIC ACID 19.81 g meglumine ioxithalamate (equiv. to 9 g iodine) SOLUTION 30 mL BOTTLE
IOXIT00000000918SOL3200055BOTTL	IOXITHALAMIC ACID 33.015 g meglumine ioxithalamate (equiv. to 15 g iodine) SOLUTION 50 mL BOTTLE
IOXIT00000000919SOL3200178BOTTL	IOXITHALAMIC ACID 66.03 g meglumine ioxithalamate (equiv. to 30 g iodine) SOLUTION 100 mL BOTTLE
IPRA10000000924MDIXX00000000000	IPRATROPIUM + FENOTEROL 20 mcg + 50 mcg x 10 mL doses METERED DOSE INHALER
IPRA10000000925RSOLU01410000000	IPRATROPIUM + FENOTEROL 250 mcg + 500 mcg per mL RESPIRATORY SOLUTION 20 ml (multidose)
IPRA10000000926RSOLU00000000000	IPRATROPIUM + FENOTEROL 500 mcg + 1.25 mg per 4 mL (unit dose) RESPIRATORY SOLUTION
IPRA20000000927RSOLU00000000000	IPRATROPIUM + SALBUTAMOL 500 mcg + 2.5 mg x 2.5 mL (unit dose) RESPIRATORY SOLUTION
IPRA20000000928MDIXX00000000000	IPRATROPIUM + SALBUTAMOL 21 mcg + 120 mcg x 200 doses x 10 mL METERED DOSE INHALER
IPRA200000001210RSOLU00000000000	IPRATROPIUM + SALBUTAMOL 500 mcg + 2.5 mg x 20 mL (unit dose) RESPIRATORY SOLUTION
IPRAT00000000771RSOLU01408000000	IPRATROPIUM 250 mcg/mL RESPIRATORY SOLUTION 1 ml (unit dose)
IPRAT00000000771RSOLU01409000000	IPRATROPIUM 250 mcg/mL RESPIRATORY SOLUTION 2 ml (unit dose)
IPRAT00000000771RSOLU01410000000	IPRATROPIUM 250 mcg/mL RESPIRATORY SOLUTION 20 ml (multidose)
IPRAT00000000920MDIXX00000000000	IPRATROPIUM 20 mcg/dose x 200 doses METERED DOSE INHALER
IPVA10000000094SOL3200000VPS01	INFLUENZA POLYVALENT VACCINE 0.5 mL SOLUTION VIAL + PRE-FILLED SYRINGE DILUENT
IPVA10000000094SUS14000000AMPUL	INFLUENZA POLYVALENT VACCINE 0.5 mL SUSPENSION AMPULE
IPVA10000000094SUS14000000PFS01	INFLUENZA POLYVALENT VACCINE 0.5 mL SUSPENSION PRE-FILLED SYRINGE
IPVA100000001367SOL3200060GVIAL	INFLUENZA POLYVALENT VACCINE 15 mcg/0.5 mL multidose SOLUTION 5 mL GLASS VIAL
IRBE100000000222TAB4900000000000	IRBESARTAN + HYDROCHLOROTHIAZIDE 150 mg + 12.5 mg TABLET
IRBE100000000372TAB4900000000000	IRBESARTAN + HYDROCHLOROTHIAZIDE 300 mg + 12.5 mg TABLET
IRBE100000001218TAB4900000000000	IRBESARTAN + HYDROCHLOROTHIAZIDE 300 mg + 25 mg TABLET
IRBES00000000028TAB4900000000000	IRBESARTAN 300 mg TABLET
IRBES00000000065TAB4900000000000	IRBESARTAN 75 mg TABLET
IRBES000000000215TAB4900000000000	IRBESARTAN 150 mg TABLET
IRE010000000489TAB4900000000000	ISONIAZID + RIFAMPICIN + ETHAMBUTOL 75 mg + 150 mg + 275 mg TABLET
IRINO00000200162CONCE000000VIALX	IRINOTECAN (as HYDROCHLORIDE) 100 mg/5 mL CONCENTRATE VIAL
IRINO00000200929CONCE000000VIALX	IRINOTECAN (as HYDROCHLORIDE) 40 mg/2 mL CONCENTRATE VIAL
IRON500000000851SOL3200361AMPUL	IRON DEXTRAN 50 mg/mL SOLUTION 2 mL AMPULE

ISOCO000000000000LOTI100195BOTTL	ISOCONAZOLE LOTION 10 mL BOTTLE
ISOCO000000000000SHAMX00178BOTTL	ISOCONAZOLE SHAMPOO 100 mL BOTTLE
ISOCO000000000000SHAMX00195SAC01	ISOCONAZOLE SHAMPOO 10 mL SACHET
ISOCO000000000000SHAMX00469BOTTL	ISOCONAZOLE SHAMPOO 60 mL BOTTLE
ISOCO000000000000SHAMX00842SAC01	ISOCONAZOLE SHAMPOO 6 mL SACHET
ISOCO00000000028OVULE00000000000	ISOCONAZOLE 300 mg OVULE
ISOCO00000000047OVULE00000000000	ISOCONAZOLE 500 mg OVULE
ISOCO00000000139OVULE00000000000	ISOCONAZOLE 100 mg OVULE
ISOCO00000000194SOL3200231BOTTL	ISOCONAZOLE 10 mg/mL SOLUTION 15 mL BOTTLE
ISOCO00000000194SOL3201375BOTTL	ISOCONAZOLE 10 mg/mL SOLUTION 25 mL BOTTLE
ISOCO00000001062CREA300227TUB01	ISOCONAZOLE 1% CREAM 15 g TUBE
ISOCO00000001062CREA301393TUB01	ISOCONAZOLE 1% CREAM 3.5 g TUBE
ISOCO00000001062CREA301398TUB01	ISOCONAZOLE 1% CREAM 450 g TUBE
ISOCO00000001314CREA300227TUB01	ISOCONAZOLE 2% CREAM 15 g TUBE
ISOCO00000001314CREA301393TUB01	ISOCONAZOLE 2% CREAM 3.5 g TUBE
ISOCO00000001314CREA301398TUB01	ISOCONAZOLE 2% CREAM 450 g TUBE
ISOFL00000000020SOL3200000BOTTL	ISOFLURANE 250 mL SOLUTION BOTTLE
ISOFL00000000178SOL3200000BOTTL	ISOFLURANE 100 mL SOLUTION BOTTLE
ISON100000000931TAB4900000000000	ISONIAZID + ETHAMBUTOL 150 mg + 400 mg TABLET
ISON100000000932TAB4900000000000	ISONIAZID + ETHAMBUTOL 200 mg + 500 mg TABLET
ISON300000000488TAB4900000000000	ISONIAZID + RIFAMPICIN 75 mg + 150 mg TABLET
ISON300000000933TAB4900000000000	ISONIAZID + RIFAMPICIN 30 mg + 60 mg (pediatric) TABLET
ISON300000000934TAB4900000000000	ISONIAZID + RIFAMPICIN 60 mg + 60 mg (pediatric) TABLET
ISON300000000935TAB4900000000000	ISONIAZID + RIFAMPICIN 150 mg + 150 mg TABLET
ISON300000000936TAB4900000000000	ISONIAZID + RIFAMPICIN 100 mg + 150 mg TABLET
ISON300000000937TAB4900000000000	ISONIAZID + RIFAMPICIN 150 mg + 300 mg TABLET
ISON300000000938TAB4900000000000	ISONIAZID + RIFAMPICIN 200 mg + 225 mg TABLET
ISON300000000939TAB4900000000000	ISONIAZID + RIFAMPICIN 300 mg + 450 mg TABLET
ISON300000000940TAB4900000000000	ISONIAZID + RIFAMPICIN 400 mg + 450 mg TABLET
ISON300000000941TAB4900000000000	ISONIAZID + RIFAMPICIN 600 mg + 400 mg TABLET
ISON300000000941TAB5200000000000	ISONIAZID + RIFAMPICIN 600 mg + 400 mg TABLET MODIFIED RELEASE
ISON400000000942TAB4900000000000	ISONIAZID + RIFAMPICIN + PYRAZINAMIDE 30 mg + 60 mg + 150 mg (pediatric) TABLET
ISON400000000943TAB4900000000000	ISONIAZID + RIFAMPICIN + PYRAZINAMIDE 75 mg + 150 mg + 400 mg TABLET
ISON400000000944TAB4900000000000	ISONIAZID + RIFAMPICIN + PYRAZINAMIDE 150 mg + 150 mg + 500 mg TABLET
ISON400000000945TAB4900000000000	ISONIAZID + RIFAMPICIN + PYRAZINAMIDE 300 mg + 450 mg + 500 mg TABLET
ISON500000000490TAB4900000000000	ISONIAZID + RIFAMPICIN + PYRAZINAMIDE + ETHAMBUTOL 75 mg + 150 mg + 400 mg + 275 mg TABLET
ISON500000000946TAB4900000000000	ISONIAZID + RIFAMPICIN + PYRAZINAMIDE + ETHAMBUTOL 60 mg + 120 mg + 300 mg + 225 mg TABLET
ISON500000000947TAB4900000000000	ISONIAZID + RIFAMPICIN + PYRAZINAMIDE + ETHAMBUTOL 200 mg + 450 mg + 500 mg + 400 mg TABLET
ISON600000000948TAB4900000000000	ISONIAZID + THIAZETAZONE 300 mg + 150 mg TABLET
ISONI00000000006SYRUP00201BOTTL	ISONIAZID 200 mg/5 mL SYRUP 120 mL BOTTLE
ISONI00000000006SYRUP00469BOTTL	ISONIAZID 200 mg/5 mL SYRUP 60 mL BOTTLE
ISONI00000000028TAB4900000000000	ISONIAZID 300 mg TABLET
ISONI00000000034TAB4900000000000	ISONIAZID 400 mg TABLET
ISONI00000000139TAB4900000000000	ISONIAZID 100 mg TABLET

ISONI0000000162SYRUP00201BOTTL	ISONIAZID 100 mg/5 mL SYRUP 120 mL BOTTLE
ISONI0000000162SYRUP00469BOTTL	ISONIAZID 100 mg/5 mL SYRUP 60 mL BOTTLE
ISOS10000000007TAB490000000000	ISOSORBIDE-5-MONONITRATE 20 mg TABLET
ISOS100000000031CAP420000000000	ISOSORBIDE-5-MONONITRATE 30 mg CAPSULE MODIFIED RELEASE
ISOS100000000031TAB520000000000	ISOSORBIDE-5-MONONITRATE 30 mg TABLET MODIFIED RELEASE
ISOS100000000036TAB490000000000	ISOSORBIDE-5-MONONITRATE 40 mg TABLET
ISOS100000000062CAP420000000000	ISOSORBIDE-5-MONONITRATE 60 mg CAPSULE MODIFIED RELEASE
ISOS100000000062TAB520000000000	ISOSORBIDE-5-MONONITRATE 60 mg TABLET MODIFIED RELEASE
ISOSO00000000007CAPSU0000000000	ISOSORBIDE DINITRATE 20 mg CAPSULE
ISOSO00000000007TAB490000000000	ISOSORBIDE DINITRATE 20 mg TABLET
ISOSO00000000007TAB520000000000	ISOSORBIDE DINITRATE 20 mg TABLET MODIFIED RELEASE
ISOSO000000000036CAP420000000000	ISOSORBIDE DINITRATE 40 mg CAPSULE MODIFIED RELEASE
ISOSO000000000036TAB520000000000	ISOSORBIDE DINITRATE 40 mg TABLET MODIFIED RELEASE
ISOSO000000000056TAB490000000000	ISOSORBIDE DINITRATE 5 mg TABLET
ISOSO000000000069TAB490000000000	ISOSORBIDE DINITRATE 10 mg TABLET
ISOSO00000000261SOL3200195AMPUL	ISOSORBIDE DINITRATE 1 mg/mL SOLUTION 10 mL AMPULE
ISOXS00000200036TAB490000000000	ISOXSUPRINE (as HYDROCHLORIDE) 40 mg TABLET
ISOXS00000200069TAB490000000000	ISOXSUPRINE (as HYDROCHLORIDE) 10 mg TABLET
ISOXS00000200459SOL3200361AMPUL	ISOXSUPRINE (as HYDROCHLORIDE) 5 mg/mL SOLUTION 2 mL AMPULE
ITRAC00000000139CAPSU0000000000	ITRACONAZOLE 100 mg CAPSULE
IVERM00000000063TAB490000000000	IVERMECTIN 6 mg TABLET
IVS0100000001233SUS1400000BOTTL	INFLUENZA VACCINE (SPLIT VIRION, INACTIVATED) 7.5 mcg/0.25 mL SUSPENSION BOTTLE
IVS0100000001234SOL3200000PFG01	INFLUENZA VACCINE (SPLIT VIRION, INACTIVATED) 0.25 mL (single dose) SOLUTION PRE-FILLED GLASS SYRIN
KANAM0006000241SOL3200000VIALX	KANAMYCIN (as SULFATE) 1 g SOLUTION VIAL
KETAM00000000851SOL3200195VIALX	KETAMINE 50 mg/mL SOLUTION 10 mL VIAL
KETOC00000000000LOTI100195BOTTL	KETOCONAZOLE LOTION 10 mL BOTTLE
KETOC00000000000SHAMX00178BOTTL	KETOCONAZOLE SHAMPOO 100 mL BOTTLE
KETOC00000000000SHAMX00195SAC01	KETOCONAZOLE SHAMPOO 10 mL SACHET
KETOC00000000000SHAMX00469BOTTL	KETOCONAZOLE SHAMPOO 60 mL BOTTLE
KETOC00000000000SHAMX00842SAC01	KETOCONAZOLE SHAMPOO 6 mL SACHET
KETOC00000000004TAB490000000000	KETOCONAZOLE 200 mg TABLET
KETOC00000000028OVULE0000000000	KETOCONAZOLE 300 mg OVULE
KETOC00000000047OVULE0000000000	KETOCONAZOLE 500 mg OVULE
KETOC00000000139OVULE0000000000	KETOCONAZOLE 100 mg OVULE
KETOC00000000194SOL3200231BOTTL	KETOCONAZOLE 10 mg/mL SOLUTION 15 mL BOTTLE
KETOC00000000194SOL3201375BOTTL	KETOCONAZOLE 10 mg/mL SOLUTION 25 mL BOTTLE
KETOC00000001062CREA300227TUB01	KETOCONAZOLE 1% CREAM 15 g TUBE
KETOC00000001062CREA301393TUB01	KETOCONAZOLE 1% CREAM 3.5 g TUBE
KETOC00000001062CREA301398TUB01	KETOCONAZOLE 1% CREAM 450 g TUBE
KETOC00000001314CREA300227TUB01	KETOCONAZOLE 2% CREAM 15 g TUBE
KETOC00000001314CREA301393TUB01	KETOCONAZOLE 2% CREAM 3.5 g TUBE
KETOC00000001314CREA301398TUB01	KETOCONAZOLE 2% CREAM 450 g TUBE
KETOP00000000139SOL3200000VIALX	KETOPROFEN 100 mg SOLUTION VIAL
KETOP00000000851SOL3200361VIALX	KETOPROFEN 50 mg/mL SOLUTION 2 mL VIAL
KETOR00000400785SOL3200516AMPUL	KETOROLAC (as TROMETHAMOL) 30 mg/mL SOLUTION 1 mL AMPULE
L-ASP00000000135SOL3200000VIALX	L-ASPARAGINASE 10,000 IU SOLUTION VIAL

LACTA0000000044SOL3200000BOTTL	LACTATED RINGER'S SOLUTION 500 mL SOLUTION BOTTLE
LACTA00000000247SOL3200000BOTTL	LACTATED RINGER'S SOLUTION 1 L SOLUTION BOTTLE
LACTU00000000951SYRUP00201BOTTL	LACTULOSE 3.3 g/5 mL SYRUP 120 mL BOTTLE
LACTU00000000951SYRUP01125BOTTL	LACTULOSE 3.3 g/5 mL SYRUP 200 mL BOTTLE
LACTU00000001324SYRUP00201BOTTL	LACTULOSE 3.35 g/5 mL SYRUP 120 mL BOTTLE
LACTU00000001324SYRUP01125BOTTL	LACTULOSE 3.35 g/5 mL SYRUP 200 mL BOTTLE
LAGUN00000000028TAB490000000000	LAGUNDI 300 mg TABLET
LAGUN00000000061TAB490000000000	LAGUNDI 600 mg TABLET
LAGUN00000000374SYRUP00201BOTTL	LAGUNDI 300 mg/5 mL SYRUP 120 mL BOTTLE
LAGUN00000000374SYRUP00469BOTTL	LAGUNDI 300 mg/5 mL SYRUP 60 mL BOTTLE
LAMIV00000000139TAB490000000000	LAMIVUDINE 100 mg TABLET
LAMIV00000000194SOL3200322BOTTL	LAMIVUDINE 10 mg/mL SOLUTION 240 mL BOTTLE
LAMIV00000000215TAB490000000000	LAMIVUDINE 150 mg TABLET
LAMOT00000000050TAB490000000000	LAMOTRIGINE 50 mg TABLET
LAMOT00000000139TAB490000000000	LAMOTRIGINE 100 mg TABLET
LANP10000000015TAB490000000000	LAPATINIB 250 mg TABLET
LANSO00000000031CAPSU00000000000	LANSOPRAZOLE 30 mg CAPSULE
LANSO00000000031TAB520000000000	LANSOPRAZOLE 30 mg TABLET MODIFIED RELEASE
LANSO00000000022CAPSU00000000000	LANSOPRAZOLE 15 mg CAPSULE
LANSO00000000022TAB520000000000	LANSOPRAZOLE 15 mg TABLET MODIFIED RELEASE
LAT0100000001198TAB490000000000	LAMIVUDINE + TENOFOVIR 300 mg/300 mg TABLET
LATAN00000000443EDS0201124BOTTL	LATANOPROST 50mcg/mL EYE DROPS SOLUTION 2.5 mL BOTTLE
LAZ0100000001195TAB490000000000	LAMIVUDINE + ZIDOVUDINE 150 mg/300 mg TABLET
LET0100000001197TAB490000000000	LAMIVUDINE + EFAVIRENZ + TENOFOVIR 300 mg/600 mg/300 mg TABLET
LEUPR00000000953SOL3200000VS001	LEUPRORELINE 3.75 mg/2 mL SOLUTION VIAL + SYRINGE
LEUPR00000100954POW2700000SDS01	LEUPRORELINE (as ACETATE) 1.88 mg POWDER SINGLE DOSE WITH SYRINGE
LEUPR00000100955POW2700000SDS01	LEUPRORELINE (as ACETATE) 3.75 mg POWDER SINGLE DOSE WITH SYRINGE
LEUPR00000100956SOL3900000VS001	LEUPRORELINE (as ACETATE) 11.25 mg DEPOT SOLUTION VIAL + SYRINGE
LEVO100000000152TAB490000000000	LEVODOPA + CARBIDOPA 100 mg + 25 mg TABLET
LEVO100000000335TAB490000000000	LEVODOPA + CARBIDOPA 250 mg + 25 mg TABLET
LEVO1000000001223TAB490000000000	LEVODOPA + CARBIDOPA 200 mg + 50 mg TABLET
LEVOF00000000015TAB490000000000	LEVOFLOXACIN 250 mg TABLET
LEVOF00000000047TAB490000000000	LEVOFLOXACIN 500 mg TABLET
LEVOF00000000064TAB490000000000	LEVOFLOXACIN 750 mg TABLET
LEVOF00000000459SOL3200178VIALX	LEVOFLOXACIN 5 mg/mL SOLUTION 100 mL VIAL
LEVOT000000000148TAB490000000000	LEVOTHYROXINE 100 mcg TABLET
LEVOT000000000203TAB490000000000	LEVOTHYROXINE 125 mcg TABLET
LEVOT000000000218TAB490000000000	LEVOTHYROXINE 150 mcg TABLET
LEVOT000000000341TAB490000000000	LEVOTHYROXINE 25 mcg TABLET
LEVOT000000000437TAB490000000000	LEVOTHYROXINE 50mcg TABLET
LEVOT00000000486TAB490000000000	LEVOTHYROXINE 75 mcg TABLET
LEVOT00000000508TAB490000000000	LEVOTHYROXINE 200 mcg TABLET
LEVOT00000001214TAB490000000000	LEVOTHYROXINE 112 mcg TABLET
LEVOT00000001215TAB490000000000	LEVOTHYROXINE 137 mcg TABLET
LEVOT00000001312TAB490000000000	LEVOTHYROXINE 80 mcg TABLET

LIDOC0000001062SOL3200055VIALX	LIDOCAINE 1% SOLUTION 50 mL VIAL
LIDOC0000200133SPRA300055BOTTL	LIDOCAINE (as HYDROCHLORIDE) 10% SPRAY 50 mL BOTTLE
LIDOC0000200142SOL3200060AMPUL	LIDOCAINE (as HYDROCHLORIDE) 100 mg/ml SOLUTION 5 mL AMPULE
LIDOC0000200264SOL3200000CAR01	LIDOCAINE (as HYDROCHLORIDE) 2%, 1.8mL w/ epinephrine SOLUTION CARPULE
LIDOC0000200414SOL3200020PRM01	LIDOCAINE (as HYDROCHLORIDE) 4 mg/mL SOLUTION 250 mL PRE-MIXED
LIDOC0000200658SOL4100379BOTTL	LIDOCAINE (as HYDROCHLORIDE) 4% TOPICAL SOLUTION 30 mL BOTTLE
LIDOC0000201062SOL3200011AMPUL	LIDOCAINE (as HYDROCHLORIDE) 1% SOLUTION 20 mL AMPULE
LIDOC0000201062SOL3200060AMPUL	LIDOCAINE (as HYDROCHLORIDE) 1% SOLUTION 5 mL AMPULE
LIDOC0000201314JELLY00367TUB01	LIDOCAINE (as HYDROCHLORIDE) 2% JELLY 30 g TUBE
LIDOC0000201314SOL3200011AMPUL	LIDOCAINE (as HYDROCHLORIDE) 2% SOLUTION 20 mL AMPULE
LIDOC0000201314SOL3200011VIALX	LIDOCAINE (as HYDROCHLORIDE) 2% SOLUTION 20 mL VIAL
LIDOC0000201314SOL3200055VIALX	LIDOCAINE (as HYDROCHLORIDE) 2% SOLUTION 50 mL VIAL
LIDOC0000201314SOL3200060AMPUL	LIDOCAINE (as HYDROCHLORIDE) 2% SOLUTION 5 mL AMPULE
LIDOC0000201314SOL3200060VIALX	LIDOCAINE (as HYDROCHLORIDE) 2% SOLUTION 5 mL VIAL
LIDOC0000201314SOL3200195VIALX	LIDOCAINE (as HYDROCHLORIDE) 2% SOLUTION 10 mL VIAL
LIDOC0000201314SOL3200361AMPUL	LIDOCAINE (as HYDROCHLORIDE) 2% SOLUTION 2 mL AMPULE
LIDOC0000201315OINTX01384TUB01	LIDOCAINE (as HYDROCHLORIDE) 5% OINTMENT 50 g TUBE
LIDOC0000201315OINTX01411TUB01	LIDOCAINE (as HYDROCHLORIDE) 5% OINTMENT 35 g TUBE
LIPID0000000133SOL3200044BOTTL	LIPIDS 10% SOLUTION 500 mL BOTTLE
LIPID0000000133SOL3200178BOTTL	LIPIDS 10% SOLUTION 100 mL BOTTLE
LIPID0000000133SOL3201125BOTTL	LIPIDS 10% SOLUTION 200 mL BOTTLE
LIPID0000000289SOL3200044BOTTL	LIPIDS 20% SOLUTION 500 mL BOTTLE
LIPID0000000289SOL3200178BOTTL	LIPIDS 20% SOLUTION 100 mL BOTTLE
LIQUI0000001315CONCE0134700000	LIQUID BICARBONATE CONCENTRATE 5% CONCENTRATE 5 L
LIQUI0000001315CONCE0134800000	LIQUID BICARBONATE CONCENTRATE 5% CONCENTRATE 10 L
LIQUI0000001462CONCE0134700000	LIQUID BICARBONATE CONCENTRATE 8.40% CONCENTRATE 5 L
LIQUI0000001462CONCE0134800000	LIQUID BICARBONATE CONCENTRATE 8.40% CONCENTRATE 10 L
LITHI0000000406TAB520000000000	LITHIUM CARBONATE 450 mg MR TABLET MODIFIED RELEASE
LIVE10000000967SOL3200000VIALX	LIVE ATTENUATED MEASLES VACCINE monodose vial + 0.5 mL diluent SOLUTION VIAL
LIVE10000000968SOL3200000VIALX	LIVE ATTENUATED MEASLES VACCINE multidose vial + 5 mL diluent SOLUTION VIAL
LIVE20000000967SOL3200000VIALX	LIVE ATTENUATED MEASLES, MUMPS, AND RUBELLA (MMR) VACCINE monodose vial + 0.5 mL diluent SOLUTION V
LIVE20000000968SOL3200000VIALX	LIVE ATTENUATED MEASLES, MUMPS, AND RUBELLA (MMR) VACCINE multidose vial + 5 mL diluent SOLUTION VI
LIVE30000000967SOL3200000VIALX	LIVE ATTENUATED MUMPS VACCINE monodose vial + 0.5 mL diluent SOLUTION VIAL
LIVE40000000967SOL3200000VIALX	LIVE ATTENUATED RUBELLA VACCINE monodose vial + 0.5 mL diluent SOLUTION VIAL
LIVE40000000968SOL3200000VIALX	LIVE ATTENUATED RUBELLA VACCINE multidose vial + 5 mL diluent SOLUTION VIAL
LIVE50000000094SOL3200000TUB01	LIVE ATTENUATED TRIVALENT ORAL POLIO VACCINE 0.5 mL SOLUTION TUBE
LIVE50000000094SOL3200000VIALX	LIVE ATTENUATED TRIVALENT ORAL POLIO VACCINE 0.5 mL SOLUTION VIAL
LIVE50000000969SOL3200000PLTo1	LIVE ATTENUATED TRIVALENT ORAL POLIO VACCINE 1 mL (10) doses SOLUTION PLASTIC TUBE
LIVE50000000969SOL3200000VIALX	LIVE ATTENUATED TRIVALENT ORAL POLIO VACCINE 1 mL (10) doses SOLUTION VIAL
LIVE50000000970SOL3200000PLTo1	LIVE ATTENUATED TRIVALENT ORAL POLIO VACCINE 2 mL (20 doses) SOLUTION PLASTIC TUBE

LIVE50000000970SOL3200000VIALX	LIVE ATTENUATED TRIVALENT ORAL POLIO VACCINE 2 mL (20 doses) SOLUTION VIAL
LIVE60000000971SOL3200000VA001	LIVE ATTENUATED VARICELLA VACCINE not less than 2000 PFU freeze-dried powder monodose vial + diluent
LIVER0000000122SOL3200000VIALX	LIVE ATTENUATED ROTAVIRUS 1.5 mL SOLUTION VIAL
LOMUS0000000028CAPSU00000CP001	LOMUSTINE 300 mg CAPSULE COMBO-PACK
LOPER0000200026CAPSU0000000000	LOPERAMIDE (as HYDROCHLORIDE) 2 mg CAPSULE
LORAT0000000069TAB490000000000	LORATADINE 10 mg TABLET
LORAT0000000454SYRUP00379BOTTL	LORATADINE 5 mg/5 mL SYRUP 30 mL BOTTLE
LORAZ0000000414SOL3200516AMPUL	LORAZEPAM 4 mg/mL SOLUTION 1 mL AMPULE
LOSA10000000152TAB490000000000	LOSARTAN + HYDROCHLOROTHIAZIDE 100 mg + 25 mg TABLET
LOSA10000000444TAB490000000000	LOSARTAN + HYDROCHLOROTHIAZIDE 50mg +12.5 mg TABLET
LOSA1000000146500000000000	LOSARTAN + HYDROCHLOROTHIAZIDE 100 mg + 12.5 mg TABLET
LOSAR0000000050TAB490000000000	LOSARTAN 50 mg TABLET
LOSAR0000000139TAB490000000000	LOSARTAN 100 mg TABLET
LRT010000001194CAPSU0000000000	LOPINAVIR + RITONAVIR 200 mg/50 mg CAPSULE
LRT010000001194TAB490000000000	LOPINAVIR + RITONAVIR 200 mg/50 mg TABLET
LYNES0000000045TAB49000000CYC01	LYNESTRENOL 500 mcg TABLET CYCLE
LZN010000001196TAB490000000000	LAMIVUDINE + ZIDOVUDINE + NEVIRAPINE 150 mg/300 mg/200 mg TABLET
MAGN10000000435SOL3200361AMPUL	MAGNESIUM SULFATE 500mg/mL SOLUTION 2 mL AMPULE
MAGN10000000972SOL3200011VIALX	MAGNESIUM SULFATE 250 mg/mL SOLUTION 20 mL VIAL
MAGN10000000972SOL3200055VIALX	MAGNESIUM SULFATE 250 mg/mL SOLUTION 50 mL VIAL
MAGN10000000972SOL3200195AMPUL	MAGNESIUM SULFATE 250 mg/mL SOLUTION 10 mL AMPULE
MAGN10000000972SOL3200195VIALX	MAGNESIUM SULFATE 250 mg/mL SOLUTION 10 mL VIAL
MAGN10000000972SOL3200361AMPUL	MAGNESIUM SULFATE 250 mg/mL SOLUTION 2 mL AMPULE
MAGN10000000973SOL3200195AMPUL	MAGNESIUM SULFATE 500 mg/mL SOLUTION 10 mL AMPULE
MANNI0000000289SOL3200020BOTTL	MANNITOL 20% SOLUTION 250 mL BOTTLE
MANNI0000000289SOL3200044BOTTL	MANNITOL 20% SOLUTION 500 mL BOTTLE
MEBEN0000000047TAB490000000000	MEBENDAZOLE 500 mg TABLET
MEBEN0000000139CAPSU0000000000	MEBENDAZOLE 100 mg CAPSULE
MEBEN0000000139TAB490000000000	MEBENDAZOLE 100 mg TABLET
MEBEN0000000162SUS1400379BOTTL	MEBENDAZOLE 100 mg/5 mL SUSPENSION 30 mL BOTTLE
MEBEN0000000162SUS1400469BOTTL	MEBENDAZOLE 100 mg/5 mL SUSPENSION 60 mL BOTTLE
MEBEN0000000851SUS1400195BOTTL	MEBENDAZOLE 50 mg/mL SUSPENSION 10 mL BOTTLE
MEBEV0000000139TAB490000000000	MEBEVERINE 100 mg TABLET
MEC010000000021TAB490000000000	MECLOZINE (MECLIZINE) 25 mg TABLET
MEC010000000974TABL900000000000	MECLOZINE (MECLIZINE) 12.5 mg TABLET CHEWABLE
MECOB0000000045TAB490000000000	MECOBALAMIN 500 mcg TABLET
MECOB0000000668SOL3200516AMPUL	MECOBALAMIN 500 mcg/mL SOLUTION 1 mL AMPULE
MEDRO0000000003TAB490000000000	MEDROXYPROGESTERONE 2.5 mg TABLET
MEDRO0000000015TAB490000000000	MEDROXYPROGESTERONE 250 mg TABLET
MEDRO0000000034TAB490000000000	MEDROXYPROGESTERONE 400 mg TABLET
MEDRO0000000047TAB490000000000	MEDROXYPROGESTERONE 500 mg TABLET
MEDRO0000000056TAB490000000000	MEDROXYPROGESTERONE 5 mg TABLET
MEDRO0000000069TAB490000000000	MEDROXYPROGESTERONE 10 mg TABLET
MEDRO0000000139TAB490000000000	MEDROXYPROGESTERONE 100 mg TABLET
MEDRO00000100851SOL3200388VS001	MEDROXYPROGESTERONE (as ACETATE) 50 mg/mL SOLUTION 3 mL VIAL + SYRINGE

MEDRO00002100224SOL3200516VIALX	MEDROXYPROGESTERONE (as ENANTHATE) 150 mg/mL SOLUTION 1 mL VIAL
MEFEN0000000015CAPSU0000000000	MEFENAMIC ACID 250 mg CAPSULE
MEFEN0000000015TAB490000000000	MEFENAMIC ACID 250 mg TABLET
MEFEN0000000047CAPSU0000000000	MEFENAMIC ACID 500 mg CAPSULE
MEFEN0000000047TAB490000000000	MEFENAMIC ACID 500 mg TABLET
MEFEN0000001228SUS1400469BOTTL	MEFENAMIC ACID 50 mg/5 mL SUSPENSION 60 mL BOTTLE
MEFLO0000000015TAB490000000000	MEFLOQUINE 250 mg TABLET
MEGES0000000036TAB490000000000	MEGESTROL 40 mg TABLET
MEGES00000000233TAB490000000000	MEGESTROL 160 mg TABLET
MELOX00000000228TAB490000000000	MELOXICAM 15 mg TABLET
MELOX0000001458TAB490000000000	MELOXICAM 7.5 mg TABLET
MELPH0000000026TAB490000000000	MELPHALAN 2 mg TABLET
MEMAN0000000007TAB490000000000	MEMANTINE HYDROCHLORIDE 20 mg TABLET
MEMAN0000000069TAB490000000000	MEMANTINE HYDROCHLORIDE 10 mg TABLET
MEN500000000975SOL3200000VIALX	MENINGOCOCCAL POLYSACCHARIDE (NEISSERIA MENINGITIDIS) VACCINE 50 mcg/0.5 mL dose (Group A + C) multi
MEN500000000976SOL3200000LPS01	MENINGOCOCCAL POLYSACCHARIDE (NEISSERIA MENINGITIDIS) VACCINE 50 mcg/0.5 mL dose(Group A + C) single
MEN500000000977SOL3200000LPV01	MENINGOCOCCAL POLYSACCHARIDE (NEISSERIA MENINGITIDIS) VACCINE 50 mcg/0.5 mL dose (Serogroup A + Sero
MERCA0000000050TAB490000000000	MERCAPTOPURINE 50 mg TABLET
MEROP0005100047POW2700000VIALX	MEROPENEM (as SODIUM) 500 mg POWDER VIAL
MEROP0005100241POW2700000VIALX	MEROPENEM (as SODIUM) 1 g POWDER VIAL
MESAL0000000015SUPP400000000000	MESALAZINE 250 mg SUPPOSITORY
MESAL0000000015TAB490000000000	MESALAZINE 250 mg TABLET
MESAL0000000047TAB490000000000	MESALAZINE 500 mg TABLET
MESAL00000000241SUPP400000000000	MESALAZINE 1 g SUPPOSITORY
MESN10000000142SOL3200060AMPUL	MESNA (SODIUM-2-MERCAPTOETHANE SULPHONATE) 100 mg/ml SOLUTION 5 mL AMPULE
MESN10000000142SOL3200195AMPUL	MESNA (SODIUM-2-MERCAPTOETHANE SULPHONATE) 100 mg/ml SOLUTION 10 mL AMPULE
MESN10000000142SOL3201386AMPUL	MESNA (SODIUM-2-MERCAPTOETHANE SULPHONATE) 100 mg/ml SOLUTION 4 mL AMPULE
MESNA00000000340000000000000000	MESNA 400 mg
METFO0000200047TAB490000000000	METFORMIN (as HYDROCHLORIDE) 500 mg TABLET
METFO0000200241TAB490000000000	METFORMIN (as HYDROCHLORIDE) 1 g TABLET
METFO0000200978TAB490000000000	METFORMIN (as HYDROCHLORIDE) 850 mg TABLET
METHO0000000003TAB490000000000	METHOTREXATE 2.5 mg TABLET
METHO0000000056TAB490000000000	METHOTREXATE 5 mg TABLET
METHO0000000069TAB490000000000	METHOTREXATE 10 mg TABLET
METHO0000000142SOL3200055VIALX	METHOTREXATE 100 mg/ml SOLUTION 50 mL VIAL
METHO0000000142SOL3200195VIALX	METHOTREXATE 100 mg/ml SOLUTION 10 mL VIAL
METHO0000000194SOL3200060VIALX	METHOTREXATE 10 mg/mL SOLUTION 5 mL VIAL
METHO0000000194SOL3200516VIALX	METHOTREXATE 10 mg/mL SOLUTION 1 mL VIAL
METHO0000000459SOL3200361VIALX	METHOTREXATE 5 mg/mL SOLUTION 2 mL VIAL
METHO0000000459SOL3200388VIALX	METHOTREXATE 5 mg/mL SOLUTION 3 mL VIAL
METHO0000000851SOL3200516VIALX	METHOTREXATE 50 mg/mL SOLUTION 1 mL VIAL
METHO0000000980SOL3200361AMPUL	METHOTREXATE 2.5 mg/ml SOLUTION 2 mL AMPULE
METHO0000000980SOL3200361VIALX	METHOTREXATE 2.5 mg/ml SOLUTION 2 mL VIAL
METHO0000000980SOL3201401VIALX	METHOTREXATE 2.5 mg/ml SOLUTION 8 mL VIAL

METHO0000000985SOL3200011VIALX	METHOTREXATE 25 mg/mL SOLUTION 20 mL VIAL
METHO0000000985SOL3200361AMPUL	METHOTREXATE 25 mg/mL SOLUTION 2 mL AMPULE
METHO0000000985SOL3200361VIALX	METHOTREXATE 25 mg/mL SOLUTION 2 mL VIAL
METHO0000000985SOL3201386VIALX	METHOTREXATE 25 mg/mL SOLUTION 4 mL VIAL
METHO0000000985SOL3201401VIALX	METHOTREXATE 25 mg/mL SOLUTION 8 mL VIAL
METHO00005101453000000000000000	METHOTREXATE (as SODIUM) 50 mg/2 mL
METHY0000000015TAB490000000000	METHYLDOPA 250 mg TABLET
METHY00000000205TAB490000000000	METHYLDOPA 125 mg TABLET
METOC0000200069TAB490000000000	METOCLOPRAMIDE (as HYDROCHLORIDE) 10 mg TABLET
METOC0000200454SYRUP00469BOTTLE	METOCLOPRAMIDE (as HYDROCHLORIDE) 5 mg/5 mL SYRUP 60 mL BOTTLE
METOC0000200459SOL3200361AMPUL	METOCLOPRAMIDE (as HYDROCHLORIDE) 5 mg/mL SOLUTION 2 mL AMPULE
METOC0000200459SOL3200388AMPUL	METOCLOPRAMIDE (as HYDROCHLORIDE) 5 mg/mL SOLUTION 3 mL AMPULE
METOP00005600050TAB490000000000	METOPROLOL (as TARTRATE) 50 mg TABLET
METOP00005600139TAB490000000000	METOPROLOL (as TARTRATE) 100 mg TABLET
METP200000000407TAB490000000000	METOPROLOL SUCCINATE 47.5 mg TABLET
METRO0000000015TAB490000000000	METRONIDAZOLE 250 mg TABLET
METRO0000000047SOL3200178VIALX	METRONIDAZOLE 500 mg SOLUTION 100 mL VIAL
METRO0000000047TAB490000000000	METRONIDAZOLE 500 mg TABLET
METRO00000000241SUPP400000000000	METRONIDAZOLE 1 g SUPPOSITORY
METRO00000000459SOL3200178VIALX	METRONIDAZOLE 5 mg/mL SOLUTION 100 mL VIAL
METRO00000800209SUS1400379BOTTLE	METRONIDAZOLE (as BENZOATE) 125 mg/5 mL SUSPENSION 30 mL BOTTLE
METRO00000800209SUS1400469BOTTLE	METRONIDAZOLE (as BENZOATE) 125 mg/5 mL SUSPENSION 60 mL BOTTLE
METY10000000194SOL3200195AMPUL	METHYLENE BLUE 10 mg/mL SOLUTION 10 mL AMPULE
METY10000000194SOL3200195VIALX	METHYLENE BLUE 10 mg/mL SOLUTION 10 mL VIAL
METY10000000194SOL3200516AMPUL	METHYLENE BLUE 10 mg/mL SOLUTION 1 mL AMPULE
METY10000000194SOL3200516VIALX	METHYLENE BLUE 10 mg/mL SOLUTION 1 mL VIAL
METY10000000989TAB490000000000	METHYLENE BLUE 1% solution 55 mg TABLET
METY10000000990TAB490000000000	METHYLENE BLUE 1% solution 65 mg TABLET
METY20003400203TAB490000000000	METHYLERGOMETRINE (as MALEATE) 125 mcg TABLET
METY20003401018SOL3200516AMPUL	METHYLERGOMETRINE (as MALEATE) 200 mcg/mL SOLUTION 1 mL AMPULE
METY40000000991TAB520000000000	METHYLPHENIDATE 18 mg TABLET MODIFIED RELEASE
METY40000000992TAB520000000000	METHYLPHENIDATE 36 mg TABLET MODIFIED RELEASE
METY50000000039SUS1400516VIALX	METHYLPREDNISOLONE 40 mg/mL SUSPENSION 1 mL VIAL
METY50000000041TAB490000000000	METHYLPREDNISOLONE 4 mg TABLET
METY50000000047LYP0100000VIALX	METHYLPREDNISOLONE 500 mg LYOPHILIZED POWDER VIAL
METY50000000235TAB490000000000	METHYLPREDNISOLONE 16 mg TABLET
METY50000000241LYP0100000VIALX	METHYLPREDNISOLONE 1 g LYOPHILIZED POWDER VIAL
METY50000000243POW2700000VIALX	METHYLPREDNISOLONE 1 g/16 mL + Diluent POWDER VIAL
METY50000000434POW2700000VIALX	METHYLPREDNISOLONE 500mg/7.7mL + Diluent POWDER VIAL
METY50000000994POW2701386VIALX	METHYLPREDNISOLONE 62.5 mg/mL POWDER 4 mL VIAL
METY50000000998POW2700000VIALX	METHYLPREDNISOLONE 500 mg/8.0 mL + diluent vial POWDER VIAL
METY50000001221POW2700000VIALX	METHYLPREDNISOLONE 40 mg (single dose) POWDER VIAL
METY50000001368POW270036100000	METHYLPREDNISOLONE 125 mg/mL + diluent POWDER 2 mL
METY50000001369POW2700361VIALX	METHYLPREDNISOLONE 62.5 mg/mL + diluent POWDER 2 mL VIAL
MICO1000000000000000000195SAC01	MICONAZOLE 10 mL SACHET

MICO100000000000LOT1100195BOTTL	MICONAZOLE LOTION 10 mL BOTTLE
MICO100000000000SHAMX00178BOTTL	MICONAZOLE SHAMPOO 100 mL BOTTLE
MICO100000000000SHAMX00469BOTTL	MICONAZOLE SHAMPOO 60 mL BOTTLE
MICO100000000000SHAMX00842SAC01	MICONAZOLE SHAMPOO 6 mL SACHET
MICO100000000028OVULE00000000000	MICONAZOLE 300 mg OVULE
MICO100000000047OVULE00000000000	MICONAZOLE 500 mg OVULE
MICO100000000139OVULE00000000000	MICONAZOLE 100 mg OVULE
MICO100000000194SOL3200231BOTTL	MICONAZOLE 10 mg/mL SOLUTION 15 mL BOTTLE
MICO100000000194SOL3201375BOTTL	MICONAZOLE 10 mg/mL SOLUTION 25 mL BOTTLE
MICO100000001062CREA300227TUB01	MICONAZOLE 1% CREAM 15 g TUBE
MICO100000001062CREA301393TUB01	MICONAZOLE 1% CREAM 3.5 g TUBE
MICO100000001062CREA301398TUB01	MICONAZOLE 1% CREAM 450 g TUBE
MICO100000001314CREA300227TUB01	MICONAZOLE 2% CREAM 15 g TUBE
MICO100000001314CREA301393TUB01	MICONAZOLE 2% CREAM 3.5 g TUBE
MICO100000001314CREA301398TUB01	MICONAZOLE 2% CREAM 450 g TUBE
MICRO00000000241POW2700000SAC01	MICRONUTRIENT 1 g POWDER SACHET
MIDAZ00000000228TAB4900000000000	MIDAZOLAM 15 mg TABLET
MIDAZ00000000261SOL3200060AMPUL	MIDAZOLAM 1 mg/mL SOLUTION 5 mL AMPULE
MIDAZ00000000261SOL3200060VIALX	MIDAZOLAM 1 mg/mL SOLUTION 5 mL VIAL
MIDAZ00000000459SOL3200060AMPUL	MIDAZOLAM 5 mg/mL SOLUTION 5 mL AMPULE
MIDAZ00000000459SOL3200195AMPUL	MIDAZOLAM 5 mg/mL SOLUTION 10 mL AMPULE
MIDAZ00000000459SOL3200361AMPUL	MIDAZOLAM 5 mg/mL SOLUTION 2 mL AMPULE
MIDAZ00000000459SOL3200388AMPUL	MIDAZOLAM 5 mg/mL SOLUTION 3 mL AMPULE
MIDAZ00000000459SOL3200516AMPUL	MIDAZOLAM 5 mg/mL SOLUTION 1 mL AMPULE
MINOC00000000050CAPSU00000000000	MINOCYCLINE 50 mg CAPSULE
MINOC00000000139CAPSU00000000000	MINOCYCLINE 100 mg CAPSULE
MITOM00000000069SOL3200000VIALX	MITOMYCIN 10 mg SOLUTION VIAL
MITOX00000000357SOL3200060VIALX	MITOXANTRONE 2 mg/ml SOLUTION 5 mL VIAL
MITOX00000000357SOL3200195VIALX	MITOXANTRONE 2 mg/ml SOLUTION 10 mL VIAL
MODIF00000000364SOL3200044BOTTL	MODIFIED FLUID GELATIN (FORMERLY POLYMERISATE OF DEGRADED SUCCINYLATED GELATIN) 3% SOLUTION 500 mL B
MODIF00000000658SOL3200044BOTTL	MODIFIED FLUID GELATIN (FORMERLY POLYMERISATE OF DEGRADED SUCCINYLATED GELATIN) 4% SOLUTION 500 mL B
MOL0100000000218SOL3200000VIALX	MOLGRAMOSTIN 150 mcg SOLUTION VIAL
MOL0100000000808SOL3200000VIALX	MOLGRAMOSTIN 400 mcg SOLUTION VIAL
MONOB00000000240SOL3201421BOTTL	MONOBASIC SODIUM PHOSPHATE, DIBASIC SODIUM PHOSPHATE 19 g/7 g SOLUTION 133 mL BOTTLE
MONOB00000000240SOL3201422BOTTL	MONOBASIC SODIUM PHOSPHATE, DIBASIC SODIUM PHOSPHATE 19 g/7 g SOLUTION 66 mL BOTTLE
MONOB00000000409SOL3201390BOTTL	MONOBASIC SODIUM PHOSPHATE, DIBASIC SODIUM PHOSPHATE 48 g/18 g SOLUTION 45 mL BOTTLE
MONTE0005100041SAC0100000000000	MONTELUKAST (as SODIUM) 4 mg SACHET
MONTE0005100041TABL9000000000000	MONTELUKAST (as SODIUM) 4 mg TABLET CHEWABLE
MONTE0005100056TABL9000000000000	MONTELUKAST (as SODIUM) 5 mg TABLET CHEWABLE
MONTE0005100069TAB4900000000000	MONTELUKAST (as SODIUM) 10 mg TABLET
MORPH0006000007CAPSU00000000000	MORPHINE (as SULFATE) 20 mg CAPSULE
MORPH0006000007TAB4900000000000	MORPHINE (as SULFATE) 20 mg TABLET
MORPH0006000031CAPSU00000000000	MORPHINE (as SULFATE) 30 mg CAPSULE
MORPH0006000031TAB4900000000000	MORPHINE (as SULFATE) 30 mg TABLET
MORPH0006000062TAB5200000000000	MORPHINE (as SULFATE) 60 mg TABLET MODIFIED RELEASE

MORPH0006000069CAPSU0000000000	MORPHINE (as SULFATE) 10 mg CAPSULE
MORPH0006000069TAB490000000000	MORPHINE (as SULFATE) 10 mg TABLET
MORPH0006000139TAB520000000000	MORPHINE (as SULFATE) 100 mg TABLET MODIFIED RELEASE
MORPH0006000194SOL3200516AMPUL	MORPHINE (as SULFATE) 10 mg/mL SOLUTION 1 mL AMPULE
MORPH00060001004SOL3200516AMPUL	MORPHINE (as SULFATE) 15 mg/mL SOLUTION 1 mL AMPULE
MULT2000000000000000000195VIALX	MULTIVITAMINS 10 mL VIAL
MULT200000000000CAPSU0000000000	MULTIVITAMINS CAPSULE
MULT200000000000EMUL500195AMPUL	MULTIVITAMINS EMULSION 10 mL AMPULE
MULT200000000000FDP0100195VIALX	MULTIVITAMINS FREEZE-DRIED POWDER + DILUENT 10 mL VIAL
MULT200000000000OD00000195BOTTL	MULTIVITAMINS ORAL DROPS 10 mL BOTTLE
MULT200000000000OD00000379BOTTL	MULTIVITAMINS ORAL DROPS 30 mL BOTTLE
MULT200000000000TAB490000000000	MULTIVITAMINS TABLET
MULT200000000516DROPS00231BOTTL	MULTIVITAMINS 1 mL DROPS 15 mL BOTTLE
MULT200000001256SYRUP00201BOTTL	MULTIVITAMINS per 5 mL SYRUP 120 mL BOTTLE
MULT200000001256SYRUP00469BOTTL	MULTIVITAMINS per 5 mL SYRUP 60 mL BOTTLE
MULT200000001370EMUL500195AMPUL	MULTIVITAMINS Pedia EMULSION 10 mL AMPULE
MUPIR00000001314CREA300227TUB01	MUPIROCIN 2% CREAM 15 g TUBE
MUPIR00000001314CREA301245SAC01	MUPIROCIN 2% CREAM 5 g SACHET
MUPIR00000001314OINTX00227TUB01	MUPIROCIN 2% OINTMENT 15 g TUBE
MUPIR00000001314OINTX01245TUB01	MUPIROCIN 2% OINTMENT 5 g TUBE
MUPOI00000001314OINTX00805TUB01	MUPIROCIN 2% OINTMENT 2.5 g TUBE
MYA010004000236TAB490000000000	MYCOPHENOLIC ACID (as MYCOPHENOLATE SODIUM) 180 mg TABLET
MYA010004000383TAB490000000000	MYCOPHENOLIC ACID (as MYCOPHENOLATE SODIUM) 360 mg TABLET
MYCOP00000000047TAB490000000000	MYCOPHENOLATE MOFETIL 500 mg TABLET
N-ACE00000000000CRY0001063BOTTL	N-ACETYL PENICILLAMINE CRYSTALS 25 g BOTTLE
NADRO00000001455000000000000000	NADROPARIN 3,800 I.U. / 0.4 mL
NADRO00001201006SOL3200516PFS01	NADROPARIN (as CALCIUM) 950 IU SOLUTION 1 mL PRE-FILLED SYRINGE
NADRO00001201006SOL3201396PFS01	NADROPARIN (as CALCIUM) 950 IU SOLUTION 0.4 mL PRE-FILLED SYRINGE
NADRO00001201006SOL3201404PFS01	NADROPARIN (as CALCIUM) 950 IU SOLUTION 0.2 mL PRE-FILLED SYRINGE
NADRO00001201006SOL3201405PFS01	NADROPARIN (as CALCIUM) 950 IU SOLUTION 0.6 mL PRE-FILLED SYRINGE
NADRO00001201006SOL3201406PFS01	NADROPARIN (as CALCIUM) 950 IU SOLUTION 0.8 mL PRE-FILLED SYRINGE
NADRO00001201006SOL3201423PFS01	NADROPARIN (as CALCIUM) 950 IU SOLUTION 0.1 mL PRE-FILLED SYRINGE
NADRO00001201006SOL3201424PFS01	NADROPARIN (as CALCIUM) 950 IU SOLUTION 0.3 mL PRE-FILLED SYRINGE
NALBU00000200194SOL3200516AMPUL	NALBUPHINE (as HYDROCHLORIDE) 10 mg/mL SOLUTION 1 mL AMPULE
NALBU00000201199SOL3200195VIALX	NALBUPHINE (as HYDROCHLORIDE) 20 mg/mL SOLUTION 10 mL VIAL
NALID00000000047TAB490000000000	NALIDIXIC ACID 500 mg TABLET
NALID00000001265SUS1400469BOTTL	NALIDIXIC ACID 250 mg/5 mL SUSPENSION 60 mL BOTTLE
NALOX00000200398SOL3200516AMPUL	NALOXONE (as HYDROCHLORIDE) 400 mcg/mL SOLUTION 1 mL AMPULE
NALOX00000200824SOL3200361AMPUL	NALOXONE (as HYDROCHLORIDE) 20 mcg/mL SOLUTION 2 mL AMPULE
NALTR00000000050TAB490000000000	NALTREXONE 50 mg TABLET
NAPRO00005101287TAB490000000000	NAPROXEN (as SODIUM) 275 mg TABLET
NAPRO00005101288TAB490000000000	NAPROXEN (as SODIUM) 550 mg TABLET

NELFI0000001013POW2701425BOTTL	NELFINAVIR 50 mg/scoop (1 g) POWDER 144 g BOTTLE
NELFI0003700015TAB240000000000	NELFINAVIR (as MESILATE) 250 mg TABLET FILM COATED
NELFI0003700015TAB490000000000	NELFINAVIR (as MESILATE) 250 mg TABLET
NEOST0000000668SOL3200516AMPUL	NEOSTIGMINE 500 mcg/mL SOLUTION 1 mL AMPULE
NEOST0001100228TAB490000000000	NEOSTIGMINE (as BROMIDE) 15 mg TABLET
NEPAF0000000261EDS0100060BOTTL	NEPAFENAC 1 mg/mL EYE DROPS SUSPENSION 5 mL BOTTLE
NETIL0000000142SOL3200122AMPUL	NETILMICIN 100 mg/ml SOLUTION 1.5 mL AMPULE
NETIL0000000142SOL3200122VIALX	NETILMICIN 100 mg/ml SOLUTION 1.5 mL VIAL
NETIL0000000142SOL3200361AMPUL	NETILMICIN 100 mg/ml SOLUTION 2 mL AMPULE
NETIL0000000142SOL3200361VIALX	NETILMICIN 100 mg/ml SOLUTION 2 mL VIAL
NETIL0000000851SOL3200361VIALX	NETILMICIN 50 mg/mL SOLUTION 2 mL VIAL
NETIL0000000985SOL3200361AMPUL	NETILMICIN 25 mg/mL SOLUTION 2 mL AMPULE
NETIL00000001016SOL3200361AMPUL	NETILMICIN 75 mg/mL SOLUTION 2 mL AMPULE
NEVIR0000000004TAB490000000000	NEVIRAPINE 200 mg TABLET
NEVIR00000001228SUS1400322BOTTL	NEVIRAPINE 50 mg/5 mL SUSPENSION 240 mL BOTTLE
NICAR0000000007TAB490000000000	NICARDIPINE 20 mg TABLET
NICAR0000000069TAB490000000000	NICARDIPINE 10 mg TABLET
NICAR0000000261SOL3200195AMPUL	NICARDIPINE 1 mg/mL SOLUTION 10 mL AMPULE
NICAR0000000261SOL3200361AMPUL	NICARDIPINE 1 mg/mL SOLUTION 2 mL AMPULE
NICB30000000050TAB490000000000	NICOTINAMIDE (VITAMIN B3) 50 mg TABLET
NICB30000000139TAB490000000000	NICOTINAMIDE (VITAMIN B3) 100 mg TABLET
NIFED0000000056CAPSU0000000000	NIFEDIPINE 5 mg CAPSULE
NIFED0000000069CAPSU0000000000	NIFEDIPINE 10 mg CAPSULE
NIFED0000000538TAB490000000000	NIFEDIPINE 30 mg MR TABLET
NIFED0000000949TAB490000000000	NIFEDIPINE 20 mg MR TABLET
NIMOD0000000031TAB490000000000	NIMODIPINE 30 mg TABLET
NIMOD0000000069SOL3200055VIALX	NIMODIPINE 10 mg SOLUTION 50 mL VIAL
NIMOD00000001018SOL3200055VIALX	NIMODIPINE 200 mcg/mL SOLUTION 50 mL VIAL
NITRO0003300050CAPSU0000000000	NITROFURANTOIN (as MACROCRYSTALS) 50 mg CAPSULE
NITRO0003300139CAPSU0000000000	NITROFURANTOIN (as MACROCRYSTALS) 100 mg CAPSULE
NORE10000000261SOL3200195AMPUL	NOREPINEPHRINE 1 mg/mL SOLUTION 10 mL AMPULE
NORE10000000261SOL3201386AMPUL	NOREPINEPHRINE 1 mg/mL SOLUTION 4 mL AMPULE
NORE10001000261SOL3200361AMPUL	NOREPINEPHRINE (as BITARTRATE) 1 mg/mL SOLUTION 2 mL AMPULE
NORET0000100056TAB490000000000	NORETHISTERONE (as ACETATE) 5 mg TABLET
NORET0002101123SOL3200516AMPUL	NORETHISTERONE (as ENANTHATE) 200 mg/mL SOLUTION 1 mL AMPULE
NORFL0000000004TAB490000000000	NORFLOXACIN 200 mg TABLET
NORFL0000000034TAB490000000000	NORFLOXACIN 400 mg TABLET
NPFA10000001014SOL3700060BOTTL	NEOMYCIN + POLYMYXIN B + FLUOCINOLONE ACETONIDE EAR DROPS 3.5 mg neomycin + 10,000 units polymyxin B
NYSTA00000001021SUS1400379BOTTL	NYSTATIN 100,000 units/mL SUSPENSION 30 mL BOTTLE
NYSTA00000001022TAB490000000000	NYSTATIN 500,000 units TABLET
NYSTA00000001023TAB530000000000	NYSTATIN 100,000 units VAGINAL TABLET
OCTRE0000000668SOL3200516AMPUL	OCTREOTIDE 500 mcg/mL SOLUTION 1 mL AMPULE
OCTRE0000000767SOL3200516AMPUL	OCTREOTIDE 100 mcg/mL SOLUTION 1 mL AMPULE
OFLOX0000000004TAB490000000000	OFLOXACIN 200 mg TABLET
OFLOX0000000034TAB490000000000	OFLOXACIN 400 mg TABLET
OFLOX0000000089EYD0100060BOTTL	OFLOXACIN 0.30% EYE DROPS 5 mL BOTTLE

OFLOX0000000089EYO0101393TUB01	OFLOXACIN 0.30% EYE OINTMENT 3.5 g TUBE
OFLOX00000000357SOL3200178VIALX	OFLOXACIN 2 mg/ml SOLUTION 100 mL VIAL
OLANZ00000000056TAB490000000000	OLANZAPINE 5 mg TABLET
OLANZ00000000069POW2700000VIALX	OLANZAPINE 10 mg POWDER VIAL
OLANZ00000000069TAB490000000000	OLANZAPINE 10 mg TABLET
OND010000200357SOL3200361AMPUL	ONDANSETRON (as HYDROCHLORIDE) 2 mg/ml SOLUTION 2 mL AMPULE
OND010000200357SOL3201386AMPUL	ONDANSETRON (as HYDROCHLORIDE) 2 mg/ml SOLUTION 4 mL AMPULE
OND010003000068TAB490000000000	ONDANSETRON (as HYDROCHLORIDE DIHYDRATE) 8 mg TABLET
ORAL200000000000POW2701273SAC01	ORAL REHYDRATION SALTS (ORS 75-REPLACEMENT) POWDER 20.5 g SACHET
ORAL200000000000POW2701279SAC01	ORAL REHYDRATION SALTS (ORS 75-REPLACEMENT) POWDER 4.1 g SACHET
ORAL200000000000POW2701323SAC01	ORAL REHYDRATION SALTS (ORS 75-REPLACEMENT) POWDER 5.125 g SACHET
ORAL200000000000POW2701426SAC01	ORAL REHYDRATION SALTS (ORS 75-REPLACEMENT) POWDER 2.17 g SACHET
ORAL200000000000SOL3200020BOTTL	ORAL REHYDRATION SALTS (ORS 75-REPLACEMENT) SOLUTION 250 mL BOTTLE
ORAL200000000483POW2700000SAC01	ORAL REHYDRATION SALTS (ORS 75-REPLACEMENT) 75 replacement POWDER SACHET
OSELT0000000065CAPSU0000000000	OSELTAMIVIR 75 mg CAPSULE
OSELT00000001026POW2400469BOTTL	OSELTAMIVIR 12 mg/mL POWDER FOR SUSPENSION 60 mL BOTTLE
OXACI0005300015SOL3200000VIALX	OXACILLIN (as SODIUM SALT) 250 mg SOLUTION VIAL
OXACI0005300047SOL3200000VIALX	OXACILLIN (as SODIUM SALT) 500 mg SOLUTION VIAL
OXALI00000000357SOL3200055VIALX	OXALIPLATIN 2 mg/ml SOLUTION 50 mL VIAL
OXALI00000000357SOL3201375VIALX	OXALIPLATIN 2 mg/ml SOLUTION 25 mL VIAL
OXALI00000000459SOLU200011VIALX	OXALIPLATIN 5 mg/mL SOLUTION CONCENTRATE 20 mL VIAL
OXALI00000000459SOLU200195VIALX	OXALIPLATIN 5 mg/mL SOLUTION CONCENTRATE 10 mL VIAL
OXALI00000000459SOLU201427VIALX	OXALIPLATIN 5 mg/mL SOLUTION CONCENTRATE 40 mL VIAL
OXALI00000001032SOL3200000VIALX	OXALIPLATIN 50 mg powder SOLUTION VIAL
OXANT00000001027TAB490000000000	OXANTEL + PYRANTEL 100 mg oxantel + 100 mg pyrantel TABLET
OXANT00000001028SUS1400195BOTTL	OXANTEL + PYRANTEL 100 mg oxantel + 100 mg pyrantel per 5 mL SUSPENSION 10 mL BOTTLE
OXCAR0000000028TAB490000000000	OXCARBAZEPINE 300 mg TABLET
OXCAR0000000061TAB490000000000	OXCARBAZEPINE 600 mg TABLET
OXCAR00000001326SUS1400020BOTTL	OXCARBAZEPINE 60 mg/mL SUSPENSION 250 mL BOTTLE
OXCAR00000001326SUS1400178BOTTL	OXCARBAZEPINE 60 mg/mL SUSPENSION 100 mL BOTTLE
OXYCO0000000007TAB490000000000	OXYCODONE 20 mg TABLET
OXYCO0000000036TAB490000000000	OXYCODONE 40 mg TABLET
OXYCO0000000056CAPSU0000000000	OXYCODONE 5 mg CAPSULE
OXYCO0000000066TAB490000000000	OXYCODONE 80 mg TABLET
OXYCO0000000069TAB490000000000	OXYCODONE 10 mg TABLET
OXYCO00000000194SOL3200361AMPUL	OXYCODONE 10 mg/mL SOLUTION 2 mL AMPULE
OXYCO00000000194SOL3200516AMPUL	OXYCODONE 10 mg/mL SOLUTION 1 mL AMPULE
OXYME0000200078NS00100195BOTTL	OXYMETAZOLINE (as HYDROCHLORIDE) 0.05% NASAL SPRAY 10 mL BOTTLE
OXYME0000200078NS00100231BOTTL	OXYMETAZOLINE (as HYDROCHLORIDE) 0.05% NASAL SPRAY 15 mL BOTTLE
OXYME0000201033NDS0100231BOTTL	OXYMETAZOLINE (as HYDROCHLORIDE) 0.03% NASAL DROPS SOLUTION 15 mL BOTTLE
OXYTO00000000124SOL3200516AMPUL	OXYTOCIN 10 IU/mL SOLUTION 1 mL AMPULE
OXYTO00000001036SOL3200516AMPUL	OXYTOCIN 5 IU/mL SOLUTION 1 mL AMPULE

PACLI00000000474SOL3200060VIALX	PACLITAXEL 6 mg/mL SOLUTION 5 mL VIAL
PACLI00000000474SOL3201375VIALX	PACLITAXEL 6 mg/mL SOLUTION 25 mL VIAL
PACLI00000000474SOL3201428VIALX	PACLITAXEL 6 mg/mL SOLUTION 16.7 mL VIAL
PACLI00000000474SOL3201429VIALX	PACLITAXEL 6 mg/mL SOLUTION 17 mL VIAL
PACLI00000000474SOL3201430VIALX	PACLITAXEL 6 mg/mL SOLUTION 43.4 mL VIAL
PACLI00000001040SOL3200231VIALX	PACLITAXEL 100 mg/15 mL SOLUTION 15 mL VIAL
PANCU00001100357SOL3200361AMPUL	PANCURONIUM (as BROMIDE) 2 mg/ml SOLUTION 2 mL AMPULE
PARAC00000000015SUPP40000000000	PARACETAMOL 250 mg SUPPOSITORY
PARAC00000000028TAB490000000000	PARACETAMOL 300 mg TABLET
PARAC00000000047TAB490000000000	PARACETAMOL 500 mg TABLET
PARAC00000000142DROPS00231BOTTL	PARACETAMOL 100 mg/ml DROPS 15 mL BOTTLE
PARAC00000000194SOL3200055VIALX	PARACETAMOL 10 mg/mL SOLUTION 50 mL VIAL
PARAC00000000194SOL3200178VIALX	PARACETAMOL 10 mg/mL SOLUTION 100 mL VIAL
PARAC00000000205SUPP40000000000	PARACETAMOL 125 mg SUPPOSITORY
PARAC00000000224SOL3200361AMPUL	PARACETAMOL 150 mg/mL SOLUTION 2 mL AMPULE
PARAC00000001042SUS1400201BOTTL	PARACETAMOL 120 mg/5 mL (125 mg/5 mL) SUSPENSION 120 mL BOTTLE
PARAC00000001042SUS1400379BOTTL	PARACETAMOL 120 mg/5 mL (125 mg/5 mL) SUSPENSION 30 mL BOTTLE
PARAC00000001042SYRUP00201BOTTL	PARACETAMOL 120 mg/5 mL (125 mg/5 mL) SYRUP 120 mL BOTTLE
PARAC00000001042SYRUP00379BOTTL	PARACETAMOL 120 mg/5 mL (125 mg/5 mL) SYRUP 30 mL BOTTLE
PARAC00000001265SUS1400201BOTTL	PARACETAMOL 250 mg/5 mL SUSPENSION 120 mL BOTTLE
PARAC00000001265SUS1400379BOTTL	PARACETAMOL 250 mg/5 mL SUSPENSION 30 mL BOTTLE
PARAC00000001265SYRUP00201BOTTL	PARACETAMOL 250 mg/5 mL SYRUP 120 mL BOTTLE
PARAC00000001265SYRUP00379BOTTL	PARACETAMOL 250 mg/5 mL SYRUP 30 mL BOTTLE
PARAC00000001265SYRUP00469BOTTL	PARACETAMOL 250 mg/5 mL SYRUP 60 mL BOTTLE
PARO10000000007TAB490000000000	PAROXETINE 20 mg TABLET
PEG2A00000001044SOL3200000PFS01	PEGINTERFERON ALFA 2A 135 mcg/0.5 mL SOLUTION PRE-FILLED SYRINGE
PEG2A00000001045SOL3200000PFS01	PEGINTERFERON ALFA 2A 180 mcg/0.5 mL SOLUTION PRE-FILLED SYRINGE
PENB100000001046SOL3200000VIALX	PENICILLIN G BENZATHINE (BENZATHINE BENZYL PENICILLIN) 1,200,000 units SOLUTION VIAL
PENB100000001047SOL3200000VIALX	PENICILLIN G BENZATHINE (BENZATHINE BENZYL PENICILLIN) 2,400,000 units SOLUTION VIAL
PENI100000000263SOL3200000VIALX	PENICILLIN G CRYSTALLINE 1 MU SOLUTION VIAL
PENI10000000046000000000000VIALX	PENICILLIN G CRYSTALLINE 5 MU VIAL
PENI100000001022SOL3200000VIALX	PENICILLIN G CRYSTALLINE 500,000 units SOLUTION VIAL
PERIT00000001051SOL3200247BOTTL	PERITONEAL DIALYSIS SOLUTION Sterile with 1.5% dextrose SOLUTION 1 L BOTTLE
PERIT00000001051SOL3201347BOTTL	PERITONEAL DIALYSIS SOLUTION Sterile with 1.5% dextrose SOLUTION 5 L BOTTLE
PERIT00000001051SOL3201431BOTTL	PERITONEAL DIALYSIS SOLUTION Sterile with 1.5% dextrose SOLUTION 1.5 L BOTTLE
PERIT00000001051SOL3201432BOTTL	PERITONEAL DIALYSIS SOLUTION Sterile with 1.5% dextrose SOLUTION 2 L BOTTLE
PERIT00000001055SOL3200247BOTTL	PERITONEAL DIALYSIS SOLUTION Sterile with 2.5% dextrose SOLUTION 1 L BOTTLE
PERIT00000001055SOL3201347BOTTL	PERITONEAL DIALYSIS SOLUTION Sterile with 2.5% dextrose SOLUTION 5 L BOTTLE
PERIT00000001055SOL3201431BOTTL	PERITONEAL DIALYSIS SOLUTION Sterile with 2.5% dextrose SOLUTION 1.5 L BOTTLE
PERIT00000001055SOL3201432BOTTL	PERITONEAL DIALYSIS SOLUTION Sterile with 2.5% dextrose SOLUTION 2 L BOTTLE
PERIT00000001059SOL3200247BOTTL	PERITONEAL DIALYSIS SOLUTION Sterile with 4.25% dextrose SOLUTION 1 L BOTTLE

PERIT0000001059SOL3201347BOTTL	PERITONEAL DIALYSIS SOLUTION Sterile with 4.25% dextrose SOLUTION 5 L BOTTLE
PERIT0000001059SOL3201431BOTTL	PERITONEAL DIALYSIS SOLUTION Sterile with 4.25% dextrose SOLUTION 1.5 L BOTTLE
PERIT0000001059SOL3201432BOTTL	PERITONEAL DIALYSIS SOLUTION Sterile with 4.25% dextrose SOLUTION 2 L BOTTLE
PERME0000001062LOT1101433BOTTL	PERMETHRIN 1% LOTION 125 mL BOTTLE
PERME0000001062SHAMX00379BOTTL	PERMETHRIN 1% SHAMPOO 30 mL BOTTLE
PERME0000001062SHAMX00469BOTTL	PERMETHRIN 1% SHAMPOO 60 mL BOTTLE
PERME0000001315LOT1100379BOTTL	PERMETHRIN 5% LOTION 30 mL BOTTLE
PERME0000001315LOT1100469BOTTL	PERMETHRIN 5% LOTION 60 mL BOTTLE
PETHI0000200851SOL3200361AMPUL	PETHIDINE (as HYDROCHLORIDE) 50 mg/mL SOLUTION 2 mL AMPULE
PETHI0000200851SOL3200379VIALX	PETHIDINE (as HYDROCHLORIDE) 50 mg/mL SOLUTION 30 mL VIAL
PETRO0000000000JELLY01063JAR01	PETROLATUM/PETROLEUM JELLY 25 g JAR
PETRO0000000000JELLY01064JAR01	PETROLATUM/PETROLEUM JELLY 100 g JAR
PETRO0000000000JELLY01065JAR01	PETROLATUM/PETROLEUM JELLY 200 g JAR
PHEN10004900015CAPSU0000000000	PHENOXYMETHYL PENICILLIN (PENICILLIN V) (as POTASSIUM SALT) 250 mg CAPSULE
PHEN10004900015TAB490000000000	PHENOXYMETHYL PENICILLIN (PENICILLIN V) (as POTASSIUM SALT) 250 mg TABLET
PHEN10004900047CAPSU0000000000	PHENOXYMETHYL PENICILLIN (PENICILLIN V) (as POTASSIUM SALT) 500 mg CAPSULE
PHEN10004900047TAB490000000000	PHENOXYMETHYL PENICILLIN (PENICILLIN V) (as POTASSIUM SALT) 500 mg TABLET
PHEN10004900209SUS1400379BOTTL	PHENOXYMETHYL PENICILLIN (PENICILLIN V) (as POTASSIUM SALT) 125 mg/5 mL SUSPENSION 30 mL BOTTLE
PHEN10004900209SUS1400469BOTTL	PHENOXYMETHYL PENICILLIN (PENICILLIN V) (as POTASSIUM SALT) 125 mg/5 mL SUSPENSION 60 mL BOTTLE
PHEN10004901265SUS1400469BOTTL	PHENOXYMETHYL PENICILLIN (PENICILLIN V) (as POTASSIUM SALT) 250 mg/5 mL SUSPENSION 60 mL BOTTLE
PHENO0000000031TAB490000000000	PHENOBARBITAL 30 mg TABLET
PHENO0000000062TAB490000000000	PHENOBARBITAL 60 mg TABLET
PHENO0000000200SOL3200516AMPUL	PHENOBARBITAL 120 mg/mL (130 mg/mL) SOLUTION 1 mL AMPULE
PHENO0000000228TAB490000000000	PHENOBARBITAL 15 mg TABLET
PHENO0000000503TAB490000000000	PHENOBARBITAL 90 mg TABLET
PHENY0000000279EDS0200060BOTTL	PHENYLEPHRINE 2.50% EYE DROPS SOLUTION 5 mL BOTTLE
PHEY40000000209SUS1400201BOTTL	PHENYTOIN 125 mg/5 mL SUSPENSION 120 mL BOTTLE
PHEY40000000851SOL3200060AMPUL	PHENYTOIN 50 mg/mL SOLUTION 5 mL AMPULE
PHEY40000000851SOL3200361AMPUL	PHENYTOIN 50 mg/mL SOLUTION 2 mL AMPULE
PHEY40005100031CAPSU0000000000	PHENYTOIN (as SODIUM) 30 mg CAPSULE
PHEY40005100139CAPSU0000000000	PHENYTOIN (as SODIUM) 100 mg CAPSULE
PHEY40005100377SUS1400201BOTTL	PHENYTOIN (as SODIUM) 30 mg/5 mL SUSPENSION 120 mL BOTTLE
PHYSO0000300261SOL3200361AMPUL	PHYSOSTIGMINE (as SALICYLATE) 1 mg/mL SOLUTION 2 mL AMPULE
PHYTO0000000194SOL3200516AMPUL	PHYTOMENADIONE 10 mg/mL SOLUTION 1 mL AMPULE
PHYTO0000001066SOL3200000AMPUL	PHYTOMENADIONE 2 mg/0.2 mL SOLUTION AMPULE
PILOC0000000658EDS0200195BOTTL	PILOCARPINE 4% EYE DROPS SOLUTION 10 mL BOTTLE
PILOC0000000658EDS0200231BOTTL	PILOCARPINE 4% EYE DROPS SOLUTION 15 mL BOTTLE
PILOC0000001314EDS0200195BOTTL	PILOCARPINE 2% EYE DROPS SOLUTION 10 mL BOTTLE
PILOC0000001314EDS0200231BOTTL	PILOCARPINE 2% EYE DROPS SOLUTION 15 mL BOTTLE
PILOC0000201062EDS0200195BOTTL	PILOCARPINE (as HYDROCHLORIDE) 1% EYE DROPS SOLUTION 10 mL BOTTLE
PILOC0000201062EDS0200231BOTTL	PILOCARPINE (as HYDROCHLORIDE) 1% EYE DROPS SOLUTION 15 mL BOTTLE
PIPER0005100346SOL3200000VIALX	PIPERACILLIN + TAZOBACTAM (as SODIUM) 2 g + 250 mg SOLUTION VIAL

PIPER0005100411SOL3200000VIALX	PIPERACILLIN + TAZOBACTAM (as SODIUM) 4 g + 500 mg SOLUTION VIAL
PIRIB0000000050TAB520000000000	PIRIBEDIL 50 mg TABLET MODIFIED RELEASE
PNE010000001070SUS1400000PFS01	PNEUMOCOCCAL CONJUGATE VACCINE 7-valent SUSPENSION PRE-FILLED SYRINGE
PNE130000000094SOL3200000VIALX	PNEUMOCOCCAL 13 VALENT CONJUGATE VACCINE 0.5 mL SOLUTION VIAL
PNEUM0000000342SOL3200094PFS01	PNEUMOCOCCAL POLYVALENT VACCINE 25 mcg/0.5 mL SOLUTION 0.5 mL PRE-FILLED SYRINGE
PNEUM0000000342SOL3200094VIALX	PNEUMOCOCCAL POLYVALENT VACCINE 25 mcg/0.5 mL SOLUTION 0.5 mL VIAL
PNEUM0000000342SOL3201434VIALX	PNEUMOCOCCAL POLYVALENT VACCINE 25 mcg/0.5 mL SOLUTION 2.5 ml (multidose) VIAL
POLY40000001202OINTX01393TUB01	POLYMIXIN B SULFATE + BACITRACIN ZINC + NEOMYCIN SULFATE 5000 units + 400 units + 5 mg/g OINTMENT
POLYG0000000366SOL4000044BOTTL	POLYGELINE 3.50% COLLOIDAL SOLUTION 500 mL BOTTLE
POTA10000000545TAB490000000000	POTASSIUM CHLORIDE 750 mg durules (equiv. to approximately 10 mEq) TABLET
POTA10000001075SYRUP00379BOTTL	POTASSIUM CHLORIDE 1 mmol/mL SYRUP 30 mL BOTTLE
POTA10000001075SYRUP00469BOTTL	POTASSIUM CHLORIDE 1 mmol/mL SYRUP 60 mL BOTTLE
POTA10000001077SOL3200011VIALX	POTASSIUM CHLORIDE 2 mEq/mL SOLUTION 20 mL VIAL
POTA10000001077SOL3200060AMPUL	POTASSIUM CHLORIDE 2 mEq/mL SOLUTION 5 mL AMPULE
POTA10000001077SOL3200195VIALX	POTASSIUM CHLORIDE 2 mEq/mL SOLUTION 10 mL VIAL
POTA10000001077SOL3200361AMPUL	POTASSIUM CHLORIDE 2 mEq/mL SOLUTION 2 mL AMPULE
POTSS0000000180TAB490000000000	POTASSIUM CITRATE 10 mEq TABLET
POVID0000000133OINTX00227TUB01	POVIDONE IODINE 10% OINTMENT 15 g TUBE
POVID0000000133OINTX00367TUB01	POVIDONE IODINE 10% OINTMENT 30 g TUBE
POVID0000000133OINTX01245TUB01	POVIDONE IODINE 10% OINTMENT 5 g TUBE
POVID0000000133PAI0100195BOTTL	POVIDONE IODINE 10% PAINT 10 mL BOTTLE
POVID0000000133SOL3200201BOTTL	POVIDONE IODINE 10% SOLUTION 120 mL BOTTLE
POVID0000000133SOL3200231BOTTL	POVIDONE IODINE 10% SOLUTION 15 mL BOTTLE
POVID0000000133SOL3200247BOTTL	POVIDONE IODINE 10% SOLUTION 1 L BOTTLE
POVID0000000133SOL3200379BOTTL	POVIDONE IODINE 10% SOLUTION 30 mL BOTTLE
POVID0000000133SOL3200469BOTTL	POVIDONE IODINE 10% SOLUTION 60 mL BOTTLE
POVID0000000133SOL3201435BOTTL	POVIDONE IODINE 10% SOLUTION 1 Gallon BOTTLE
POVID0000000322SOL3200000BOTTL	POVIDONE IODINE 240 mL SOLUTION BOTTLE
POVID0000001062OA00100201BOTTL	POVIDONE IODINE 1% ORAL ANTISEPTIC 120 mL BOTTLE
POVID0000001062OA00100322BOTTL	POVIDONE IODINE 1% ORAL ANTISEPTIC 240 mL BOTTLE
POVID0000001062OA00100469BOTTL	POVIDONE IODINE 1% ORAL ANTISEPTIC 60 mL BOTTLE
POVID0000001083SSC0100201BOTTL	POVIDONE IODINE 7.50% SURGICAL SKIN CLEANSER 120 mL BOTTLE
POVID0000001083SSC0100247BOTTL	POVIDONE IODINE 7.50% SURGICAL SKIN CLEANSER 1 L BOTTLE
POVID0000001083SSC0100469BOTTL	POVIDONE IODINE 7.50% SURGICAL SKIN CLEANSER 60 mL BOTTLE
POVID0000001083SSC0101435BOTTL	POVIDONE IODINE 7.50% SURGICAL SKIN CLEANSER 1 Gallon BOTTLE
POVID0000001083SSC0101436BOTTL	POVIDONE IODINE 7.50% SURGICAL SKIN CLEANSER 480 mL BOTTLE
POVID0000001314EDS0200000BOTTL	POVIDONE IODINE 2% EYE DROPS SOLUTION BOTTLE
POVID0000001314SOL3200000BOTTL	POVIDONE IODINE 2% SOLUTION BOTTLE
POVID0000001315EDS0200000BOTTL	POVIDONE IODINE 5% EYE DROPS SOLUTION BOTTLE
PRALI0000000851SOL3200011VIALX	PRALIDOXIME CHLORIDE 50 mg/mL SOLUTION 20 mL VIAL
PRAZI0000000061TAB490000000000	PRAZQUANTEL 600 mg TABLET
PRED20000000007TAB490000000000	PREDNISONE 20 mg TABLET
PRED20000000031TAB240000000000	PREDNISONE 30 mg TABLET FILM COATED
PRED20000000056TAB490000000000	PREDNISONE 5 mg TABLET

PRED20000000069TAB490000000000	PREDNISONE 10 mg TABLET
PRED20000000184SUS1400379BOTTL	PREDNISONE 10 mg/5 mL SUSPENSION 30 mL BOTTLE
PRED20000000184SUS1400469BOTTL	PREDNISONE 10 mg/5 mL SUSPENSION 60 mL BOTTLE
PREDN0000000056TAB490000000000	PREDNISOLONE 5 mg TABLET
PREDN0005200007TAB490000000000	PREDNISOLONE (as SODIUM PHOSPHATE) 20 mg TABLET
PREDN0005200314SYRUP000469BOTTL	PREDNISOLONE (as SODIUM PHOSPHATE) 20 mg/5 mL SYRUP 60 mL BOTTLE
PREDN0005200706SYRUP00011BOTTL	PREDNISOLONE (as SODIUM PHOSPHATE) 15 mg/5 mL SYRUP 20 mL BOTTLE
PREDN0005200706SYRUP00379BOTTL	PREDNISOLONE (as SODIUM PHOSPHATE) 15 mg/5 mL SYRUP 30 mL BOTTLE
PREDN0005200706SYRUP000469BOTTL	PREDNISOLONE (as SODIUM PHOSPHATE) 15 mg/5 mL SYRUP 60 mL BOTTLE
PREN10000000649EYD0100060BOTTL	PREDNISOLONE 0.50% EYE DROPS 5 mL BOTTLE
PREN10000001062EYD0100060BOTTL	PREDNISOLONE 1% EYE DROPS 5 mL BOTTLE
PRIMA0000001090TAB490000000000	PRIMAQUINE 26.3 mg TABLET
PROCA0000200050CAPSU0000000000	PROCARBAZINE (as HYDROCHLORIDE) 50 mg CAPSULE
PROPA0000200649SOL3200231BOTTL	PROPARACAINE (as HYDROCHLORIDE) 0.50% SOLUTION 15 mL BOTTLE
PROPO0000000194SOL3200011AMPUL	PROPOFOL 10 mg/mL SOLUTION 20 mL AMPULE
PROPO0000000194SOL3200011VIALX	PROPOFOL 10 mg/mL SOLUTION 20 mL VIAL
PROPO0000000194SOL3200055PFS01	PROPOFOL 10 mg/mL SOLUTION 50 mL PRE-FILLED SYRINGE
PROPO0000000194SOL3200055VIALX	PROPOFOL 10 mg/mL SOLUTION 50 mL VIAL
PROPO00000001199SOL3200055VIALX	PROPOFOL 20 mg/mL SOLUTION 50 mL VIAL
PROPR0000200036CAP420000000000	PROPRANOLOL (as HYDROCHLORIDE) 40 mg CAPSULE MODIFIED RELEASE
PROPR0000200036TAB490000000000	PROPRANOLOL (as HYDROCHLORIDE) 40 mg TABLET
PROPR0000200066CAP420000000000	PROPRANOLOL (as HYDROCHLORIDE) 80 mg CAPSULE MODIFIED RELEASE
PROPR0000200066TAB490000000000	PROPRANOLOL (as HYDROCHLORIDE) 80 mg TABLET
PROPR0000200069TAB490000000000	PROPRANOLOL (as HYDROCHLORIDE) 10 mg TABLET
PROPY0000000050TAB490000000000	PROPYLTHIOURACIL 50 mg TABLET
PROTA0000000194SOL3200060AMPUL	PROTAMINE SULFATE 10 mg/mL SOLUTION 5 mL AMPULE
PROTA0000000194SOL3201375AMPUL	PROTAMINE SULFATE 10 mg/mL SOLUTION 25 mL AMPULE
PROX10000000649EDS0200060BOTTL	PROXYMETACAINE (PROPARACAINE) 0.50% EYE DROPS SOLUTION 5 mL BOTTLE
PROX10000000649EDS0200231BOTTL	PROXYMETACAINE (PROPARACAINE) 0.50% EYE DROPS SOLUTION 15 mL BOTTLE
PYRAZ0000000047TAB490000000000	PYRAZINAMIDE 500 mg TABLET
PYRAZ0000000645SUS1400000BOTTL	PYRAZINAMIDE 500 mg/5 mL SUSPENSION BOTTLE
PYRAZ00000001265SUS1400201BOTTL	PYRAZINAMIDE 250 mg/5 mL SUSPENSION 120 mL BOTTLE
PYRAZ00000001265SUS1400469BOTTL	PYRAZINAMIDE 250 mg/5 mL SUSPENSION 60 mL BOTTLE
PYRI10000000021TAB490000000000	PYRIDOXINE (VIT. B6) 25 mg TABLET
PYRI10000200050TAB490000000000	PYRIDOXINE (VIT. B6) (as HYDROCHLORIDE) 50 mg TABLET
PYRI10000200142SOL3200195AMPUL	PYRIDOXINE (VIT. B6) (as HYDROCHLORIDE) 100 mg/ml SOLUTION 10 mL AMPULE
PYRID0001100062TAB490000000000	PYRIDOSTIGMINE (as BROMIDE) 60 mg TABLET
PYRIM0000000021TAB490000000000	PYRIMETHAMINE 25 mg TABLET
QUADV0000000094SUS1400000GPS01	HUMAN PAPILLOMAVIRUS QUADRIVALENT RECOMBINANT VACCINE (TYPE 6, 11, 16 AND 18) 0.5 mL SUSPENSION GLA
QUE010002300004TAB490000000000	QUETIAPINE (as FUMARATE) 200 mg TABLET
QUE010002300021TAB490000000000	QUETIAPINE (as FUMARATE) 25 mg TABLET
QUE010002300028TAB490000000000	QUETIAPINE (as FUMARATE) 300 mg TABLET
QUE010002300139TAB490000000000	QUETIAPINE (as FUMARATE) 100 mg TABLET

QUIN10001800030SOL3200516AMPUL	QUININE (as DIHYDROCHLORIDE) 300 mg/mL SOLUTION 1 mL AMPULE
QUIN10001800972SOL3200516AMPUL	QUININE (as DIHYDROCHLORIDE) 250 mg/mL SOLUTION 1 mL AMPULE
QUIN10006000551TAB4900000000000	QUININE (as SULFATE) 325 mg TABLET
RABEP0005100007TAB4900000000000	RABEPRAZOLE (as SODIUM) 20 mg TABLET
RABEP0005100069TAB4900000000000	RABEPRAZOLE (as SODIUM) 10 mg TABLET
RABIE00000001095SOL3200060AMPUL	RABIES IMMUNOGLOBULIN (HUMAN) 150 IU/mL SOLUTION 5 mL AMPULE
RABIE00000001095SOL3200060VIALX	RABIES IMMUNOGLOBULIN (HUMAN) 150 IU/mL SOLUTION 5 mL VIAL
RABIE00000001095SOL3200195AMPUL	RABIES IMMUNOGLOBULIN (HUMAN) 150 IU/mL SOLUTION 10 mL AMPULE
RABIE00000001095SOL3200361AMPUL	RABIES IMMUNOGLOBULIN (HUMAN) 150 IU/mL SOLUTION 2 mL AMPULE
RABIE00000001095SOL3200361VIALX	RABIES IMMUNOGLOBULIN (HUMAN) 150 IU/mL SOLUTION 2 mL VIAL
RALOX0000200062TAB4900000000000	RALOXIFENE (as HYDROCHLORIDE) 60 mg TABLET
RANIT0000000065TAB4900000000000	RANITIDINE 75 mg TABLET
RANIT0000200028TAB4900000000000	RANITIDINE (as HYDROCHLORIDE) 300 mg TABLET
RANIT0000200215TAB4900000000000	RANITIDINE (as HYDROCHLORIDE) 150 mg TABLET
RANIT0000200491SYRUP00469BOTTL	RANITIDINE (as HYDROCHLORIDE) 75 mg/5 mL SYRUP 60 mL BOTTLE
RANIT0000200491SYRUP01374BOTTL	RANITIDINE (as HYDROCHLORIDE) 75 mg/5 mL SYRUP 150 mL BOTTLE
RANIT0000200985SOL3200361AMPUL	RANITIDINE (as HYDROCHLORIDE) 25 mg/mL SOLUTION 2 mL AMPULE
RANIT0000200985SOL3200361VIALX	RANITIDINE (as HYDROCHLORIDE) 25 mg/mL SOLUTION 2 mL VIAL
RECOM0000000050POW2700000VIALX	RECOMBINANT TISSUE PLASMINOGEN ACTIVATOR (RTPA) 50 mg POWDER VIAL
RETA10000000135CAP3100000000000	RETINOL (VITAMIN A) 10,000 IU CAPSULE SOFT GELATIN
RETA10000000294CAP3100000000000	RETINOL (VITAMIN A) 200,000 IU CAPSULE SOFT GELATIN
RETA10000000747SOL3200231BOTTL	RETINOL (VITAMIN A) 10,000 IU/mL SOLUTION 15 mL BOTTLE
RETA10000000747SOL3200379BOTTL	RETINOL (VITAMIN A) 10,000 IU/mL SOLUTION 30 mL BOTTLE
RETA100000001099CAP3100000000000	RETINOL (VITAMIN A) 25,000 IU CAPSULE SOFT GELATIN
RETA100000001100CAP3100000000000	RETINOL (VITAMIN A) 50,000IU CAPSULE SOFT GELATIN
RETA100000001103CAP3100000000000	RETINOL (VITAMIN A) 100,000 IU CAPSULE SOFT GELATIN
RETIN00000001100CAPSU00000000000	RETINOL 50,000IU CAPSULE
RIBOF0000000050TAB4900000000000	RIBOFLAVIN (VITAMIN B2) 50 mg TABLET
RIEH10000000489TAB4900000000000	RIFAMPICIN + ISONIAZID + ETHAMBUTOL HYDROCHLORDE 75 mg + 150 mg + 275 mg TABLET
RIFAB00000000215CAPSU00000000000	RIFABUTIN 150 mg CAPSULE
RIFAM0000000006SUS1400201BOTTL	RIFAMPICIN 200 mg/5 mL SUSPENSION 120 mL BOTTLE
RIFAM0000000006SUS1400379BOTTL	RIFAMPICIN 200 mg/5 mL SUSPENSION 30 mL BOTTLE
RIFAM0000000006SUS1400469BOTTL	RIFAMPICIN 200 mg/5 mL SUSPENSION 60 mL BOTTLE
RIFAM00000000028CAPSU00000000000	RIFAMPICIN 300 mg CAPSULE
RIFAM00000000028TAB4900000000000	RIFAMPICIN 300 mg TABLET
RIFAM00000000040CAPSU00000000000	RIFAMPICIN 450 mg CAPSULE
RIFAM00000000040TAB4900000000000	RIFAMPICIN 450 mg TABLET
RIFAM00000000061CAPSU00000000000	RIFAMPICIN 600 mg CAPSULE
RIFAM00000000061TAB4900000000000	RIFAMPICIN 600 mg TABLET
RIFAM00000000162SUS1400201BOTTL	RIFAMPICIN 100 mg/5 mL SUSPENSION 120 mL BOTTLE
RIFAM00000000162SUS1400379BOTTL	RIFAMPICIN 100 mg/5 mL SUSPENSION 30 mL BOTTLE
RIFAM00000000162SUS1400469BOTTL	RIFAMPICIN 100 mg/5 mL SUSPENSION 60 mL BOTTLE

RIFAM0000000215CAPSU0000000000	RIFAMPICIN 150 mg CAPSULE
RIFAM0000000215TAB490000000000	RIFAMPICIN 150 mg TABLET
RIH010000000761SOL3200195VIALX	REGULAR INSULIN (RECOMBINANT DNA, HUMAN) 100 IU/mL SOLUTION 10 mL VIAL
RIH010000000761SOL3200388PFS01	REGULAR INSULIN (RECOMBINANT DNA, HUMAN) 100 IU/mL SOLUTION 3 mL PRE-FILLED SYRINGE
RIH010000001339SOL3200060VIALX	REGULAR INSULIN (RECOMBINANT DNA, HUMAN) 100 IU/mL (Equiv. to 3.571 mg/mL) SOLUTION 5 mL VIAL
RISPE0000000026TAB490000000000	RISPERIDONE 2 mg TABLET
RISPE0000000033TAB490000000000	RISPERIDONE 3 mg TABLET
RISPE0000000041TAB490000000000	RISPERIDONE 4 mg TABLET
RISPE00000000252TAB490000000000	RISPERIDONE 1 mg TABLET
RISPE00000000261OS00100178BOTTL	RISPERIDONE 1 mg/mL ORAL SOLUTION 100 mL BOTTLE
RISPE00000001225POW2700000VIALX	RISPERIDONE 50 mg prolonged-release powder + 2 mL diluent POWDER VIAL
RISPE00000001371POW2400000VPS02	RISPERIDONE 25 mg + 2 mL diluent POWDER FOR SUSPENSION VIAL + PRE-FILLED SYRINGE
RISPE00000001372POW2400000VPS02	RISPERIDONE 37.5 mg + 2 mL diluent POWDER FOR SUSPENSION VIAL + PRE-FILLED SYRINGE
RITON00000000139CAPSU0000000000	RITONAVIR 100 mg CAPSULE
RITUX00000000194SOL3200055VIALX	RITUXIMAB 10 mg/mL SOLUTION 50 mL VIAL
RITUX00000000194SOL3200195VIALX	RITUXIMAB 10 mg/mL SOLUTION 10 mL VIAL
RIVAS0003100033CAPSU0000000000	RIVASTIGMINE (as HYDROGEN TARTRATE) 3 mg CAPSULE
RIVAS0003100063CAPSU0000000000	RIVASTIGMINE (as HYDROGEN TARTRATE) 6 mg CAPSULE
RIVAS0003100121CAPSU0000000000	RIVASTIGMINE (as HYDROGEN TARTRATE) 1.5 mg CAPSULE
RIVAS0003100392CAPSU0000000000	RIVASTIGMINE (as HYDROGEN TARTRATE) 4.5 mg CAPSULE
ROCUR0001100194SOL3200060AMPUL	ROCURONIUM (as BROMIDE) 10 mg/mL SOLUTION 5 mL AMPULE
ROCUR0001100194SOL3200060VIALX	ROCURONIUM (as BROMIDE) 10 mg/mL SOLUTION 5 mL VIAL
ROCUR0001100194SOL3200195AMPUL	ROCURONIUM (as BROMIDE) 10 mg/mL SOLUTION 10 mL AMPULE
ROCUR0001100194SOL3201124VIALX	ROCURONIUM (as BROMIDE) 10 mg/mL SOLUTION 2.5 mL VIAL
ROPIV0000200194SOL3200195AMPUL	ROPIVACAINE (as HYDROCHLORIDE) 10 mg/mL SOLUTION 10 mL AMPULE
ROSE100000001062SOL3200094BOTTL	ROSE BENGAL 1% SOLUTION 0.5 mL BOTTLE
ROSE100000001316STRI10000000000	ROSE BENGAL 1.3 mg STRIPS
ROSEC0000000007TAB490000000000	ROSUVASTATIN (AS CALCIUM) 20 mg TABLET
ROSEC00000000069TAB490000000000	ROSUVASTATIN (AS CALCIUM) 10 mg TABLET
ROSIG00000000041TAB490000000000	ROSIGLITAZONE 4 mg TABLET
ROSIG00000000068TAB490000000000	ROSIGLITAZONE 8 mg TABLET
SALBU00000000261RSOLU0001100000	SALBUTAMOL 1 mg/mL RESPIRATORY SOLUTION 20 mL
SALBU00000000261RSOLU000379BOTTL	SALBUTAMOL 1 mg/mL RESPIRATORY SOLUTION 30 mL BOTTLE
SALBU00000000261RSOLU01412NEB01	SALBUTAMOL 1 mg/mL RESPIRATORY SOLUTION 2.5 ml (unit dose) NEBULE
SALBU0006000026TAB490000000000	SALBUTAMOL (as SULFATE) 2 mg TABLET
SALBU0006000041TAB520000000000	SALBUTAMOL (as SULFATE) 4 mg TABLET MODIFIED RELEASE
SALBU0006000068TAB520000000000	SALBUTAMOL (as SULFATE) 8 mg TABLET MODIFIED RELEASE
SALBU0006000349SYRUP000379BOTTL	SALBUTAMOL (as SULFATE) 2 mg/5 mL SYRUP 30 mL BOTTLE
SALBU0006000350SYRUP000201BOTTL	SALBUTAMOL (as SULFATE) 2 mg/5mL SYRUP 120 mL BOTTLE
SALBU0006000350SYRUP000469BOTTL	SALBUTAMOL (as SULFATE) 2 mg/5mL SYRUP 60 mL BOTTLE
SALBU0006000357SOL3201412NEB01	SALBUTAMOL (as SULFATE) 2 mg/ml SOLUTION 2.5 ml (unit dose) NEBULE
SALBU0006000459RSOLU01410NEB01	SALBUTAMOL (as SULFATE) 5 mg/mL RESPIRATORY SOLUTION 20 ml (multidose) NEBULE
SALBU0006000459RSOLU01413NEB01	SALBUTAMOL (as SULFATE) 5 mg/mL RESPIRATORY SOLUTION 10 ml (multidose) NEBULE

SALBU0006000564DPI010000000000	SALBUTAMOL (as SULFATE) 200 mcg/dose DRY POWDER
SALBU0006000565DPI010000000000	SALBUTAMOL (as SULFATE) 400 mcg/dose DRY POWDER
SALBU0006000569BAM010000000000	SALBUTAMOL (as SULFATE) 100 mcg/dose x 200 doses BREATH ACTUATED MDI
SALBU0006000569MDIXX0000000000	SALBUTAMOL (as SULFATE) 100 mcg/dose x 200 doses METERED DOSE INHALER
SALBU0006001110MDIXX0000000000	SALBUTAMOL (as SULFATE) 100 mcg/dose x 300 doses METERED DOSE INHALER
SALBU0006001111BAM010000000000	SALBUTAMOL (as SULFATE) 100 mcg/dose x 400 doses BREATH ACTUATED MDI
SALIC0000000133SOL3200201BOTTL	SALICYLIC ACID 10% SOLUTION 120 mL BOTTLE
SALIC0000000133SOL3200231BOTTL	SALICYLIC ACID 10% SOLUTION 15 mL BOTTLE
SALIC0000000133SOL3200379BOTTL	SALICYLIC ACID 10% SOLUTION 30 mL BOTTLE
SALIC0000000289SOL3200000BOTTL	SALICYLIC ACID 20% SOLUTION BOTTLE
SALIC00000001117SOL3201414BOTTL	SALICYLIC ACID 17% SOLUTION 13.3 mL BOTTLE
SALIC00000001315SOL3200379BOTTL	SALICYLIC ACID 5% SOLUTION 30 mL BOTTLE
SALIC00000001315SOL3200469BOTTL	SALICYLIC ACID 5% SOLUTION 60 mL BOTTLE
SALM1000000014590000000000000000	SALMETEROL (AS XINAFOATE) + FLUTICASONE PROPIONATE 25 mcg + 50 mcg metered dose inhaler + dose Count
SALM1000000014600000000000000000	SALMETEROL (AS XINAFOATE) + FLUTICASONE PROPIONATE 25 mcg + 125 mcg metered dose inhaler + dose Coun
SALM1000000014610000000000000000	SALMETEROL (AS XINAFOATE) + FLUTICASONE PROPIONATE 25 mcg + 250 mcg metered dose inhaler + dose Coun
SAMBO0000000015TAB4900000000000	SAMBONG 250 mg TABLET
SAMBO0000000047TAB4900000000000	SAMBONG 500 mg TABLET
SAQUI0000000004CAPSU0000000000	SAQUINAVIR 200 mg CAPSULE
SBSC100000014500000000000POUC1	SODIUM BICARBONATE WITH SODIUM CHLORIDE 845 g POUCH
SELEG0000200056TAB490000000000	SELEGILINE (as HYDROCHLORIDE) 5 mg TABLET
SELEN0000000279LOT100178BOTTL	SELENIUM SULFIDE 2.50% LOTION 100 mL BOTTLE
SELEN00000001062SHAMX00020BOTTL	SELENIUM SULFIDE 1% SHAMPOO 250 mL BOTTLE
SELEN00000001062SHAMX000201BOTTL	SELENIUM SULFIDE 1% SHAMPOO 120 mL BOTTLE
SELEN00000001062SHAMX000379BOTTL	SELENIUM SULFIDE 1% SHAMPOO 30 mL BOTTLE
SELEN00000001062SHAMX000469BOTTL	SELENIUM SULFIDE 1% SHAMPOO 60 mL BOTTLE
SERTR0000200050TAB490000000000	SERTRALINE (as HYDROCHLORIDE) 50 mg TABLET
SEVOF00000000000RSOLU00020BOTTL	SEVOFLURANE RESPIRATORY SOLUTION 250 mL BOTTLE
SILV10000001062CREA300227TUB01	SILVER SULFADIAZINE 1% CREAM 15 g TUBE
SILV10000001062CREA300367TUB01	SILVER SULFADIAZINE 1% CREAM 30 g TUBE
SILV10000001062CREA301063TUB01	SILVER SULFADIAZINE 1% CREAM 25 g TUBE
SILV10000001062CREA301245TUB01	SILVER SULFADIAZINE 1% CREAM 5 g TUBE
SILV10000001062CREA301334JAR01	SILVER SULFADIAZINE 1% CREAM 400 g JAR
SILV10000001062CREA301380TUB01	SILVER SULFADIAZINE 1% CREAM 10 g TUBE
SILV10000001062CREA301381JAR01	SILVER SULFADIAZINE 1% CREAM 500 g JAR
SILV10000001062CREA301384TUB01	SILVER SULFADIAZINE 1% CREAM 50 g TUBE
SILV10000001062CREA301398JAR01	SILVER SULFADIAZINE 1% CREAM 450 g JAR
SILV10000001062CREA301415TUB01	SILVER SULFADIAZINE 1% CREAM 20 g TUBE
SIMVA0000000007TAB490000000000	SIMVASTATIN 20 mg TABLET
SIMVA0000000036TAB490000000000	SIMVASTATIN 40 mg TABLET
SIMVA0000000069TAB490000000000	SIMVASTATIN 10 mg TABLET
SIROL0000000026TAB490000000000	SIROLIMUS 2 mg TABLET
SIROL00000000252TAB490000000000	SIROLIMUS 1 mg TABLET
SIROL00000000261SOL3200469BOTTL	SIROLIMUS 1 mg/mL SOLUTION 60 mL BOTTLE

SODo10000001457EYDo10000000000	SODIUM DIHYDROAZAPENTACENE POLYSULFONATE 150 mcg/mL x 15 mL EYE DROPS
SODCC0000000649SOL3200231BOTTL	SODIUM CARBOXYMETHYLCELLULOSE 0.50% SOLUTION 15 mL BOTTLE
SODCC0000000649SOL3201386BOTTL	SODIUM CARBOXYMETHYLCELLULOSE 0.50% SOLUTION 4 mL BOTTLE
SODCC0005100649EYDo100231BOTTL	SODIUM CARBOXYMETHYLCELLULOSE (as SODIUM) 0.50% EYE DROPS 15 mL BOTTLE
SODHY0000000194OPS0101416BOTTL	SODIUM HYALURONATE 10 mg/mL OPTHALMIC SOLUTION 0.85 mL BOTTLE
SODI10000001123SOL3200060AMPUL	SODIUM CALCIUM EDETATE 200 mg/mL SOLUTION 5 mL AMPULE
SODI20000000283SOL3200011VIALX	SODIUM CHLORIDE 2.5 mEq/mL SOLUTION 20 mL VIAL
SODI20000000283SOL3200055VIALX	SODIUM CHLORIDE 2.5 mEq/mL SOLUTION 50 mL VIAL
SODI40000000974TAB490000000000	SODIUM DICHLOROCYANURATE 12.5 mg TABLET
SODI60000000649SOL3200000BOTTL	SODIUM HYPOCHLORITE 0.50% SOLUTION BOTTLE
SODI60000001333SOL3200000BOTTL	SODIUM HYPOCHLORITE 1.25% SOLUTION BOTTLE
SODI70000001131CAPSU0000000000	SODIUM IODIDE 1311 0.8 to 100 mCi CAPSULE
SODI70000001132CAPSU0000000000	SODIUM IODIDE 1311 1.0 to 250 mCi CAPSULE
SODI70000001133SOL3200000VIALX	SODIUM IODIDE 1311 3.5 to 150 mCi SOLUTION VIAL
SODI80000000785SOL3200195AMPUL	SODIUM NITRITE 30 mg/mL SOLUTION 10 mL AMPULE
SODI80000000785SOL3200195VIALX	SODIUM NITRITE 30 mg/mL SOLUTION 10 mL VIAL
SODI90000000050POW2700000AMPUL	SODIUM NITROPRUSSIDE 50 mg POWDER AMPULE
SODIC0000000472TAB490000000000	SODIUM DICHLOROISOCYANURATE 67 mg TABLET
SODIC00000001126TAB490000000000	SODIUM DICHLOROISOCYANURATE 3.5 mg TABLET
SODIC00000001127TAB490000000000	SODIUM DICHLOROISOCYANURATE 8.68 mg TABLET
SODIC00000001129TAB490000000000	SODIUM DICHLOROISOCYANURATE 17 mg TABLET
SODIT0000000279SOL3200000BOTTL	SODIUM THIOSULFATE 2.50% SOLUTION BOTTLE
SODIT00000001184SOL3200000AMPUL	SODIUM THIOSULFATE 250 mg/50 mL SOLUTION AMPULE
SODIT00000001315SOL3200000BOTTL	SODIUM THIOSULFATE 5% SOLUTION BOTTLE
SODIU0000000061TAB490000000000	SODIUM BICARBONATE 600 mg TABLET
SODIU0000000251SOL3200055AMPUL	SODIUM BICARBONATE 1 meq/mL SOLUTION 50 mL AMPULE
SODIU0000000251SOL3200055VIALX	SODIUM BICARBONATE 1 meq/mL SOLUTION 50 mL VIAL
SODIU0000000251SOL3200178AMPUL	SODIUM BICARBONATE 1 meq/mL SOLUTION 100 mL AMPULE
SODIU0000000251SOL3200178VIALX	SODIUM BICARBONATE 1 meq/mL SOLUTION 100 mL VIAL
SODIU0000000251SOL3200195VIALX	SODIUM BICARBONATE 1 meq/mL SOLUTION 10 mL VIAL
SODIU0000000551TAB490000000000	SODIUM BICARBONATE 325 mg TABLET
SODIU00000001337TAB490000000000	SODIUM BICARBONATE 650 mg TABLET
SODIU000000014480000000000POUC1	SODIUM BICARBONATE 840 g POUCH
SODIU00000001449POW2700000CAR10	SODIUM BICARBONATE 440 g POWDER CARTRIDGE
SOMAT0000000013POW2700000AMPUL	SOMATOSTATIN 250 mcg POWDER AMPULE
SOMAT0000000013POW2700000VIALX	SOMATOSTATIN 250 mcg POWDER VIAL
SOMAT0000000033POW2700000AMPUL	SOMATOSTATIN 3 mg POWDER AMPULE
SOMAT0000000033POW2700000VIALX	SOMATOSTATIN 3 mg POWDER VIAL
SOTA10000000233TAB490000000000	SOTALOL 160 mg TABLET
SPECT0000000025POW2700000VIALX	SPECTINOMYCIN 2 g POWDER VIAL
SPIRO0000000021TAB490000000000	SPIRONOLACTONE 25 mg TABLET
SPIRO0000000050TAB490000000000	SPIRONOLACTONE 50 mg TABLET
SPIRO00000000139TAB490000000000	SPIRONOLACTONE 100 mg TABLET
STAND0000000237TAB490000000000	STANDARD SENNA CONCENTRATE 187 mg TABLET
STAND0000000386TAB490000000000	STANDARD SENNA CONCENTRATE 374 mg TABLET

STAND00000001135GRA1200367BOTTL	STANDARD SENNA CONCENTRATE 337 mcg/3 g GRANULE 30 g BOTTLE
STAVU00000000007CAPSU00000000000	STAVUDINE 20 mg CAPSULE
STAVU000000000031CAPSU00000000000	STAVUDINE 30 mg CAPSULE
STAVU000000000036CAPSU00000000000	STAVUDINE 40 mg CAPSULE
STER100000000000SOL3200011AMPUL	STERILE WATER FOR INJECTION SOLUTION 20 mL AMPULE
STER100000000000SOL3200044BAGXX	STERILE WATER FOR INJECTION SOLUTION 500 mL BAG
STER100000000000SOL3200044BOTTL	STERILE WATER FOR INJECTION SOLUTION 500 mL BOTTLE
STER100000000000SOL3200044GB001	STERILE WATER FOR INJECTION SOLUTION 500 mL GLASS BOTTLE
STER100000000000SOL3200055BAGXX	STERILE WATER FOR INJECTION SOLUTION 50 mL BAG
STER100000000000SOL3200055BOTTL	STERILE WATER FOR INJECTION SOLUTION 50 mL BOTTLE
STER100000000000SOL3200060AMPUL	STERILE WATER FOR INJECTION SOLUTION 5 mL AMPULE
STER100000000000SOL3200178BAGXX	STERILE WATER FOR INJECTION SOLUTION 100 mL BAG
STER100000000000SOL3200178BOTTL	STERILE WATER FOR INJECTION SOLUTION 100 mL BOTTLE
STER100000000000SOL3200195AMPUL	STERILE WATER FOR INJECTION SOLUTION 10 mL AMPULE
STER100000000000SOL3200247BAGXX	STERILE WATER FOR INJECTION SOLUTION 1 L BAG
STER100000000000SOL3200247BOTTL	STERILE WATER FOR INJECTION SOLUTION 1 L BOTTLE
STER100000000000SOL3200247GB001	STERILE WATER FOR INJECTION SOLUTION 1 L GLASS BOTTLE
STER100000000000SOL3200361AMPUL	STERILE WATER FOR INJECTION SOLUTION 2 mL AMPULE
STRE10006000241POW2700000VIALX	STREPTOMYCIN (as SULFATE) 1 g POWDER VIAL
STREP00000000117POW2700000VIALX	STREPTOKINASE 1,500,000 IU POWDER VIAL
STREP00000001136POW2700000VIALX	STREPTOKINASE 750,000 IU POWDER VIAL
SUCCI00000000004CAPSU00000000000	SUCCIMER 200 mg CAPSULE
SUCCI00000000139CAPSU00000000000	SUCCIMER 100 mg CAPSULE
SUCRA00000000047TAB4900000000000	SUCRALFATE 500 mg TABLET
SUCRA00000000241GELXX00060SAC01	SUCRALFATE 1 g GEL 5 mL SACHET
SUCRA00000000241TAB4900000000000	SUCRALFATE 1 g TABLET
SUF010005300133EDS0200060BOTTL	SULFACETAMIDE (as SODIUM SALT) 10% EYE DROPS SOLUTION 5 mL BOTTLE
SUF010005300133EDS0200195BOTTL	SULFACETAMIDE (as SODIUM SALT) 10% EYE DROPS SOLUTION 10 mL BOTTLE
SUF010005300133EDS0200231BOTTL	SULFACETAMIDE (as SODIUM SALT) 10% EYE DROPS SOLUTION 15 mL BOTTLE
SUFAP00000001140EDS0100060BOTTL	SULFACETAMIDE + PREDNISOLONE 10% + 0.25% EYE DROPS SUSPENSION 5 mL BOTTLE
SUG0100000001420000000361VIALX	SUGAMMADEX 100 mg/ml 2 mL VIAL
SULF300000001141TAB4900000000000	SULFADOXINE + PYRIMETHAMINE 500 mg + 25 mg TABLET
SULFU00000000133OINTX00367TUB01	SULFUR 10% OINTMENT 30 g TUBE
SULFU00000001315CREA300227TUB01	SULFUR 5% CREAM 15 g TUBE
SULFU00000001315CREA300367TUB01	SULFUR 5% CREAM 30 g TUBE
SULFU00000001315OINTX00227TUB01	SULFUR 5% OINTMENT 15 g TUBE
SULFU00000001315OINTX00367TUB01	SULFUR 5% OINTMENT 30 g TUBE
SUMAT00000000050TAB4900000000000	SUMATRIPTAN 50 mg TABLET
SUX010001301199SOL3200195VIALX	SUXAMETHONIUM (SUCCINYLCHOLINE) (as CHLORIDE) 20 mg/mL SOLUTION 10 mL VIAL
TACRO00000000056CAPSU00000000000	TACROLIMUS 5 mg CAPSULE
TACRO00000000252CAPSU00000000000	TACROLIMUS 1 mg CAPSULE
TACRO00000000459SOL3200516AMPUL	TACROLIMUS 5 mg/mL SOLUTION 1 mL AMPULE
TAMOX0000100007TAB4900000000000	TAMOXIFEN (as ACETATE) 20 mg TABLET
TAMOX0000100031TAB4900000000000	TAMOXIFEN (as ACETATE) 30 mg TABLET
TAMOX0000100036TAB4900000000000	TAMOXIFEN (as ACETATE) 40 mg TABLET

TAMOX0000100069TAB490000000000	TAMOXIFEN (as ACETATE) 10 mg TABLET
TAMSU0000200508CAPSU0000000000	TAMSULOSIN (as HYDROCHLORIDE) 200 mcg CAPSULE
TAMSU0000200508TAB490000000000	TAMSULOSIN (as HYDROCHLORIDE) 200 mcg TABLET
TAMSU0000200808TAB490000000000	TAMSULOSIN (as HYDROCHLORIDE) 400 mcg TABLET
TDF010000000028TAB490000000000	TENOFOVIR DISOPROXIL FUMARATE 300 mg TABLET
TEGAF0000001143CAPSU0000000000	TEGAFUR + URACIL 100 mg + 224 mg CAPSULE
TELA20000000067TAB490000000000	TELMISARTAN + HYDROCHLOROTHIAZIDE 80 mg + 12.5 mg TABLET
TELA200000000403TAB490000000000	TELMISARTAN + HYDROCHLOROTHIAZIDE 40 mg + 12.5 mg TABLET
TELAM0000000036TAB490000000000	TELMISARTAN 40 mg TABLET
TELAM0000000066TAB490000000000	TELMISARTAN 80 mg TABLET
TERAZ0000000026TAB490000000000	TERAZOSIN 2 mg TABLET
TERAZ0000000056TAB490000000000	TERAZOSIN 5 mg TABLET
TERAZ00000000252TAB490000000000	TERAZOSIN 1 mg TABLET
TERBI0006001062CREA300227TUB01	TERBINAFINE (as SULFATE) 1% CREAM 15 g TUBE
TERBI0006001062CREA301245TUB01	TERBINAFINE (as SULFATE) 1% CREAM 5 g TUBE
TERBI0006001062CREA301380TUB01	TERBINAFINE (as SULFATE) 1% CREAM 10 g TUBE
TERBI0006001062CREA301407TUB01	TERBINAFINE (as SULFATE) 1% CREAM 3 g TUBE
TERBI0006001062SOL3200379BOTTL	TERBINAFINE (as SULFATE) 1% SOLUTION 30 mL BOTTLE
TERBU00000000668SOL3200516AMPUL	TERBUTALINE 500 mcg/mL SOLUTION 1 mL AMPULE
TERBU0000001147DPI010000000000	TERBUTALINE 500 mcg/dose x 100 doses DRY POWDER
TERBU0006000003TAB490000000000	TERBUTALINE (as SULFATE) 2.5 mg TABLET
TERBU0006000026TAB490000000000	TERBUTALINE (as SULFATE) 2 mg TABLET
TERBU0006000056TAB490000000000	TERBUTALINE (as SULFATE) 5 mg TABLET
TERBU0006000980SOL3201409NEB01	TERBUTALINE (as SULFATE) 2.5 mg/ml SOLUTION 2 ml (unit dose) NEBULE
TERBU0006000980SOL3201417NEB01	TERBUTALINE (as SULFATE) 2.5 mg/ml SOLUTION 15 ml (multi dose) NEBULE
TERBU0006000980SOL3201418NEB01	TERBUTALINE (as SULFATE) 2.5 mg/ml SOLUTION 30 ml (multi dose) NEBULE
TERBU0006001145SYRUP00201BOTTL	TERBUTALINE (as SULFATE) 1.5 mg/5 mL SYRUP 120 mL BOTTLE
TERBU0006001145SYRUP00469BOTTL	TERBUTALINE (as SULFATE) 1.5 mg/5 mL SYRUP 60 mL BOTTLE
TERIZ0000000015CAPSU0000000000	TERIZODONE 250 mg CAPSULE
TESTO0000000036CAPSU0000000000	TESTOSTERONE 40 mg CAPSULE
TETA10000000060SOL3200000VIALX	TETANUS TOXOID 5 mL SOLUTION VIAL
TETA10000000094SOL3200000AMPUL	TETANUS TOXOID 0.5 mL SOLUTION AMPULE
TETA100000000195SOL3200000VIALX	TETANUS TOXOID 10 mL SOLUTION VIAL
TETC100000000261SOL3200516AMPUL	TETRACOSACTIDE (CONSYNTROPIN) 1 mg/mL SOLUTION 1 mL AMPULE
TETC100000000771SOL3200516AMPUL	TETRACOSACTIDE (CONSYNTROPIN) 250 mcg/mL SOLUTION 1 mL AMPULE
TETR20000000015CAPSU0000000000	TETRACYCLINE 250 mg CAPSULE
TETR20000000047CAPSU0000000000	TETRACYCLINE 500 mg CAPSULE
TETRA0000000007POW2700000AMPUL	TETRACAINE 20 mg POWDER AMPULE
THEO10000000004TAB490000000000	THEOPHYLLINE (ANHYDROUS) 200 mg TABLET
THEO10000000028TAB490000000000	THEOPHYLLINE (ANHYDROUS) 300 mg TABLET
THEO100000000205TAB490000000000	THEOPHYLLINE (ANHYDROUS) 125 mg TABLET
THEO100000000215TAB490000000000	THEOPHYLLINE (ANHYDROUS) 150 mg TABLET
THEO100000000548SYRUP00201BOTTL	THEOPHYLLINE (ANHYDROUS) 25 mg/5 mL SYRUP 120 mL BOTTLE
THEO100000000548SYRUP00469BOTTL	THEOPHYLLINE (ANHYDROUS) 25 mg/5 mL SYRUP 60 mL BOTTLE
THEO100000000627CAPSU0000000000	THEOPHYLLINE (ANHYDROUS) 400 mg MR CAPSULE

THEO10000000627TAB490000000000	THEOPHYLLINE (ANHYDROUS) 400 mg MR TABLET
THEO10000000715TAB490000000000	THEOPHYLLINE (ANHYDROUS) 250 mg MR TABLET
THIA10000000028TAB490000000000	THIAMINE (HYDROCHLORIDE) 300 mg TABLET
THIA10000000050TAB490000000000	THIAMINE (HYDROCHLORIDE) 50 mg TABLET
THIA100000000139TAB490000000000	THIAMINE (HYDROCHLORIDE) 100 mg TABLET
THIA100000000142SOL3200195VIALX	THIAMINE (HYDROCHLORIDE) 100 mg/ml SOLUTION 10 mL VIAL
THIA100000000142SOL3200516AMPUL	THIAMINE (HYDROCHLORIDE) 100 mg/ml SOLUTION 1 mL AMPULE
THIA100000000142SOL3200516VIALX	THIAMINE (HYDROCHLORIDE) 100 mg/ml SOLUTION 1 mL VIAL
THIA60000000069TAB490000000000	THIAMINE HYDROCHLORIDE 10 mg TABLET
THIAC00000000215TAB490000000000	THIACETAZONE 150 mg TABLET
THIAM00000000056TAB490000000000	THIAMAZOLE (METHIMAZOLE) 5 mg TABLET
THIAM00000000069TAB490000000000	THIAMAZOLE (METHIMAZOLE) 10 mg TABLET
THIOP00000000047POW2700000VIALX	THIOPENTAL SODIUM 500 mg POWDER VIAL
THIOP00000001373SOL3200000VIALX	THIOPENTAL SODIUM 1 g + 50 mL diluent SOLUTION VIAL
TIHU100000000145SOL3200122VIALX	TETANUS IMMUNOGLOBULIN (HUMAN) 1000 IU/mL SOLUTION 1.5 mL VIAL
TIHU100000000216SOL3200516AMPUL	TETANUS IMMUNOGLOBULIN (HUMAN) 1500 IU/mL SOLUTION 1 mL AMPULE
TIHU1000000001152SOL3200361AMPUL	TETANUS IMMUNOGLOBULIN (HUMAN) 250 IU/mL SOLUTION 2 mL AMPULE
TIHU1000000001152SOL3200516AMPUL	TETANUS IMMUNOGLOBULIN (HUMAN) 250 IU/mL SOLUTION 1 mL AMPULE
TIHU1000000001152SOL3200516PFS01	TETANUS IMMUNOGLOBULIN (HUMAN) 250 IU/mL SOLUTION 1 mL PRE-FILLED SYRINGE
TIHU1000000001152SOL3201386AMPUL	TETANUS IMMUNOGLOBULIN (HUMAN) 250 IU/mL SOLUTION 4 mL AMPULE
TIHU1000000001154SOL3200361VIALX	TETANUS IMMUNOGLOBULIN (HUMAN) 250 units/mL SOLUTION 2 mL VIAL
TIHU1000000001154SOL3200516VIALX	TETANUS IMMUNOGLOBULIN (HUMAN) 250 units/mL SOLUTION 1 mL VIAL
TIMOL00000000085EDS0200060BOTTL	TIMOLOL 0.25% EYE DROPS SOLUTION 5 mL BOTTLE
TIMOL000000000649EDS0200060BOTTL	TIMOLOL 0.50% EYE DROPS SOLUTION 5 mL BOTTLE
TIN010005101163SOL3200361VIALX	TINZAPARIN (as SODIUM) 10,000 anti-XA IU/mL SOLUTION 2 mL VIAL
TIN010005101163SOL3201419PFS01	TINZAPARIN (as SODIUM) 10,000 anti-XA IU/mL SOLUTION 0.35 mL PRE-FILLED SYRINGE
TIN010005101163SOL3201420PFS01	TINZAPARIN (as SODIUM) 10,000 anti-XA IU/mL SOLUTION 0.45 ml PRE-FILLED SYRINGE
TIOCO00000000000LOTH100195BOTTL	TIOCONAZOLE LOTION 10 mL BOTTLE
TIOCO00000000000SHAMX00178BOTTL	TIOCONAZOLE SHAMPOO 100 mL BOTTLE
TIOCO00000000000SHAMX00195SAC01	TIOCONAZOLE SHAMPOO 10 mL SACHET
TIOCO00000000000SHAMX00469BOTTL	TIOCONAZOLE SHAMPOO 60 mL BOTTLE
TIOCO00000000028OVULE00000000000	TIOCONAZOLE 300 mg OVULE
TIOCO00000000047OVULE00000000000	TIOCONAZOLE 500 mg OVULE
TIOCO000000000139OVULE00000000000	TIOCONAZOLE 100 mg OVULE
TIOCO000000000194SOL3200231BOTTL	TIOCONAZOLE 10 mg/mL SOLUTION 15 mL BOTTLE
TIOCO000000000194SOL3201375BOTTL	TIOCONAZOLE 10 mg/mL SOLUTION 25 mL BOTTLE
TIOCO000000000842SHAMX00842SAC01	TIOCONAZOLE 6 mL SHAMPOO 6 mL SACHET
TIOCO000000001062CREA300227TUB01	TIOCONAZOLE 1% CREAM 15 g TUBE
TIOCO000000001062CREA301393TUB01	TIOCONAZOLE 1% CREAM 3.5 g TUBE
TIOCO000000001062CREA301398TUB01	TIOCONAZOLE 1% CREAM 450 g TUBE
TIOCO000000001314CREA300227TUB01	TIOCONAZOLE 2% CREAM 15 g TUBE
TIOCO000000001314CREA301393TUB01	TIOCONAZOLE 2% CREAM 3.5 g TUBE
TIOCO000000001314CREA301398TUB01	TIOCONAZOLE 2% CREAM 450 g TUBE

TIOTR0001100238DPI0100000INH01	TIOTROPIUM (as BROMIDE) 18 mcg/dose DRY POWDER INHALER
TOBR100000000088EDS0100060BOTTL	TOBRAMYCIN + DEXAMETHASONE 0.3% + 0.1% EYE DROPS SUSPENSION 5 mL BOTTLE
TOBR100000000088EYO0101393TUB01	TOBRAMYCIN + DEXAMETHASONE 0.3% + 0.1% EYE OINTMENT 3.5 g TUBE
TOBRA00000000089EDS0200060BOTTL	TOBRAMYCIN 0.30% EYE DROPS SOLUTION 5 mL BOTTLE
TOBRA00000000089EYO0101393TUB01	TOBRAMYCIN 0.30% EYE OINTMENT 3.5 g TUBE
TOLNA000000CREA30008200227TUB01	TOLNAFTATE 0.10% CREAM 15 g TUBE
TOPIR00000000004TAB4900000000000	TOPIRAMATE 200 mg TABLET
TOPIR00000000021CAPSU00000000000	TOPIRAMATE 25 mg CAPSULE
TOPIR00000000021TAB4900000000000	TOPIRAMATE 25 mg TABLET
TOPIR00000000050TAB4900000000000	TOPIRAMATE 50 mg TABLET
TOPIR00000000139TAB4900000000000	TOPIRAMATE 100 mg TABLET
TOPIR00000000228CAPSU00000000000	TOPIRAMATE 15 mg CAPSULE
TRA0100000000215LYP0100000VIALX	TRASTUZUMAB 150 mg LYOPHILIZED POWDER VIAL
TRAMA00000000004TAB5200000000000	TRAMADOL 200 mg TABLET MODIFIED RELEASE
TRAMA00000000050CAPSU00000000000	TRAMADOL 50 mg CAPSULE
TRAMA00000000050DISPE00000000000	TRAMADOL 50 mg DISPERSIBLE TABLET
TRAMA00000000139TAB5200000000000	TRAMADOL 100 mg TABLET MODIFIED RELEASE
TRAMA00000000215TAB5200000000000	TRAMADOL 150 mg TABLET MODIFIED RELEASE
TRAMA00000000851SOL3200361AMPUL	TRAMADOL 50 mg/mL SOLUTION 2 mL AMPULE
TRAMA00000000851SOL3200516AMPUL	TRAMADOL 50 mg/mL SOLUTION 1 mL AMPULE
TRANE00000000015CAPSU00000000000	TRANEXAMIC ACID 250 mg CAPSULE
TRANE00000000047CAPSU00000000000	TRANEXAMIC ACID 500 mg CAPSULE
TRANE00000000047TAB4900000000000	TRANEXAMIC ACID 500 mg TABLET
TRANE00000000142SOL3200060AMPUL	TRANEXAMIC ACID 100 mg/ml SOLUTION 5 mL AMPULE
TRANE00000000142SOL3201124AMPUL	TRANEXAMIC ACID 100 mg/ml SOLUTION 2.5 mL AMPULE
TRAVO00000000073SOL3201124BOTTL	TRAVOPROST 0.004% SOLUTION 2.5 mL BOTTLE
TRIAM0000500039SOL3200516AMPUL	TRIAMCINOLONE (as ACETONIDE) 40 mg/mL SOLUTION 1 mL AMPULE
TRIAM0000500194SOL3200060VIALX	TRIAMCINOLONE (as ACETONIDE) 10 mg/mL SOLUTION 5 mL VIAL
TRIAM0000500194SOL3200516VIALX	TRIAMCINOLONE (as ACETONIDE) 10 mg/mL SOLUTION 1 mL VIAL
TRIME00000000007TAB4900000000000	TRIMETAZIDINE 20 mg TABLET
TRIME00000000382TAB4900000000000	TRIMETAZIDINE 35 mg TABLET
TROEY00000000649EDS0200060BOTTL	TROPICAMIDE 0.50% EYE DROPS SOLUTION 5 mL BOTTLE
TSAAN00000000015TAB4900000000000	TSAANG GUBAT 250 mg TABLET
TUBER00000001165SOL3200060VIALX	TUBERCULIN, PURIFIED PROTEIN DERIVATIVE (PPD) 2 TU/0.1 mL SOLUTION 5 mL VIAL
TUBER00000001165SOL3200516VIALX	TUBERCULIN, PURIFIED PROTEIN DERIVATIVE (PPD) 2 TU/0.1 mL SOLUTION 1 mL VIAL
TUBER00000001166SOL3200000VA001	TUBERCULIN, PURIFIED PROTEIN DERIVATIVE (PPD) 5 TU/0.1 mL + 2 mL diluent SOLUTION VIAL + AMPULE
TYPHO00000001167TAB4900000000000	TYPHOID VACCINE live-attenuated S. typhi (not less than 109) viable strain TABLET
TYPHO0000000116800000000094PFS01	TYPHOID VACCINE Vi-capsular polysaccharide S. typhi 25 mcg 0.5 mL PRE-FILLED SYRINGE
URSOD00000000004TAB4900000000000	URSODEOXYCHOLIC ACID 200 mg TABLET
URSOD00000000015CAPSU00000000000	URSODEOXYCHOLIC ACID 250 mg CAPSULE
URSOD00000000139TAB4900000000000	URSODEOXYCHOLIC ACID 100 mg TABLET
VALAC0000200047TAB4900000000000	VALACICLOVIR (as HYDROCHLORIDE) 500 mg TABLET
VALGA00000000040TAB4900000000000	VALGANCICLOVIR 450 mg TABLET
VALP100000000015TAB4900000000000	VALPROIC ACID 250 mg TABLET

VALP10000000645SOL3200060VIALX	VALPROIC ACID 500 mg/5 mL SOLUTION 5 mL VIAL
VALP10000001265SYRUP00201BOTTL	VALPROIC ACID 250 mg/5 mL SYRUP 120 mL BOTTLE
VALPR0000000015TAB490000000000	VALPROATE DISODIUM/VALPROIC ACID 250 mg TABLET
VALS10000000067TAB490000000000	VALSARTAN + HYDROCHLOROTHIAZIDE 80 mg + 12.5 mg TABLET
VALS100000000234TAB490000000000	VALSARTAN + HYDROCHLOROTHIAZIDE 160 mg + 12.5 mg TABLET
VALSA0000000066TAB490000000000	VALSARTAN 80 mg TABLET
VALSA00000000233TAB490000000000	VALSARTAN 160 mg TABLET
VANCO00000200047POW2700000VIALX	VANCOMYCIN (as HYDROCHLORIDE) 500 mg POWDER VIAL
VANCO00000200241POW2700000VIALX	VANCOMYCIN (as HYDROCHLORIDE) 1 g POWDER VIAL
VARIC00000001169SOL3200000VIALX	VARICELLA ZOSTER IMMUNOGLOBULIN (VZIG) 125 units/1.25 mL SOLUTION VIAL
VECUR00000001171FDP0100000AMPUL	VECURONIUM 4 mg/mL + 1 mL solvent FREEZE-DRIED POWDER + DILUENT AMPULE
VERAP0000000036TAB490000000000	VERAPAMIL 40 mg TABLET
VERAP0000000066TAB490000000000	VERAPAMIL 80 mg TABLET
VERAP00000000199CAPSU00000000000	VERAPAMIL 120 mg CAPSULE
VERAP00000000236TAB520000000000	VERAPAMIL 180 mg TABLET MODIFIED RELEASE
VERAP00000000321CAPSU00000000000	VERAPAMIL 240 mg CAPSULE
VERAP00000000321TAB520000000000	VERAPAMIL 240 mg TABLET MODIFIED RELEASE
VERAP00000000507SOL3200361AMPUL	VERAPAMIL 2.5mg/ml SOLUTION 2 mL AMPULE
VERO100000001174LYP0100000VIALX	VERO CELL (PURIFIED) 2.5 IU/0.5 mL + diluent LYOPHILIZED POWDER VIAL
VERO100000001175SUS1400516VIALX	VERO CELL (PURIFIED) 2.5 IU/mL SUSPENSION 1 mL VIAL
VINBL0006000069POW2700000VIALX	VINBLASTINE (as SULFATE) 10 mg POWDER VIAL
VINBL0006000261SOL3200195VIALX	VINBLASTINE (as SULFATE) 1 mg/mL SOLUTION 10 mL VIAL
VINCR0006000261SOL3200060VIALX	VINCRISTINE (as SULFATE) 1 mg/mL SOLUTION 5 mL VIAL
VINCR0006000261SOL3200195VIALX	VINCRISTINE (as SULFATE) 1 mg/mL SOLUTION 10 mL VIAL
VINCR0006000261SOL3200361VIALX	VINCRISTINE (as SULFATE) 1 mg/mL SOLUTION 2 mL VIAL
VINCR0006000261SOL3200516VIALX	VINCRISTINE (as SULFATE) 1 mg/mL SOLUTION 1 mL VIAL
VINCR0006001176SOL3200000VIALX	VINCRISTINE (as SULFATE) 1 mg + diluent SOLUTION VIAL
VINCR0006001177SOL3200000VIALX	VINCRISTINE (as SULFATE) 2 mg + diluent SOLUTION VIAL
VIT160000000150SOL3200388AMPUL	VITAMIN B1 B6 B12 100 mg + 100 mg + 1 mg SOLUTION 3 mL AMPULE
VIT160000000153CAPSU00000000000	VITAMIN B1 B6 B12 100 mg + 5 mg + 50 mcg CAPSULE
VIT160000000153TAB490000000000	VITAMIN B1 B6 B12 100 mg + 5 mg + 50 mcg TABLET
VIT1600000001179DROPS0023100000	VITAMIN B1 B6 B12 10 mg + 5 mg + 5 mcg/0.6 mL DROPS 15 mL
VIT1600000001181SOL3200000VIALX	VITAMIN B1 B6 B12 100 mg + 100 mg + 1 mg/10 mL SOLUTION VIAL
VIT1600000001331TAB490000000000	VITAMIN B1 B6 B12 250 mg + 250 mg + 1000 mcg TABLET
VITAA00000000294CAP310000000000	VITAMIN A 200,000 IU CAPSULE SOFT GELATIN
VITAT00000000195SOL3200000AMPUL	VITAMINS IV - TRACE ELEMENTS 10 mL SOLUTION AMPULE
VORIC0000000004LPS0300379VIALX	VORICONAZOLE 200 mg LYOPHILIZED POWDER + SOLUTION 30 mL VIAL
VORIC0000000004TAB490000000000	VORICONAZOLE 200 mg TABLET
VSIN100000001170SOL3200516AMPUL	VASOPRESSIN 20 pressor units/mL SOLUTION 1 mL AMPULE
WARFA0000000003TAB490000000000	WARFARIN 2.5 mg TABLET
WARFA0000000056TAB490000000000	WARFARIN 5 mg TABLET
WARFA00000000252TAB490000000000	WARFARIN 1 mg TABLET
YELLO00000001182SOL3200000VS001	YELLOW FEVER VACCINE 1000 DL 50 mouse min (attenuated) + 0.5 mL solvent SOLUTION VIAL + SYRINGE
YERBA0000000015TAB490000000000	YERBA BUENA 250 mg TABLET
YERBA0000000047TAB490000000000	YERBA BUENA 500 mg TABLET

ZALCI00000001183TAB490000000000	ZALCITABINE 375 mcg TABLET
ZIDOV00000000028TAB490000000000	ZIDOVUDINE 300 mg TABLET
ZIDOV00000000139CAPSU0000000000	ZIDOVUDINE 100 mg CAPSULE
ZIDOV00000001228SUS1400322BOTTL	ZIDOVUDINE 50 mg/5 mL SUSPENSION 240 mL BOTTLE
ZINC200000001238DISPE0000000000	ZINC SULFATE 88 mg DISPERSIBLE TABLET
ZINCX00000000031TAB490000000000	ZINC 30 mg TABLET
ZINCX00000000069TABL900000000000	ZINC 10 mg TABLET CHEWABLE
ZINCX00000001335OD000000231BOTTL	ZINC 27.5 mg/mL (Equiv. to 10 mg Elemental Zinc) ORAL DROPS 15 mL BOTTLE
ZINCX00000001336SYRUP00469BOTTL	ZINC 55 mg/5 mL (Equiv. to 20 mg Elemental Zinc) SYRUP 60 mL BOTTLE
ZINCX00000001344SYRUP00201BOTTL	ZINC 70 mg/5 mL (Equiv. 10 mg Elemental Zinc) SYRUP 120 mL BOTTLE
ZINCX00000001344SYRUP00469BOTTL	ZINC 70 mg/5 mL (Equiv. 10 mg Elemental Zinc) SYRUP 60 mL BOTTLE
ZOLMI00000000003TAB490000000000	ZOLMITRIPTAN 2.5 mg TABLET
ZOLPI00000000069TAB490000000000	ZOLPIDEM 10 mg TABLET
NOMED00000000000000000000000000	DRUGS AND MEDICINES NOT NEEDED DURING THIS PARTICULAR EPISODE OF CARE