

Validating XML File using Netbeans

1. Copy the Konsulta XML file to a folder
2. Copy the DTD_KonsultaData_v1.14.dtd file (from the Konsulta dev kit) to the same folder
3. Launch NetBeans
4. In NetBeans, click the File>Open menu and select the Konsulta XML file
5. In the NetBeans editor, insert the following lines at the top of the displayed content of the eclaims XML

```
<?xml version="1.0"?>  
<!DOCTYPE Konsulta PUBLIC "-//PHIC-ITMD//DTD Konsulta File //EN"  
"DTD_KonsultaData_v1.14.dtd" >
```

6. Right click on any of the whitespace of displayed XML text.
7. On the context menu that appeared, click the Validate XML menu.
8. Netbeans will display the validation results on the window below the editor.